Standard Hot Cold Pack

WARNING: Read and follow all directions carefully to avoid possible injury.

Sunbeam

HOW TO USE FOR COLD THERAPY:

1. Place the gel pack in the freezer for at least 1 hour before use.
2. Once frozen, place gel pack into the mesh pocket of the wrap and secure it with the hook and loop. Then apply the wrap onto the injured area with the gel pack pointing towards the body.

Note: Always use product with the provided cover. Never apply directly to skin.

HOW TO USE FOR HEAT THERAPY:

BE EXTREMELY CAREFUL WHEN HANDLING PACK FROM MICROWAVE TO AVOID BURNS. HEAT ONLY IN A MICROWAVE. HEATING TIME MAY VARY DEPENDING ON MICROWAVE RATED OUTPUT.

Heating Instructions:

From Room Temperature:

1. Heat in microwave on HIGH in 30 second intervals until a safe and desired heat level is reached, with total heating time not to exceed 90 seconds. Massage and turn gel pack over in microwave between intervals to avoid hot spots.
2. Monitor gel pack in microwave while heating. If any swelling of the gel pack occurs, stop heating immediately and wait for the swelling to subside before removing.
3. Before removing gel pack from the microwave, carefully check for any signs of possible rupture or leakage. If any leakage is found, do not use the gel pack.

From Frozen:

1. Heat in microwave on HIGH for 60 seconds, then in intervals of 30 seconds until a safe and desired heat level is reached, with total heating time not to exceed 120 seconds. Massage and turn gel pack over in microwave between intervals to avoid hot spots.
2. Monitor gel pack in microwave while heating. If any swelling of the gel pack occurs, stop heating immediately and wait for the swelling to subside before removing.
3. Before removing gel pack from the microwave, carefully check for any signs of possible rupture or leakage. If any leakage is found, do not use the gel pack.

**PLACE GEL PACK INSIDE THE MESH POCKET OF THE WRAP PRIOR TO USE.
NEVER APPLY DIRECTLY TO SKIN.**

CAUTION

1. Consult a physician prior to use.
2. People with poor circulation, diabetes, nerve damage, paralysis or sensitive skin should not use this product.
3. Although non-toxic, discard gel pack if product leaks or is punctured.
4. Adult supervision is required for use on invalids or children. Do not use on infants or children under the age of 2 years.
5. Do not use the gel pack in conjunction with liniments, sports creams, balms or lotions.
6. Do not use on broken skin.
7. If gel pack is too hot, allow to cool before use.
8. Do not fall asleep with the gel pack.

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