



Technical Publications

Direction 5506903-1EN




Rev. 3



GE Healthcare

DaTQUANT™

Quantitative Software for DaTscan™ Images



DaTQUANT™

Quantitative Software for DaTscan™ Images

User Installation & Activation Guide

© 2020 by GE. All rights reserved.



Legal Notice

All Material Copyrighted ©2020 by GE. All rights reserved.

The software and related documentation are confidential and proprietary information of GE, protected by copyright laws of the United States, other countries, and by international treaty provisions. Only licensees of GE have a right to use the information contained herein. Only licensees specifically granted copy and/or transfer rights have the right to copy and/or transfer the information. Any unauthorized use, disclosure, assignment, transfer or reproduction of this confidential information will be prosecuted to the full extent of the law.

This software and documentation is provided "as is." GE makes no representations or warranties, express or implied, including but not limited to, warranties of merchantability or fitness for any particular purpose or that the use of the software or documentation will not infringe any third party patents, copyrights, trademarks or other rights. GE will not be liable for any direct, indirect, special, or consequential damages arising out of the user of the software or documentation. GE will not be liable nor obligated in any manner in respect of bodily injury and/or property damage arising from the use of this software if such use is not in strict compliance with instructions and safety precautions contained in the relevant operating manuals and in all supplements thereto, in all product labels, and according to all terms of warranty and sale of this software, nor if any change not authorized by GE is made to the software contained herein. In addition, user provided programs or scripts are NOT validated nor warranted by GE. The use of data obtained using such user provided programs or scripts is the sole responsibility of the party using such programs or protocols.

Trademark Notice

Xeleris™ – Xeleris is a trademark of General Electric Company.

All third party trademarks are property of their respective owners.



CAUTION

Federal law restricts this device to sale by or on the order of a physician.



Council Directive 93/42/EEC concerning medical devices:

The "CE 0459" label affixed to the product testifies compliance to the Directive.



Name and Address of Authorized Representative in EU:

GE Medical Systems SCS
283, Rue de la Minière
78530 BUC France



The original manual was written in English.

Revision History

Revision	Date	Description of Changes	Chapter/Pages
Rev. 1	January 2015	Initial Release	All
Rev. 2	June 2016	Various edits throughout. Added new sections: ■ Running in Trial Mode ■ Network Addressing Selection	All 2-14 2-21
Rev. 3	March 2020	■ Updating minimal system requirements ■ Updating DaTQUANT installation screens to SW version 2 ■ InSite Express not supported: relevant section removed and License activation updated.	2-1 Throughout Ch.2, Installation Chapter 2, Chapter 5



Topics

Chapter 1: Overview

- Safety Notice Legends 1-1
- System Documentation Set and Online Access. 1-2
- Manual Overview 1-2
- Printing this Document 1-3
- Conventions in This Document 1-4

Chapter 2: Installation

- Minimal System Requirements 2-1
- Installing DaTQUANT 2-3
- Monitor Quality Test. 2-7
- System Registration and License Activation. 2-13
- Running in Trial Mode 2-14
 - System Registration. 2-16
 - License Activation 2-18
- End User License Agreement (EULA). 2-19
- Network Addressing Selection. 2-21



Configure DICOM Server 2-22

Uninstalling DaTQUANT 2-25

Computer Replacement 2-26

Chapter 3: Managing Patient Data and Connectivity

Using DaTQUANT Data Manager 3-1

Import Patient Data 3-3

 Import a DICOM Folder 3-3

 Import a DICOM Study. 3-4

Save or Delete Patient Data. 3-6

 Saving or Deleting Patient Study 3-6

 Saving or Deleting Selected Patient Series 3-6

 Viewing Patient Data 3-7

Chapter 4: Backup and Restore

Backup. 4-1

Restore 4-2

Chapter 5: Troubleshooting

System ID From Activation Form. 5-1

DICOM Connectivity 5-2



Creating a Problem Report. 5-3



Chapter 1: Overview

Safety Notice Legends



WARNING

Potentially hazardous situation, which if not avoided, could result in death or serious injury.



CAUTION

Potentially hazardous situation, which if not avoided, may result in minor or moderate injury.

NOTE

Non-hazardous situation, which if not avoided, could result in equipment damage, lost time, or reduced image quality.

IMPORTANT

Refer to *DaTQUANT Operator Guide* for important safety and regulatory information that pertains to the DaTQUANT application.



System Documentation Set and Online Access

The DaTQUANT documentation set includes:

- *DaTQUANT User Installation & Activation Guide* (this document)
- *DaTQUANT Operator Guide*
- *eIFU and User Information booklet*

The DaTQUANT documentation is available in the DaTQUANT Help folder. Once the application is installed and activated, you can access the manuals from within the application: Select **How To...** to access the **How To** page that links to this guide, the *DaTQUANT Operator Guide*, and the *eIFU and User Information booklet*.

Manual Overview




This guide contains and describes the following chapters and topics:

- Overview (this chapter)
- [Chapter 2: Installation](#), which includes minimal requirements, installing DaTQUANT, Monitor check, registration and license activation, End User License Agreement (EULA), DICOM configuration and computer replacement.
- [Chapter 3: Managing Patient Data and Connectivity](#)
- [Chapter 4: Backup and Restore](#)
- [Chapter 5: Troubleshooting](#)



Printing this Document

This document is created using A5 sheet size. Use the following guidelines when printing:

Print Dialog	Property ^a	Two pages per sheet (A4 or Letter)	One page per sheet (A5 or A4/Letter)	
			A5 is recommended for a compact book and paper saving A4/Letter is recommended when large format and print are needed	
				
Main	Page Scaling	Fit to Printable Area (removes extra white margins)		N/A
Printer Properties (Advanced)	Orientation	Landscape		
	Double-sided	Double-sided (or Print on Both Sides)		
	Page direction / Binding	Open to Side (or Open to Left)	Open to Top (or Flip pages Up)	Open to Side (or Open to Left)
	Number of pages per sheet	2	1	1
	Page Borders	Print Page Borders (adds visual separation between the two pages)	N/A	N/A

^a Names of properties and options can differ, depending on your specific printer driver

Conventions in This Document

The following conventions are used throughout the manual:

Description	Example
Keys on the operator keyboard, hand-held controller and gantry	<SET> , <Ctrl>
Software Interface Buttons	[OK] , [Apply] , [Cancel]
Names of items in the graphical interface including: - Names of dialog boxes, windows, tabs, areas and lists - Menu items - Field and icon labels	Configuration tab; To Do List File menu Gantry icon; Properties field
System messages	Press Y to continue.
System parameters whose actual values must be defined by the user	Type-in the <i>Patient ID</i>
Hyperlinks	Figure 3-1
Paths	D:\Einstein\
References to other documents	<i>Operator Manual</i>
End of a procedure	◆



Chapter 2: Installation

Topics covered in this chapter:

- Minimal System Requirements
- Installing DaTQUANT
- Monitor Quality Test
- System Registration and License Activation
- End User License Agreement (EULA)
- Uninstalling DaTQUANT
- Computer Replacement

Minimal System Requirements

- Minimal memory – 1024 MB
- CPU speed – 1GHz
- Minimal supported screen resolution – height 768, width 1024 (a typical 17" monitor)
Recommended resolutions – 1366x768, 1280x1024, 1920x1080

IMPORTANT

Minimal resolution is a mandatory requirement, and DaTQUANT will not be launched if the system does not meet this requirement.



- Supported Windows versions:
 - Windows 10 32/64 bit
 - Windows 7 32/64 bit (SP1 mandatory)
- For DICOM connectivity purposes, the following ports need to be open:
 - For fetching 4242
 - For sending 104

For details, refer to [Configure DICOM Server, p.2-22](#).

IMPORTANT

If DaTQUANT is installed on a computer that already has an installed DaTQUANT (a previous or current version): Use the **Backup** tool to back up the DaTQUANT configuration and licenses, and all user-created normal databases and patient exams (refer to [page 4-1](#)).

Keep a copy of the DaTQUANT license activator file in a separate location (for example, a removable drive or network location). The patient exams, the user-created normal databases and all customization settings are retained from the previous DaTQUANT version.

The DaTQUANT license activator file used in version 1.0 is valid for version 2.0. No new license activator file is needed when upgrading an active license.



Installing DaTQUANT

The installation sequence is as follows: software installation, monitor quality test, registration, activation and activation run. The necessary steps are described in this chapter.

1. Run the installation package. The package will be extracted to a temporary folder.
2. In the **Welcome to DaTQUANT Setup Wizard**, click **[Next]**.




3. In the following **DaTQUANT Installation Wizard**:

- a. You can change the location of installed files. Click **[Browse]** and set the new location.
- b. Clear the **Sample Studies** checkbox if you do not want to copy them during installation.

NOTE

In case of an upgrade, the sample studies are added to the existing sample studies and do not replace them.


- c. Click **[Next]** to continue.




4. In the following **DaTQUANT Installation Wizard**, click **[Install]**.

Installation begins and the progress indicator appears.

5. The **Monitor Quality Test** utility is launched.



6. The following **DaTQUANT Installation Wizard** appears when the installation and the **Monitor Quality Test** are completed. Click **[Finish]**. A computer restart will complete the process.



Monitor Quality Test



WARNING

Monitor quality test must be performed on every monitor that will be used for clinical evaluations. A failure in one of the tests will **prevent** that monitor from being used for clinical reading. In such cases, you should replace the monitor to enable its use for clinical reading. The Monitor quality test must be performed in the following cases:

- During installation.
- When a monitor is replaced.


Make sure to perform this procedure on every monitor that will be used for running DaTQUANT.

1. The **Monitor Quality Test** utility is automatically launched after installation (refer to [page 2-5](#)).

After replacing the monitor, you must initiate this utility manually:


From the **Start** menu, select

DaTQUANT → Tools → Monitor Quality Test.




2. The **Monitor Quality Test** screen appears:


- Optional: Select language (default **English**)
 B. The selected language will be the DaTQUANT user interface language (that is, the language used in the application).
- If your configuration includes more than one monitor, you need to select the monitor to be used in the DaTQUANT application. Click on the relevant monitor **A**.
- Optional: Click **[Close]** to cancel the procedure.



3. Read **Introduction** text **C** in the **Quality Test** window.
4. Click **[Start]** **D** to begin the test.
(To exit the procedure, click **[Close]** **F**.)
5. Perform all steps **E** as they appear on screen.




A ✓ marks each step that was completed successfully. A Ø marks each step that failed.





6. When all steps are completed, the **Quality Test Summary** window appears.

Optional: If Microsoft Outlook is defined on the computer, click **[Email Summary]**  to send the summary.



The following window appears when the monitor does not meet minimal requirements.



7. Click **[Close]** to exit the utility.




System Registration and License Activation

This procedure is performed after installing DaTQUANT and performing the **Monitor Quality Test** utility.

IMPORTANT

The following instructions are mandatory in order to activate the application.



Running in Trial Mode


If the license is missing or has expired, the application can only be run in **Trial Mode**.

After installation, double-click the DaTQUANT icon on the desktop to launch the application. The **About** screen is displayed.

NOTE

Once the active license is installed, you can select **Do not show this dialog again** ③.

Click **[OK]** to launch in **Trial Mode** or click **Registration Form** ① to acquire an active license and register the software.





In Trial mode, you can process and review the sample studies provided as part of the installation, and the UI indicates that you are in **Trial Mode** .

A tooltip is displayed when you hover over the text **Trial Version** .

NOTE

If an active license is not found, only demo studies are displayed.

If an active license was previously available, all studies already saved in the local database are available for processing and review when the license is re-activated.



System Registration

1. The first time the application is run, when you double-click the **DaTQUANT** icon on the desktop, the **End-user License Agreement (EULA)** screen opens. To register the product, open the **Product Registration Form** from the EULA screen or select **Start → DaTQUANT → Tools**.
2. Optional: Select **Language** (default **English**) **E**. The selected language will be the **DaTQUANT** user interface language (that is, the language used in the application).
3. For **Site Information**, fill in all the mandatory fields **A** marked with an asterisk (*).
4. Optional: **Yes** or **No** **B** to indicate if the site address is your billing address.
5. Optional: Under **Product Information** **C**, enter the **Order Number** (from the GE invoice mail).
6. Click **[Save]** **D**.

DaTQUANT Tools

Registration Activation

Product Registration Form

Site Information

* Name:

* Institution:

Address:

City/State: Country:

Phone #: Zip code:

* E-mail:

B Should the site address be used also as your billing address? ☒ Yes ☐ No

C **Product Information**

Order Number (Optional):

Unique ID:

PSI Code:

Install Date:

Version:

D

E Language: **English**

* Indicates mandatory fields.



The following text then appears in the **Registration Form**:

1. To get the Activation Information, [click here](#). An email will automatically be created containing the registration information.
2. Once you receive the license information from GE, fill in the Activation form in the [Activation tab](#).

7. To obtain license activation information, you need to send the Site and Product information that you provided in the **Registration** form by using the [click here](#) link. One of the following options will occur:
 - If Microsoft Outlook is installed on the computer, an email message is automatically created containing the registration information that you filled in. Send the email message to standalone_sw.registration@ge.com
 - If Microsoft Outlook is **not** installed on the computer, a notepad file opens instructing you to create an email message containing the necessary registration information. Follow the instructions and send the information to standalone_sw.registration@ge.com.

The GE registration team will provide you with the license and relevant system information within one business day.

8. Once you receive a reply email with the active license attached, proceed to [License Activation](#).



License Activation

Perform the following after DaTQUANT is installed on your system and an active license file is received via email:

1. After receiving a reply from the GE registration team, double-click the **DaTQUANT** icon on the desktop and then click **Activation** in the **Registration** form.
2. Copy the information received from the GE registration team into the appropriate fields **A**.
3. Save the **License File** (.bin file) attached to the email message in a designated location on the computer. **[Browse]** to that location **B**.
4. Click **[Save]** **C**.
5. It is recommended that you back up your configuration ([Backup, p.4-1](#)).
6. The application is now activated. You can click **[Run DaTQUANT]** **D** to launch the application.

The screenshot shows the 'Product Activation Form' in the 'DaTQUANT Tools' application. The form has a 'Registration' tab and an 'Activation' tab. The 'Activation' tab is active, showing fields for 'GE Service Phone R:' (972-4-8563799) and 'GE System ID:' (test4). A yellow box labeled 'A' highlights these fields. Below these fields is a 'License File:' field with a 'Browse...' button, labeled 'B'. The 'License Status:' field shows 'OK' and 'Expiration date: 2 Feb 2022'. A 'Save' button is labeled 'C'. At the bottom, there is a 'Run DaTQUANT' button labeled 'D' and a 'Close' button. A note at the bottom left states '* Indicates mandatory fields.' and the language is set to 'English'.




End User License Agreement (EULA)

Once the license is activated and DaTQUANT is launched for the first time, the **End User License Agreement (EULA)** window appears.

The purpose of this stage is to verify that:

- The following requirements are met:
CPU speed, physical memory, minimal resolution.
- Adobe Reader is installed to enable using Help files.
- Monitor quality test was performed.
- The license is valid.

1. EULA will appear each time you launch DaTQUANT unless you select **Do not show this dialog again** **A**.



2. The Unique Device Identifier (UDI) is displayed in the same screen. This is a unique identification mark per medical device. The UDI code can be viewed from within the **DaTQUANT Patient Selector** by clicking **Applications → About**.

NOTE

The UDI in this screen is an example and may vary per DaTQUANT version.

IMPORTANT

Even if you select **Do not show this dialog again**, the system will display EULA in the following cases:

- Two months prior to license expiration date, EULA will show expiration date. If the license is **not** renewed by expiration date, the **Registration** form automatically appears when you try to run the application.
- When CPU and/or physical memory and/or monitor quality test do not meet minimal requirements. Requirements that are not met will show in red in EULA.

3. Optional: Select **Language**  (default **English**). The selected language will be the **DaTQUANT** user interface language (that is, the language used in the application).

NOTE



You can also select the **DaTQUANT** language in the **Registration** and **Activation** forms, and from within **DaTQUANT** in the **Data Manager** window.



Network Addressing Selection

You can determine the IP version used as follows:

1. Select **Applications** → **Configuration** **A**.
2. In the **Configuration** dialog, choose **IPv4** (default) or **IPv6** **B**.
3. Restart DaTQUANT to apply the configuration change.



Configure DICOM Server



Data can be transferred to any DICOM station, such as PACS or workstation. In order to use a DICOM station for importing or exporting data, the DICOM station needs to be configured.

Open port for DICOM retrieve (fetch):

To enable the retrieval of DICOM files, add a rule to the Windows firewall to allow port 4242. If an alternative firewall (not Windows) is used, manually add an exception rule for port 4242.



1. Click  (A) in the **DaTQUANT Data Manager** window on [page 3-2](#)) to open the **Connect to server** window:
2. To define a new DICOM server, click [**New connection...**] (F).
3. A dialog appears requesting DICOM server details from your IT department:
 - If you have the necessary details, proceed to step 4.
 - If you **do not** have the remote host settings, contact your IT department. The **Local Host Settings** (C): **Application Entity Title**, **Network Address** and **Port Number**, are needed for the configuration on the DICOM remote station.
4. When the **Remote Host Settings** (B) are received, fill in all fields and click [**OK**]. The new server is added to the list of servers (A).



5. Optional: In the **Connect to server** window, select a DICOM server and click **[Edit]**  to open the **DICOM Server connection details** dialog and edit the DICOM server properties.


Additional DICOM options in the Connect to server window :

- Click **[Delete]** to delete a configured DICOM server and confirm the deletion.
- Click **[Export Connections]** to save the configured DICOM servers as a file on system.
- Click **[Import Connections]** to import a previously saved DICOM servers' file.



Uninstalling DaTQUANT

1. From the **Start** menu, select
DaTQUANT → Tools → Uninstall DaTQUANT.
2. Click **[Yes]** and follow the screen instructions.





NOTE

Uninstalling DaTQUANT does **not** automatically remove the following:

- Local database
- Application logs
- User preferences

During uninstall, you are prompted to approve removal of these items from the hard disk.

Uninstall does not wipe-erase files.



Computer Replacement

NOTE

If the replaced computer is damaged. proceed to step 3.

1. On your old computer, create a backup ([Backup, p.4-1](#)).
2. Uninstall DaTQUANT from the computer you are replacing before proceeding ([Uninstalling DaTQUANT](#)).
3. Install the application on the new computer ([Installing DaTQUANT, p.2-3](#)).
4. Perform a monitor quality test ([Monitor Quality Test, p.2-7](#))
5. Restore the backup ([Restore, p.4-2](#)).
6. Request a new license and activate it ([System Registration and License Activation, p.2-13](#)).



Chapter 3: Managing Patient Data and Connectivity

This chapter includes:

- Using DaTQUANT Data Manager
- Import Patient Data
- Save or Delete Patient Data

Using DaTQUANT Data Manager

- ◆ Click the **DaTQUANT** icon on your desktop; The **DaTQUANT Data Manager** window lists the available patient studies, and includes various options for managing patient data, as described in this section.

The following sections describe the main operations that you can perform in this window.



The screenshot shows the DaTQUANT Data Manager application window. The interface includes a menu bar (File, Applications), a toolbar with icons for file operations, and a main data table. The table is divided into two sections: '1/48 examinations' and '1/4 series'. The '1/48 examinations' section contains a table with columns for Patient's Name, Study Date, Study ID, Study Description, Patient ID, Patient's Birth Date, Patient's Sex, Modalities, Accession Number, Performing Physic..., Study Time, and Referring Physicia... The '1/4 series' section contains a table with columns for Series Number, Series Date, Series Type, Series Description, Images, and Modality. A preview window on the right shows a brain scan image. The bottom of the window features a progress bar and buttons for 'Open', 'Start', and 'Customize'.

1/48 examinations

Patient's Name	Study Date	Study ID	Study Description	Patient ID	Patient's Birth Date	Patient's Sex	Modalities	Accession Number	Performing Physic...	Study Time	Referring Physicia...
Dat DaT phantom	Jul 6, 2012	NUH DaT	NUH DaT	060712		F	NM			11:21:20 AM	
DQ	Feb 17, 2010	Em+ATTMAP+Hyb...	User&Crane&Cra...	024443++1+	Jul 28, 1952	F	CT NM	6030825-1		11:49:14 AM	zzzzzz
DQ	Jul 27, 2004	Brain Tomo I123	Factory&Brain&Br...	0916573++2+		O	CT NM			1:22:59 PM	
DQ	Jul 27, 2004	Brain Tomo I123	Factory&Brain&Br...	0916573++2+	Apr 8, 1970	O	NM			1:22:59 PM	
DGJ, St.Jean X3_...	Jan 5, 2011	Crane LEHR Neuro	HMPRT	5247	Jun 5, 1955	M	CT NM	6037324-1		12:43:23 PM	HOFMAN, FrankD3
DGJ, St.Jean X3_...	Jan 6, 2011	Crane Neuro Diam	HMPRT	5247	Jun 5, 1955	M	CT NM	6037370-1		12:16:53 PM	HOFMAN, FrankD3
DATSCAN	Dec 29, 2010	DATSCAN	DAT SCAN	DATSCAN	Sep 1, 1955	F	NM			12:58:07 PM	
DATQUANT, exper...	Feb 19, 2013	Brain Tomo I123	Brain Tomo I123	D670		M	NM			12:55:34 PM	
D670 BRAIN EANM...	May 25, 2011	Tc99m DIAG DOSE	NEUROLITE	HYPOPERFUSION+...	Apr 29, 1922	F	CT NM			10:20:07 AM	
Brain Scan, Head	Feb 9, 2011	Brain Tomo I123	Brain Tomo I123	09032011		M	CT NM			3:40:23 PM	

1/4 series

Series Number	Series Date	Series Type	Series Description	Images	Modality
35	Nov 11, 2014		DaTQUANT RESUL...	3	NM
3	Jul 27, 2004		ATT MAP TOMO-1	1	NM
1	Jul 27, 2004		CT TOMO-1	34	CT
1	Jul 27, 2004		TOMO-1	1	NM

Preview Images

Ex: Brain Tomo I123
Se: 1
Im: 01
Fr: 01

Open: Start Customize

34%

Import Patient Data

NOTE

The **Import Patient Data** feature is not available in **Trial Mode**.

To import patient studies, click the relevant icon in **A** (alternatively, select an option from the **File** menu).

Import a DICOM Folder

- Click  to import files from a DICOM folder, and then select a folder in the file browser.
- Click  to import a DICOM folder from a removable media (DICOMDIR), and then select the removable media.

Import a DICOM Study

1. Click  to open the **Connect to server** window.
If required, refer to [Network Addressing Selection, p.2-21](#) for network configuration options.
2. Select the DICOM server that you want to import from **A** and click **[Connect]** **B**. The **DICOM Server Query** dialog appears
3. Fill in the patient information **C** and click **[OK]**.



4. In the **Browse repository** dialog, select the data that you want to import **D**.

D

Save or Delete Patient Data

Saving or Deleting Patient Study

To save or delete a selected patient study and its series, click the relevant icon in **1** in the **DaTQUANT Data Manager** window ([page 3-2](#)):

- Click  to send study to DICOM server. The **Connect to server** window appears ([page 3-4](#)).
- Click  to save study to system folder. Select location in the file browser.
- Click  to delete a study. Use the CTRL key for multiple selection. Confirm deletion.

Saving or Deleting Selected Patient Series


To save or delete a selected patient series within a study, click the relevant icon in **2** the **DaTQUANT Data Manager** window ([page 3-2](#)):

- Click  to send selected patient series to DICOM server. The **Connect to server** window appears ([page 3-4](#)).
- Click  to save selected patient series to system folder. Select location in the file browser.
- Click  to delete a selected patient series. Use the CTRL key for multiple selection. Confirm deletion.

Viewing Patient Data

In the **DaTQUANT Data Manager** window (page 3-2):

1. Select a patient study **C**; its series are listed **E**.
2. Search the studies according to one of the study attributes (**B**): select an attribute from the drop-down list **L** and type a value in the **Search...** field **K**.
3. A preview of the selected series raw data is shown in the Preview window **H**. Use the scroll bar to view different frames.
4. Click **Images I** (page 3-2) to view the series content, instead of the preview.
5. Click **N** to view DICOM attributes of the selected image. In the **Meta data** window **M**, you can search for a particular attribute.



Additional buttons and indicators in the DaTQUANT Data Manager window (page 3-2):

- Click **[Customize]** **N** to customize DaTQUANT parameters. Each parameter is described in the relevant customization tab.
- Click **[Start]** **M** to launch DaTQUANT.
- Free computer space indicator is displayed in **F**.
- When there are background processes, click **G** to view the processes running in the background.



Chapter 4: Backup and Restore

NOTE

The patient database is not included in the backup.

Backup

The following process should be used after installation and activation of the software on a new system.

1. From the **Start** menu, select **DaTQUANT → Backup/Restore**. The **Backup** window appears.

2. Select the folders to be backed up **A**.


NOTE

It is important to select the licenses for backup.

3. Click **[Browse]** **B** to change the backup file location.
4. Click **[Backup]** **C**.

NOTE

It is recommended to keep the backup files on an external device, such as USB or remote storage device.



Restore


IMPORTANT

- Restoring a user normal database will override any existing user normal database **B**.
- Do not restore backup files from the network.

1. From the **Start** menu, select **DaTQUANT → Backup/Restore**. The **Backup** window appears.
2. Click the **Restore** tab **A** to display the **Restore** window.
3. Browse to the folder of the required backup file and select the file from the list **C**. The meta data of the backup file is shown in **F**.
4. Click **[Restore]** **D**.

Before performing restore, the system verifies that the backup file:

- contains DaTQUANT data and not data from another application.
- and that the data is from a compatible DaTQUANT version.



NOTE

If you click **[Cancel]**  during the restore, the system deletes the partial restore that was performed.



Chapter 5: Troubleshooting

This chapter contains the following topics:

- [System ID From Activation Form](#)
- [DICOM Connectivity](#)
- [Creating a Problem Report](#)

System ID From Activation Form

You are required to provide your system ID when connecting to GE for support. The system ID can be obtained from the **Activation Form**:

1. From the **Start** menu, select **DaTQUANT → Tools → Registration Form**.
2. Select the **Activation** tab.
3. The System ID located in the **Activation Form** is the number requested by the GE engineer.



DICOM Connectivity

If DICOM query/retrieve does not work, verify that:

- The remote station AE, IP and port are configured correctly. Notice that AE is case sensitive.
- The local application (AE, IP and port) is configured on the remote system.
- If a firewall is running on the local system, an exception rule is configured for the DICOM port.
- Try to temporary disable the firewall to check if it is blocking communication.
- The network address version is defined correctly (see [Network Addressing Selection, p.2-21](#)).



Creating a Problem Report

1. From the **Start** menu, select **DaTQUANT → Tools → Problem Report**. The **Create Report** window appears.
2. Optional: Select the study **A**.
3. Enter a problem description **C**.
4. **Important:** To avoid capturing patient information, make sure that **Anonymize Data B** is selected.
5. Click **[Browse] D** and browse to the location in which the report zip file will be saved.
6. Click **[Create Report] E**. A message confirms that the report was created **F**.

