



Instruction guide

DaTQUANT™ 2.0



Before you upgrade to DaTQUANT for PC version 2.0 ...

WHAT YOU WILL NEED



- 1 Flexera portal user email address and password from Version 1, to download the installer for the new version.

Note:

If you don't remember the Flexera portal password, you can request a reset on the login page on the Flexera portal.

If you don't remember the e-mail address associated with the DaTQUANT for PC initial purchase and installation, or that user is no longer with the organization, you can contact GE Healthcare customer support at datquant.calendarBS@ge.com for assistance.

- 2 Windows® administrator rights on the computer onto which DaTQUANT for PC version 2.0 will be installed (if in doubt, check with your IT department).
- 3 Ensure that you know the technical support contact information and your assigned GE System ID for DaTQUANT for PC version 1.0 in case you need to contact the technical support team during the upgrade process.


BEFORE YOU BEGIN



Provided that DaTQUANT for PC version 2.0 is installed on the same PC where version 1.0 is already installed, the installation and upgrade process should retain the license activator file, patient database, user-created normal database and user-created normal databases during the upgrade process. However, it is still prudent to first make a backup of those items onto a separate location using the built-in tool.

To transfer customizations and settings to the new version, the provided Backup and Restore Utility must be used in version 1.0 as follows:

- 1 Login to Windows as an Administrator-Type user.
- 2 In Windows 10, use Start -> DaTQUANT -> Backup Restore



- 3 Save the backup file on a portable drive or a remote location

DOWNLOADING THE DaTQUANT FOR PC UPGRADE INSTALLER



1 Login to the Flexera Portal

<https://gehealthcare.flexnetoperations.com/flexnet/operationsportal/logon.do>

2 In the “Recent Files” section, three items should be visible; the software installer, the Instructions for Use/User Information Booklet and the Privacy and Security Manual.

Click on each item
to download

Recent Files
Description
DaTQUANT software for your PC Privacy and Security Manual
DaTQUANT Software for your PC V2 0
DaTQUANT software for your PC eIFU and User Information Booklet

The Booklet contains information on how to access the Installation Manual and the Operator Manual for DaTQUANT for PC version 2.0.


Note: Access to the Flexera portal and download of the files can be done from any PC with appropriate internet access, then the downloaded files can be transferred to the PC running DaTQUANT using a removable drive or shared storage location. The files do not have to be downloaded directly to the PC running DaTQUANT.

Note: Problems accessing the Flexera portal or downloading the files are issues related to the internet connection settings at your location. Any such problems should be reported to and resolved by the user's local IT department. Or you can try using a different internet connection (e.g. a hot spot from a mobile phone) to connect to the Flexera portal. GE Healthcare technical support cannot assist users with issues related to logging into or downloading files from the Flexera portal.

INSTALLING THE UPGRADE TO VERSION 2.0




- 1 Copy the installer file to a location on the hard drive of the PC running DaTQUANT 1.0.
- 2 Close all open applications.
- 3 Either:
 - (1) login to Windows as an Administrator-type user and double-click the installer file
 - or
 - (2) right-click on the installer file, select “Run as administrator” and provide the username and password for an Administrator-type user
- 4 In the User Account Control warning message about allowing the software to make changes, click “Yes”.
- 5 The installer will extract the files required and then show a prompt with a list of actions. Select “Install” and click “Next”.



Allow the installation of files to complete.

The installer will place a shortcut and icon for DaTQUANT for PC on the desktop.

6 Restart the computer when prompted to at the end of installation



7 Once the PC has completed restarting, login as a standard user type and launch DaTQUANT for PC. In the upper right corner of the patient data manager window, verify that the version is listed as 2.0000.

If DaTQUANT opens with the Registration screen or if the previously processed patients are not visible, then:

- 1 Close and Exit DaTQUANT.
- 2 Login to Windows as an Administrator-Type user.
- 3 In Windows 10, use Start -> DaTQUANT -> Backup Restore.
- 4 Select Restore and locate the previously saved backup file.



- 8 To expand the list of customization options and to restore the previously saved customization settings:

- 1 Click “Customize” in the lower right corner of the main window.
- 2 Click “Reset...” in the upper right corner of the pop-up window.



- 3 In the pop-up window, select “Parameter Customization” and click “Ok”.
- 4 Click “Ok” to close the Customize window.
- 5 Repeat sub-steps 1-4 from step 7 to restore the previously saved backup file (even if previously done, it needs to be restored again).

GE Healthcare Limited, Pollards Wood, Nightingales Lane,
Chalfont St Giles, Buckinghamshire, England HP8 4SP
www.gehealthcare.com

General Electric Company reserves the right to make changes in specifications and features shown herein, or discontinue the product described at any time without notice or obligation. Contact your GE Representative for the most current information.

© 2020 General Electric Company.

GE, the GE Monogram and DaTQUANT are trademarks of General Electric Company.

Windows is a registered trademark of Microsoft Corporation in the United States.

06-2020 JB8509/PRT/OS INT'L ENGLISH

DaTQUANT is CE marked

