GE HealthCare

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1. Suzumura EA, Amato MB, Cavalcanti AB: Understanding recruitment maneuvers. Intensive Care Med 2016; 42: 908-11
2. Futier E, Marret E, Jaber S: Perioperative positive pressure ventilation: an integrated approach to improve pulmonary care. Anesthesiology 2014; 121: 400-8
3. Duggan M, Kavanagh BP: Pulmonary atelectasis: a pathogenic perioperative entity. Anesthesiology 2005; 102: 838-54
4. Lundquist H, Hedenstierna G, Strandberg A, Tokics L, Brismar B: CT-assessment of dependent lung densities in man during general anaesthesia. Acta Radiol 1995; 36: 626-32
5. Ladha K, Vidal Melo MF, McLean DJ, Wanderer JP, Grabitz SD, Kurth T, Eikermann M: Intraoperative protective mechanical ventilation and risk of postoperative respiratory complications: hospital based registry study. BMJ 2015; 351: h3646
6. Pelosi P, Brusasco C, Gama de Abreu M: Mechanical ventilation during general anesthesia, Principles and practice of mechanical ventilation, 3rd edition. Edited by Tobin MJ. New York, McGraw Hill Medical, pp 597-627
7. Guldner A, Kiss T, Serpa Neto A, Hemmes SN, Canet J, Spieth PM, Rocco PR, Schultz MJ, Pelosi P, Gama de Abreu M: Intraoperative Protective Mechanical Ventilation for Prevention of Postoperative Pulmonary Complications: A Comprehensive Review of the Role of Tidal Volume, Positive End-expiratory Pressure, and Lung Recruitment Maneuvers. Anesthesiology 2015; 123: 692-713
8. Slutsky AS, Ranieri VM: Ventilator-induced lung injury. N Engl J Med 2013; 369: 2126-36
9. Jabaudon M, Hamroun N, Roszyk L, Guerin R, Bazin JE, Sapin V, Pereira B, Constantin JM: Effects of a recruitment maneuver on plasma levels of soluble RAGE in patients with diffuse acute respiratory distress syndrome: a prospective randomized crossover study. Intensive Care Med 2015; 41: 846-55
10. Unzueta C, Tusman G, Suarez-Sipmann F, Bohm S, Moral V: Alveolar recruitment improves ventilation during thoracic surgery: a randomized controlled trial. Br J Anaesth 2012; 108: 517-24
11. Ferrando C, Mugarra A, Gutierrez A, Carbonell JA, Garcia M, Soro M, Tusman G, Belda FJ: Setting individualized positive end-expiratory Pressure level with a positive end-expiratory Pressure decrement trial after a recruitment maneuver improves oxygenation and lung mechanics during one-lung ventilation. Anesth Analg 2014; 118: 657-65
12. Futier E, Constantin JM, Polisi P, Chanques G, Kwiatkowski F, Jaber S, Bazin JE: Intraoperative recruitment maneuver reverses detrimental pneumoperitoneum-induced respiratory effects in healthy weight and obese patients undergoing laparoscopy. Anesthesiology 2010; 113: 1310-9
13. Futier E, Constantin JM, Paugam-Burtz C, Pascal J, Eurin M, Neuschwander A, Marret E, Beaussier M, Gutton C, Lefrant JY, Allaouchiche B, Verzilli D, Leone M, De Jong A, Bazin JE, Pereira B, Jaber S: A trial of intraoperative low-tidal-volume ventilation in abdominal surgery. N Engl J Med 2013; 369: 428-37
14. Severgnini P, Selmo G, Lanza C, Chiesa A, Frigerio A, Bucuzzi A, Dionigi G, Novario R, Gregoret C, de Abreu MG, Schultz MJ, Jaber S, Futier E, Chiaranda M, Pelosi P: Protective mechanical ventilation during general anesthesia for open abdominal surgery improves postoperative pulmonary function. Anesthesiology 2013; 118: 1307-21
15. Hemmes SN, Gama de Abreu M, Pelosi P, Schultz MJ: High versus low positive end-expiratory pressure during general anaesthesia for open abdominal surgery (PROVHILO trial): a multicentre randomised controlled trial. Lancet 2014; 384: 495-503
16. Yang D, Grant MC, Stone A, Wu CL, Wick EC: A Meta-analysis of Intraoperative Ventilation Strategies to Prevent Pulmonary Complications: Is Low Tidal Volume Alone Sufficient to Protect Healthy Lungs? Ann Surg 2016; 263: 881-7
17. 中山 禎人: 麻酔中の肺胞リクルートメント手技 基礎と実践. 臨床麻酔 2016; 40: 21-29
18. Rothen HU, Neumann P, Berglund JE, Valtysson J, Magnusson A, Hedenstierna G: Dynamics of re-expansion of atelectasis during general anaesthesia. Br J Anaesth 1999; 82: 551-6

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© 2024 GE HealthCare. GE is a trademark of General Electric Company used under trademark license.