

# Air-Purifying Respirator Types & Protection Levels

## Recommended by CDC for COVID-19



Brief Issue Date: **March 24, 2020**

By now most people are familiar with the disposable half-mask respirators commonly called “N95s” in the healthcare setting. But a much broader array of respiratory products are available that meet United States Centers for Disease Control and Prevention (CDC) guidance for protection against COVID-19.

The CDC has specified that reusable National Institute for Occupational Safety and Health (NIOSH)-approved elastomeric respirators are a viable option for use by healthcare workers and the first responder community. This document outlines MSA Safety’s NIOSH-approved respirators and filter configurations that meet CDC’s guidance.

Two key considerations when selecting an air-purifying respirator (APR) for protection against COVID-19 are respirator type and protection level. While there are many types of respiratory protection, CDC recommends NIOSH-approved masks with a protection factor of N95 or higher for certain healthcare workers who may be exposed to COVID-19. Below you will find a comparison of respirator types, including supplied-air respirators; filter classifications regulated by NIOSH that meet CDC recommendations; and in Table 1, a detailed comparison of all respirator types, including supplied-air respirators.

### Types of Respirators



#### Half Mask, Reusable

An elastomeric, reusable mask that covers nose and mouth. The mask includes an inlet receptacle(s) for particulate filters and a one-way exhalation valve. The filters are disposable. The reusable mask can be disinfected.

*Filter Classifications: All; N95 and P100 most common*



#### Full Facepiece, Reusable

An elastomeric, reusable mask that covers eyes, nose, and mouth. The mask includes an inlet receptacle(s) for particulate filters and a one-way exhalation valve. The filters are disposable. The reusable mask can be disinfected.

*Filter Classifications: All; N95 and P100 most common*



#### Powered Air-Purifying Respirator (PAPR), Mask-Mounted

An elastomeric, reusable mask that covers eyes, nose and mouth. The mask includes an inlet receptacle for a motor/blower with attached particulate filter and a one-way exhalation valve. The motor/blower actively draws air through the particulate filter and delivers breathable air to the respirator inlet. The motor blower is typically powered by a belt-mounted battery. The filters are disposable. The reusable mask/PAPR assembly can be disinfected.

*Filter Classifications: All; HEPA/Type HE (similar particulate efficiency to P100)*



#### Powered Air-Purifying Respirator (PAPR), Belt-Mounted

A belt-mounted, battery-powered motor/blower that provides filtered air to a full facepiece or hood. The filters are disposable. The reusable mask/PAPR assembly can be disinfected.

*Filter Classifications: All; HEPA/Type HE (similar particulate efficiency to P100)*

### Filter Classifications (*all meet CDC recommendations for COVID-19*)

Filter Types	95	99	100
N (NON-OIL PROOF)	N95*	N99	N100
R (OIL-RESISTANT)	R95	R99	R100
P (OIL PROOF)	P95	P99	P100*

\* N95 and P100 are the most common and most available filter types. Higher levels of particulate protection are inclusive of all lower levels.

Term	Definition
N	Effective against particulate aerosols free of oil; time use restrictions may apply
R	Effective against all particulate aerosols; time use restrictions may apply
P	Effective against all particulate aerosols
95	95% efficient against 0.3 micron (most penetrating particle size)
99	99% efficient against 0.3 micron (most penetrating particle size)
100	99.97% efficient against 0.3 micron (most penetrating particle size)

## Filter Replacement Guidance for Reusable Respiratory Protection

Disposable filters for use in reusable respiratory protection should be discarded and replaced at least daily or more frequently as circumstances warrant according to CDC guidance. CDC guidance on COVID-19 is rapidly evolving. Always consult the CDC website for the latest recommendations.

Table 1: Comparative Analysis of Various Forms of Respiratory Protection Provided by MSA Safety

Respirator Type	Elastomeric Half-Mask Respirator	Full-Facepiece Respirator	Powered Air-Purifying Respirator, Tight-Fitting Facepiece	Powered Air-Purifying Respirator, Loose-Fitting Hood	Self-Contained Breathing Apparatus
<b>RESPIRATOR NAME</b>	<ul style="list-style-type: none"> <li>• Advantage® 420</li> <li>• Advantage 200LS</li> <li>• Comfo® Classic</li> </ul>	<ul style="list-style-type: none"> <li>• Advantage 1000</li> <li>• Advantage 3200</li> <li>• Advantage 4100</li> <li>• Advantage 4200</li> <li>• Ultra-Elite</li> <li>• Ultra-Twin/Ultravue</li> <li>• Millennium</li> <li>• G1 Facepiece with APR Adapter</li> </ul>	<ul style="list-style-type: none"> <li>• OptimAir Mask- Mounted</li> <li>• Optimair TL Tight-Fitting Facepiece</li> <li>• Responder PAPR</li> </ul>	<ul style="list-style-type: none"> <li>• OptimAir TL</li> </ul>	<ul style="list-style-type: none"> <li>• G1 SCBA</li> <li>• G1 Industrial</li> </ul>
<b>COMPATIBLE PARTICULATE FILTER TYPES</b>	<ul style="list-style-type: none"> <li>• N95 Flexi Filter or prefilter with chemical cartridge</li> <li>• R95 Flexi Filter or prefilter with chemical cartridge</li> <li>• P100</li> </ul>	<ul style="list-style-type: none"> <li>• N95 Flexi Filter or prefilter with chemical cartridge</li> <li>• R95 Flexi Filter or prefilter with chemical cartridge</li> <li>• P100</li> <li>• Riot control</li> <li>• CBRN</li> </ul>	<ul style="list-style-type: none"> <li>• HEPA</li> <li>• CBRN for Responder PAPR</li> </ul>	<ul style="list-style-type: none"> <li>• HEPA</li> </ul>	<ul style="list-style-type: none"> <li>• Supplied Air</li> </ul>
<b>MEETS CDC GUIDANCE FOR COVID-19</b>	✓	✓	✓	✓	✓
<b>Type and Scope of Protection Provided</b>					
<b>DISPOSABLE VS. REUSABLE</b>	Reusable after disinfection, and replacement of filtering element	Reusable after disinfection, and replacement of filtering element	Reusable after disinfection, and replacement of filtering element	Reusable after disinfection, and replacement of filtering element	Reusable after disinfection
<b>ORGANS PROTECTED</b>	<ul style="list-style-type: none"> <li>• Respiratory</li> </ul>	<ul style="list-style-type: none"> <li>• Respiratory</li> <li>• Eye</li> <li>• Face</li> </ul>	<ul style="list-style-type: none"> <li>• Respiratory</li> <li>• Eye</li> <li>• Face</li> </ul>	<ul style="list-style-type: none"> <li>• Respiratory</li> <li>• Eye</li> <li>• Face</li> </ul>	<ul style="list-style-type: none"> <li>• Respiratory</li> <li>• Eye</li> <li>• Face</li> </ul>
<b>COMPATIBLE WITH OTHER FORMS OF PPE</b>	Yes	Yes	Yes	Yes	Yes
<b>DURATION OF USE</b>	Dependent upon use conditions	Dependent upon use conditions	Dependent upon use conditions	Dependent upon use conditions	Dependent upon use conditions
<b>LEVEL OF TRAINING REQUIRED FOR PROPER USE</b>	Moderate	Moderate	High	High	Highest

Table 1: Comparative Analysis of Various Forms of Respiratory Protection Provided by MSA Safety

Continued

Respirator Type	Elastomeric Half-Mask Respirator	Full-Facepiece Respirator	Powered Air-Purifying Respirator, Tight-Fitting Facepiece	Powered Air-Purifying Respirator, Loose-Fitting Hood	Self-Contained Breathing Apparatus
<b>RESPIRATOR NAME</b>	<ul style="list-style-type: none"> <li>• Advantage® 420</li> <li>• Advantage 200LS</li> <li>• Comfo® Classic</li> </ul>	<ul style="list-style-type: none"> <li>• Advantage 1000</li> <li>• Advantage 3200</li> <li>• Advantage 4100</li> <li>• Advantage 4200</li> <li>• Ultra-Elite</li> <li>• Ultra-Twin/Ultravue</li> <li>• Millennium</li> <li>• G1 Facepiece with APR Adapter</li> </ul>	<ul style="list-style-type: none"> <li>• OptimAir Mask- Mounted</li> <li>• Optimair TL Tight-Fitting Facepiece</li> <li>• Responder PAPR</li> </ul>	<ul style="list-style-type: none"> <li>• OptimAir TL</li> </ul>	<ul style="list-style-type: none"> <li>• G1 SCBA</li> <li>• G1 Industrial</li> </ul>
<b>COMPATIBLE PARTICULATE FILTER TYPES</b>	<ul style="list-style-type: none"> <li>• N95 Flexi Filter or prefilter with chemical cartridge</li> <li>• R95 Flexi Filter or prefilter with chemical cartridge</li> <li>• P100</li> </ul>	<ul style="list-style-type: none"> <li>• N95 Flexi Filter or prefilter with chemical cartridge</li> <li>• R95 Flexi Filter or prefilter with chemical cartridge</li> <li>• P100</li> <li>• Riot control</li> <li>• CBRN</li> </ul>	<ul style="list-style-type: none"> <li>• HEPA</li> <li>• CBRN for Responder PAPR</li> </ul>	<ul style="list-style-type: none"> <li>• HEPA</li> </ul>	<ul style="list-style-type: none"> <li>• Supplied Air</li> </ul>
<b>MEETS CDC GUIDANCE FOR COVID-19</b>	✓	✓	✓	✓	✓
<b>Government Certification and Performance</b>					
<b>GOVERNMENT CERTIFICATION(S)</b>	NIOSH 42 CFR Part 84	NIOSH 42 CFR Part 84	NIOSH 42 CFR Part 84	NIOSH 42 CFR Part 84	NIOSH 42 CFR Part 84, NFPA 1918 (2018) for structural firefighting applications
<b>BREATHING RESISTANCE</b>	Moderate	Moderate	Low	Low	Supplied air
<b>ADAPTABLE TO A RANGE OF ADULT FACE SIZES</b>	<ul style="list-style-type: none"> <li>• Small</li> <li>• Medium</li> <li>• Large</li> </ul>	<ul style="list-style-type: none"> <li>• Small</li> <li>• Medium</li> <li>• Large</li> </ul>	Adaptable to a range of adult face sizes	Adaptable to a range of adult face sizes	Adaptable to a range of adult face sizes
<b>ASSIGNED PROTECTION FACTOR (OSHA)</b>	10	50	1000	1000 SWPF Available	10,000
<b>FIT TESTING REQUIRED</b>	Yes	Yes	Yes	—	Yes

Table 1: Comparative Analysis of Various Forms of Respiratory Protection Provided by MSA Safety

Continued

Respirator Type	Elastomeric Half-Mask Respirator	Full-Facepiece Respirator	Powered Air-Purifying Respirator, Tight-Fitting Facepiece	Powered Air-Purifying Respirator, Loose-Fitting Hood	Self-Contained Breathing Apparatus
<b>RESPIRATOR NAME</b>	<ul style="list-style-type: none"> <li>• Advantage® 420</li> <li>• Advantage 200LS</li> <li>• Comfo® Classic</li> </ul>	<ul style="list-style-type: none"> <li>• Advantage 1000</li> <li>• Advantage 3200</li> <li>• Advantage 4100</li> <li>• Advantage 4200</li> <li>• Ultra-Elite</li> <li>• Ultra-Twin/Ultravue</li> <li>• Millennium</li> <li>• G1 Facepiece with APR Adapter</li> </ul>	<ul style="list-style-type: none"> <li>• OptimAir Mask- Mounted</li> <li>• Optimair TL Tight-Fitting Facepiece</li> <li>• Responder PAPR</li> </ul>	<ul style="list-style-type: none"> <li>• OptimAir TL</li> </ul>	<ul style="list-style-type: none"> <li>• G1 SCBA</li> <li>• G1 Industrial</li> </ul>
<b>COMPATIBLE PARTICULATE FILTER TYPES</b>	<ul style="list-style-type: none"> <li>• N95 Flexi Filter or prefilter with chemical cartridge</li> <li>• R95 Flexi Filter or prefilter with chemical cartridge</li> <li>• P100</li> </ul>	<ul style="list-style-type: none"> <li>• N95 Flexi Filter or prefilter with chemical cartridge</li> <li>• R95 Flexi Filter or prefilter with chemical cartridge</li> <li>• P100</li> <li>• Riot control</li> <li>• CBRN</li> </ul>	<ul style="list-style-type: none"> <li>• HEPA</li> <li>• CBRN for Responder PAPR</li> </ul>	<ul style="list-style-type: none"> <li>• HEPA</li> </ul>	<ul style="list-style-type: none"> <li>• Supplied Air</li> </ul>
<b>MEETS CDC GUIDANCE FOR COVID-19</b>	✓	✓	✓	✓	✓
<b>Design Related Characteristics and Attributes</b>					
<b>COMPOSITION</b>	Synthetic, metallic, and natural materials	Synthetic, metallic, and natural materials	Synthetic, metallic, and natural materials	Synthetic, metallic, and natural materials	Synthetic, metallic, and natural materials
<b>EXTERNAL FORMS OF POWER REQUIRED</b>	No	No	Yes, battery	Yes, battery	Yes, cylinder of breathable air, battery on some models
<b>COST</b>	Moderate	Moderate	High	High	Highest

**The situation with COVID-19 is rapidly evolving. Organizations like the CDC are providing updated information and guidance as it becomes available.**

For the most-up-to date information on COVID-19, you should regularly consult guidance being published by national and international organizations, such as the CDC, the National Institutes of Health (NIH), the World Health Organization (WHO), the European Centre for Disease Prevention and Control (ECDC), and/or your local health authority. Guidance on COVID-19, including information on actions needed to prevent, control, and manage contact with the virus, is available at the following websites:

<b>US CDC</b>	<a href="https://www.cdc.gov/coronavirus/index.html">https://www.cdc.gov/coronavirus/index.html</a>
<b>US NIH</b>	<a href="https://www.nih.gov/news-events/news-releases/nih-officials-discuss-novel-coronavirus-recently-emerged-china">https://www.nih.gov/news-events/news-releases/nih-officials-discuss-novel-coronavirus-recently-emerged-china</a>
<b>WHO</b>	<a href="https://www.who.int/health-topics/coronavirus">https://www.who.int/health-topics/coronavirus</a>
<b>ECDC</b>	<a href="https://www.ecdc.europa.eu/en/coronavirus">https://www.ecdc.europa.eu/en/coronavirus</a>

Note: This Bulletin contains only a general description of the products shown. While product uses and performance capabilities are generally described, the products shall not, under any circumstances, be used by untrained or unqualified individuals. The products shall not be used until the product instructions/user manual, which contains detailed information concerning the proper use and care of the products, including any warnings or cautions, have been thoroughly read and understood. Specifications are subject to change without prior notice.

MSA operates in over 40 countries worldwide. To find an MSA office near you, please visit [MSAsafety.com/offices](https://www.MSA.com/offices).