3M Science. Applied to Life.™

**3M™ Emphaze™ AEX Hybrid Purifier
Scale-Up Capsules**

**3M™ Emphaze™ AEX Hybrid Purifier
Scale-Up Capsulen**

Cápsulas filtrantes a escala

3M™ Emphaze™ AEX Hybrid

**Capsule 3M™ Emphaze™ AEX Hybrid
Purifier Scale-Up**

**Capsules Scale-Up 3M™ Emphaze™ AEX
Hybrid Purifier**


**3M™ Emphaze™ AEX Hybrid Purifier
Scale-Up-capsules**

BV60R (EMP503AEX020R)

BV120R (EMP513AEX020R)

BV360R (EMP533AEX020R)

EN	Safety Information & Installation and Operation Instructions	1-11
DE	Sicherheitshinweise	13-24
ES	Información de seguridad	26-37
IT	Informazioni sulla sicurezza e istruzioni di installazione e funzionamento	39-50
FR	Consignes de sécurité	52-63
NL	Veiligheidsinformatie	65-76

EN English Instructions Are the Original Instructions. All Other Languages Are Translations of the Original Instructions.

Failure to follow installation, operation, and maintenance instructions may result in system failure and will void warranty.

Safety Information

Read, understand, and follow all safety information contained in these instructions prior to installation and use of the 3M™ Emphaze™ AEX Hybrid Purifier Scale-Up Capsules. Retain these instructions for future reference.

Intended Use:

Single use processing of aqueous based biological pharmaceuticals (drugs) and vaccines strictly following the product operating instructions and cGMP requirements, where applicable. Customers must determine whether the 3M product is suitable for a specific application based on a risk assessment that considers the product leachable characteristics and its impact on drug safety.

Prohibited Use:

Do not use as a component in a medical device that is regulated by any agency, and/or globally exemplary agencies, including but not limited to: a) FDA, b) European Medical Device Directive (MDD), c) Japan Pharmaceuticals and Medical Devices Agency (PMDA). Do not use in applications involving permanent implantation into the body, life-sustaining medical applications, applications requiring Food Contact compliance.

Explanation of Signal Word Consequences	
△ WARNING:	Indicates a hazardous situation which, if not avoided, could result in serious injury or death.
△ CAUTION:	Indicates a hazardous situation which, if not avoided, could result in minor or moderate injury and/or property damage.
NOTICE:	Indicates a situation which, if not avoided, could result in property damage.
Explanation of Safety and Related Symbols	
	WARNING: Explosion hazard
	Single use
	Wear PPE
	Storage conditions

 **WARNING**

To reduce the risk of injury from possible end user exposure associated with a) biological contamination and b) degradation products released from the purification media; and

To reduce the risk of product loss associated with a) chemical release from purification media, and b) nullification of Supplier's USP<88> Class VI compliance of filter capsule components:

- Use the installation procedure specified in Operation Instructions.
- The capsule must be subjected to a preconditioning flush prior to use. Use the procedure specified in Operation Instructions.
- At least 54 L/m² preconditioning flush volume
- Carry out preconditioning flush procedures of 3M™ Zeta Plus™ capsules and 3M™ Emphaze™ AEX Hybrid Purifier capsules separately.
- After pre-use NaOH sanitization, flush the system with appropriate aqueous buffer until the pH range of 5 to 9 has been attained and monitor the pH of the capsule effluent to determine when the baseline has been reached (e.g., ± 0.5 pH units from the pH of the buffer used to flush the media).
- Capsules are designed for single use only. Do not reuse.
- Do not store capsules under conditions outside of the recommended storage temperature range.
- Do not use past the expiration date.
- Do not use with organic solvents or flammable liquids.
- Use only with aqueous based process fluids in the pH range of 5 to 9.
- Use only alkaline resistant capsules.
- Do not use with or expose this product to hot liquids >40°C (104°F).
- Do not use the product on unlevelled surfaces.
- Do not gamma or ethylene oxide sterilize the capsules.
- Do not steam-in-place sterilize the capsules.
- Only expose the capsule to 1 sterilization or 1 sanitization cycle prior to use.
- Do not autoclave and base sanitize the same capsule prior to use.
- Do not autoclave or base sanitize capsules with the following part numbers BV120 (EMP513AEX020A), and BV360 (EMP533AEX020A).
- During pre-use autoclave sterilization do not exceed a temperature of 121°C or a duration longer than 30 minutes for numbers BV60R (EMP503AEX020R), BV120R (EMP513AEX020R), and 40 minutes for BV360R (EMP533AEX020R).
- Do not use with or expose this product to process liquids until capsule temperature goes <40°C (104°F).
- Only use NaOH as sanitizing chemical and do not exceed 1M NaOH concentration or 1-hour duration for sanitization.
- Do not NaOH sanitize 3M Emphaze AEX Hybrid Purifier and 3M Zeta Plus Encapsulated capsules together in series.
- Use only 70% ethyl and/or isopropyl alcohol to wipe the exterior of the capsules.
- Do not use with or expose this product to pressurized steam or sanitizing chemicals other than NaOH.

 **WARNING**

To reduce the risks associated with chemical exposure, impact, crush, tipping hazards, property damage from product loss, bypass or cross contamination due to pressure barrier failure, connector failure, purification media or capsule failure:

- Use the installation procedure specified in Operation Instructions.
- Wear Personal Protective Equipment (PPE) to prevent chemical splash.
- Hardware Qualification testing is recommended prior to use.
- The capsule must be subjected to a preconditioning flush prior to use. Use the procedure specified in Operation Instructions.
- After pre-use NaOH sanitization, flush the system with appropriate aqueous buffer until the pH range of 5 to 9 has been attained and monitor the pH of the capsule effluent to determine when the baseline has been reached (e.g., ± 0.5 pH units from the pH of the buffer used to flush the media).
- Capsules are designed for single use only. Do not reuse.
- Do not modify or alter the capsules.
- Do not put this product into service if any damage is occurred to the capsule or any of the connecting parts.
- Do not install where fluid pressure exceeds 3.1 bar (45 psig).
- Do not use in organic solvent based processing.
- Use only with aqueous based process fluids in the pH range of 5 to 9.
- Do not use with or expose this product to hot liquids $>40^{\circ}\text{C}$ (104°F).
- Do not use the product on unlevelled surfaces.
- Do not use the product for more than 8 hours.
- Do not use this product for continuous service with compressed gasses. Use of compressed gas is permissible for post-use integrity testing and blow-down purposes only. Limit post-use integrity test and blow-down gas pressures to less than 1.4 bar (20 psig) and temperature less than 25°C for no longer than 30 minutes.
- Do not hang other items on this product.
- Do not gamma or ethylene oxide sterilize the capsules.
- Do not steam-in-place sterilize the capsules.
- Only expose the capsule to 1 sterilization or 1 sanitization cycle prior to use.
- Do not autoclave and base sanitize the same capsule prior to use.
- Do not autoclave or base sanitize capsules with the following part numbers BV120 (EMP513AEX020A), and BV360 (EMP533AEX020A).
- During pre-use autoclave sterilization do not exceed a temperature of 121°C or a duration longer than 30 minutes for numbers BV60R (EMP503AEX020R), BV120R (EMP513AEX020R), and 40 minutes for BV360R (EMP533AEX020R).
- Do not use with or expose this product to process liquids until capsule temperature goes $<40^{\circ}\text{C}$ (104°F).
- Only use NaOH as sanitizing chemical and do not exceed 1M NaOH concentration or 1-hour duration for sanitization.

⚠ WARNING

To reduce the risks associated with chemical exposure, impact, crush, tipping hazards, property damage from product loss, bypass or cross contamination due to pressure barrier failure, connector failure, purification media or capsule failure:

- Do not expose the purification capsule to base sanitization chemicals for more than 60 minutes.
- Do not pressurize the capsule more than 5 psi during introduction of NaOH solution and static soak during sanitization.
- Use only 70% ethyl and/or isopropyl alcohol to wipe the exterior of the capsules.
- Do not use with or expose this product to pressurized steam or sanitizing chemicals other than NaOH.

⚠ WARNING

To reduce the risks associated with chemical exposure due to pressure barrier failure, connector failure, purification media or capsule failure from use of disinfectants within the product:

- Use up to 1 M Sodium Hydroxide (NaOH) for pre-use and post-use sanitization and up to 5% Bleach (NaClO) for post-use sanitization
- Do not expose the purification capsule to base sanitization chemicals more than 60 minutes.
- Do not pressurize the capsule more than 5 psi during introduction of NaOH solution and static soak during sanitization.
- Wear personal protective equipment (PPE) to prevent chemical splash.

⚠ CAUTION

To reduce the risk of burn or exposure injuries associated with autoclave sterilization.

- Wear personal protective equipment (PPE) to handle hot capsules.

⚠ CAUTION

To reduce the risk of fluid exposure:

- Evaluate the risk of your process liquids and determine an appropriate protocol to dispose of spent capsules and all other waste.
- Dispose of capsules and all other waste in accordance with federal, state, and local regulations or your local country's laws and regulations.
- The protective caps on the capsule connections are designed for the protection of capsule parts. The protective caps are not intended to contain process liquids during the disposal of spent capsules.

Table of Contents

1. Operating Parameters and Product Specifications page 5
2. Single Stage Use Instructions. page 7
3. Capsule Use with Pre-use Autoclave Sterilization page 8
4. Capsule Use with Pre-use NaOH Sanitization. page 9
5. Disposal page 10

1. Operating Parameters and Product Specifications

Table 1: Operating Parameters

Liquid	
Maximum Inlet Pressure	3.1 bar (45 psig)
Maximum Differential Pressure	2.4 bar (35 psid)
Maximum Operating Temperature	40°C (104°F)
Post-Use Compressed Gas for Blow Down Only ¹	
Maximum Inlet Pressure ¹	1.4 bar (20 psig)
Maximum Operating Temperature	25°C (77°F)
Required Preconditioning Flush	
Preconditioning Flush Volume	54 L/m ² (15.3 mL flush solution per mL of bed volume)
Recommended Flow Rate	210 L/m ² /h (1 mL/min per mL of bed volume)

Product Specifications			
	BV60R	BV120R	BV360R
	EMP503AEX020R	EMP513AEX020R	EMP533AEX020R
Dimensions			
Height x Diameter	4" x 9 1/2" (10.1 cm x 24.1 cm)		6" x 9 1/2" (15.2 cm x 24.1 cm)
Bed Volume	60 mL	120 mL	360 mL
Filter Surface Area	170 cm ²	340 cm ²	1020 cm ²
Approximate Weight			
Dry	1.0 kg	1.1 kg	1.6 kg
Wet Post Blow down	1.1 kg	1.2 kg	2.1 kg
Materials of Construction			
Filter Media	Quaternary Ammonium Functional Polymer Grafted to a Polypropylene Nonwoven		
Membrane	Polyamide		
Membrane Support	Polypropylene Nonwoven		
Capsule Shell	Polysulfone		
Separator	Polypropylene	Filled Polypropylene	
Flow Inhibitor	Polypropylene	N/A	
End Cap, Edge Seal	Filled Polypropylene		
O-ring Support, Vent Cap	Polypropylene		
O-ring	Fluorocarbon		
Approximate Volume			
Capsule Fill Volume ²	550 mL	550 mL	1.4 L
Holdup Volume Post Blow Down ³	100 mL	150 mL	460 mL

Product Specifications			
	BV60R	BV120R	BV360R
	EMP503AEX020R	EMP513AEX020R	EMP533AEX020R
Pre-use Sterilization and Sanitization			
Autoclave Sterilization (Only 1 autoclave cycle)	121°C for 30 minutes using Pre-vac cycle		121°C for 40 minutes using Pre-vac cycle
Base Sanitization (Only 1 base sanitization cycle)	1M Sodium Hydroxide (NaOH) for 60 minutes		
Post-use Sterilization and Sanitization			
Autoclave Sterilization	Post use 1 autoclave cycle at 126°C for 30 minutes		
Sanitization	1M Sodium Hydroxide (NaOH) or 5% Bleach (NaClO) for post-use sanitization		
Required Preconditioning Flush Volume⁴			
Flush Volume, No Pre-treatment	0.9 L	1.8 L	5.5 L
Flush Volume, Post-autoclave	0.9 L	1.8 L	5.5 L
Flush Volume, Post Base Sanitization - 50 mM phosphate buffered saline at pH 7.4 ⁵	4.5 L	5.4 L	10.2 L
Miscellaneous			
Connector	¾" Sanitary Connectors		
Recommended Use Flow Rate	60 mL/min	120 mL/min	360 mL/min
Maximum Flow Rate ⁶	180 mL/min	360 mL/min	1080 mL/min
Storage Conditions ⁷	Controlled indoor temperature: 0 - 30°C (32 - 86°F) in original packaging		
Shelf-Life ⁸	2 years from date of manufacture @ 30°C maximum storage		

⚠ WARNING

¹ Do not use this product for continuous service with compressed gasses. The use of compressed gas is permissible for integrity testing and blow-down purposes. Limit integrity test and blow-down gas pressures to less than 1.4 bar (20 psig) and temperature less than 25°C no more than 30 min.

² Capsule Fill Volume is defined as the volume of liquid that is required to fill the capsule.

³ Post Blow down Hold up Volume is defined as the volume of the residual liquid after gas (air or nitrogen) blow down.

⁴ The Required Preconditioning Flush is a saline solution, having a conductivity in the same range as that of 25 to 150 mM NaCl (typical 3-16 mS/cm). DO NOT use water to flush the capsule.

⁵ Review Regulatory Support File for recommended alternate buffer selection.

⁶ Do not allow the liquid differential pressure across the capsule to exceed 2.4 bar (35 psid). Maximum flow rate may be limited by system differential pressure.

⁷ Do not allow product to be exposed to freezing or extreme heat conditions.

⁸ USP <88> Class VI compliance and high adsorption capacity for up to 2 years from date of manufacture if stored as indicated. Product labeling includes the expiration date.

The Scale-Up Capsule must be flushed prior to use to remove the glycerin that is present as a media stabilizer. Flush the capsule with 54 L/m² (15.3 mL/bed volume) of a saline solution or buffer having a conductivity in the same range as that of 25 to 150 mM NaCl (3 - 16 mS/cm) at a flow rate no higher than 210 L/m²/h (LMH). DO NOT use water as the flush solution. Flush in the forward flow direction only. Intermittently bleed air off using the Luer Capped Air Vent on the capsule. Drain fluid from the capsule, if required, by opening the vent and outlet.

IMPORTANT NOTICE: Precondition flush the capsule prior to use.

 **WARNING**

To reduce the risk of Injury from possible end user exposure associated with degradation products released from the purification media; and To reduce the risk of product loss associated with a) chemical release from purification media and b) nullification of Supplier's USP<88> Class VI compliance of filter capsule components.

- The capsule must be subjected to a preconditioning flush prior to use.

2. Single Stage Use Instructions

2.A. Required Preconditioning Flush Procedure

1. Remove the sanitary connection protective covers.
2. Connect the capsule inlet and outlet with the appropriate tubing and process equipment.
3. Connect the capsule inlet to the reservoir containing the pre-conditioning flush solution.
4. Open the capsule vent.
5. Close the outlet valve and open the inlet valve slowly to begin the flow of pre-conditioning flush solution through the capsule inlet.
6. When fluid emerges from the vent, close the vent while simultaneously opening the outlet valve slowly.
7. Pump the appropriate amount of flush solution per capsule size to complete the required preconditioning flush.
8. Close the inlet valve and stop the flow of preconditioning flush solution.

IMPORTANT NOTICE: Do not blow down 3M™ Emphaze™ AEX Hybrid Purifier capsules after pre-conditioning flush procedure.

2.B. Single-stage Use Instructions

1. Connect the inlet and outlet of the capsule to the upstream and downstream process equipment, respectively.
2. Connect the capsule inlet to the reservoir containing the bioprocess fluid.
3. Start the flow of bioprocess fluid and collect purified permeate.

IMPORTANT NOTICE: For optimal performance, run the filtration process at a flux of 160 - 630 LMH (0.75-3.0 Bed Volume per minute).

4. Monitor the efficiency of the bioprocess purification experiment by taking aliquots and analyzing the performance.
5. Stop the flow of bioprocess fluid at the appropriate volume based by on the laboratory study design.
6. The bioprocess fluid remaining on the upstream side of the capsule may then be processed through the system:
 - a. Open the vent to reduce system pressure to 0 bar (0 psig).
 - b. Close the inlet valve and open the outlet to drain.
 - c. Apply 0.3-0.7 bar (5-10 psig) air or nitrogen gas pressure to vent. Using a non-surgng pump or clean air supply to deliver regulated low-pressure air to push the process fluid through the capsule. Do not allow the inlet pressure to pulse or exceed 1.4 bar (20 psig).

- d. Monitor the capsule outlet. When liquid flow is reduced to a trickle, reduce the system pressure to 0 bar (0 psig).

2.C. Instructions To Use Multiple Scale-Up Capsules In Series Or Parallel Flow Path

3M™ Emphaze™ AEX Hybrid Purifier Scale-Up capsules can be used in parallel to add capacity to the process stream. 3M Emphaze AEX Hybrid Purifier Scale-Up Capsules can be used in series flow path with 3M™ Zeta Plus™ Scale-Up Capsules. When using combinations of 3M Zeta Plus and 3M Emphaze Scale-Up capsules, the 3M Emphaze AEX Hybrid Purifier capsules should be placed downstream from the 3M Zeta Plus capsules. The capsules can be stacked on top of each other.


⚠ WARNING

To reduce the risk of Injury from possible end user exposure associated with degradation products released from the purification media; and To reduce the risk of product loss associated with a) chemical release from purification media and b) nullification of Supplier's USP<88> Class VI compliance of filter capsule components.

- The capsule must be subjected to a preconditioning flush prior to use.

Figure 1. Series Operation for Two Scale-Up Capsules

For Series Operation for Two Scale-Up Capsules:

IMPORTANT NOTE: Precondition flush each capsule individually prior to use.

1. Follow the preconditioning flush procedure given in II. Single Stage Use Instructions for each capsule used (See IMPORTANT NOTE).
2. Configure the system as shown in Figure 2.
3. To process fluids, close valves V1, V2, and V3.

4. Open both vents V4 and V5 if closed.
5. Open V1 slowly to begin fluid flow.
6. When fluid exits vent V5, close V5, and simultaneously open V2.
7. When fluid exits vent V4, close V4 and simultaneously open V3.
8. The fluid remaining on the upstream side may then be processed through the system:
 - a. Allow the system pressure to go to 0 bar (0 psig).
 - b. Close V1 and open V5 and reduce the system pressure to 0 bar (0 psig).
 - c. Apply 0.3 – 0.7 bar (5 – 10 psi) air or nitrogen gas pressure to V5. Using a non-surging pump or regulated, low pressure, clean air supply, push process fluid through the capsule. Do not allow the inlet pressure to pulse or exceed 1.4 bar (20 psig).
 - d. Monitor the system outlet. When liquid flow is reduced to a trickle from the top capsule, reduce the system pressure to 0 bar (0 psig) and then close V5 and V2.
 - e. Open V4 and apply 0.3 – 0.7 bar (5 – 10 psi) air or nitrogen gas pressure.
 - f. Monitor the system outlet. When liquid flow is reduced to a trickle, reduce the system pressure to 0 bar (0 psig) and close V4.

3. Capsule Use with Pre-use Autoclave Sterilization

3.A. Recommended Pre-use Autoclave Sterilization of Capsules

1. Remove the sanitary connection protective covers and vent cap from the capsule.

2. Wrap the sanitary connections and vent with sterilization wrap.
3. Place the capsule in the autoclave.
4. Set-up the autoclave cycle:
 - a. Select pre-vac cycle for all capsule sizes.
 - b. Set the sterilization temperature at 121°C.
 - c. Set the autoclave duration.
 - i. BV60R and BV120R capsules for 30 minutes.
 - ii. BV360R capsules for 40 minutes.
5. Start the autoclave cycle.
6. At the end of autoclave cycle, remove the capsules from the autoclave and let the capsule cool down at ambient temperature for at least 60 minutes.
7. Remove the Sterilization Wrap.

⚠ CAUTION

To reduce the risk of burn or exposure injuries associated with autoclave sterilization:

- Wear personal protective equipment (PPE) to handle hot capsules.

3.B. Required Preconditioning Flush Procedure for the Pre-use Autoclaved Capsules

⚠ WARNING

To reduce the risk of injury from possible end user exposure associated with a) biological contamination and b) degradation products released from the purification media; and To reduce the risk of product loss associated with a) chemical release from purification media, and b) nullification of Supplier's USP<88> Class VI compliance of filter capsule components:

To reduce the risks associated with chemical exposure, impact, crush, tipping hazards, property damage from product loss, bypass or cross contamination due to pressure barrier failure, connector failure, purification media or capsule failure:

Do not use with or expose this product to process liquids until capsule temperature goes <40°C (104°F).

Perform Required Preconditioning Flush procedure as described in Section 2.A.

IMPORTANT NOTICE: Do not blow down 3M™ Emphaze™ AEX Hybrid Purifier capsules after pre-conditioning flush procedure.

3.C. Single-stage Use

Perform Single-Stage Use as described in Section 2.B.

3.D. Two-stage Use

Perform Two-Stage Use as described in Section 2.C.

4. Capsule Use with Pre-use NaOH Sanitization

4.A. Pre-use NaOH Sanitization of Capsules

Recommended maximum conditions: Only use NaOH as sanitizing chemical and do not exceed 1M NaOH concentration or 1 hour duration for sanitization.

⚠ WARNING

To reduce the risk of product loss associated with a) chemical release from purification media, and b) nullification of Supplier's USP<88> Class VI compliance of filter capsule components:

Do not NaOH sanitize 3M™ Emphaze™ AEX Hybrid Purifier and 3M™ Zeta Plus™ Encapsulated capsules together in series.

1. Remove the sanitary connection protective covers and vent cap from the capsule.

2. Connect the capsule inlet to the reservoir containing the 1 M NaOH solution.
3. Open the capsule vent.

⚠ WARNING

To reduce the risks associated with chemical exposure due to pressure barrier failure, connector failure, purification media or capsule failure from use of disinfectants within the product: Do not pressurize the capsule more than 5 psi during introduction of NaOH solution and static soak during sanitization.

4. When fluid emerges from the vent, close the vent while simultaneously opening the outlet valve slowly.
5. When fluid emerges from the capsule outlet, stop the flow and close the capsule outlet valve.
6. Allow the capsule to incubate with the 1 M NaOH solution. Incubation time should not exceed 1 hour.
7. After the incubation is complete, drain the NaOH solution from the capsule.

IMPORTANT NOTICE: Do not blow down the 3M™ Emphaze™ AEX Hybrid Purifier capsules after NaOH sanitization procedure.

4.B. Required Preconditioning Flush Procedure for the Pre-use NaOH Sanitized Capsules

⚠ WARNING

To reduce the risk of product loss associated with a) chemical release from purification media, and b) nullification of Supplier's USP<88> Class VI compliance of filter capsule components:
After pre-use NaOH sanitization, flush the system with appropriate aqueous buffer until the pH range of 5 to 9 has been attained and monitor the pH of the capsule effluent to determine when the baseline has been reached (e.g., ± 0.5 pH units from the pH of the buffer used to flush the media).

Perform Required Preconditioning Flush procedure as described in Section 2.A.

IMPORTANT NOTICE: Do not blow down with a gas the 3M Emphaze AEX Hybrid Purifier capsules after preconditioning flush procedure.

4.C. Single-stage Use

Perform Single-Stage Use as described in Section 2.B.

4.D. Two-stage Use

Perform Two-Stage Use as described in Section 2.C.

5. Disposal

1. After fluid processing is complete, drain the fluid from the capsule until the internal pressure on the capsule is reduced to 0 bar (0 psig).
2. Drain any remaining fluid from the system and remove the external connections from the inlet, outlet, and vent.

⚠ CAUTION

To reduce the risk of fluid exposure:

- Evaluate the risk of your process liquids and determine an appropriate protocol to dispose of spent capsules and all other waste.
- Dispose of capsules and all other waste in accordance with federal, state, and local regulations or your local country's laws and regulations.

Intended Use:

Single use processing of aqueous based biological pharmaceuticals (drugs) and vaccines strictly following the product operating instructions and cGMP requirements, where applicable. Customers must determine whether the 3M product is suitable for a specific application based on a risk assessment that considers the product leachable characteristics and its impact on drug safety.

Prohibited Use:

Do not use as a component in a medical device that is regulated by any agency, and/or globally exemplary agencies, including but not limited to: a) FDA, b) European Medical Device Directive (MDD), c) Japan Pharmaceuticals and Medical Devices Agency (PMDA). Do not use in applications involving permanent implantation into the body, life-sustaining medical applications, applications requiring Food Contact compliance.

Product Selection and Use:

Many factors beyond 3M's control and uniquely within user's knowledge and control can affect the use and performance of a 3M product in a particular application. As a result, customer is solely responsible for evaluating the product and determining whether it is appropriate and suitable for customer's application, including conducting a workplace hazard assessment and reviewing all applicable regulations and standards (e.g., OSHA, ANSI, etc.). Failure to properly evaluate, select, and use a 3M product and appropriate safety products, or to meet all applicable safety regulations, may result in injury, sickness, death, and/or harm to property.

Warranty, Limited Remedy, and Disclaimer:

Unless a different warranty is specifically stated on the applicable 3M product packaging or product literature (in which case such warranty governs), 3M warrants that each 3M product meets the applicable 3M product specification at the time 3M ships the product. 3M MAKES NO OTHER WARRANTIES OR CONDITIONS, EXPRESS OR IMPLIED, INCLUDING, BUT NOT LIMITED TO, ANY IMPLIED WARRANTY OR CONDITION OF MERCHANTABILITY, FITNESS FOR A PARTICULAR PURPOSE, OR ARISING OUT OF A COURSE OF DEALING, CUSTOM, OR USAGE OF TRADE. If a 3M product does not conform to this warranty, then the sole and exclusive remedy is, at 3M's option, replacement of the 3M product or refund of the purchase price.

Limitation of Liability

Except for the limited remedy stated above, and except to the extent prohibited by law, 3M will not be liable for any loss or damage arising from or related to the 3M product, whether direct, indirect, special, incidental, or consequential (including, but not limited to, lost profits or business opportunity), regardless of the legal or equitable theory asserted, including, but not limited to, warranty, contract, negligence, or strict liability.

3M Purification Inc.

400 Research Parkway
Meriden, CT 06450
U.S.A.
(800) 243-6894
(203) 237-5541
Fax (203) 630-4530
www.3Mpurification.com

3M, Emphaze and Zeta Plus are
trademarks of 3M Company.
Please recycle. Printed in U.S.A.
© 2018 3M Company. All rights reserved.
34-8723-4281-0

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