

- | | | | | | | | | | | | |
|---|---|------------------|-----------------------------------|---------------------|---|------------------|-------------------|--------------|---|--|--|
|  SGS
1300 Winetka Ave. N
Brooklyn Park, MN 55445
(763) 488-5700 |  | CUST ID # | | PRINT METHOD | | SGS # | | ROUND | | | |
| | | | | | | 6191839-2 | | 2 | | | |
| Production Art | Description: 3M™ Ranger™
Blood/Fluid Warming
Standard Flow Set | | Item Spec#: 34-8725-9950-0 | | DATE | | CHANGE LOG | | INT | | |
| | Base File Name:
PA: 34872599500.indd
GA:
Software: InDesign 2020 | | Supersedes#: | | 12.09.20 Update page count and RXOnly
12.15.20 Final Release | | | | kat
kat | | |
| | Structure#: | | Cat/Product#: 24355 | | | | | | | | |
| | Reference: 24355 | | Requester: Frank Boeshart | | | | | | | | |
| Die # / Doc. Size: 6.5" x 11" | | | | | | | | | | | |
| INKS | 
PROCESS BLOCK | | MATCH | | | | SPECIAL | |  | | |
| <small> LIMITATION OF LIABILITY – The supplied materials should be thoroughly inspected immediately upon receipt by the recipient for all content, placement and copy accuracy. It is also the responsibility of the recipient printer to examine and proofread all supplied materials before producing image carriers and printing final product. The recipient materials are found to be incorrect please contact the SGS representative immediately. The liability of SGS shall be limited to correcting the copy error. SGS is not responsible for any loss or damage including but not limited to the cost of product, print, image carriers, press delays, or other ancillary cost. </small> | | | | | | | | | | | |
| <small> UPC NOTE – Because of the difference in scanning equipment, SGS cannot guarantee the scanability of UPC barcodes. The Uniform Code Council recommends a UPC symbol be printed black on white at 100% with no truncation. All UPC codes are test scanned and verified before leaving our facility. </small> | | | | | | | | | | | |

en Indications for Use

The Ranger blood/fluid warmer is intended to warm blood, blood products, and liquids.

Patient Population and Settings

Adult and pediatric patients being treated in operating rooms, emergency trauma settings, or other areas where blood/fluids are infused.

Explanation of Signal Word Consequences
WARNING: Indicates a hazardous situation which, if not avoided, could result in death or serious injury.
CAUTION: Indicates a hazardous situation which, if not avoided, could result in minor or moderate injury.

WARNING: To reduce the risks associated with air embolism:

- Never infuse fluids if air bubbles are present in the fluid line
- Ensure all luer connections are tightened
- Do not mount the Ranger blood fluid warming system higher than patient level, as air embolism may result

WARNING: To reduce the risks of biohazard exposure, cross-contamination or infection:

- Do not reuse. This product is designed for single patient use only.
- Do not use if caps are not in place
- Ensure that all luer connections are securely tightened prior to priming set
- Do not exceed 300 mmHg pressure; ensure inline hand pressure pumps or other pumps do not exceed 300 mm/Hg pressure.

WARNING: To reduce the risk of leaking and loss of therapy:

- Do not reuse. This product is designed for single patient use only.
- Ensure that all luer connections are securely tightened prior to priming set.

WARNING: To reduce the risks associated with potential blood loss:

- Do not use in combination with an extracorporeal circuit.

WARNING: To reduce the risk associated with patient exposure to leakage current that exceeds Type BF equipment limits:

- Do not use in combination with other medical equipment without verifying that the total leakage current from the combined equipment does not exceed the safety limits for Type BF equipment.

CAUTION: To reduce the risk of leaking and loss of therapy:

- Do not exceed 300 mmHg pressure; ensure inline hand pressure pumps or other pumps do not exceed 300 mm/Hg pressure.

- Do not remove warming set from warming unit while in use.

CAUTION: To reduce the risks associated with warming platelets, cryoprecipitate, or granulocyte suspensions:

- Follow the AABB Guidelines for the use of blood warming devices which caution against warming when administering platelets, cryoprecipitate, or granulocyte suspensions.

CAUTION: To reduce the risks associated with potential material toxicity:

- Do not use in combination with an extracorporeal circuit as the Ranger blood/fluid warming system is not designed for use with a circulating blood system; therefore, biocompatibility testing has not been completed to assess genotoxicity.

Instructions for Use

Do not use if package has been previously opened or damaged. Refer to the 3M™ Ranger™ blood/fluid warming unit model 245 Operator's Manual for set up and use of the warming system.


Figure 1: Illustration of Disposable Set

- | | |
|-----------------------------|-----------------------------|
| A. Bag spikes | H. White outlet pinch clamp |
| B. Luer connections (4) | I. Roller Clamp |
| C. Drip chamber with filter | J. Injection port |
| D. Roller clamp | K. White outlet clamp |
| E. Blue inlet pinch clamp | L. Patient connection |
| F. Fluid warming cassette | M. Bubble trap holder |
| G. Auto-venting bubble trap | |

NOTE: Change the filter when it becomes clogged or according to institutional protocol. Additional filters may be obtained from 3M Company.

Figure 2


Peel To Open

Peel To Open

Figure 3

Instructions

Set-up and Priming the warming set

1. Adjust warming unit height to patient level or lower.
2. Open warming set pouch where indicated by the arrows (Fig. 2)
3. Before priming, slide the fluid warming cassette (Fig. 1, F) into the slot in the warming unit (Fig. 3, #1).
4. Place the auto-venting bubble trap (Fig. 1, G) into the holder (Fig. 1, M) on the warming unit.
5. Tighten luer connections (Fig. 1, B).
6. Close the clamp below the auto-venting bubble trap (Fig. 1, H)
7. Close the 2 clamps below the spikes
8. Insert the Ranger set spikes (Fig. 1, A) into the I.V. bag ports. Open the clamp (Fig. 3, #2).
9. Invert the drip chamber (Fig. 1, C) and fill 2/3 full. Turn upright.
10. Allow the auto-venting bubble trap (Fig. 1, G) to fill. Flow will stop when the chamber is properly filled. Gently tap the auto-venting bubble trap against the warming unit to expel any trapped air, if necessary.
11. Open the clamps below the auto-venting bubble trap and continue to prime the patient line.
12. Close any clamp to stop flow.
13. Turn the warming unit on (Fig. 3, #3).

The fluid warming set is now ready for use.

Removing Air from the Auto-venting bubble Trap

1. The auto-venting bubble trap is designed to continuously vent air from the patient line. Monitor fluid lines to ensure that they are air-free.
2. Gently tap the auto-venting bubble trap against the warming unit to expel any trapped air.

Removing the Disposable Set

1. Close the inlet clamp (Fig. 1, E) proximal to the fluid warming cassette and open all clamps distal to the cassette.
2. Allow fluid to flow into the patient (this may take 2 to 3 seconds) and close outlet clamp (Fig. 1, H).
3. Remove the fluid warming cassette from the warming unit and discard according to institutional protocol.

Please report a serious incident occurring in relation to the device to 3M and the local competent authority (EU) or local regulatory authority.

Symbol Glossary

Symbol Title	Symbol	Description and Reference
Manufacturer		Indicates the medical device manufacturer as defined in Medical Device Regulation (EU) 2017/745 formerly EU Directive 93/42/EEC. Source: ISO 15223, 5.1.1
Authorized Representative in European Community		Indicates the authorized representative in the European Community. Source: ISO 15223, 5.1.2, 2014/35/EU, and/or 2014/30/EU
Date of Manufacture		Indicates the date when the medical device was manufactured. Source: ISO 15223, 5.1.3
Use-by date		Indicates the date after which the medical device is not to be used. Source: ISO 15223, 5.1.4
Batch code		Indicates the manufacturer's batch code so that the batch or lot can be identified. Source : ISO 15223, 5.1.5
Catalogue number		Indicates the manufacturer's catalogue number so that the medical device can be identified. Source : ISO 15223, 5.1.6

Sterilized using ethylene oxide		Indicates a medical device that has been sterilized using ethylene oxide. Source: ISO 15223, 5.2.3
Do not use if package is damaged or open		Indicates a medical device that should not be used if the package has been damaged or opened. Source: ISO 15223, 5.2.8
Sterile fluid path		Indicates the presence of a sterile fluid path within the medical device in cases when other parts of the medical device, including the exterior, might not be supplied sterile. Source: ISO 15223, 5.2.9
Do not re-use		Indicates a medical device that is intended for one use or for use on a single patient during a single procedure. Source: ISO 15223, 5.4.2
Caution		Indicates the need for the user to consult the instructions for use for important cautionary information such as warnings and precautions that cannot, for a variety of reasons, be presented on the medical device itself. Source: ISO 15223, 5.4.4
Natural rubber latex is not present		Indicates natural rubber or dry natural rubber latex is not present as a material of construction within the medical device or the packaging of a medical device. Source: ISO 15223, 5.4.5 and Annex B
Medical Device		Indicates the item is a medical device.
Unique device identifier		Indicates bar code to scan product information into patient electronic health record
Importer		Indicates the entity importing the medical device into the EU.
Single sterile barrier system		To identify a single layer that makes up the sterile barrier system (SBS) and differentiate it from layers of protective packaging designed to prevent damage to the sterile barrier system and its contents
CE Mark		Indicates conformity to all applicable European Union Regulations and Directives with notified body involvement.
Rx Only		Indicates that U.S. Federal Law restricts this device to sale by or on the order of a healthcare professional. 21 Code of Federal Regulations (CFR) sec. 801.109(b)(1)

For more information see, [HCBRegulatory.3M.com](https://www.3m.com/HCBRegulatory)

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