3M

Bair Hugger™

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3M™ Bair Hugger™ Temperature Monitoring Sensor Model 360

Product Description

The 3M Bair Hugger Temperature Monitoring Sensor Model 360 is a non-invasive, non-sterile, single use, disposable temperature sensor that is used exclusively with the 3M Bair Hugger Temperature Monitoring System to measure core temperature of the patient. The sensor connects to the Model 370 control unit through the sensor cable.

Indications for Use

Measure, monitor, and trend body temperature of adult and pediatric patients.

This product is intended to be used by trained medical professionals in a clinical/surgical setting.

Patient Population and Settings

Adult and pediatric patients being treated in operating rooms, emergency departments and other departments in the hospital setting where patient temperature management is required.

Explanation of signal word consequences

Warning: Indicates a hazardous situation which, if not avoided, could result in death or serious injury.

Caution: Indicates a hazardous situation which, if not avoided, could result in minor or moderate injury.

Notice: Indicate a situation which, if not avoided, could result in property damage only.

WARNING: To reduce the risks associated with access to critical patient information or equipment:

- The 3M™ Bair Hugger™ control unit is to be attached to other equipment by authorized service personnel only.
- Do not re-install or re-locate the control unit.

WARNING: To reduce the risks associated with hazardous voltage, fire, and thermal energy hazards:

- Use only with 3M™ Bair Hugger™ temperature monitoring system components (sensors, cables, and power supply). Do not substitute other devices for the control unit, sensor, cables, or power supply.
- Do not use the temperature monitoring system (sensor, control unit, cables, or cords) in an MRI environment.
- Do not connect the monitor output of the device to a non-CF rated monitor input.

WARNING: To reduce the risks associated with hazardous voltage and fire:

- Keep power supply visible and accessible at all times. The plug on the power supply serves as the disconnect device. The power outlet shall be as close as practical and shall be easily accessible.
- Use only the power supply specified for this product and certified for the country of use.
- Use only a properly grounded power outlet; do not use extension cords or multiple portable socket outlets.
- Do not allow the power supply to get wet.

- Do not use the control unit when it appears the unit or system components are damaged. Contact your local 3M technical support representative.
- Do not service or modify the power supply, control unit, cables, sensor, or any part of the temperature monitoring system. There are no user serviceable parts.

WARNING: To reduce the risks associated with fire:

- The temperature monitoring system is not suitable for use in the presence of a flammable anesthetic mixture with air nitrous oxide.

WARNING: To reduce the risks associated with pressure, and crushed or entangled cords and connectors:

- Do not allow the patient to lie on any sensor cable or connector.
- Always position cables and cords away from the patient's body.
- Do not use a headband or other device to secure the sensor to the patient.

WARNING: To reduce the risks associated with incorrect system or sensor use:

- Limit use of the sensor to 24 hours. Extended use may compromise skin, cause degradation of material or performance.
- Avoid exposing the temperature monitoring sensor, sensor cable, power supply, and control unit to surgical skin surface preparation solutions or other fluids.
- Use an additional independent thermometer to measure body temperature during intentional hyperthermia or hypothermia therapy.
- Confirm unanticipated temperature reading with an independent thermometer, if needed.
- Do not reposition the sensor; repositioning may weaken the sensor adhesive, damage the sensor, or compromise the device performance.
- Avoid placing the sensor in the center of the forehead as this may affect sensor accuracy.

WARNING: To reduce the risks associated with exposure to biohazards:

- Follow facilities policies and procedures for disposal of contaminated materials.
- Always perform the decontamination procedure prior to returning the temperature monitoring system for service and prior to disposal.

WARNING: To reduce the risks associated with entanglement:

- Do not leave pediatric patients unattended while using the temperature monitoring system.

CAUTION: To reduce the risks associated with skin maceration:

- Do not use the sensor on damaged or compromised skin.

CAUTION: To reduce the risks associated with cross-contamination:

- Clean the sensor cable before connecting to a new sensor.

CAUTION: To reduce the risks associated with impact and facility medical device damage:

- Do not use the control unit as a handle to transport or to move the device to which it is attached.
- Do not initiate temperature monitoring unless the control unit is safely placed on a hard, flat surface or is securely mounted.

CAUTION: To reduce the risks associated with environmental contamination:

- Follow applicable regulations when disposing of this device or any of its electronic components.

NOTICE

1. The temperature monitoring system meets medical electronic interference requirements. If radio frequency interference with other equipment should occur, connect the unit to a different power source.
2. Protection against electrostatic discharge (ESD) is important when working with the temperature monitoring system. Always be mindful of the immediate work area and any potentially hazardous ESD situations that could impact temperature monitoring system performance.
3. To avoid temperature monitoring system damage that may impact performance:
 - Do not store the control unit and system components in a wet or damp place.
 - Do not spray cleaning solutions onto the control unit or into the sensor cable connector.
 - Do not immerse the control unit or system components in any liquid, use a dripping wet cloth for cleaning, or subject them to any sterilization process.
 - Do not use solvents such as acetone or thinner to clean the control unit; avoid abrasive cleaners.
 - Cleaning should be performed in accordance with hospital practices for cleaning OR equipment. After every use wipe the device surfaces, cords and cables. Avoid getting liquid into electronic ports. Use a damp soft cloth and a hospital approved mild detergent, germicidal disposable wipes, disinfecting towelettes, or antimicrobial spray. The following active ingredients are acceptable for cleaning:
 - Oxidizers (e.g. 10% Bleach)
 - Quaternary Ammonium Compounds (e.g. 3M™ Quat Disinfectant Cleaner)
 - Phenolics (e.g. 3M™ Phenolic Disinfectant Cleaner)
 - Alcohols (e.g. 70% Isopropyl Alcohol)
4. The sensor is not made with natural rubber latex.
5. To the full extent permitted by law, the manufacturer and/or importer declines all responsibility for injury resulting from the unit being used in conjunction with unapproved system components.

Warranties/3M Liability**Model 360 Product Warranty****WARRANTY INFORMATION**

3M warrants only that the Bair Hugger sensor (the "Product") will meet 3M's applicable physical product specification at the time of shipment to the customer. The terms of this Warranty do not affect or prejudice the statutory rights of a purchaser acquiring the Product otherwise than in the normal course of a business.

DISCLAIMER

3M DOES NOT OFFER ANY OTHER WARRANTY AND DOES NOT WARRANT THE PERFORMANCE, SAFETY OR SUCH OTHER

CHARACTERISTICS OF THE PRODUCT IN COMBINATION WITH OTHER MATERIALS. 3M SPECIFICALLY DOES NOT WARRANTY THE PRODUCT FOR ANY INTENDED OR UNINTENDED USES (WHETHER OR NOT FORESEEABLE FOR COMPATIBILITY OR SUITABILITY WITH OTHER COMPONENTS OR COMPATIBILITY WITH ANY METHODS OF MANUFACTURE OR CONVERSION). THE FOREGOING WARRANTY IS MADE IN LIEU OF ALL OTHER WARRANTIES, EXPRESSED OR IMPLIED, INCLUDING THE IMPLIED WARRANTIES OF MERCHANTABILITY, FITNESS FOR A PARTICULAR PURPOSE AND FREEDOM FROM NON-INFRINGEMENT.

LIMITATION OF REMEDY

If Product is proven to not meet 3M's specifications, the sole and exclusive remedy available and 3M's only obligation shall be, at 3M's option, to replace such quantity of Product which are proven out of specification or to refund the purchase price paid for the Product.

LIMITATIONS OF LIABILITIES

THE REMEDIES PROVIDED HEREIN ARE EXCLUSIVE REMEDIES AGAINST 3M FOR ANY ALLEGED OR ACTUAL NONCONFORMANCE TO SPECIFICATIONS OR DEFECT OR OTHER FAILURE IN PRODUCTS OR FOR 3M'S PERFORMANCE OF ITS SUPPLY OBLIGATIONS. UNDER NO CIRCUMSTANCES IS 3M LIABLE FOR ANY DIRECT, INDIRECT, INCIDENTAL, SPECIAL OR CONSEQUENTIAL DAMAGES (INCLUDING LOST PROFITS) IN ANY WAY RELATED TO THE PRODUCT OR TO PERFORMANCE OF THIS AGREEMENT UNDER ANY THEORY OF LAW INCLUDING, BUT NOT LIMITED TO, NEGLIGENCE AND STRICT LIABILITY.

Directions for Use

1. Ensure the power cord is connected to the back of the control unit and the power cord is plugged into an appropriate outlet. If desired, ensure the provided monitor cable is connected to the back of the control unit and to the YSI-400 type input on the patient monitor.


Figure 1

2. Use an alcohol wipe to clean and disinfect the patient's forehead above the orbital ridge. Allow the patient's forehead to completely dry.
3. Connect the sensor to the sensor cable connector. Ensure the sensor tab is properly connected and is fully inserted into the sensor cable connector (Figure 1).

Figure 2

4. Remove backing liner on the sensor to expose adhesive (Figure 2).
5. Position the sensor cable away from the patient's face and apply the sensor to the patient's forehead above the orbital ridge, on either side, avoiding the frontal sinus. (Figure 3).

Figure 3

- Gently press the edge of the sensor to ensure good adhesion with the skin.
Optional: Support the sensor cable with the provided sensor cable clip.
- The control unit will sense patient attachment after several seconds.
- The patient's temperature will display on the patient monitor once system equilibrium is reached.

Removing the sensor

- Disconnect the sensor cable from the sensor.
Disconnect the sensor cable clip from the bed or sheet, if used.
- Starting at the edge of the sensor, gently remove the sensor from the patient (Figure 4).


Figure 4

Note: Do not use the sensor tab to remove the sensor from the patient (Figure 4).

Note: If needed, use an alcohol swab along the edge of the sensor to help remove the sensor from the patient.

- Discard the sensor and clean the sensor cable according to institutional protocol.

Storage/Shelf Life/Disposal

Sensor Shelf Life: 36 months

Storage Conditions:

-20°C to 60°C

Store all components at room temperature and in a dry place when not in use.

Dispose of sensor according to institutional protocol.

Please report a serious incident occurring in relation to the device to 3M and the local competent authority (EU) or local regulatory authority.

For further information, please contact your local 3M representative.

Explanation of Symbols:

Symbol Title	Symbol	Description and Reference
Authorized Representative in European Community / European Union	EC REP	Indicates the authorized representative in the European Community / European Union. Source: ISO 15223, 5.1.2, 2014/35/EU, and/or 2014/30/EU
Catalogue number	REF	Indicates the manufacturer's catalogue number so that the medical device can be identified. Source: ISO 15223, 5.1.6
Caution	⚠	Indicates that caution is necessary when operating the device or control close to where the symbol is placed, or that the current situation needs operator awareness or operator action in order to avoid undesirable consequences. Source: ISO 15223, 5.4.4

CE Mark		Indicates conformity to all applicable European Union Medical Device Regulations and Directives.
Date of manufacture		Indicates the date when the medical device was manufactured. Source: ISO 15223, 5.1.3
Green Dot		Indicates a financial contribution to national packaging recovery company per European Directive No. 94/62 and corresponding national law. Packaging Recovery Organization Europe.
Importer		Indicates the entity importing the medical device into the locale. Source: ISO 15223, 5.1.8
Keep dry		Indicates a medical device that needs to be protected from moisture. Source: ISO 15223, 5.3.4
Manufacturer		Indicates the medical device manufacturer. Source: ISO 15223, 5.1.1
Medical Device		Indicates the item is a medical device. Source: ISO 15223, 5.7.7
Consult instructions for use or consult electronic instructions for use		Indicates the need for the user to consult the instructions for use. Source: ISO 15223, 5.4.3
Recycle electronic equipment		DO NOT throw this unit into a municipal trash bin when this unit has reached the end of its lifetime. Please recycle. Source: Directive 2012/19/EC on waste electrical and electronic equipment (WEEE)
Temperature limit		Indicates the temperature limits to which the medical device can be safely exposed. Source: ISO 15223, 5.3.7
Rx Only		Indicates that U.S. Federal Law restricts this device to sale by or on the order of a healthcare professional. 21 Code of Federal Regulations (CFR) sec. 801.109(b)(1)
Unique device identifier		Indicates a carrier that contains Unique Device Identifier information. Source: ISO 15223, 5.7.10
Do not re-use		Indicates a medical device that is intended for one single use only. Source: ISO 15223, 5.4.2
Batch code		Indicates the manufacturer's batch code so that the batch or lot can be identified. Source: ISO 15223, 5.1.5
Use-by date		Indicates the date after which the medical device is not to be used. Source: ISO 15223, 5.1.4

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Made in USA with Globally
Sourced Materials

2510 Conway Ave.
St. Paul, MN 55144 USA
1-800-228-3957 (USA Only)

3M Deutschland GmbH
Health Care Business
Carl-Schurz-Str. 1
41453 Neuss, Germany

3M, Bair Hugger and the Bair Hugger logo are
trademarks of 3M. Used under license in Canada.
© 2020, 3M. Unauthorized use prohibited.
All rights reserved.

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Match Colors:

Scale:  1 Inch

This artwork has been created as requested by 3M. 3M is responsible for the artwork AS APPROVED and assumes full responsibility for its correctness.

PROOF REQUIRED

Please refer to 3M Corporate General Spec. 205 for details