


Warming Blanket


Full Body
REF 30000
84 × 36 in
(213 × 91 cm)


Pediatric
REF 31000
60 × 36 in
(152 × 91 cm)


Multi-Access
REF 31500
84 × 36 in
(213 × 91 cm)


Lower Body
REF 52500
60 × 36 in
(152 × 91 cm)


Small Lower Body
REF 53700
35 × 24 in
(89 × 61 cm)

Dual Port Torso
REF 54200
42 × 36 in
(107 × 91 cm)

Surgical Access
REF 57000
84 × 36 in
(213 × 91 cm)

Full Body Surgical
REF 61000
72 × 36 in
(183 × 91 cm)

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

3M Deutschland GmbH
Health Care Business
Carl-Schurz-Str. 1
41453 Neuss, Germany

Internal Use

Internal Use

9300 Winnetka Ave. N Brooklyn Park, MN 55445 (763) 488-5700			CUST ID#	PRINT METHOD	SGS#	ROUND		
					6444659-2	1		
Production Art	Description: 3M™ Bair Hugger™ Combined Over-Body Warming Blankets - IFU	Item Spec#: 34-8726-6465-0	DATE				CHANGE LOG	INT
	Base File Name: PA: 34872664650.indd	Supersedes#:	1.04.22	New Build		MY		
	GA:	Structure#:	05.04.22	Final release		MY		
	Software: InDesign CC 2020	Cat/Product#: multiple						
	SGS Contact:	Reference:						
	Email:	Requester: Justin Haas						
	Printer:	Die # / Doc. Size: 12" x 11" Flat, 6" x 11" Folded						
	Printer Location:							
	Supplier:							
	INKS	PROCESS BLACK	MATCH	SPECIAL				
LIMITATION OF LIABILITY – The supplied materials should be thoroughly inspected immediately upon receipt by the recipient for all content, placement and copy accuracy. It is also the responsibility of the recipient printer to examine and proofread all supplied materials before producing image carriers and printing final product. Should the supplied materials be found to be incorrect please contact an SGS representative immediately. The liability of SGS shall be limited to correcting its own product. SGS will not be responsible for any loss or damage including but not limited to the cost of product, print, image carriers, press delays, or other ancillary cost. UPC NOTE – Because of the differences in scanning equipment, SGS cannot guarantee the scanability of UPC barcodes. The Uniform Code Council recommends a UPC symbol be printed black on white at 100% with no truncation. All UPC codes are test scanned and verified before leaving our facility.								

Indications for Use

The Bair Hugger family of temperature management systems are indicated for hypothermic patients or normothermic patients for whom induced hypothermia or localized temperature therapy is clinically indicated. In addition, the temperature management systems can be used to provide patient thermal comfort when conditions exist that may cause patients to become too warm or too cold. The temperature management systems can be used with adult and pediatric patients.

The Bair Hugger temperature management system should only be used by trained medical professionals.

Patient Population and Settings

Adult and pediatric patients being treated in operating rooms, emergency departments and other departments in the hospital setting where patient temperature management is required.

Contraindications, Warnings, and Cautions

WARNING: Indicates a hazardous situation which, if not avoided, could result in death or serious injury

CAUTION: Indicates a hazardous situation which, if not avoided, could result in minor or moderate injury

CONTRAINDICATION: To reduce the risk of thermal injury:

- Do not apply heat to lower extremities during aortic cross-clamping. Thermal injury may occur if heat is applied to ischemic limbs.

WARNING: To reduce the potential for thermal injury:

- Do not treat patients with the Bair Hugger warming unit hose alone. Always attach the hose to a Bair Hugger warming blanket before providing warming therapy.
- Do not allow the patient to lie on the warming unit hose.
- Do not allow the warming unit hose to directly contact the patient's skin during warming therapy.
- Do not leave neonates, infants, children and other vulnerable patient populations unattended during warming therapy.
- Do not leave patients with poor perfusion unmonitored during prolonged warming therapy.
- Do not place the non-perforated side of the warming blanket on the patient. Always place the perforated side (with the small holes) directly on top of the patient in contact with the patient's skin.
- In the operating room, do not use this warming blanket with any device other than a Bair Hugger 500 series, 700 series, or 675 warming unit.
- Do not use a Bair Hugger 200 series warming unit in the operating room.
- Do not use a Bair Hugger 800 series patient adjustable warming unit with any Bair Hugger warming blanket.
- Do not continue warming therapy if the red Over-temp indicator light illuminates and the alarm sounds. Unplug the warming unit and contact a qualified service technician.
- Do not place patient securement device (i.e. safety strap or tape) over the warming blanket.
- Do not place the warming blanket directly over a dispersive electrode pad.

WARNING: To reduce the risk of patient injury or death due to altered drug delivery:

- Do not use a warming blanket over transdermal medication patches.

WARNING: To reduce the potential for injury due to interference with ventilation:

- Do not allow the warming blanket or head drape to cover the patient's head or airway when the patient is not mechanically vented.

WARNING: To reduce the potential for injury due to patient falls:

- Do not use a warming blanket to transfer or move the patient.

CAUTION: To reduce the risk of cross-contamination:

- This warming blanket is not sterile and is intended for single patient use ONLY. Placing a sheet between the warming blanket and the patient does not prevent contamination of the product.

CAUTION: To reduce the risk of fire:

- This product is classified as Class I Normal Flammability as defined by the Consumer product Safety Commission's flammable fabric regulation, 16 CFR 1610. Follow standard safety protocols when using high intensity heat sources.

CAUTION: To reduce the risk of thermal injury:

- Do not use if primary packaging has been previously opened or is damaged

CAUTION: To reduce the risk of thermal injury, hyperthermia, or hypothermia:

- 3M recommends continuously monitoring core temperature. In the absence of continuous monitoring, monitor the temperature of patients who are incapable of reacting, communicating and/or who cannot sense temperature a minimum of every 15 minutes or according to institutional protocol.
- Monitor cutaneous responses of patients who are incapable of reacting, communicating and/or who cannot sense temperature a minimum of every 15 minutes or according to institutional protocol.


- Adjust air temperature or discontinue therapy when the therapeutic goal is reached, if elevated temperatures are recorded or if there is an adverse cutaneous response under the warming blanket.

Intraoperative Warming


Follow your institution's protocols and adhere to guidelines for sterile technique, including good surgical draping practice.

Directions for Use

- Place the perforated side of the 3M™ Bair Hugger™ warming blanket (the side with small holes) directly on top of the patient in contact with the patient's skin (see Figure A).

- Where applicable, to help secure the warming blanket, remove the backing from the adhesive tape strip and adhere the warming blanket to the patient (see Figure B). The warming blanket should be placed with the adhesive strip towards the surgical site.


Warning: Do not place patient securement device (i.e. safety strap or tape) over the warming blanket.

Warning: Do not place the warming blanket directly over a dispersive electrode pad.

Warning: do not treat patients with the Bair Hugger warming unit hose alone. Always attach the hose to a Bair Hugger warming blanket before providing warming therapy.

- Insert the end of the Bair Hugger warming unit hose in the hose port (see Figure C). Use a twisting motion to ensure a snug fit. A visual marker is located around the mid-section of the hose end to guide the depth of hose insertion. Support hose to ensure secure attachment.

NOTE: See special considerations for Bair Hugger blankets shown below.

- Select the desired temperature setting on the warming unit to initiate warming therapy. (See the Operator's Manual for your specific warming unit model)


Caution: Patient Monitoring Recommendations:

- 3M recommends continuously monitoring core temperature. In the absence of continuous monitoring, monitor the temperature of patients who are incapable of reacting, communicating and/or who cannot sense temperature a minimum of every 15 minutes or according to institutional protocol.
 - Monitor cutaneous responses of patients who are incapable of reacting, communicating and/or who cannot sense temperature a minimum of every 15 minutes or according to institutional protocol.
 - Adjust air temperature or discontinue therapy when the therapeutic goal is reached, if elevated temperatures are recorded or if there is an adverse cutaneous response under the warming blanket.
- Based on the warming unit model utilized, turn the unit off or to standby mode to discontinue warming therapy. Disconnect the warming unit hose from the blanket and discard the warming blanket per hospital policy.

Special Considerations:

Dual Port Torso 54200 and Surgical Access 57000 Warming Blankets


Two hose ports are provided for clinician preference. Place the removable hose port card in the hose port that is not being used during warming therapy (see Figure D).

Head Drape for Dual Port Torso 54200 and Surgical Access 57000 Warming Blankets


If the patient is intubated and ventilated, lay the head drape over the patient's head and neck (see figure e) otherwise, tuck the drape between the channels of the warming blanket, away from the patient's head.

Warning: Do not allow the warming blanket or head drape to cover the patient's head or airway when the patient is not mechanically ventilated.

Multi-Access 31500 and Full Body 61000 Warming Blankets

To use the access panels, tear the uncut tab at the edge of the warming blanket. Remove the backing from the tape strip on the center of the blanket. Fold the access panel back and press against the exposed tape. Pull the panel away from the tape to release (see Figure F).


Please report serious incident occurring in relation to 3M and the local competent authority (EU) or local regulatory authority.

Symbol Glossary

Symbol Title	Symbol	Description and Reference
Manufacturer		Indicates the medical device manufacturer. Source: ISO 15223, 5.1.1
Authorized Representative in European Community / European Union		Indicates the authorized representative in the European Community / European Union. Source: ISO 15223, 5.1.2, 2014/35/EU, and/or 2014/30/EU
Date of manufacture		Indicates the date when the medical device was manufactured. Source: ISO 15223, 5.1.3
Batch code		Indicates the manufacturer's batch code so that the batch or lot can be identified. Source : ISO 15223, 5.1.5
Catalogue number		Indicates the manufacturer's catalogue number so that the medical device can be identified. Source : ISO 15223, 5.1.6
Non Sterile		Indicates a medical device that has not been subjected to a sterilization process. Source: ISO 15223, 5.2.7
Do not use if package is damaged and consult instructions for use		Indicates a medical device that should not be used if the package has been damaged or opened and that the user should consult the instructions for use for additional information. Source: ISO 15223, 5.2.8
Do not re-use		Indicates a medical device that is intended for one single use only. Source: ISO 15223, 5.4.2
Consult instructions for use or consult electronic instructions for use		Indicates the need for the user to consult the instructions for use. Source: ISO 15223, 5.4.3
Natural rubber latex is not present		Indicates natural rubber or dry natural rubber latex is not present as a material of construction within the medical device or the packaging of a medical device. Source: ISO 15223, 5.4.5 and Annex B
Medical Device		Indicates the item is a medical device. Source: ISO 15223, 5.7.7
Unique device identifier		Indicates a carrier that contains Unique Device Identifier information. Source: ISO 15223, 5.7.10
Importer		Indicates the entity importing the medical device into the locale. Source: ISO 15223, 5.1.8
UKCA Mark		Indicates conformity to all applicable regulations and/or directives in the United Kingdom (UK), for products placed on the market in Great Britain (GB).
Swiss Authorised Representative		Indicates the authorised representative in Switzerland. Source: Swissmedic.ch
CE Mark		Indicates conformity to all applicable European Union Medical Device Regulations and Directives.
Rx Only		Indicates that U.S. Federal Law restricts this device to sale by or on the order of a healthcare professional. 21 Code of Federal Regulations (CFR) sec. 801.109(b)(1)
Green Dot		Indicates a financial contribution to national packaging recovery company per European Directive No. 94/62 and corresponding national law. Packaging Recovery Organization Europe.








For more information see, HCBGregulatory.3M.com

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3M, Bair Hugger and the Bair Hugger logo are trademarks of 3M.
© 2021, 3M. Unauthorized use prohibited. All rights reserved.