
# 3M Bair Hugger™

(en) Warming Blanket	(sk) Ohrievacia deka
(fr) Couverture de réchauffement	(sl) Grelna odeja
(de) Wärmedecke	(et) Soojendav tekk
(it) Coperta termica	(lv) Sasilšanas sega
(es) Manta térmica	(lt) Atšilimo antklodė
(nl) Verwarmingsdeken	(ro) Pătură încălzitoare
(sv) Värmetäcke	(ru) Обогревающее одеяло
(da) Varmetæppe	(hr) Pokrivač za zagrijavanje
(no) Varmeteppe	(bg) Затоплящо одеяло
(fi) Lämpöpeite	(sr) Deka za grejanje
(pt) Manta de aquecimento	(tr) Isıtma Battanisi
(el) Κουβέρτα θέρμανσης	(zh) 加温毯
(pl) Koc ogrzewający	(ar) بطانية تدفئة
(hu) Melegítőtakaró	(sq) Batanije për ngrohje
(cs) Vyhřívání přikrývka	(mk) ќебе за затоплување

## (en) Large Pediatric Underbody Warming Blanket 55000 Pediatric Underbody Warming Blanket 55501

### Indications for Use

The Bair Hugger family of temperature management systems are indicated for hypothermic patients or normothermic patients for whom induced hypothermia or localized temperature therapy is clinically indicated. In addition, the temperature management systems can be used to provide patient thermal comfort when conditions exist that may cause patients to become too warm or too cold. The temperature management systems can be used with adult and pediatric patients.

The Bair Hugger temperature management system should only be used by trained medical professionals.

### Patient Population and Settings

Adult and pediatric patients being treated in operating rooms, emergency departments and other departments in the hospital setting where patient temperature management is required.

### Contraindication: To reduce the risk of thermal injury:

- Do not apply heat to lower extremities during aortic cross-clamping. Thermal injury may occur if heat is applied to ischemic limbs.

**Warning:** indicates a hazard which, if not avoided, could result in death or serious injury.

### WARNING: To reduce the risk of thermal injury:

- Do not treat patients with the Bair Hugger warming unit hose alone. Always attach the hose to a Bair Hugger warming blanket before providing warming therapy.
- Do not allow the patient to lie on the warming unit hose.
- Do not allow the patient's arms or legs to rest on either hose port.
- Do not allow the warming unit hose to directly contact the patient's skin during warming therapy.
- Do not leave neonates, infants, children and other vulnerable patient populations unattended during warming therapy.
- Do not leave patients with poor perfusion unmonitored during prolonged warming therapy.
- Do not place the non-perforated side of the warming blanket on the patient. Always place the perforated side (with the small holes) directly in contact with the patient's skin.
- In the operating room, do not use this warming blanket with any device other than a Bair Hugger 500 series, 700 series, or 675 warming unit.
- Do not use a Bair Hugger 200 series warming unit in the operating room.
- Do not use a Bair Hugger 800 series patient adjustable warming unit with any Bair Hugger warming blanket.
- Do not continue warming therapy if the red Over-temp indicator light illuminates and the alarm sounds. Unplug the warming unit and contact a qualified service technician.
- If a securement device (i.e. safety strap, tape) is used, ensure the warming channels are not occluded.
- Do not place the warming blanket directly over a dispersive electrode pad.

**WARNING: To reduce the risk of patient injury or death due to altered drug delivery:**

- Do not use a warming blanket over transdermal medication patches.

**WARNING: To reduce the risk of injury due to interference with ventilation:**

- Do not allow the warming blanket or head drape to cover the patient's head or airway when the patient is not mechanically ventilated.

**WARNING: To reduce the risk of injury due to patient falls:**

- Do not use a warming blanket to transfer or move the patient.

**Caution:** Indicates a hazardous situation which, if not avoided, could result in minor or moderate injury.

**CAUTION: To reduce the risk of cross-contamination:**

- This warming blanket is not sterile and is intended for single patient use ONLY. Placing a sheet between the warming blanket and the patient does not prevent contamination of the product.

**CAUTION: To reduce the risk of fire:**

- This product is classified as Class I Normal Flammability as defined by the Consumer Product Safety Commission's flammable fabric regulation, 16 CFR 1610. Follow standard safety protocols when using high intensity heat sources.

**CAUTION: To reduce the risk of thermal injury:**

- Do not use if primary packaging has been previously opened or is damaged.

**CAUTION: To reduce the risk of thermal injury, hyperthermia or hypothermia:**

- 3M recommends continuously monitoring core temperature. In the absence of continuous monitoring, monitor the temperature of patients who are incapable of reacting, communicating and/or who cannot sense temperature a minimum of every 15 minutes or according to institutional protocol.
- Monitor cutaneous responses of patients who are incapable of reacting, communicating and/or who cannot sense temperature a minimum of every 15 minutes or according to institutional protocol.
- Adjust air temperature or discontinue warming therapy when the therapeutic goal is reached, if elevated temperatures are recorded or if there is an adverse cutaneous response in the warmed area.

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

EC REP 3M Deutschland GmbH  
Health Care Business  
Carl-Schurz-Str. 1  
41453 Neuss, Germany

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Internal Use


### Intraoperative Warming

Follow your institution's protocols and adhere to guidelines for sterile technique, including good surgical draping practice.

### Directions for Use

**Note:** In cases where fluid runoff is expected, place an absorbent cloth beneath the Bair Hugger warming blanket.

- Using the adhesive tape strips on the underside of the warming blanket, secure the blanket to the table to prevent it from slipping (Figures A and B). The side flaps (located on underside of 55000) can be tucked under the mattress or table to provide extra stability.


- Place the patient on top of the blanket while it is in the non-inflated state. Do not position patient's arms out across blanket; airflow may be restricted.

**WARNING:** If a securement device (i.e. safety strap, tape) is used, ensure the warming channels are not occluded.

- Apply the clear drapes as desired for the surgical procedure. The head drape should be used only if the patient is intubated and ventilated. Use adhesive tape strips to adhere drape(s) in place along the upper and/or lower edge of the prep site (Figures C and D).

**WARNING:** Do not allow the head drape to cover the patient's head or airway when the patient is not mechanically ventilated.


Large Pediatric Underbody Warming Blanket 55000

Pediatric Underbody Warming Blanket 55501

- Hose ports are provided at either end of the blanket for clinician preference (Figures E and F). A removable card is located in the lower hose port (Figure G). Place the card into the hose port that is not being used during warming.

**WARNING:** Do not allow the patient's arms or legs to rest on either hose port.


Model 55000  
Large Pediatric Underbody Warming Blanket 55000

Model 55501  
Pediatric Underbody Warming Blanket 55501

- Insert the end of the Bair Hugger warming unit hose in the hose port (Figure H). Use a twisting motion to ensure a snug fit. A visual marker is located around the mid-section of the hose end to guide the depth of hose insertion. Support hose to insure secure attachment.

**WARNING:** Do not treat patients with the Bair Hugger hose alone. Always attach the hose to a Bair Hugger blanket before providing warming therapy.

- Select the desired temperature setting on the warming unit to initiate warming therapy. (See the Operator Manual for your specific Warming Unit Model)

**CAUTION:** Patient Monitoring Recommendations:

- 3M recommends continuously monitoring core temperature. In the absence of continuous monitoring, monitor the temperature of patients who are incapable of reacting, communicating and/or who cannot sense temperature a minimum of every 15 minutes or according to institutional protocol.

- Monitor cutaneous responses of patients who are incapable of reacting, communicating and/or who cannot sense temperature a minimum of every 15 minutes or according to institutional protocol.

- Adjust air temperature or discontinue warming therapy when the therapeutic goal is reached, if elevated temperatures are recorded or if there is an adverse cutaneous response in the warmed area.

- Based on the warming unit model utilized, turn the unit off or to standby mode to discontinue warming therapy. Disconnect the hose from the warming blanket and discard the blanket per hospital policy.

**Please report a serious incident occurring in relation to the device to 3M and the local competent authority (EU) or local regulatory authority**

### Symbol Glossary

Symbol Title	Symbol	Description and Reference
Manufacturer		Indicates the medical device manufacturer. Source: ISO 15223, 5.1.1
Authorized Representative in European Community / European Union		Indicates the authorized representative in the European Community / European Union. Source: ISO 15223, 5.1.2, 2014/35/EU, and/or 2014/30/EU
Date of manufacture		Indicates the date when the medical device was manufactured. Source: ISO 15223, 5.1.3
Use-by date		Indicates the date after which the medical device is not to be used. Source: ISO 15223, 5.1.4
Batch code		Indicates the manufacturer's batch code so that the batch or lot can be identified. Source: ISO 15223, 5.1.5
Catalogue number		Indicates the manufacturer's catalogue number so that the medical device can be identified. Source: ISO 15223, 5.1.6
Non Sterile		Indicates a medical device that has not been subjected to a sterilization process. Source: ISO 15223, 5.2.7
Do not use if package is damaged and consult instructions for use		Indicates a medical device that should not be used if the package has been damaged or opened and that the user should consult the instructions for use for additional information. Source: ISO 15223, 5.2.8
Do not re-use		Indicates a medical device that is intended for one single use only. Source: ISO 15223, 5.4.2
Consult instructions for use or consult electronic instructions for use		Indicates the need for the user to consult the instructions for use. Source: ISO 15223, 5.4.3
Natural rubber latex is not present		Indicates natural rubber or dry natural rubber latex is not present as a material of construction within the medical device or the packaging of a medical device. Source: ISO 15223, 5.4.5 and Annex B
Medical Device		Indicates the item is a medical device. Source: ISO 15223, 5.7.7
Unique device identifier		Indicates a carrier that contains Unique Device Identifier information. Source: ISO 15223, 5.7.10
Importer		Indicates the entity importing the medical device into the locale. Source: ISO 15223, 5.1.8
UKCA Mark		Indicates conformity to all applicable regulations and/or directives in the United Kingdom (UK), for products placed on the market in Great Britain (GB).
Swiss Authorised Representative		Indicates the authorised representative in Switzerland. Source: Swissmedic.ch
CE Mark		Indicates conformity to all applicable European Union Medical Device Regulations and Directives.
Rx Only		Indicates that U.S. Federal Law restricts this device to sale by or on the order of a healthcare professional. 21 Code of Federal Regulations (CFR) sec. 801.109(b)(1)
Green Dot		Indicates a financial contribution to national packaging recovery company per European Directive No. 94/62 and corresponding national law. Packaging Recovery Organization Europe.

For more information see, [HCBGRegulatory.3M.com](http://HCBGRegulatory.3M.com)

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<p>9300 Winnetka Ave. N Brooklyn Park, MN 55445 (763) 488-5700</p>		CUST ID#	PRINT METHOD	SGS#	ROUND
				6448277-2	2
<p><b>Production Art</b></p> <p>Description: 3M™ Bair Hugger™ Pediatric Underbody Warming Blanket - IFU</p> <p>Base File Name: PA: 34872664684.indd</p> <p>GA: Software: InDesign CC 2020</p> <p>SGS Contact: Email: Printer: Printer Location: Supplier:</p>	<p>Item Spec#: 34-8726-6468-4</p> <p>Supersedes#:</p> <p>Structure#:</p> <p>Cat/Product#: 54500</p> <p>Reference:</p> <p>Requester: Justin Haas</p> <p>Die # / Doc. Size: 12" x 11" Flat, 6" x 11" Folded</p>	<p>DATE CHANGE LOG INT</p> <p>04.03.22 Update as per annotation SM</p> <p>04.12.22 Final release MY</p>			
		<p>INKS PROCESS BLACK</p> <p>MATCH</p> <p>SPECIAL</p>			
<p><small>LIMITATION OF LIABILITY - The supplied materials should be thoroughly inspected immediately upon receipt by the recipient for all content, placement and copy accuracy. It is also the responsibility of the recipient printer to examine and proofread all supplied materials before producing image carriers and printing final product. Should the supplied materials be found to be incorrect please contact an SGS representative immediately. The liability of SGS shall be limited to correcting its own product. SGS will not be responsible for any loss or damage including but not limited to the cost of product, print, image carriers, press delays, or other ancillary cost.</small></p> <p><small>UPC NOTE - Because of the differences in scanning equipment, SGS cannot guarantee the scanability of UPC barcodes. The Uniform Code Council recommends a UPC symbol be printed black on white at 100% with no function. All UPC codes are first scanned and verified before leaving our facility.</small></p>					

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3M, Bair Hugger and the Bair Hugger logo are trademarks of 3M.
All rights reserved.