
3M

Bair Hugger™

- en

Warming Blanket
- fr

Couverture de réchauffement
- de

Wärmedecke
- it

Coperta riscaldante
- es

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Verwarmingsdeken
- sv

Värmetäcke
- da

Varmetæppe
- no

Varmeteppe
- fi

Lämpöpeite
- pt

Manta de aquecimento
- el

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Lithotomy Underbody Warming Blanket 58501

Indications for Use

The Bair Hugger family of temperature management systems are indicated for hypothermic patients or normothermic patients for whom induced hypothermia or localized temperature therapy is clinically indicated. In addition, the temperature management systems can be used to provide patient thermal comfort when conditions exist that may cause patients to become too warm or too cold. The temperature management systems can be used with adult and pediatric patients.

The Bair Hugger temperature management system should only be used by trained medical professionals.

Patient Population and Settings

Adult and pediatric patients being treated in operating rooms, emergency departments and other departments in the hospital setting where patient temperature management is required.

CONTRAINDICATION: To reduce the risk of thermal injury:

- Do not apply heat to lower extremities during aortic cross-clamping. Thermal injury may occur if heat is applied to ischemic limbs.


Warnings: Indicates a hazardous situation which, if not avoided, could result in death or serious injury.

WARNING: To reduce the risk of thermal injury:


- Do not treat patients with the Bair Hugger warming unit hose alone. Always attach the hose to a Bair Hugger warming blanket before providing warming therapy.
- Do not allow the patient to lie on the warming unit hose.
- Do not allow the warming unit hose to directly contact the patient's skin during warming therapy.
- Do not leave neonates, infants, children and other vulnerable patient populations unattended during warming therapy.
- Do not leave patients with poor perfusion unmonitored during prolonged warming therapy.
- Do not place the non-perforated side of the warming blanket on the patient. Always place the perforated side (with the small holes) directly in contact with the patient's skin.
- In the operating room, do not use this warming blanket with any device other than a Bair Hugger 500 series, 700 series, or 675 warming unit.
- Do not use a Bair Hugger 200 series warming unit in the operating room
- Do not use a Bair Hugger 800 series patient adjustable warming unit with any Bair Hugger warming blanket.
- Do not continue warming therapy if the red Over-temp indicator light illuminates and the alarm sounds. Unplug the warming unit and contact a qualified service technician.
- If a securement device (i.e. safety strap, tape) is used, ensure the warming channels are not occluded.
- Do not place the warming blanket directly over a dispersive electrode pad.

WARNING: To reduce the risk of patient injury or death due to altered drug delivery:

- Do not use a warming blanket over transdermal medication patches.

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
3M Company
2510 Conway Ave.
St. Paul, MN 55144 USA
1-800-228-3957 (USA Only)

3M, Bair Hugger and the Bair Hugger logo are trademarks of 3M.
Used under license in Canada.
© 2021, 3M. Unauthorized use prohibited.
All rights reserved.

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Issue Date: 2021-10, rev 4
Date d'émission: 2021-10, rev. 4
34-8726-6469-2

Note: In cases where fluid runoff is expected, place an absorbent cloth beneath the model 58501 warming blanket.


-
- A schematic diagram labeled 'A' showing a cross-section of a device. It features a central cavity with a rounded bottom and two side chambers. The top of the device has a flange with a central rectangular slot. The bottom of the device has a wide, flat base. The central cavity contains several small, dark, oval-shaped components. The side chambers also contain similar components. The entire device is shown in a symmetrical, cross-sectional view.

-
- B**

[illegible]

-

-

-

-

Symbol Title	Symbol	Description and Reference
Manufacturer		Indicates the medical device manufacturer. Source: ISO 15223, 5.1.1
Authorized Representative in European Community / European Union		Indicates the authorized representative in the European Community / European Union. Source: ISO 15223, 5.1.2, 2014/35/EU, and/or 2014/30/EU
Date of manufacture		Indicates the date when the medical device was manufactured. Source: ISO 15223, 5.1.3
Use-by date		Indicates the date after which the medical device is not to be used. Source: ISO 15223, 5.1.4
Batch code		Indicates the manufacturer's batch code so that the batch or lot can be identified. Source : ISO 15223, 5.1.5
Catalogue number		Indicates the manufacturer's catalogue number so that the medical device can be identified. Source : ISO 15223, 5.1.6

Unique device identifier		Indicates a carrier that contains Unique Device Identifier information. Source: ISO 15223, 5.7.10
Importer		Indicates the entity importing the medical device into the locale. Source: ISO 15223, 5.1.8
UKCA Mark		Indicates conformity to all applicable regulations and/or directives in the United Kingdom (UK), for products placed on the market in Great Britain (GB).
Swiss Authorised Representative		Indicates the authorised representative in Switzerland. Source: Swissmedic.ch
CE Mark		Indicates conformity to all applicable European Union Medical Device Regulations and Directives.
Rx Only		Indicates that U.S. Federal Law restricts this device to sale by or on the order of a healthcare professional. 21 Code of Federal Regulations (CFR) sec. 801.109(b)(1)
Green Dot		Indicates a financial contribution to national packaging recovery company per European Directive No. 94/62 and corresponding national law. Packaging Recovery Organization Europe.

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Production Art	SGS 1300 Winnetka Ave. N Brooklyn Park, MN 55445 (763) 488-5700	3M	CUST ID #		PRINT METHOD		SGS #	ROUND	
							6448363-3	2	
			Item Spec#: 34-8726-6469-2		DATE		CHANGE LOG		INT
			Supersedes#:		04.02.22 Update as per amend				GA
			Structure#:		04.12.22 Final release				MY
			Cat/Product#: 58501						
			Reference:						
			Requester: Justin Haas						
			Die # / Doc. Size: 12" x 11" Flat, 6" x 11" Folded						
INKS	PROCESS BLACK	MATCH						SPECIAL	
Limitation of Liability: The supplied materials should be thoroughly inspected immediately upon receipt by the recipient for all content, placement and copy accuracy. It is also the responsibility of the recipient printer to examine and proofread all supplied materials before producing image carriers and printing final product. The liability of SGS shall be limited to the cost of the product. SGS will not be responsible for any loss or damage including but not limited to the cost of product, print image carriers, press delays, or other ancillary costs. UPC NOTE: Because of the differences in scanning equipment, SGS cannot guarantee the scannability of UPC barcodes. The Uniform Code Council recommends a UPC symbol be printed black on white at 100% with no truncation. All UPC codes are test scanned and verified before leaving our facility.									

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