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
Spinal Underbody Warming Blanket 57501

Indications for Use

The Bair Hugger family of temperature management systems are indicated for hypothermic patients or normothermic patients for whom induced hypothermia or localized temperature therapy is clinically indicated. In addition, the temperature management systems can be used to provide patient thermal comfort when conditions exist that may cause patients to become too warm or too cold. The temperature management systems can be used with adult and pediatric patients.

The Bair Hugger temperature management system should only be used by trained medical professionals.

Patient Population and Settings

Adult and pediatric patients being treated in operating rooms, emergency departments and other departments in the hospital setting where patient temperature management is required.


Contraindication: To reduce the risk of thermal injury:

- Do not apply heat to lower extremities during aortic cross-clamping. Thermal injury may occur if heat is applied to ischemic limbs.


Warning: Indicates a hazardous situation which, if not avoided, could result in death or serious injury.

WARNING: To reduce the risk of thermal injury:

- Do not treat patients with the Bair Hugger warming unit hose alone. Always attach the hose to a Bair Hugger warming blanket before providing warming therapy.
- Do not allow the patient to lie on the warming unit hose.
- Do not allow the warming unit hose to directly contact the patient's skin during warming therapy.
- Do not leave neonates, infants, children and other vulnerable patient populations unattended during warming therapy.
- Do not leave patients with poor perfusion unmonitored during prolonged warming therapy.
- Do not place the non-perforated side of the warming blanket on the patient. Always place the perforated side (with the small holes) directly on top of the patient in contact with the patient's skin.
- In the operating room, do not use this warming blanket with any device other than a Bair Hugger 500 series, 700 series, or 675 warming unit.
- Do not use a Bair Hugger 200 series warming unit in the operating room
- Do not use a Bair Hugger 800 series patient adjustable warming unit with any Bair Hugger warming blanket.
- Do not continue warming therapy if the red Over-temp indicator light illuminates and the alarm sounds. Unplug the warming unit and contact a qualified service technician.
- If a securement device (i.e. safety strap, tape) is used, ensure the warming channels are not occluded.
- Do not place the warming blanket directly over a dispersive electrode pad.

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

Issue Date: 2021-10, rev. 4
Date d'émission : 2021-10, rev. 4
34-8726-6470-0

3M Deutschland GmbH
Health Care Business
Carl-Schurz-Str. 1
41453 Neuss, Germany


1. Open and unfold warming blanket to locate the white product label on the blanket. Unfold the blanket so the product label faces the patient.
2. Position the warming blanket onto the surface of the patient footboard or onto the surface of the sling assembly. Using the adhesive tape strip on the underside of the blanket, secure the blanket to the table to prevent it from slipping (Figure A).


3. Pass the warming blanket between the patient footboard (➔) and the patient top thigh set (●) through the frame to the underside of the frame of the spinal surgery table (Figure B). Unfold the blanket under the top thigh pad set (★), top hip pad set (●), patient chest pad (▲), and under the patient head support plate (■).
4. Draw the end of the warming blanket to the end of the spinal surgery table and secure it to the end of the table frame using the tie strips located at the end of the blanket (Figure B).


5. Open the four perforations on both sides of the warming blanket and position the blanket around the edges of the spinal surgery frame (Figure C).

- Before the patient has been positioned, draw tie strips across the top of the frame of the spinal surgery table and around each side of the patient chest pad (Figure D, Option 1) or directly to the frame (Figure D, Option 2).


7. Pillows and other positioning modalities used for the lower extremities should be positioned under the blanket.
8. **Optional:** Tear out the circular (face) section at the head of the warming blanket to visualize the patient's eyes and airway from under the table (Figure E).

9. Insert the end of the Bair Hugger warming unit hose into the hose port (Figures F and G). Use a twisting motion to ensure a snug fit. A visual marker is located around the mid-section of the hose end to guide the depth of hose insertion. Support hose to insure secure attachment.

Warning: Do not treat patients with the Bair Hugger hose alone. Always attach the hose to a Bair Hugger blanket before providing therapy.


10. Select the desired temperature setting on the warming unit to initiate warming therapy. (See the Operator Manual for your specific Warming Unit Model)

- 3M recommends continuously monitoring core temperature. In the absence of continuous monitoring, monitor the temperature of patients who are incapable of reacting, communicating and/or who cannot sense temperature a minimum of every 15 minutes or according to institutional protocol.
- Monitor cutaneous responses of patients who are incapable of reacting, communicating and/or who cannot sense temperature a minimum of every 15 minutes or according to institutional protocol.
- Adjust air temperature or discontinue therapy when the therapeutic goal is reached, if elevated temperatures are recorded or if there is an adverse cutaneous response in the warmed area.

11. Based on the warming unit model utilized, turn the unit off or to standby mode to discontinue warming therapy. Disconnect the warming unit hose from the warming blanket and discard the blanket per hospital policy.

Please report a serious incident occurring in relation to the device to 3M and the local competent authority (EU) or local regulatory authority.

Symbol Title	Symbol	Description and Reference
Manufacturer		Indicates the medical device manufacturer. Source: ISO 15223, 5.1.1
Authorized Representative in European Community / European Union		Indicates the authorized representative in the European Community / European Union. Source: ISO 15223, 5.1.2, 2014/35/EU, and/or 2014/30/EU
Date of manufacture		Indicates the date when the medical device was manufactured. Source: ISO 15223, 5.1.3
Use-by date		Indicates the date after which the medical device is not to be used. Source: ISO 15223, 5.1.4
Batch code		Indicates the manufacturer's batch code so that the batch or lot can be identified. Source: ISO 15223, 5.1.5
Catalogue number		Indicates the manufacturer's catalogue number so that the medical device can be identified. Source: ISO 15223, 5.1.6
Non Sterile		Indicates a medical device that has not been subjected to a sterilization process. Source: ISO 15223, 5.2.7
Do not use if package is damaged and consult instructions for use		Indicates a medical device that should not be used if the package has been damaged or opened and that the user should consult the instructions for use for additional information. Source: ISO 15223, 5.2.8
Do not re-use		Indicates a medical device that is intended for one single use only. Source: ISO 15223, 5.4.2
Consult instructions for use or consult electronic instructions for use		Indicates the need for the user to consult the instructions for use. Source: ISO 15223, 5.4.3

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Swiss Authorised Representative		Indicates the authorised representative in Switzerland. Source: Swissmedic.ch
CE Mark		Indicates conformity to all applicable European Union Medical Device Regulations and Directives.
Rx Only		Indicates that U.S. Federal Law restricts this device to sale by or on the order of a healthcare professional. 21 Code of Federal Regulations (CFR) sec. 801.109(b)(1)
Green Dot		Indicates a financial contribution to national packaging recovery company per European Directive No. 94/62 and corresponding national law. Packaging Recovery Organization Europe.

For more information see, HCBGRegulatory.3M.com

 SGS 1300 Winnetka Ave. N Brooklyn Park, MN 55445 (763) 488-5700		CUST ID #		PRINT METHOD		SGS #	ROUND				
						6448345-2	2				
Production Art	Description: 3M™ Bair Hugger™ Spinal Underbody Warming Blanket - IFU		Item Spec#: 34-8726-6470-0		DATE			CHANGE LOG		INT	
	Base File Name: PA: 34872664700.indd		Supersedes#:		04.03.22 Update as per annotations 04.07.22 Final File Release					SM GA	
	GA:		Structure#:								
	Software: InDesign CC 2020		Cat/Product#: 57501								
	SGS Contact:		Reference:								
	Email:		Requester: Justin Haas								
	Printer:		Die # / Doc. Size: 12" x 11" Flat, 6" x 11" Folded								
	Printer Location:										
	Supplier:										
INKS PROCESS BLACK	MATCH				SPECIAL						
LIMITATION OF LIABILITY: The supplied materials should be thoroughly inspected immediately upon receipt by the recipient for all content, placement and copy accuracy. It is also the responsibility of the recipient printer to examine and proofread all supplied materials before producing image carriers and printing final product. The recipient printer is bound to be incorrect please contact an SGS representative immediately. The liability of SGS shall be limited to the cost of the product, print image carriers, press delays, or other ancillary costs. UPC NOTE: Because of the differences in scanning equipment, SGS cannot guarantee the scannability of UPC symbols. The Uniform Code Council recommends a UPC symbol be printed back on white at 100% with no truncation. All UPC codes are test scanned and verified before leaving our facility.											

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3M, Bair Hugger and the Bair Hugger logo are trademarks of 3M.
All rights reserved.