
3M Bair Hugger™

- en Warming Unit Model 675 Operator's Manual
- fr Manuel d'utilisation de l'unité de réchauffement modèle 675
- de Benutzerhandbuch zur Wärmeeinheit, Modell 675
- it Unità riscaldante modello 675 - Manuale d'uso
- es Manual del operador de la unidad de calentamiento modelo 675
- pt Manual do Operador da Unidade de Aquecimento Modelo 675
- el Εγχειρίδιο Χρήστη Μονάδας Θέρμανσης, Μοντέλο 675
- pl Aparat do terapii grzewczej Model 675 — instrukcja obsługi
- cs Ohřívací jednotka, model 675 – Návod k obsluze
- lt 675 modelio šildymo įrenginio operatoriaus vadovas
- ru Руководство по эксплуатации устройства для обогрева, модель 675
- tr Isıtma Ünitesi Model 675 Kullanım Kılavuzu
- zh 675 型加温器操作手册
- ar دليل المشغل للطراز 675 لوحدة التدفئة

MD CE 2797

Made in USA with Globally Sourced Materials

3M Company
2510 Conway Ave.,
St. Paul, MN 55144 USA
1-800-228-3957 (USA Only)

3M, BAIR HUGGER and BAIR HUGGER logo are trademarks of 3M. Used under license in Canada.
© 2021, 3M. Unauthorized use prohibited. All rights reserved.

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
Issue Date: 2021-10, rev. 8
34-8727-8204-9

EC REP 3M Deutschland GmbH
Health Care Business
Carl-Schurz-Str. 1
41453 Neuss, Germany

3M Bair Hugger™


- en English 1
- fr Français 15
- de Deutsch 29
- it Italiano 45
- es Español 59
- pt Português 73
- el Ελληνικά 87
- pl Polski 103
- cs Čeština 117
- lt Lietuviškai 131
- ru РУССКИЙ 145
- tr Türkçe 161
- zh 简体中文 175
- ar العربية 187

<p>9300 Winnetka Ave. N Brooklyn Park, MN 55445 (763) 488-5700</p>		CUST ID#	PRINT METHOD	SGS#	ROUND																			
				6595581-1	6																			
<p>Description: Bair Hugger Warming Unit - IFU</p> <p>Base File Name: PA: 34872782049.indd GA: Software: InDesign CC SGS Contact: Email: Printer: Printer Location: Supplier:</p>	<p>Item Spec#: 34-8727-8204-9</p> <p>Supersedes#:</p> <p>Structure#:</p> <p>Cat/Product#: MD675</p> <p>Reference:</p> <p>Requester: Frank Boeshart</p> <p>Die # / Doc. Size: 8.5" x 11"</p>	CHANGE LOG																						
		<table border="1" style="width: 100%; border-collapse: collapse;"> <thead> <tr> <th style="width: 10%;">DATE</th> <th style="width: 80%;">UPDATE</th> <th style="width: 10%;">INT</th> </tr> </thead> <tbody> <tr> <td>09.11.21</td> <td>Update as per annotations & symbol chart</td> <td>SM</td> </tr> <tr> <td>29.11.21</td> <td>Update as per Brief</td> <td>SR</td> </tr> <tr> <td>07.12.21</td> <td>Update as per annotation</td> <td>SM</td> </tr> <tr> <td>09.12.21</td> <td>Rx content only update for all languages</td> <td>SM</td> </tr> <tr> <td>12.24.21</td> <td>Update as per annotation</td> <td>GA</td> </tr> <tr> <td>12.28.21</td> <td>Final Release</td> <td>kat</td> </tr> </tbody> </table>				DATE	UPDATE	INT	09.11.21	Update as per annotations & symbol chart	SM	29.11.21	Update as per Brief	SR	07.12.21	Update as per annotation	SM	09.12.21	Rx content only update for all languages	SM	12.24.21	Update as per annotation	GA	12.28.21
DATE	UPDATE	INT																						
09.11.21	Update as per annotations & symbol chart	SM																						
29.11.21	Update as per Brief	SR																						
07.12.21	Update as per annotation	SM																						
09.12.21	Rx content only update for all languages	SM																						
12.24.21	Update as per annotation	GA																						
12.28.21	Final Release	kat																						
<p>INKS</p> <table style="display: inline-table; border: none;"> <tr> <td style="width: 20px; height: 20px; background-color: black; margin-right: 5px;"></td> <td style="font-size: 8px; margin-right: 5px;">PROCESS BLACK</td> <td style="width: 20px; height: 20px; background-color: blue; margin-right: 5px;"></td> <td style="font-size: 8px; margin-right: 5px;">PANTONE 294C</td> <td style="width: 20px; height: 20px; background-color: red; margin-right: 5px;"></td> <td style="font-size: 8px;">3M RED</td> </tr> </table>			PROCESS BLACK		PANTONE 294C		3M RED	MATCH	SPECIAL															
	PROCESS BLACK		PANTONE 294C		3M RED																			
<p>LIMITATION OF LIABILITY – The supplied materials should be thoroughly inspected immediately upon receipt by the recipient for all content, placement and copy accuracy. It is also the responsibility of the recipient printer to examine and proofread all supplied materials before producing image carriers and printing final product. Should the supplied materials be found to be incorrect please contact an SGS representative immediately. The liability of SGS shall be limited to correcting its own product. SGS will not be responsible for any loss or damage including but not limited to the cost of product, print, image carriers, press delays, or other ancillary cost.</p> <p>UPC NOTE – Because of the differences in scanning equipment, SGS cannot guarantee the scanability of UPC barcodes. The Uniform Code Council recommends a UPC symbol be printed black on white at 100% with no truncation. All UPC codes are test scanned and verified before leaving our facility.</p>																								

en **Table of Contents**

Product Description of the 3M™ Bair Hugger™ Temperature Management System	3
Indications for Use	3
Patient Population and Settings	3
Contraindications, Warnings, Cautions and Notices.....	3
Explanation of Signal Word Consequences.....	3
CONTRAINDICATIONS.....	3
WARNINGS	3
CAUTIONS.....	4
Notices	4
Using the Model 675 Bair Hugger Warming Unit	5
Overview of Model 675 Warming Unit.....	5
Overview of Control Panel	5
Important Alarm Symbols and Indicators.....	5
Filter time limit has been reached:.....	6
Instructions for Use.....	6
Caution:	6
Connecting the Warming Unit Hose and Power Cord	7
Connecting the Warming Unit Hose to the Warming Unit.....	7
Connecting the Power Cord to the Warming Unit.....	7
Warming Unit Mounting Instructions.....	7
Mounting on IV Pole	7
Mounting onto Rolling Cart Accessory	7
Cleaning, Storage, and Maintenance	7
Cleaning Instructions	7
Storage	7
Maintenance: Filter.....	7
Proper Use and Maintenance.....	8
Diagnostics	8
Read Before Servicing Equipment.....	8
In-Warranty Repair and Exchange.....	8
USA:	8
Outside of the USA:	8
Symbol Glossary.....	8
Specifications.....	10
Physical Characteristics.....	10
Recommended operating environment	10
Operating temperatures	10
Safety System	10
Alarm/Fault Codes.....	10
Electrical Characteristics.....	11
Classifications.....	13
Frequently Asked Questions (FAQ) on Filtration	14
What filter is used in the 3M™ Bair Hugger™ Model 675 warming unit?	14
What is a MERV rating?	14
What is the MERV rating recommended for operating rooms?.....	14
What is the MERV rating recommended for patient warming devices in the operating room?	14
Is a HEPA filter necessary?	14
What is the ASHRAE 52.2-2012 standard?	14
References:.....	14

Operator's Manual

Product Description of the 3M™ Bair Hugger™ Temperature Management System

The Bair Hugger temperature management system consists of a Model 675 warming unit (with optional rolling cart and other accessories) along with a 3M warming blanket or a Bair Hugger warming gown.

The Bair Hugger warming unit is attached to the warming blanket or warming gown by means of a flexible hose. Warm air is generated in the unit and flows through the warming unit hose into the blanket or gown. Depending on the model, the 3M blanket or gown is placed either over, around or underneath the patient. Small perforations on the blanket or gown allow the warm air to be dispersed over the patient's skin.

This manual includes operating instructions and warming unit specifications for the Model 675 warming unit. You can use the Model 675 warming unit in clinical settings where the patient may feel too warm or too cold; including the operating room to provide patient temperature management. For information about using 3M Bair Hugger warming blankets or gowns with the Model 675 warming unit, refer to the respective Instructions for Use.

Indications for Use

The Bair Hugger temperature management system is intended to prevent and treat hypothermia. In addition, the temperature management system can be used to provide patient thermal comfort when conditions exist that may cause patients to feel too warm or too cold. The temperature management system can be used with adult and pediatric patients.

The Bair Hugger temperature management system should only be used by trained medical professionals.

Patient Population and Settings

Adult and pediatric patients being treated in operating rooms, emergency departments and other departments in the hospital setting where patient temperature management is required.

Contraindications, Warnings, Cautions and Notices

Explanation of Signal Word Consequences

WARNING: Indicates a hazardous situation which, if not avoided, could result in death or serious injury.

CAUTION: Indicates a hazardous situation which, if not avoided, could result in minor or moderate injury.

NOTICE: Indicates a situation which, if not avoided, could result in property damage only.

CONTRAINDICATIONS

CONTRAINDICATION: To reduce the risk of thermal injury:

- Do not apply heat to lower extremities during aortic cross-clamping. Thermal injury may occur if heat is applied to ischemic limbs.

WARNINGS

WARNING: To reduce the risk of thermal injury:

- The Bair Hugger warming unit Model 675 has been designed to operate safely ONLY with 3M disposable warming blankets and gowns. Use with other products may cause thermal injury. To the full extent permitted by law, the manufacturer and/or importer declines all responsibility for thermal injury resulting from the warming unit being used in conjunction with non-3M warming blankets and gowns.
- Do not treat patients with the warming unit hose alone. Always attach the hose to a 3M warming blanket/gown before providing warming therapy.
- Do not place the non-perforated side of the warming blanket/gown on the patient. Always place the perforated side (with the small holes) directly on top of the patient in contact with the patient's skin.
- Do not place the warming blanket/gown directly over a dispersive electrode pad.
- For the Bair Hugger flex warming gown; ensure that the blood pressure cuff, ECG, IV or other lines or cords are not between the upper sleeve insert and the gown prior to deploying the upper sleeve warming insert(s), as this could result in tearing of the insert during deployment.
- For over-body warming blankets and gowns; do not place patient securement device (i.e. safety strap or tape) over the warming blanket/gown.
- For underbody or side channel warming blankets; if a securement device (i.e. safety strap, tape) is used, ensure the warming channels are not occluded.
- Do not allow the patient to lie on the warming unit hose.
- Do not allow the warming unit hose to directly contact the patient's skin during warming therapy.
- Do not connect a torn or damaged warming blanket/gown to the warming unit.
- Do not leave patients with poor perfusion unmonitored during prolonged warming therapy.
- Do not leave neonates, infants, children, and other vulnerable patient populations unattended during warming therapy.
- Do not continue warming therapy if the red Over-temp symbol illuminates and the alarm sounds. Unplug the warming unit and contact a qualified service technician.
- Do not perform the over-temperature detection system test while the warming unit is being used for warming therapy.

WARNING: To reduce the risk of patient injury or death due to altered drug delivery:

- Do not use a warming blanket/gown over transdermal medication patches.

WARNING: To reduce the risk of injury due to interference with ventilation

- Do not allow the warming blanket/gown or head drape to cover the patient's head or airway when the patient is not mechanically ventilated.

WARNING: To reduce the risk of injury due to patient falls

- Do not use a warming blanket/gown to transfer or move the patient.

WARNING: To reduce the risks associated with hazardous voltage and fire

- Keep power cord visible and accessible at all times. The plug on the power cord serves as the disconnect device.
- Only connect to outlets marked "Hospital only," "Hospital Grade," or a reliable grounded outlet.

- Use only the 3M-supplied power cord intended for this product and certified for the country of use.
- Do not allow the power cord to get wet.
- Do not use the warming unit when it appears the warming unit, power cord or any component is damaged. Contact 3M Technical Support.
- Do not disassemble the temperature management unit unless you are a qualified service technician. There are electrically live parts within the unit when it is connected to a power source.
- Connect each warming unit being tested to a separate power source.
- Do not attempt to replace the filter while the warming unit is in use.

CAUTIONS

CAUTION: To reduce the risk of cross-contamination:

- Except for specific Bair Hugger warming blanket models, 3M warming blankets/gowns are not sterile. Each warming blanket/gown is intended for single patient use ONLY. Placing a sheet between the warming blanket/gown and the patient does not prevent contamination of the product.
- Follow applicable regulations when disposing of this warming unit or any of its electrical components.
- Clean the warming unit and the outside of the warming unit hose after each patient use. See "Cleaning Instructions" in the Cleaning Instructions section. Do not attempt to clean the air filter as it may be contaminated from use. Discard the filter in a manner consistent with institutional protocol
- Do not operate the Bair Hugger warming unit with the hose detached from the 3M warming blanket/gown. The Bair Hugger warming unit is intended to be used with the hose properly attached to its corresponding warming blanket/gown and in accordance with good practices for operating room sterile technique.

CAUTION: To reduce the risk of fire:

- 3M warming blankets and gowns are classified as Class I Normal Flammability as defined by the Consumer Product Safety Commissions flammable fabric regulation, 16 CFR 1610. Follow standard safety protocols when using high intensity heat sources.

CAUTION: To reduce the risk of patient or caregiver injury:

- If mounted on an IV pole, the distance from the bottom of the warming unit to the floor must be less than 44" (112 cm) and the IV pole wheelbase diameter must be at least 28" (71 cm) to prevent tipping.
- Do not place the warming unit and rolling cart on an inclined surface of more than 7 degrees with rear casters facing downward angle.
- Do not attempt to move the rolling cart when the front casters are locked.

CAUTION: To reduce the risk of thermal injury, hyperthermia or hypothermia:

- 3M recommends continuously monitoring core temperature. In the absence of continuous monitoring, monitor the temperature of patients who are incapable of reacting, communicating and/or who cannot sense temperature a minimum of every 15 minutes or according to institutional protocol.
- Monitor cutaneous responses of patients who are incapable of reacting, communicating and/or who cannot sense temperature a minimum of every 15 minutes or according to institutional protocol.
- Adjust air temperature or discontinue therapy when the therapeutic goal is reached, if elevated temperatures are recorded or if there is an adverse cutaneous response in the warmed area.

- Do not place the warming unit on a soft uneven surface, such as a bed, or a visibly wet surface as the air intake may become blocked and cause the warming unit to overheat, compromising the warming unit's performance.
- Perform all temperature testing of the warming unit with a 3M Model 22110 temperature test unit. The Model 22110 temperature test unit itself requires no calibration.
- Only use 3M-supplied power cord with the warming unit.

CAUTION: To reduce the risk associated with electromagnetic interference (EMI) due to portable and mobile radio frequency (RF) communications equipment:


- The 3M Bair Hugger temperature management system has been tested to be resistant to both EMI and electrostatic discharge (ESD).
- Install and put into service the 3M Bair Hugger temperature management system according to the electromagnetic compatibility (EMC) information provided in the Guidance and Manufacturer's Declaration.
- Should interference occur, move away from the portable or mobile RF communications equipment.

Notices

1. The Bair Hugger warming unit meets medical electronic interference requirements. If radio frequency interference with other equipment should occur, connect the warming unit to a different power source.
2. To avoid warming unit damage:
 - Use proper Electrostatic Discharge (ESD) procedures when performing maintenance.
 - Do not modify this equipment without authorization from the manufacturer.
 - Do not immerse the warming unit, warming unit parts or accessories in any liquid or subject them to any sterilization process.
 - Do not use cleaning solutions with greater than 80% alcohol or solvents, including acetone and thinner, to clean the warming unit or hose. Solvents may damage the labels and other plastic parts.

Using the Model 675 Bair Hugger Warming Unit

Overview of Model 675 Warming Unit

Overview of Control Panel

Display Screen (LCD):

- Shows temperature of air coming out of the hose

Temperature Settings:


- (Ambient)
- 32° C (89.6° F)
- 38° C (100.4° F)
- 43° C (109.4° F)

Green LED:

- Shows selected temperature setting



Power:

- Green = On
- Amber = Standby


Important Alarm Symbols and Indicators


(The Model 675 unit's visual and audible alarms are classified as *Low Priority* per IEC 80601-2-35):

Service required


Red over-temperature

Blue under-temperature

- **Stop use:** Service required
- **Audible alarm will sound:** To silence the alarm press the Power button.
- **Do not re-start**

Filter time limit has been reached:

- Continue current use
- After current use, filter should be changed by service personnel
- Filter change symbol will stay illuminated, even when the warming unit is in standby mode


Warning: Do not attempt to replace the filter while the warming unit is in use.

Caution: Do not attempt to clean the air filter as it may be contaminated from use. Discard the filter in a manner consistent with institutional protocol.

Instructions for Use

Read and comply with all CONTRAINDICATIONS, WARNINGS and CAUTIONS, and NOTICES included in the Operator's Manual. The instructions below describe how to operate the Model 675 warming unit. For information about using 3M warming blankets/gowns, refer to the Instructions for Use included with the blanket/gown.


PREPARE	<ol style="list-style-type: none"> Place the unit on a flat, hard, dry surface (e.g., on the rolling cart) or mount it on an IV pole (See mounting instructions on page 7.) <p>Caution: If mounted on an IV pole, the distance from the bottom of the warming unit to the floor must be less than 44" (112 cm) and the IV pole wheelbase diameter must be at least 28" (71 cm) to prevent tipping.</p> <p>Caution: Do not place the warming unit on a soft uneven surface, such as a bed, or a visibly wet surface as the air intake may become blocked and cause the warming unit to overheat, compromising the warming unit's performance.</p>
CONNECT	<ol style="list-style-type: none"> Connect the warming unit to a properly grounded power source. Power button will illuminate amber showing the warming unit has power. Insert the end of the Bair Hugger warming unit hose into the warming blanket or gown hose port up to the stop collar. <p>Warning: Do not treat patients with the warming unit hose alone.</p> <p>Warning: Do not allow the patient to lie on the warming unit hose or the hose to have direct contact with the patient's skin during warming.</p> <p>Caution: Do not operate the Bair Hugger warming unit with the hose detached from the 3M warming blanket/gown.</p> Support the hose using the sheet clip to ensure secure attachment.
START	<ol style="list-style-type: none"> Select the appropriate temperature setting to start the warming unit. <ul style="list-style-type: none"> • Green LED shows selected temperature setting • The LCD display shows the temperature of the air exiting the hose.
MONITOR	<ol style="list-style-type: none"> Monitor the patient: <p>Caution:</p> <ul style="list-style-type: none"> • 3M recommends continuously monitoring core temperature. In the absence of continuous monitoring, monitor the temperature of patients who are incapable of reacting, communicating and/or who cannot sense temperature a minimum of every 15 minutes or according to institutional protocol. • Monitor cutaneous responses of patients who are incapable of reacting, communicating and/or who cannot sense temperature a minimum of every 15 minutes or according to institutional protocol. • Adjust air temperature or discontinue therapy when the therapeutic goal is reached, if elevated temperatures are recorded or if there is an adverse cutaneous response in the warmed area.
DISCONNECT	<ol style="list-style-type: none"> Complete warming therapy: <ul style="list-style-type: none"> • Press the Power button to turn the warming unit off • Disconnect the warming unit hose from the blanket/gown and discard the blanket/gown per hospital policy.
CLEAN	<ol style="list-style-type: none"> Clean the warming unit and the outside of the warming unit hose (See Cleaning Instructions, page 7). Properly store the unit and hose after cleaning.

Connecting the Warming Unit Hose and Power Cord

Connecting the Warming Unit Hose to the Warming Unit

To connect the warming unit hose to the warming unit, first ensure the O-ring is present on the hose end. Align the triangle on top of the warming unit hose connection with the notched groove on the hose end. Push the hose in and turn clockwise until it snaps into place. To disconnect turn hose counter-clockwise and pull.

Connecting the Power Cord to the Warming Unit


If the warming unit does not come with the power cord connected, follow the Power Cord Instructions for Use sheet included with the Model 675 warming unit.

Caution: Only use 3M-supplied power cord with the warming unit.


Warming Unit Mounting Instructions

Mounting on IV Pole

To mount the Model 675 warming unit to an IV pole, turn the clamp handle clockwise to tighten the clamp on the IV pole and counterclockwise to release it.


Caution: If mounted on an IV pole, the distance from the bottom of the warming unit to the floor must be less than 44" (112 cm) and the IV pole wheelbase diameter must be at least 28" (71 cm) to prevent tipping.

Mounting onto Rolling Cart Accessory

See Instructions included with the Model 96750 Rolling Cart Accessory.

- Phenolics (e.g. 3M™ Phenolic Disinfectant Cleaner)
- Aldehydes (e.g. 2% Activated Gluteraldehyde Solution)
- Alcohols (e.g. 70% Isopropyl Alcohol)

3. Let air dry.

Notice:

- Do not use cleaning solutions with greater than 80% alcohol or solvents, including acetone and thinner, to clean the warming unit or hose. Solvents may damage the labels and other plastic parts.
- Do not immerse the Bair Hugger warming unit, warming unit parts, or any of its accessories in any liquid or subject it to any sterilization process as this may damage the warming unit or its parts.

Storage

- Store all components in a cool, dry place when not in use.

Maintenance: Filter

When you see this symbol

- Continue current use.
- After current use, the filter should be changed by service personnel.
- Recommended filter change: Every 12 months or 500 hours of use.

Warning: Do not attempt to replace the filter while the warming unit is in use.

Caution: Do not attempt to clean the air filter as it may be contaminated from use. Discard the filter in a manner consistent with institutional protocol.

Cleaning, Storage, and Maintenance

Cleaning Instructions

- Disconnect the warming unit from the power source before cleaning.
- Cleaning should be performed in accordance with hospital practices for cleaning OR equipment. After every use; wipe the warming unit, the outside of the warming unit hose, and any other surfaces that may have been touched. Use a damp, soft cloth and a hospital approved mild detergent, germicidal disposable wipes, disinfecting towelettes, or antimicrobial spray. The following active ingredients are acceptable for use in cleaning the warming unit:
 - Oxidizers (e.g. 10% Bleach)
 - Quaternary Ammonium Compounds (e.g. 3M™ Quat Disinfectant Cleaner)

Proper Use and Maintenance

3M assumes no responsibility for the reliability, performance, or safety of the temperature warming unit if the following events occur:

- Modifications or repairs are not performed by a qualified, medical equipment service technician who is familiar with good practice for medical device repair.
- The unit is used in a manner other than that described in the Operator's or Preventive Maintenance Manual.
- The unit is installed in an environment that does not provide grounded electrical outlets.
- The warming unit is not maintained in accordance with the procedures described in the Preventive Maintenance Manual

Diagnostics

A qualified service technician can perform over-temperature detection system testing, temperature output testing, operating temperature calibration, and error/fault code troubleshooting.

Read Before Servicing Equipment

All repair, calibration, and servicing of the temperature warming unit require the skill of a qualified, medical equipment service technician who is familiar with good practice for medical device repair. (No further 3M warming unit-specific special training, licensing or certification required for servicing of Model 675 warming unit.) If service does not require the manufacturer's attention, the Model 675 Preventive Maintenance Manual provides the technical information needed to service the unit. Perform all repairs and maintenance in accordance with the instructions in the Model 675 Preventive Maintenance Manual For additional service information please contact 3M Technical Support.

In-Warranty Repair and Exchange

USA:

- Call 3M at 1-800-228-3957 if your Model 675 warming unit requires factory service.
- Provide the serial number located on the back of the warming unit.
- Technical Support will give you a Return Authorization (RA) number. Use this RA number on all correspondence. Your Technical Support representative will also send a shipping carton to you at no charge, if needed.

Outside of the USA:

Contact your local distributor concerning in-warranty repair and exchange. Provide the serial number located on the back of the Bair Hugger warming unit when you call.

PLEASE REPORT A SERIOUS INCIDENT OCCURRING IN RELATION TO THE DEVICE TO 3M AND THE LOCAL COMPETENT AUTHORITY (EU) OR LOCAL REGULATORY AUTHORITY.

Symbol Glossary

Symbol Title	Symbol	Description and Reference
Manufacturer		Indicates the medical device manufacturer. Source: ISO 15223, 5.1.1
Authorized Representative in European Community / European Union		Indicates the authorized representative in the European Community / European Union. Source: ISO 15223, 5.1.2, 2014/35/EU, and/or 2014/30/EU
Date of manufacture		Indicates the date when the medical device was manufactured. Source: ISO 15223, 5.1.3
Catalogue number		Indicates the manufacturer's catalogue number so that the medical device can be identified. Source : ISO 15223, 5.1.6
Serial number		Indicates the manufacturer's serial number so that a specific medical device can be identified. Source: ISO 15223, 5.1.7
Keep dry		Indicates a medical device that needs to be protected from moisture. Source: ISO 15223, 5.3.4
Consult instructions for use or consult electronic instructions for use		Indicates the need for the user to consult the instructions for use. Source: ISO 15223, 5.4.3
Caution		Indicates that caution is necessary when operating the device or control close to where the symbol is placed, or that the current situation needs operator awareness or operator action in order to avoid undesirable consequences. Source: ISO 15223, 5.4.4
Medical Device		Indicates the item is a medical device. Source: ISO 15223, 5.7.7
Unique device identifier		Indicates a carrier that contains unique device identifier information. Source: ISO 15223, 5.7.10
Importer		Indicates the entity importing the medical device into the locale. Source: ISO 15223, 5.1.8
No Free Hosing		Safety sign indicating no free hosing
CE Mark		Indicates conformity to all applicable European Union Regulations and Directives with notified body involvement.

Rx Only		Indicates that U.S. Federal Law restricts this device to sale by or on the order of a healthcare professional. 21 Code of Federal Regulations (CFR) sec. 801.109(b)(1)
Maximum safe working load		Indicates the maximum safe working load in less than the number reported.
UL Classified		Indicates product was evaluated and Listed by UL for the USA and Canada.
Green Dot		Indicates a financial contribution to national packaging recovery company per European Directive No. 94/62 and corresponding national law. Packaging Recovery Organization Europe.
Refer to instruction manual/booklet		To signify that the instruction manual/booklet must be read. ISO 7010-M002
Protective earth (ground)		To identify any terminal which is intended for connection to an external conductor for protection against electric shock in case of a fault, or the terminal of a protective earth (ground) electrode. Source: IEC 60417, 5019
Equipotentiality		To identify the terminals which, when connected together, bring the various parts of an equipment or of a system to the same potential, not necessarily being the earth (ground) potential. Source: IEC 60417-5021
Fuse		Indicates a replaceable fuse
Recycle electronic equipment		DO NOT throw this unit into a municipal trash bin when this unit has reached the end of its lifetime. Please recycle. Source: Directive 2012/19/EC on waste electrical and electronic equipment (WEEE)
IP Code		Indicates vertically falling drips shall have no harmful effects. Source: IEC 60529+AMD1:1999+AMD2:2013CSV/COR2:2015
Type BF applied part		To identify a type BF applied part complying with IEC 60601-1. Source: IEC 60417-5333
Use only 3M warming units and 3M warming blankets or gowns together		Safety sign indicating the Bair Hugger warming unit Model 675 has been designed to operate safely ONLY with 3M disposable warming blankets and gowns. Use with other products may cause thermal injury.

Over-temperature - Stop use: Service required		Indicates over-temperature set point has been reached. Stop use: Service required
Under-temperature - Stop use: Service required		The under-temperature set point has been reached. Stop use: Service required
Filter time limit has been reached		Indicates that after current use, the filter should be changed by service personnel
Power button (power indicator)		Color indicates the power status.
Ambient air (Room temperature air)		Indicates ambient air (Room temperature air) set point
Stop use: Service required		Indicates to stop use and service is required

For more information see, HCBRegulatory.3M.com

Specifications

Physical Characteristics

Dimensions	14.5" high x 10.5" deep x 11" wide (36.2 cm high x 26.7 cm deep x 28 cm wide)
Weight	9.93 lb (4.5 kg)
Relative noise level	50 dBA
Hose	Flexible Length: 80" (203 cm) Weight: 1.51 lb (0.686 kg) Diameter: 2" (5.08 cm)
Unit Intake Filter	MERV 14*
Recommended Filter Change	Every 12 months or 500 hours of use.

* See Frequently Asked Questions on Filtration

Recommended operating environment

Temperature	15°C-25°C (59°-77°F)
Humidity	Up to 80% relative humidity (non-condensing)
Altitude	Max 2,000 m

Operating temperatures

Temperature control	Electronically controlled
Heat generated	1600 BTU/hr (average), 470 W (average)
Time to reach operating temperature	2 - 5 minutes (dependent on blanket model) Time required for the contact surface temperature to heat up from 23 ± 2°C to 37°C (73 ± 2°F to 99°F)
Storage/transport temperature	-20°C to 45°C (-4°F to 113°F) Store all components in a cool, dry place when not in use.

Safety System

Thermostat	Independent electronic circuit; thermal cutoff shuts the heater OFF to ensure hose end air remains below 56°C (53°C ± 3°C typical).
Alarm system	Over-temperature (≤56°C, 53°C ± 3°C typical): red <i>Over-Temp</i> indicator symbol will flash alternately with the "Wrench" symbol, alarm sounds, heater and blower shut down, operating indicator lights turn OFF, and control panel becomes unresponsive. Remove from service. Under-temperature (<3°C of selected temperature/setting for 3 minutes): blue <i>Under-Temp</i> indicator symbol will flash alternately with the "Wrench" symbol, alarm sounds, heater and blower shut down, operating indicator lights turn OFF, and control panel becomes unresponsive. Remove from service.
Fault	"Wrench" symbol displays and alarm sounds. Heater and blower shut down, operating indicator light turns OFF and control panel becomes unresponsive. Service required.
Overcurrent protection	Dual input fused lines. Fuses are to be accessed and/or replaced only by service personnel.

Alarm/Fault Codes

Alarm conditions	General fault Over-temperature Under-temperature
Fault code indicator	Yes on LCD screen
Filter change indicator	Yes: Triggered by 500 hours of use

Electrical Characteristics

Heating element	1190 W Resistive
Heater control	3 stage heater flicker free
Temperature settings	Ambient 32° C (89.6° F) 38° C (100.4° F) 43° C (109.4° F)
Temperature accuracy	+/- 3° C
Temperature display resolution	0.1° C
Leakage current	Meets IEC 60601-1 requirements
Blower motor	Non-adjustable, one speed Airflow: up to 44 cfm or 20.77 L/s Operating speed: 4,700 RPM
Power consumption	Peak: 1255 W Average: 650 W

Power cord	15 ft., SJT, 3 cond., 13 A 15 ft., SJT, 3 cond., 15 A 4.6 m, HAR, 3 cond., 10 A
Device ratings	110-120 VAC, 50/60 Hz, 11 A, or 220-240 VAC, 50/60 Hz, 6 A, or 100 VAC, 50/60 Hz, 12 A
Fuses	Type: Fast acting ceramic fuses, 250 VAC AMP rating: 12 A (110 - 120 VAC) 8 A (220 - 240 VAC) 15 A (100 VAC) Operating speed: F (Fast Acting) Breaking capacity: 15 A, 12 A: 750 A @ 250 Vac 15 A, 12 A: 10,000 A @ 125 Vac 8 A: 200 A @ 250 Vac 8 A: 10,000 A @ 125 Vac

Guidance and manufacturer's declaration – electromagnetic emissions

The model 675 is intended for use in the electromagnetic environment specified below. The customer or the user of the model 675 should assure that it is used in such an environment.

Emissions test	Compliance	Electromagnetic environment - guidance
RF emissions, CISPR 11	Group 1	The model 675 uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF emissions, CISPR 11	Class A	NOTE: The EMISSIONS characteristics of this equipment make it suitable for use in industrial areas and hospitals (CISPR 11 class A). If it is used in a residential environment (for which CISPR 11 class B is normally required) this equipment might not offer adequate protection to radio-frequency communication services. The user might need to take mitigation measures, such as relocating or re-orienting the equipment.
Harmonic emissions, IEC 61000-3-2	Class A	
Voltage fluctuations/ Flicker emissions, IEC 61000-3-3	Complies	

Guidance and manufacturer's declaration – electromagnetic immunity

The model 675 is intended for use in the electromagnetic environment specified below. The customer or the user of the model 675 should assure that it is used in such an environment.

These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

IEC 60601 test level, Compliance level	Electromagnetic environment - guidance
<i>Electrostatic discharge (ESD), IEC 61000-4-2</i>	
±8 kV contact ±15 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.
<i>Electrical fast transient/burst, IEC 61000-4-4</i>	
±2 kV for power supply lines ±1 kV for input/output lines	Mains power quality should be that of a typical commercial or hospital environment.
<i>Surge, IEC 61000-4-5</i>	
±1 kV line to line ±2 kV line(s) to earth	Mains power quality should be that of a typical commercial or hospital environment.
<i>Voltage dips, short interruptions and voltage variations on power supply input lines, IEC 61000-4-11</i>	

<5% U_T (>95% dip in U_T) for 0,5 cycle 40% U_T (60% dip in U_T) for 5 cycles 70% U_T (30% dip in U_T) for 25 cycles <5% U_T (>95% dip in U_T) for 5 sec U_T is the a.c. mains voltage prior to application of the test level.	Mains power quality should be that of a typical commercial or hospital environment. If the user of the model 675 requires continued operation during power mains interruptions, it is recommended that the model 675 be powered from an uninterruptible power supply or a battery.
<i>Power frequency (50/60 Hz) magnetic field, IEC 61000-4-8</i>	
30 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.
IEC 60601 test level	Electromagnetic environment - guidance
<i>Conducted RF, IEC 61000-4-6, and Radiated RF, IEC 61000-4-3</i>	
6 Vrms, 150 kHz to 80 MHz 10 V/m, 80 MHz to 2,5 GHz	<p>Portable and mobile RF communications equipment should be used no closer to any part of the model 675, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter.</p> <p>Recommended separation distance</p> $d = 1.2\sqrt{P}$ $d = 1.2\sqrt{P} \text{ 80 MHz to 800 MHz}$ $d = 2.3\sqrt{P} \text{ 800 MHz to 2,5 GHz}$ <p>where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in meters (m).</p> <p>At 80 MHz and 800 MHz, the higher frequency range applies.</p> <p>Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey,^a should be less than the compliance level in each frequency range.^b</p> <p>Interference may occur in the vicinity of equipment marked with the following symbol: </p>
Notes	
<p>^a Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the model 675 is used exceeds the applicable RF compliance level above, the model 675 should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as re-orienting or relocating the model 675.</p> <p>^b Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V/m.</p>	

Recommended separation distances between portable and mobile RF communications equipment and the model 675

The model 675 is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the model 675 can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the model 675 as recommended below, according to the maximum output power of the communications equipment.

Rated maximum output power of transmitter W	Separation distance according to frequency of transmitter (m)		
	150 kHz to 80 MHz $d = 1.2\sqrt{P}$	80 MHz to 800MHz $d = 1.2\sqrt{P}$	800 MHz to 2,5 GHz $d = 2.3\sqrt{P}$
0,01	0.12	0.12	0.23
0,1	0.38	0.38	0.73
1	1.2	1.2	2.3
10	3.8	3.8	7.3
100	12	12	23

For transmitters rated at a maximum output power not listed above, the recommended separation distance d in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

Notes

At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.

These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

Classifications

- Protection against electric shock: Class I Medical Electrical Equipment with Type BF applied part.
- Protection against ingress of water: IPX1
- Mode of operation: Continuous operation.

Not suitable for use in the presence of flammable anesthetic mixtures with air or with oxygen or nitrous oxide.

MEDICAL – GENERAL MEDICAL EQUIPMENT AS TO ELECTRIC SHOCK, FIRE AND MECHANICAL HAZARDS ONLY IN ACCORDANCE WITH ANSI/AAMI ES 60601-1 (2005) + AMD (2012), CAN/CSA-C22.2 No. 60601-1 (2008) + (2014), IEC 80601-2-35 (2009) + AMD (2016), IEC 60601-1-6:2010 (Third Edition) + A1:2013, IEC 60601-1-8:2008 (Second Edition) + A1:2012; Control No. 4HZ8

Frequently Asked Questions (FAQ) on Filtration

What filter is used in the 3M™ Bair Hugger™ Model 675 warming unit?

The Model 675 uses the same filter media as has for many years been used in prior Bair Hugger 700-series models. That filter media is classified on an industry-standard scale as MERV 14 at 44 CFM (cubic feet per minute) under ASHRAE Standard 52.2-2012: Method of Testing General Air Ventilation Air-Cleaning Devices for Removal Efficiency by Particle Size.^{1,2}

What is a MERV rating?

The Minimum Efficiency Reporting Value, or MERV, is a measurement scale to rate the efficiency of an air filter as a function of particle size based on the ASHRAE Standard 52.2-2012.^{1,2}

What is the MERV rating recommended for operating rooms?

ASHRAE Standard 170-2013, Ventilation of Health Care Facilities, provides standards and guidance on the design requirements for ventilation systems utilized in health care facilities.³ The standard recommends that operating rooms utilize a filter that is equivalent to MERV 14.³ This MERV rating is specific to the ventilation system in the hospital's operating room.

What is the MERV rating recommended for patient warming devices in the operating room?

There is no filtration standard for patient warming devices utilized within the operating room. The ASHRAE Standard is a standard used to assess filter efficiency in health care ventilation systems.

Is a HEPA filter necessary?

No. We are not aware of any evidence or claim demonstrating that a HEPA filter improves patient safety.

What is the ASHRAE 52.2-2012 standard?

The American Society of Heating, Refrigerating and Air-Conditioning Engineers (ASHRAE) is a global organization focused on research, education and standards for heating, ventilation, air conditioning and refrigeration. The ASHRAE standard 52.2-2012^{1,2} provides test methods for air intake filters to ensure standardization in classification and efficiency. For more information on ASHRAE, please visit www.ashrae.org.

References:

1. ANSI/ASHRAE Standard 52.2-2012: Method of Testing General Air Ventilation Air-Cleaning Devices for Removal Efficiency by Particle Size. The American Society of Heating, Refrigerating and Air-Conditioning Engineers, Inc. © 2012 ASHRAE. ISSN 1041-2336. Available at <http://www.ashrae.org>. Accessed on 6/7/2016.
2. ANSI/ASHRAE Addenda a, b, and d to ANSI/ASHRAE Standard 52.2-2012. Method of Testing General Air Ventilation Air-Cleaning Devices for Removal Efficiency by Particle Size. The American Society of Heating, Refrigerating and Air-Conditioning Engineers, Inc. © 2015 ASHRAE. ISSN 1041-2336. Available at <http://www.ashrae.org>. Accessed on 6/7/2016.
3. ANSI/ASHRAE/ASHE Standard 170-2013: Ventilation of Health Care Facilities. The American Society of Heating, Refrigerating and Air-Conditioning Engineers, Inc. © 2013 ASHRAE. ISSN 1041-2336. Available at <http://www.ashrae.org>. Accessed on 6/7/2016.

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Références :

1. Norme ANSI/ASHRAE 52.2-2012 : Method of Testing General Air Ventilation Air-Cleaning Devices for Removal Efficiency by Particle Size. The American Society of Heating, Refrigerating and Air-Conditioning Engineers, Inc. © 2012 ASHRAE. ISSN 1041-2336. Disponible à l'adresse <http://www.ashrae.org>. Consultation : 07/06/2016.
2. Addenda ANSI/ASHRAE a, b et d à la norme ANSI/ASHRAE 52.2-2012. Method of Testing General Air Ventilation Air-Cleaning Devices for Removal Efficiency by Particle Size. The American Society of Heating, Refrigerating and Air-Conditioning Engineers, Inc. © 2015 ASHRAE. ISSN 1041-2336. Disponible à l'adresse <http://www.ashrae.org>. Consultation : 07/06/2016.
3. Norme ANSI/ASHRAE/ASHE 170-2013 : Ventilation of Health Care Facilities. The American Society of Heating, Refrigerating and Air-Conditioning Engineers, Inc. © 2013 ASHRAE. ISSN 1041-2336. Disponible à l'adresse <http://www.ashrae.org>. Consultation : 07/06/2016.

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Literaturangaben:

1. ANSI/ASHRAE Standard 52.2-2012: Method of Testing General Air Ventilation Air-Cleaning Devices for Removal Efficiency by Particle Size. The American Society of Heating, Refrigerating and Air-Conditioning Engineers, Inc. © 2012 ASHRAE. ISSN 1041-2336. Verfügbar unter <http://www.ashrae.org>. Zugriff am 7.06.2016.
2. ANSI/ASHRAE Addenda a, b, and d to ANSI/ASHRAE Standard 52.2-2012. Method of Testing General Air Ventilation Air-Cleaning Devices for Removal Efficiency by Particle Size. The American Society of Heating, Refrigerating and Air-Conditioning Engineers, Inc. © 2015 ASHRAE. ISSN 1041-2336. Verfügbar unter <http://www.ashrae.org>. Zugriff am 7.06.2016.
3. ANSI/ASHRAE/ASHE Standard 170-2013: Ventilation of Health Care Facilities. The American Society of Heating, Refrigerating and Air-Conditioning Engineers, Inc. © 2013 ASHRAE. ISSN 1041-2336. Verfügbar unter <http://www.ashrae.org>. Zugriff am 7.06.2016.

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Riferimenti:

1. ANSI/ASHRAE Standard 52.2-2012: Method of Testing General Air Ventilation Air-Cleaning Devices for Removal Efficiency by Particle Size. The American Society of Heating, Refrigerating and Air-Conditioning Engineers, Inc. © 2012 ASHRAE. ISSN 1041-2336. Disponibile su <http://www.ashrae.org>. Ultimo accesso 6/7/2016.
2. ANSI/ASHRAE Addenda a, b, and d to ANSI/ASHRAE Standard 52.2-2012. Method of Testing General Air Ventilation Air-Cleaning Devices for Removal Efficiency by Particle Size. The American Society of Heating, Refrigerating and Air-Conditioning Engineers, Inc. © 2015 ASHRAE. ISSN 1041-2336. Disponibile su <http://www.ashrae.org>. Ultimo accesso 6/7/2016.
3. ANSI/ASHRAE/ASHE Standard 170-2013: Ventilation of Health Care Facilities. The American Society of Heating, Refrigerating and Air-Conditioning Engineers, Inc. © 2013 ASHRAE. ISSN 1041-2336. Disponibile su <http://www.ashrae.org>. Ultimo accesso 6/7/2016.

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1. ANSI/ASHRAE Standard 52.2-2012: Method of Testing General Air Ventilation Air-Cleaning Devices for Removal Efficiency by Particle Size. The American Society of Heating, Refrigerating and Air-Conditioning Engineers, Inc. © 2012 ASHRAE. ISSN 1041-2336. Опубликовано на веб-сайте <http://www.ashrae.org>. По состоянию на 06.07.2016.
2. ANSI/ASHRAE Addenda a, b, and d to ANSI/ASHRAE Standard 52.2-2012. Method of Testing General Air Ventilation Air-Cleaning Devices for Removal Efficiency by Particle Size. The American Society of Heating, Refrigerating and Air-Conditioning Engineers, Inc. © 2015 ASHRAE. ISSN 1041-2336. Опубликовано на веб-сайте <http://www.ashrae.org>. По состоянию на 06.07.2016.
3. ANSI/ASHRAE/ASHE Standard 170-2013: Ventilation of Health Care Facilities. The American Society of Heating, Refrigerating and Air-Conditioning Engineers, Inc. © 2013 ASHRAE. ISSN 1041-2336. Опубликовано на веб-сайте <http://www.ashrae.org>. По состоянию на 06.07.2016.

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Referanslar:

- ANSI/ASHRAE Standard 52.2-2012: Method of Testing General Air Ventilation Air-Cleaning Devices for Removal Efficiency by Particle Size. The American Society of Heating, Refrigerating and Air-Conditioning Engineers, Inc. © 2012 ASHRAE. ISSN 1041-2336. <http://www.ashrae.org> adresinde mevcuttur. Erişim tarihi: 6/7/2016.
- ANSI/ASHRAE Addenda a, b, and d to ANSI/ASHRAE Standard 52.2-2012. Method of Testing General Air Ventilation Air-Cleaning Devices for Removal Efficiency by Particle Size. The American Society of Heating, Refrigerating and Air-Conditioning Engineers, Inc. © 2015 ASHRAE. ISSN 1041-2336. <http://www.ashrae.org> adresinde mevcuttur. Erişim tarihi: 6/7/2016.
- ANSI/ASHRAE/ASHE Standard 170-2013: Ventilation of Health Care Facilities. The American Society of Heating, Refrigerating and Air-Conditioning Engineers, Inc. © 2013 ASHRAE. ISSN 1041-2336. <http://www.ashrae.org> adresinde mevcuttur. Erişim tarihi: 6/7/2016.

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1. ANSI/ASHRAE Standard 52.2-2012: Method of Testing General Air Ventilation Air-Cleaning Devices for Removal Efficiency by Particle Size. The American Society of Heating, Refrigerating and Air-Conditioning Engineers, Inc. © 2012 ASHRAE. ISSN 1041-2336. Available at http://www.ashrae.org. Accessed on 6/7/2016.
2. ANSI/ASHRAE Addenda a, b, and d to ANSI/ASHRAE Standard 52.2-2012. Method of Testing General Air Ventilation Air-Cleaning Devices for Removal Efficiency by Particle Size. The American Society of Heating, Refrigerating and Air-Conditioning Engineers, Inc. © 2015 ASHRAE. ISSN 1041-2336. Available at http://www.ashrae.org. Accessed on 6/7/2016.
3. ANSI/ASHRAE/ASHE Standard 170-2013: Ventilation of Health Care Facilities. The American Society of Heating, Refrigerating and Air-Conditioning Engineers, Inc. © 2013 ASHRAE. ISSN 1041-2336. Available at http://www.ashrae.org. Accessed on 6/7/2016.

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