
3M

Ranger™

en	Pressure Infusor Model 145 Operator's Manual	02
fr	Manuel d'utilisation du dispositif de perfusion sous pression modèle 145	12
de	Druckinfusor, Modell 145 – Bedienungsanleitung	24
it	Manuale per l'operatore dell'infusore a pressione modello 145	35
es	Manual del operador del infusor de presión modelo 145	46
nl	Pressure Infusor Model 145 Bedieningshandleiding	57
sv	Tryckinfusor modell 145 användarhandbok	68
da	Brugsanvisning til trykinfusionsenhed model 145	79
no	Trykkinfusjonsenhet modell 145, Brukerhåndbok	90
fi	Paineinfusorin malli 145 käyttöopas	101
pt	Manual do Operador do Infusor de Pressão, Modelo 145	112
el	Εγχειρίδιο Χειριστή Εγχυτήρα Πίεσης Μοντέλο 145	123
pl	Infuzor ciśnieniowy, model 145 Instrukcja obsługi	135
hu	Infúziós pumpa, 145-ös típus – Kezelői kézikönyv	147
cs	Tlakový infuzor, model 145 – návod k obsluze	158
sk	Príručka pre obsluhu pre tlakový infúzor, model 145	169
sl	Navodila za uporabo tlačnega infuzorja, model 145	180
et	Rõhkinfuusori mudeli 145 kasutusjuhend	191
lv	Spiediena infūzijas ierīce, modelis 145, lietošanas rokasgrāmata	202
lt	Slėgio infuzoriaus modelio 145 naudotojo vadovas	213
ro	Perfuzor sub presiune Model 145 Manual de utilizare	224
ru	Руководство по эксплуатации компрессионного инфузора, модель 145	236
hr	Korisnički priručnik za infuzijsku pumpu, model 145	249
bg	Ръководство за оператора на компресионен инфузор модел	260
tr	Basınçlı İnfüzör Model 145 Kullanım Kılavuzu	272
zh	145 型压力输液器操作手册	283
ar	دليل المشغل للطراز 145 لجهاز الحقن بالضغط	293

MDCE2797

Made in USA with Globally Sourced Materials
3M Company
2510 Conway Ave.
St. Paul, MN 55144 USA
1-800-228-3957 (USA Only)

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EC REP

3M Deutschland GmbH
Health Care Business
Carl-Schurz-Str. 1
41453 Neuss, Germany

Issue Date 2021-10, rev. 3
34-8727-8209-8

Section 1: Technical Service and Order Placement

3

Technical service and order placement

3

USA

3

Proper use and maintenance

3

When you call for technical support

3

Servicing

3

Section 2: Introduction

3

Product description

3

Indications for use

3

Patient Population and Settings

3

Explanation of signal word consequences

3

WARNING:

3

CAUTION:

4

NOTICE:

4

Overview and Operation

4

Pressure infusor panel

5

Section 3: Instructions for Use

6

Attaching the pressure infusor to the I.V. pole

6

Load and pressurize the infusors

6

Changing a fluid bag

6

Section 4: Troubleshooting

7

Standby/ON mode

7

Pressure infusor

7

Section 5: Maintenance and Storage

8

General Maintenance and Storage

8

Cleaning Instructions

8

Storage

8

Servicing

8

Symbol Glossary

8

Section 6: Specifications

9

Physical characteristics

10

Electrical characteristics

11

Storage and transport conditions

11

Performance characteristics

11

Technical service and order placement

USA: TEL: 1-800-228-3957 (USA Only)

Outside of the USA: Contact your local 3M representative.

Proper use and maintenance

- 3M assumes no responsibility for the reliability, performance, or safety of the Ranger pressure infusor system if any of the following events occur:
- Modifications or repairs are not performed by a qualified, medical equipment service technician who is familiar with good practice for medical device repair.
 - The unit is used in a manner other than that described in the Operator's or Preventive Maintenance Manual.
 - The unit is installed in an environment that does not provide grounded electrical outlets.
 - The warming unit is not maintained in accordance with the procedures described in the Preventive Maintenance Manual.

When you call for technical support

We will need to know the serial number of your unit when you call us. The serial number label is located on the back of the pressure infusor system.

Servicing

All service must be performed by 3M or an authorized service technician. Call 3M at 1-800 228-3957 (USA only) for service information. Outside of the USA, contact your local 3M representative.

Section 2: Introduction

Product description

The Ranger pressure infusor consists of the pressure infusor device, which is combined with user supplied intravenous (IV) fluid bags which are pressurized in the infusor's chambers. The Ranger pressure infusor also requires user supplied disposable, sterile patient administration sets that can deliver the bag fluid to the patient under pressures up to 300 mmHg. The pressure infusor accepts solution bags ranging from 250 mL to 1000 mL. Fluids for use with the pressure infusor include, but are not limited to blood, saline, sterile water, and irrigation solution. The pressure infusor is intended to be used only with fluid bags that meet the standards of the American Association of Blood Banks. The pressure infusor is not intended for use with fluid bags and administration sets that do not meet the above specifications.

Indications for use

The 3M™ Ranger™ pressure infusor is intended to provide pressure to I.V. solution bags when rapid infusion of liquids is required.

Patient Population and Settings

Adult and pediatric patients being treated in operating rooms, emergency trauma settings, or other areas when rapid infusion of liquids is required. The liquids infused may interact with any part of the body as determined by the medical professional.

Explanation of signal word consequences

- WARNING:** Indicates a hazardous situation which, if not avoided, could result in death or serious injury.
- CAUTION:** Indicates a hazardous situation which, if not avoided, could result in minor or moderate injury.
- NOTICE:** Indicates a situation which, if not avoided, could result in property damage only.

WARNING:

- To reduce the risks associated with hazardous voltage and fire:
 - Connect power cord to receptacles marked "Hospital Only," "Hospital Grade," or a reliably grounded outlet. Do not use extension cords or multiple portable power socket outlets.
 - Do not position the equipment where unplugging is difficult. The plug serves as the disconnect device.
 - Use only the 3M-supplied power cord intended for this product and certified for the country of use.
 - Examine the pressure infusor for physical damage before each use. Never operate the equipment if the pressure infusor housing, power cord, or plug is visibly damaged. In the USA, contact 3M at 1-800-228-3957 (USA only). Outside of the USA contact your local 3M representative.
 - Do not allow the power cord to get wet.
 - Do not attempt to open the equipment or service the unit yourself. In the USA, contact 3M at 1-800-228-3957 (USA only). Outside of the USA contact your local 3M representative.
 - Do not modify any part of the pressure infusor.
 - Do not pinch the power cord of the pressure infusor when attaching other devices to the I.V. pole.
- To reduce the risks associated with exposure to biohazards:
 - Always perform the decontamination procedure prior to returning the pressure infusor for service and prior to disposal.
- To reduce the risks associated with air embolism and incorrect routing of fluids:
 - Never infuse fluids if air bubbles are present in the fluid line.
 - Ensure all luer connections are tightened.

SGS

9300 Winnetka Ave. N
Brooklyn Park, MN 55445
(763) 488-5700

Production Art

Description:
Ranger Pressure Infusor - IFU

Base File Name:
34872782098.indd

GA:

Software: InDesign CC 2020

SGS Contact:

Email:

Printer:

Printer Location:

Supplier:

Item Spec#: 34-8727-8209-8

Supersedes#:

Structure#: SS-87879

Cat/Product#: MD145

Reference:

Requester: Frank Boeshart

Die # / Doc. Size: 8.5" x 11"

CUST ID#

PRINT METHOD

SGS#

ROUND

DATE

CHANGE LOG

INT

11.09.21

New Build

SM

12.08.21

Translation Content updated

SM

12.08.21

Rx content only updated in all languages

SM

17.12.21

Amend as per annotation

GA

12.01.22

Amend as per annotation

SV

13.01.22

Final file release

GA

24.01.22

Final file release

MY

INKS

PROCESS
CYAN

PROCESS
MAGENTA

PROCESS
YELLOW

PROCESS
BLACK

PMS
294 C

MATCH

SPECIAL

LIMITATION OF LIABILITY – The supplied materials should be thoroughly inspected immediately upon receipt by the recipient for all content, placement and copy accuracy. It is also the responsibility of the recipient printer to examine and proofread all supplied materials before producing image carriers and printing final product. Should the supplied materials be found to be incorrect please contact an SGS representative immediately. The liability of SGS shall be limited to correcting its own product. SGS will not be responsible for any loss or damage including but not limited to the cost of product, print, image carriers, press delays, or other ancillary cost.

UPC NOTE – Because of the differences in scanning equipment, SGS cannot guarantee the scanability of UPC barcodes. The Uniform Code Council recommends a UPC symbol be printed black on white at 100% with no truncation. All UPC codes are test scanned and verified before leaving our facility.

2

CAUTION:

- 1. To reduce the risks associated with instability, impact and facility medical device damage:
 - Mount the Ranger pressure infusor Model 145 only on a 3M model 90068/90124 pressure infusor I.V. pole/base.
 - Do not mount this unit more than 56" (142 cm) from the floor to the base of the pressure infusor unit.
 - Do not use the power cord to transport or move the device.
 - Ensure the power cord is free from the castors during transport of the device.
 - Do not push on surfaces identified with the pushing prohibited symbol.
 - Do not pull the I.V. pole using the pressure infusor power cord.
- 2. To reduce the risks associated with environmental contamination:
 - Follow applicable regulations when disposing of this device or any of its electronic components.
- 3. This product is designed for pressure infusion only.

NOTICE:

- 1. The Ranger pressure infusor meets medical electronic interference requirements. If radio frequency interference with other equipment should occur, connect the pressure infusor to a different power source.
- 2. Do not use cleaning solutions with greater than 80% alcohol or solvents, including acetone and thinner, to clean the warming unit or hose. Solvents may damage the labels and other plastic parts.
- 3. Do not immerse the Ranger unit or accessories in any liquid or subject them to any sterilization process.

Overview and Operation

Overview and Operation

The pressure infusor system is attached to a custom I.V. pole with base. The power cord is secured allowing the system to be moved around the hospital/ healthcare facility to various locations where needed. The system is designed for frequent use whenever pressure infusion is necessary.


For operation, the Ranger pressure infusor is mounted on the custom I.V. pole and base. A fluid bag is placed in one of the pressure infusor chambers, the infusor's main power switch is turned ON, and the infusor chambers are activated using the user interface at the front of the device. Upon activation, the air is routed to the pressure infusor bladder, the bladder begins to inflate and pressurize the fluid bag. The user interface indicates when the fluid bag pressure is **In Range** and ready for use.

Note: The solution pressure of the Ranger pressure infusor is dependent on surface area and volume of the solution bag. To verify the pressure, refer to the preventive maintenance manual.

The Ranger pressure infusor has no user adjustable controls. The user slides an I.V. solution bag behind the metal fingers and against the inflation bladder located inside the pressure infusor. The Ranger pressure infusor should only be used in healthcare facilities by trained medical professionals and can be used within the patient environment.


When the infusor is attached to an external power source and the main power switch is turned ON, pushing the pressure infusor start/stop button ON inflates the inflation bladder and maintains pressure on blood and solution bags.

Turning the pressure infusor start/stop button OFF deflates the bladder. Turn the main power switch OFF when the pressure infusor is not in use

Pressure infusor panel

The pressure infusor panel displays the status of the pressure infusors. When the pressure infusor is initially turned ON, the indicators will illuminate to show operation. The pressure infusor power indicator illuminates yellow (standby) when the main power switch is turned ON and the pressure infusors are able to be turned ON. A green LED indicates that the infusor is ON. To pressurize/ depressurize the pressure infusor, verify the pressure infusor door is closed and latched, then press the pressure infusor start/stop button. Each pressure infusor is controlled independently.

1 Pressure Infusor Power

no power to unit

The indicator on the start/stop button notifies the user of the status of each pressure infusor. No light indicates that the unit is either not plugged in, the main power switch is not turned ON, or there is a system fault. See "Section 4: Troubleshooting" on 7 for more information.


standby

A yellow indicator notifies the user that the pressure infusor is in standby mode and is ready to be turned ON.

ON

A green indicator notifies the user that the infusor is pressurized.

2

High

Visual and audible indicator: The High yellow indicator illuminates and an audible indicator notifies the user when the pressure infusor bladder is above 330 mmHg. The visual and audible indicator will continue as long as the pressure remains above 330 mmHg. If the High condition is observed, the infusor chamber should be turned OFF by using the pressure infusor start/stop button. Use of the infusor chamber should be discontinued immediately, and 3M contacted for repair and servicing.

3

In Range

Visual only: The In Range green indicator flashes as the pressure is increasing in the pressure infusor. Once the pressure is within the target range of 230 330 mmHg, the indicator will be at a consistent green.

4

Low

Visual and audible indicator: The Low yellow indicator illuminates and an audible indicator notifies the user when the pressure infusor bladder has not reached 230 mmHg within approximately 30 seconds or when the pressure drops below 230 mmHg during use.

Section 3: Instructions for Use

NOTE: Assemble the model 90068/90124 pressure infusor I.V. pole/base according to the instructions for use that accompanies the I.V. pole base.

NOTE: Assembly of the I.V. pole/base and attachment of the pressure infusor to the I.V. pole should only be performed by a qualified, medical equipment service technician.


Attaching the pressure infusor to the I.V. pole

- 1. Mount the Model 145 pressure infusor onto a Model 90068/90124 pressure infusor custom I.V. pole and base.

CAUTION:

To reduce the risks associated with instability and facility medical device damage:

- Do not mount this unit more than 56" (142 cm) from the floor to the base of the pressure infusor.
- 2. Securely fasten the clamps at the back of the infusor and tighten the knob screws until the infusor is stable.
 - 3. Using the provided hook and loop strap, secure the power cord to the lower portion of the I.V. pole.

Load and pressurize the infusors

- 1. Plug the power cord into an appropriately grounded outlet.
- 2. Using the main power switch located under the pressure infusor, turn the unit ON.
- 3. If administration set is used, spike bags and prime administration set tubing, ensuring that all air is removed from the tubing set.
- 4. Slide the cassette into the slot in the 3M™ Ranger™ Warming Unit, Model 245 or 247. The cassette can only fit into the device one way
- 5. Connect warming set and continue priming, ensuring that all air is removed from the warming set. If no administration set is used, spike bags and prime warming set ensuring that all air is removed from the tubing set. For more information about priming the set, refer to instructions provided with the warming sets.
- 6. Open pressure infusor door.
- 7. Slide fluid bag to the bottom of the pressure infusor ensuring bag is completely inside the metal fingers.

Note: Ensure solution bag port and spike hang below the pressure infusor fingers.

- 8. Securely close and latch the pressure infusor door.
- 9. Close pinch clamps on tubing
- 10. Press the pressure infusor start/stop button on the pressure infusor control panel to turn the corresponding pressure chamber ON.

Note: A pressure infusor can only be turned ON when the start/stop button display status indicator is yellow. A green display status indicator notifies the infusor chamber is ON.

- 11. After pressing the start/stop button, the infusor indicator LED should be flashing green in the In Range status. When the indicator LED is solid green, open clamps to begin flow.
- 12. To depressurize, press the pressure infusor power button on the pressure infusor control panel to turn the corresponding pressure chamber off.

Changing a fluid bag

- 1. Press the pressure infusor button on the pressure infusor control panel to turn the corresponding pressure chamber OFF.
- 2. Close pinch clamps on tubing.
- 3. Open pressure infusor door and remove fluid bag.
- 4. Remove the spike from the used fluid bag.
- 5. Insert spike into the new I.V. bag port.
- 6. Push the pressure infusor's bladder to expel the remaining air and slide solution bag to the bottom of the pressure infusor, ensuring bag is completely inside the metal fingers.

Note: Ensure fluid bag port and spike hang below the pressure infusor fingers.

- 7. Securely close and latch the pressure infusor door.
- 8. Prime the warming set, ensuring all air is removed from the tubing. For more information about priming the set, refer to instructions provided with the warming sets.
- 9. Press the pressure infusor start/stop button on the pressure infusor control panel to turn the corresponding pressure chamber ON.

Note: A pressure infusor can only be turned on when the start/stop button display status indicator is yellow. A green display status indicator notifies the infusor chamber is ON.

- 10. Once the pressure infusor is *In Range*, open clamps to resume flow from the new bag of fluid.
- 11. Discard fluid bags and warming sets according to institutional protocol.

Section 4: Troubleshooting

All repair, calibration, and servicing of the Ranger pressure infusor requires the skill of a qualified, medical equipment service technician who is familiar with good practice for medical unit repair. All repairs and maintenance should be in accordance with manufacturer's instructions. Service the Ranger pressure infusor every six months or whenever service is needed. For replacement of the Ranger pressure infusor door latch, door, bladder, fingers, or power cord, contact a biomedical technician. For additional technical support refer to the preventive maintenance manual or contact 3M.

Standby/ON mode

Condition	Cause	Solution
Nothing illuminates on the pressure infusor control panel when the main power switch is turned ON.	Power cord is not plugged into the power entry module, or power cord is not plugged into an appropriately grounded outlet.	Make sure the power cord is plugged into the power entry module of the pressure infusor. Make sure the pressure infusor is plugged into a properly grounded outlet.
	Burned out LED light(s).	Contact a biomedical technician.
	Unit failure.	Contact a biomedical technician.
	Blown fuse.	Contact a biomedical technician.
Power LED status lights do not illuminate.	Power cord is not plugged into the power entry module, or power cord is not plugged into an appropriately grounded outlet.	Make sure the power cord is plugged into the power entry module of the pressure infusor. Make sure the pressure infusor is plugged into a properly grounded outlet.
	Unit is not turned ON.	Using the main power switch located under the pressure infusor turn the unit ON.
	Burned out LED light.	Push the pressure infusor start/stop button, if unit functions properly, continue use. Contact a biomedical technician after use to replace LED.
	Unit failure.	Contact a biomedical technician.
The status indicators (Low, In Range, and/or High) do not illuminate when the pressure infusor start/stop button is pushed.	Power cord is not plugged into the power entry module, or power cord is not plugged into an appropriately grounded outlet.	Make sure the power cord is plugged into the power entry module of the pressure infusor. Make sure the pressure infusor is plugged into a properly grounded outlet.
	Unit is not turned ON.	Using the main power switch located under the pressure infusor turn the unit ON.
	Burned out LED light(s).	Contact a biomedical technician.
	Unit failure.	Contact a biomedical technician.
	Blown fuse.	Contact a biomedical technician.

Pressure infusor

Condition	Cause	Solution
Pressure infusor is not working.	Power cord is not plugged into the power entry module, or power cord is not plugged into an appropriately grounded outlet.	Make sure the power cord is plugged into the power entry module of the pressure infusor. Make sure the pressure infusor is plugged into a properly grounded outlet.
	Unit is not turned ON.	Using the main power switch located under the pressure infusor turn the unit ON.
	Unit fault.	Discontinue use of unit. Contact a biomedical technician.
	Blown fuse.	Contact a biomedical technician.
Low indicator (solid yellow visual with audible indicator).	Pressure infusor bladder is loose or has become unattached.	Discontinue use of pressure infusor chamber. Use the other side of the pressure infusor. Reattach bladder by using your thumbs to fit one side of the bladder port on the bladder retaining collar and stretch into position.
	Pressure infusor door may not be closed and securely latched.	Securely close and latch the pressure infusor door.
	Detected pressure has fallen below 230 mmHg.	Continue infusion or use the other side of the pressure infusor. Contact a biomedical technician after use.
High indicator (solid yellow visual with audible indicator).	Pressure is above 330 mmHg.	Discontinue use of pressure infusor chamber. Use the other side of the pressure infusor. Contact a biomedical technician after use.
Leakage of fluid.	Bag is not spiked securely.	Secure spike in bag.
Bladder does not deflate after pressure is discontinued.	Unit fault.	Contact a biomedical technician after use.

Section 5: Maintenance and Storage

General Maintenance and Storage

All repairs and maintenance should be in accordance with manufacturer's instructions. Service the Ranger pressure infusor every six months or whenever service is needed. For replacement of the Ranger pressure infusor door latch, door, bladder, fingers, or power cord, contact a biomedical technician. For additional technical support refer to the preventive maintenance manual or contact 3M.

Cleaning Instructions

- Disconnect the Ranger unit from the power source before cleaning.
- Cleaning should be performed in accordance with hospital practices for cleaning OR equipment. After every use; wipe the warming unit, the outside of the warming unit hose, and any other surfaces that may have been touched. Use a damp, soft cloth and a hospital approved mild detergent, germicidal disposable wipes, disinfecting towelettes, or antimicrobial spray. The following active ingredients are acceptable for use in cleaning the warming unit:
 - Oxidizers (e.g. 10% Bleach)
 - Quaternary Ammonium Compounds (e.g. 3M Quat Disinfectant Cleaner)
 - Phenolics (e.g. 3M™ Phenolic Disinfectant Cleaner)
 - Alcohols (e.g. 70% Isopropyl Alcohol)

NOTICE:

- Do not use cleaning solutions with greater than 80% alcohol or solvents, including acetone and thinner, to clean the warming unit or hose. Solvents may damage the labels and other plastic parts.
- Do not immerse the Ranger unit or accessories in any liquid or subject them to any sterilization process.

Storage

Cover and store all components in a cool, dry place when not in use. Take care not to drop or jar the unit.

Servicing

All repairs and maintenance should be in accordance with manufacturer's instructions. Service the Ranger pressure infusor every six months or whenever service is needed. For replacement of the Ranger pressure infusor door latch, door, bladder, fingers, or power cord, contact a biomedical technician. For additional technical support refer to the preventive maintenance manual or contact 3M.

Please report a serious incident occurring in relation to the device to 3M and the local competent authority (EU) or local regulatory authority.

Symbol Glossary

The following symbols may appear on the product's labeling or exterior packaging.

Symbol Title	Symbol	Description and Reference
Manufacturer		Indicates the medical device manufacturer. Source: ISO 15223, 5.1.1
Authorized Representative in European Community / European Union		Indicates the authorized representative in the European Community / European Union. Source: ISO 15223, 5.1.2, 2014/35/EU, and/or 2014/30/EU
Date of manufacture		Indicates the date when the medical device was manufactured. Source: ISO 15223, 5.1.3
Catalogue number		Indicates the manufacturer's catalogue number so that the medical device can be identified. Source : ISO 15223, 5.1.6
Serial number		Indicates the manufacturer's serial number so that a specific medical device can be identified. Source: ISO 15223, 5.1 .7
Keep dry		Indicates a medical device that needs to be protected from moisture. Source: ISO 15223, 5.3.4
Caution		Indicates that caution is necessary when operating the device or control close to where the symbol is placed, or that the current situation needs operator awareness or operator action in order to avoid undesirable consequences. Source: ISO 15223, 5.4.4
Medical Device		Indicates the item is a medical device. Source: ISO 15223, 5.7.7
Unique device identifier		Indicates a carrier that contains unique device identifier information. Source: ISO 15223, 5.7.10
Importer		Indicates the entity importing the medical device into the locale. Source: ISO 15223, 5.1.8
No pushing		To prohibit pushing against an object. Source: ISO 7010-P017
CE Mark		Indicates conformity to all applicable European Union Regulations and Directives with notified body involvement.
Rx Only		Indicates that U.S. Federal Law restricts this device to sale by or on the order of a healthcare professional. 21 Code of Federal Regulations (CFR) sec. 801.109(b)(1)
Maximum safe working load		Indicates the maximum safe working load in less than the number reported.
UL Classified		Indicates product was evaluated and Listed by UL for the USA and Canada.
Refer to instruction manual/booklet		To signify that the instruction manual/booklet must be read. ISO 7010-M002

Defibrillation-proof Type BF applied part		Indicates the device applied part is Defibrillation-Proof Type BF. Source: EN 60601-1, Table D.1, Symbol 26
“OFF” (power)		To indicate disconnection from the mains, at least for main switches, or their positions, and all those cases where safety is involved. Source: IEC 60417-5008
“ON” (power)		To indicate connection to the mains, at least for mains switches, or their positions, and all those cases where safety is involved. Source: IEC 60417-5007
Protective earth (ground)		To identify any terminal which is intended for connection to an external conductor for protection against electric shock in case of a fault, or the terminal of a protective earth (ground) electrode. Source: IEC 60417, 5019
Equipotentiality		To identify the terminals which, when connected together, bring the various parts of an equipment or of a system to the same potential, not necessarily being the earth (ground) potential. Source: IEC 60417-5021
Fuse		Indicates a replaceable fuse
Recycle electronic equipment		DO NOT throw this unit into a municipal trash bin when this unit has reached the end of its lifetime. Please recycle. Source: Directive 2012/19/EC on waste electrical and electronic equipment (WEEE)
IP Code	IPX1	Indicates vertically falling drips shall have no harmful effects. Source: IEC 60529+AMD1:1999+AMD2:2013CSV/COR2:2015

For more information see, [HCBG](#)[regulatory.3M.com](#)

Section 6: Specifications

Guidance and manufacturer's declaration – electromagnetic emissions		
The model 145 is intended for use in the electromagnetic environment specified below. The customer or the user of the model 145 should assure that it is used in such an environment.		
Emissions test	Compliance	Electromagnetic environment - guidance
RF emissions CISPR 11	Group 1	The model 145 uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF emissions CISPR 11	Class B	The model 145 is suitable for use in all establishments, including domestic establishments and those directly connected to the public low voltage power supply network that supplies buildings used for domestic purposes.
Harmonic emissions IEC 61000-3-2	Class A	
Voltage fluctuations/Flicker emissions IEC 61000-3-3	Complies	

Guidance and manufacturer's declaration – electromagnetic immunity			
The model 145 is intended for use in the electromagnetic environment specified below. The customer or the user of the model 145 should assure that it is used in such an environment.			
Immunity Test	IEC 60601 test level	Compliance level	Electromagnetic environment - guidance
Electrostatic discharge (ESD) IEC 61000-4-2	±8 kV contact ±15 kV air	±8 kV contact ±15 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.
Electrical fast transient/burst IEC 61000-4-4	±2 kV power lines	±2 kV power lines	Mains power quality should be that of a typical commercial or hospital environment.
Surge IEC 61000-4-5	±1 kV line to line ±2 kV line to gnd	±1 kV line to line ±2 kV line to gnd	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	<5% U_T (>95% dip in U_T) for 0,5 cycle 40% U_T (60% dip in U_T) for 6 cycles 70% U_T (30% dip in U_T) for 30 cycles <5% U_T (>95% dip in U_T) for 5 sec	<5% U_T (>95% dip in U_T) for 0,5 cycle 40% U_T (60% dip in U_T) for 6 cycles 70% U_T (30% dip in U_T) for 30 cycles <5% U_T (>95% dip in U_T) for 5 sec	Mains power quality should be that of a typical commercial or hospital environment. If the user of the model 145 requires continued operation during power mains interruptions, it is recommended that the model 145 be powered from an uninterruptible power supply or a battery.
Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	3 A/m	3 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.

NOTE U_T is the a.c. mains voltage prior to application of the test level.

Guidance and manufacturer's declaration – electromagnetic immunity			
The model 145 is intended for use in the electromagnetic environment specified below. The customer or the user of the model 145 should assure that it is used in such an environment.			
Immunity Test	IEC 60601 test level	Compliance level	Electromagnetic environment - guidance
Conducted RF IEC 61000-4-6	3 Vrms 150 kHz to 80 MHz	3 Vrms	Portable and mobile RF communications equipment should be used no closer to any part of the model 145, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter. Recommended separation distance $d = 1.2 \sqrt{P}$ $d = 1.2 \sqrt{P}$ 80 MHz to 800 MHz $d = 2.3 \sqrt{P}$ 800 MHz to 2,5 GHz where <i>P</i> is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in meters (m). Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey^a, should be less than the compliance level in each frequency range^b. Interference may occur in the vicinity of equipment marked with the following symbol: 
Radiated RF IEC 61000-4-3	3 V/m 80 MHz to 2,5 GHz	3 V/m	
NOTE 1	At 80 MHz and 800 MHz, the higher frequency range applies.		
NOTE 2	These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.		

- ^a Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the model 145 is used exceeds the applicable RF compliance level above, the model 145 should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as re-orienting or relocating the model 145.
- ^b Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V/m.

Recommended separation distances between portable and mobile RF communications equipment and the model 145			
The model 145 is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the model 145 can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the model 145 as recommended below, according to the maximum output power of the communications equipment.			
Rated maximum output power of transmitter W	Separation distance according to frequency of transmitter m		
	150 kHz to 80 MHz $d = 1.2 \sqrt{P}$	80 MHz to 800 MHz $d = 1.2 \sqrt{P}$	800 MHz to 2,5 GHz $d = 2.3 \sqrt{P}$
0,01	0.12	0.12	0.23
0,1	0.37	0.37	0.74
1	1.17	1.17	2.30
10	3.69	3.69	7.37
100	11.67	11.67	23.30

For transmitters rated at a maximum output power not listed above, the recommended separation distance d in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

NOTE 1 At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.

NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

Physical characteristics

Dimensions

15.75 in (40 cm) high
20 in (51 cm) wide, 7.75 in (20 cm) deep

Weight

17 lb. (7.7 kg)

Mounting

Dual clamp

Classification

Protection against electric shock:

- Class I Medical Electrical Equipment with Type BF defibrillation-proof applied parts
- Protection against ingress of water: IPX1
- Mode of operation: Continuous operation.

MEDICAL – GENERAL MEDICAL EQUIPMENT AS TO ELECTRIC SHOCK, FIRE AND MECHANICAL HAZARDS ONLY IN ACCORDANCE WITH ANSI/AAMI ES 60601-1 (2005) + AMD (2012), CAN/CSA-C22.2 No. 60601-1 (2008) + (2014), and IEC 60601-1-6:2010 (Third Edition) + A1:2013; Control No. 4HZ8

Electrical characteristics

Leakage current

Meets leakage current requirements in accordance with IEC 60601-1.

Power cord

15 feet (4.6 m)

Device rating

110-120 VAC, 50/60 Hz, 1 Amp
220-240 VAC, 50/60 Hz, 0.8 Amp

Fuse

2 x F1A-H, rated 250V, for 110-120 VAC unit
2 x F0.8A-H, rated 250V, for 220-240 VAC unit

Storage and transport conditions

Storage/transport temperature

-20 to 60°C (-4°F to 140°F)

Store all components at room temperature and in a dry place when in use.

Operating humidity

Up to 90% RH, noncondensing

Atmospheric pressure range

Altitude up to 2000m or 80 kPa

Performance characteristics

Operating pressure

300 ± 10 mmHg setpoint

Note:

- Pressure system is In Range when the pressure infusor bladders are inflated to between 230 mmHg (low) and 330 mmHg (high). If pressure falls below 230 mmHg for more than approximately 30 seconds the Low yellow indicator will illuminate and an audible indicator will sound. The High yellow and audible indicator notifies the user when the pressure infusor bladder is above 330 mmHg.
- The outlet pressure of the fluid may vary with the surface area and volume of the fluid bag.

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Assistance technique et commandes

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3M and Ranger are trademarks of 3M.
All rights reserved.