3M™ Attest™ Steam Chemical Integrators 1243A and 1243B

Product Description

3M™ Attest™ Steam Chemical Integrators 1243A and 1243B are chemical indicators consisting of a paper wick and steam and temperature sensitive chemical pellet contained in a paper/film/foil laminate. The chemical pellet melts and migrates as a dark color along the paper wick. The migration is visible through a green window marked ACCEPT or red window marked REJECT; the extent of migration depends on steam and temperature. 3M™ Attest™ Steam Chemical Integrators are Type 5 (Category 15) Integrating Indicators as categorized by ISO 11410-1:2014.

Stated Values

See Stated Values listed in the table on the front of the package. The table lists the time in minutes at the corresponding temperature for chemical indicators in the lot contained in the package to reach their endpoint when tested in a steam resistometer under highly controlled conditions with saturated steam. The Stated Values are determined for each lot at the time of manufacture and are reproducible only under the exact conditions in which they were determined. The Stated Values compare with a 1 × 106 spore population of *Geobacillus stearothermophilus* with a D-value of 1.6 minutes at 121°C and a z-value of 10.0°C.

This Product is Not Made With Lead.

Indications for Use

Outside the United States

Use 3M™ Attest™ Steam Chemical Integrators for pack control monitoring of all 121-135°C (250-275°F) steam sterilization cycles.

United States

The 3M™ Attest™ Steam Chemical Integrators are designed to respond to all critical parameters over a specified range of steam sterilization cycles. The integrating indicator is intended to be placed in each pack, pouch, container, tray or other containment device to function as an independent monitor of critical parameters for the following sterilization cycles:

Cycle Type	Temperature	Exposure Time
Gravity	250°F/121°C	30 minutes
Gravity	270°F/132°C	4 minutes
Gravity	270°F/132°C	10 minutes
Gravity	270°F/132°C	15 minutes
Gravity	270°F/132°C	25 minutes
Gravity	275°F/135°C	3 minutes
Gravity	275°F/135°C	10 minutes
Dynamic Air Removal	250°F/121°C	30 minutes
Dynamic Air Removal	270°F/132°C	4 minutes
Dynamic Air Removal	270°F/132°C	10 minutes
Dynamic Air Removal	273°F/134°C	3 minutes
Dynamic Air Removal	273°F/134°C	4 minutes
Dynamic Air Removal	275°F/135°C	3 minutes

Minimum Stated Values for 3M™ Attest™ Steam Chemical Integrators as determined in a resistometer:

25.0°F/121°C	270°F/132°C	273°F/134°C	275°F/135°C
16.5 Minutes	2.0 Minutes	1.4 Minutes	1.2 Minutes

Contraindications

- None

Precautions

Do not use 3M™ Attest™ Steam Chemical Integrators to monitor dry heat, ethylene oxide, hydrogen peroxide, or other low temperature sterilization processes.

Directions for Use

1. Carefully tear or use scissors at the tear notch at the top of the foil package to make the initial opening. Remove only the number of 3M™ Attest™ Steam Chemical Integrators 1243A/1243B needed. Reseal the package.

2. Place a 3M™ Attest™ Steam Chemical Integrator 1243A/1243B in each pack, peel pouch, container system or tray to be steam sterilized in the area determined to be the least accessible to steam penetration.

3. Process the load according to established procedures.

4. After processing, the dark color should have entered the green ACCEPT window of the 3M™ Attest™ Steam Chemical Integrator 1243A/1243B. If the dark color has not entered the green ACCEPT window, a red REJECT result is indicated and the items in the pack, peel pouch, container system, or tray were not exposed to sufficient steam sterilization conditions. These items should be returned for reprocessing.

Storage/ Shelf Life/ Disposal

• Best stored at normal room temperature conditions, 15-30°C (59-86°F) and 40-60% relative humidity. Store away from direct sunlight. Do not store near strong alkaline or acidic products such as cleaning/sterilization agents.

• After use, the indicator will not change visually within 6 months when stored at above conditions.

Symbol Glossary

Symbol Title	Symbol	Description and Reference
Manufacturer		Indicates the medical device manufacturer as defined in Medical Device Regulation (EU) 2017/745 formerly EU Directive 93/42/EEC. Source: ISO 15223.5.1.1
Date of Manufacture		Indicates the date when the medical device was manufactured. Source: ISO 15223.5.1.3
Use-by date		Indicates the date after which the medical device is not to be used. Source: ISO 15223.5.1.4
Batch code	LOT	Indicates the manufacturer's batch code so that the batch or lot can be identified. Source: ISO 15223.5.1.5
Catalogue number	REF	Indicates the manufacturer's catalogue number so that the medical device can be identified. Source: ISO 15223.5.1.6
Do not re-use		Indicates a medical device that is intended for one use or for use on a single patient during a single procedure. Source: ISO 15223.5.4.2
Caution		Indicates the need for the user to consult the instructions for use for important cautionary information such as warnings and precautions that cannot, for a variety of reasons, be presented on the medical device itself. Source: ISO 15223.5.4.4
Steam indicator		Indicates product is designed for use with steam sterilization processes. Source: ISO 11410-1, 5.6

For more information see, HCBRegulatory.3M.com

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