3M™ Attest™ Steam Chemical Integrators 1243RE and 1243RES

Product Description
3M™ Attest™ Steam Chemical Integrators 1243RE and 1243RES are chemical indicators consisting of a paper wick and steam and temperature sensitive chemical pellet contained in a paper/polyester foil laminate. The chemical pellet melts and migrates as a dark color along the paper wick. The migration is visible through a green window marked ACCEPT or a red window marked REJECT; the extent of migration depends on steam, time, and temperature. 3M™ Attest™ Steam Chemical Integrators are Type 5 (Category 15) Integrating Indicators as categorized by ISO 11440-1:2014. This product comes with an extensor strip affixed to one end of a 3M™ Attest™ Steam Chemical Integrator. The affixed 1243RE extender is a 17.8 cm (7 in.) long by 1.8 cm (0.7 in.) wide strip. The affixed 1243RES extender is a 5.1 cm (2 in.) long by 1.8 cm (0.7 in.) wide rigid strip. The extensor allows a single integrator to retrieve processed integrators from inner packs.

Stated Values

See Stated Values listed in the table on the front of the package. The table sets the time in minutes at the corresponding temperature for chemical indicators in the lot contained in the package to reach their endpoint when tested in a steam resistorium under tightly controlled conditions with saturated steam. The Stated Values are determined for each lot at the time of manufacture and are reproducible only under the exact conditions in which they were determined. The Stated Values compare with a 1 × 10⁶ spore population of *Geobacillus stearothermophilus* with a D-value of 1.6 minutes at 121°C and a z-value of 10.0°C.

This Product Is Not Made with Lead.

Indications for Use

Use 3M™ Attest™ Steam Chemical Integrators for pack control monitoring of ISO 121-135°C (250–275°F) steam sterilization cycles. United States

The 3M™ Attest™ Steam Chemical Integrators are designed to respond to all critical parameters over a specified range of steam sterilization cycles. The integrating indicator is intended to be placed on each pack of a steam sterilization cycle. The integrator device to function as an independent monitor of critical parameters for the following sterilization cycles:

Cycle Type	Temperature	Exposure Time
Gravity	250°F/121°C	30 minutes
Gravity	270°F/132°C	3 minutes
Gravity	270°F/132°C	10 minutes
Gravity	270°F/132°C	15 minutes
Gravity	270°F/132°C	25 minutes
Gravity	275°F/135°C	3 minutes
Gravity	275°F/135°C	10 minutes
Dynamic Air Removal	250°F/121°C	30 minutes
Dynamic Air Removal	270°F/132°C	4 minutes
Dynamic Air Removal	270°F/132°C	10 minutes
Dynamic Air Removal	273°F/134°C	3 minutes
Dynamic Air Removal	273°F/134°C	4 minutes
Dynamic Air Removal	275°F/135°C	3 minutes

Minimum Stated Values for 3M™ Attest™ Steam Chemical Integrators as determined in a resistorium:

250°F/121°C	270°F/132°C	273°F/134°C	275°F/135°C
16.5 Minutes	2.0 minutes	1.4 minutes	1.2 minutes

Contraindications

Precautions
Do not use 3M™ Attest™ Steam Chemical Integrators 1243RE and 1243RES to monitor dry heat, ethylene oxide, hydrogen peroxide, or other low temperature sterilization processes.

Directions for Use

- Carefully tear or use scissors at the tear notch at the top of the foil package to make the initial opening. Remove the number of 3M™ Attest™ Steam Chemical Integrators 1243RE/1243RES needed. Reseal the package.
- Place a 3M™ Attest™ Steam Chemical Integrator 1243RE/1243RES in each pack, peel pouch, container system or tray to be steam sterilized in the area determined to be the least accessible to steam penetration. Position the unattached end of the extender so that it extends slightly beyond the top edge of the pack. This will permit integrator retrieval without touching the pack contents.
- Process the load according to established procedures.
- After processing, the dark color should have entered the green ACCEPT window of the 3M™ Attest™ Steam Chemical Integrator 1243RE/1243RES. If the dark color has not entered the green ACCEPT window, a red REJECT result is indicated and the items in the pack, peel pouch, container system, or tray were not exposed to sufficient steam sterilization conditions. These items should be returned for reprocessing.

Storage/ Shelf Life/ Disposal
• Best stored at normal room temperature conditions, 15–30°C (59–86°F) and 40–60% relative humidity. Store away from direct sunlight. Do not store near strong alkaline or acidic products such as cleaning/disinfecting agents.
• After use, the indicator will not change visually within 6 months when stored at above conditions.

For further information, please contact your local 3M representative or contact us at 3M.com and select your country.

Symbol Glossary

Symbol Title	Symbol	Description and Reference
Manufacturer		Indicates the medical device manufacturer as defined in Medical Device Regulation (EU) 2017/745 formerly EU Directive 93/42/EEC. Source: ISO 15223, 5.11
Date of Manufacture		Indicates the date when the medical device was manufactured. Source: ISO 15223, 5.13
Use-by date		Indicates the date after which the medical device is not to be used. Source: ISO 15223, 5.14
Batch code		Indicates the manufacturer's batch code so that the batch or lot can be identified. Source: ISO 15223, 5.15
Catalogue number		Indicates the manufacturer's catalogue number so that the medical device can be identified. Source: ISO 15223, 5.16
Do not re-use		Indicates a medical device that is intended for one use or for use on a single patient during a single procedure. Source: ISO 15223, 5.4.2
Caution		Indicates the need for the user to consult the instructions for use for important cautionary information such as warnings and precautions that cannot, for a variety of reasons, be presented on the medical device itself. Source: ISO 15223, 5.4.4
Steam indicator		Indicates product is designed for use with steam sterilization processes. Source: ISO 11440-1, 5.6

For more information see, HCBRegulatory.3M.com

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Item Stock:	34-8725-990-0
Supermarket:	
Structure:	SS-33391
SGS Contact:	1243RE-1243RES
Reference:	
Requester:	Kathryn Early
Part#: Die / Doc. Size:	21" x 21" FLAT 3" x 3" FOLDED

DESCRIPTION	3M™ Attest™ Steam Chemical Integrators 1243RE and 1243RES
Base File Name:	34-847259900.indd
GA Product:	1243RE-1243RES
Software:	Idesign CC 2017
SGS Contact:	
Printer:	
Print Location:	
Supplier:	

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3M and Attest are trademarks of 3M.
© 2019, 3M. All rights reserved.