







## 3M™ V.A.C. Veraflo™ Dressing Selection Guide

This guide may be used in deciding which wound dressing to use in conjunction with 3M™ Veraflo™ Therapy.

	3M™ V.A.C. Veraflo™ Dressing - Medium 	3M™ V.A.C. Veraflo Cleanse Choice™ Dressing 
Wound Types		
Wound Characteristics	Open wounds, including wounds with shallow undermining or tunnel areas where the distal aspect is visible	Wounds with thick fibrinous exudate, slough, infectious material and other wound bioburden
Key Goal(s) of Therapy	<ul style="list-style-type: none"><li>• When used in conjunction with Veraflo Therapy, to help facilitate the removal of infectious material and other wound bioburden</li><li>• Generation of robust granulation tissue</li></ul>	<ul style="list-style-type: none"><li>• When used in conjunction with Veraflo Therapy, to help facilitate the removal of infectious material and other wound bioburden</li><li>• Generation of robust granulation tissue</li></ul>
SKU	<b>Medium - ULTVFL05MD</b>	<b>Medium - ULTVCC05MD</b>
Shape	Spiral-cut foam	Block foam pre-slit into three layers
Application Characteristics	Easy application: <ul style="list-style-type: none"><li>• Size without scissors</li><li>• Precut area for pad application when used for bridging</li><li>• Single pad application</li></ul>	Easy application with flexibility: <ul style="list-style-type: none"><li>• Thin layers for improved conformability</li><li>• Multiple layers provide application options for wounds with varying depths</li><li>• Ideal for dirty wounds needing active therapy</li><li>• Single or duo pad application</li></ul>
Dressing Characteristics	<b>Medium:</b> 17.4 x 14.7 x 1.8cm per piece (2 pieces spiral) <ul style="list-style-type: none"><li>• Each piece 80cm uncoiled</li><li>• Each piece 74cm long bridge</li></ul> - 3 sheets 3M™ V.A.C.® Advanced Drape - 1 3M™ V.A.C. VeraT.R.A.C.™ Pad - 2pk 3M™ Cavilon™ No Sting Barrier Film - 1 disposable ruler/foam quantity label <b>FILL ASSIST SET POINT RANGE</b> <b>Medium:</b> 38-240mL (1 piece) 80-450mL (2 pieces)	<b>Medium:</b> 18.0cm x 12.5cm x Varying Thicknesses (0.8cm – 1.6cm) - 4 sheets 3M™ V.A.C.® Advanced Drape - 1 3M™ V.A.C. VeraT.R.A.C.™ Pad - 4pk 3M™ Cavilon™ No Sting Barrier Film - 1 disposable ruler/foam quantity label <b>FILL VOLUME START POINTS</b> <b>Medium:</b> 85mL (1.6cm Cover Layer); 42mL (0.8cm Cover Layer); 24mL (0.8cm Wound Contact Layer)