

3M™ Littmann® Stethoscopes

The Convenience of Tunable Technology


A quick reference guide



3M Health Care


Using a 3M™ Littmann® Stethoscope with tunable technology in a single-sided chestpiece design.

Littmann stethoscope's tunable technology provides convenience, versatility and enhanced acoustics. Its dual frequency diaphragm combines both bell (low-frequency) and diaphragm (high-frequency) modes in the convenience of a single-sided chestpiece. You control bell and diaphragm modes by pressure on the chestpiece rather than turning it over.



Bell Mode (low-frequency)

Light contact with chestpiece allows you to hear low-frequency sounds.



Diaphragm Mode (high-frequency)


Without turning the chestpiece over, use firm pressure to hear high-frequency sounds.

Available Models

- 3M™ Littmann® Master Cardiology Stethoscope
- 3M™ Littmann® Cardiology S.T.C. Stethoscope
- 3M™ Littmann® Master Classic II Stethoscope
- 3M™ Littmann® Select Stethoscope


Using tunable technology in a traditional chestpiece design.

The benefits of tunable technology are also available in a traditional design. These versatile stethoscopes include an adult tunable diaphragm and pediatric tunable diaphragm or traditional bell for small site auscultation.



Bell Mode (low-frequency)

Light contact with tunable chestpiece side allows you to hear low-frequency sounds.



Diaphragm Mode (high-frequency)


Without turning the chestpiece over, use firm pressure to hear high-frequency sounds.

Available Models

- 3M™ Littmann® Cardiology III Stethoscope
- 3M™ Littmann® Classic II S.E. Stethoscope
- 3M™ Littmann® Lightweight II S.E. Stethoscope


Using a traditional, combination-style chestpiece.

Traditional chestpieces have a bell side to hear low frequencies and a diaphragm side to hear high-frequency sounds. You must turn the chestpiece over to listen to the different sounds.



Bell Mode (low-frequency)

Use light contact on the bell side to hear low-frequency sounds.



Diaphragm Mode (high-frequency)

Turn the chestpiece over, index to the opposite side and use firm pressure to listen to high-frequency sounds.








Available Models

- 3M™ Littmann® Classic II Pediatric Stethoscope
- 3M™ Littmann® Classic II Infant Stethoscope













3M™ Littmann® Stethoscope

Ordering Information

Catalog #	Length	Tube Color	
Littmann Electronic Stethoscope Model 4000			
4000	27 in. (69 cm)	Black	
Littmann Electronic Stethoscope Model 2000			
2000	27 in. (69 cm)	Black	
Littmann Master Cardiology Stethoscope			
2159	22 in. (55.8 cm)	Black	
2160	27 in. (69 cm)	Black	
2163	27 in. (69 cm)	Burgundy	
2164	27 in. (69 cm)	Navy Blue	
2165	27 in. (69 cm)	Hunter Green	
Black Edition—Black-coated chestpiece and eartubes			
2161	27 in. (69 cm)	Black	
Littmann Cardiology S.T.C. Stethoscope			
4470	22 in. (55.8 cm)	Black	
4471	27 in. (69 cm)	Black	
4472	27 in. (69 cm)	Burgundy	
4473	27 in. (69 cm)	Navy Blue	
4474	27 in. (69 cm)	Hunter Green	
Littmann Cardiology III Stethoscope			
3127	22 in. (55.8 cm)	Black	
3128	27 in. (69 cm)	Black	
3129	27 in. (69 cm)	Burgundy	
3130	27 in. (69 cm)	Navy Blue	
3134	27 in. (69 cm)	Hunter Green	
Littmann Master Classic II Stethoscope			
2144L	27 in. (69 cm)	Black	
2143	27 in. (69 cm)	Purple	
2146	27 in. (69 cm)	Burgundy	
2147	27 in. (69 cm)	Navy Blue	
2630	27 in. (69 cm)	Caribbean Blue	
2632	27 in. (69 cm)	Hunter Green	
Black Edition—Black-coated chestpiece and eartubes			
2141	27 in. (69 cm)	Black	
Gold Edition—23 Karat gold-plated chestpiece and eartubes			
2142G	27 in. (69 cm)	Black	

The Littmann stethoscopes listed do not contain natural rubber latex or dry natural rubber as components in the product or packaging.

Catalog #	Length	Tube Color	
Littmann Classic II S.E. Stethoscope			
2201	28 in. (71.1 cm)	Black	
2203	28 in. (71.1 cm)	Gray	
2205	28 in. (71.1 cm)	Navy Blue	
2206	28 in. (71.1 cm)	Caribbean Blue	
2208	28 in. (71.1 cm)	Hunter Green	
2209	28 in. (71.1 cm)	Purple	
2210	28 in. (71.1 cm)	Raspberry	
2211	28 in. (71.1 cm)	Burgundy	
Littmann Classic II Pediatric Stethoscope			
2113	28 in. (71.1 cm)	Black	
2113R	28 in. (71.1 cm)	Red	
2119	28 in. (71.1 cm)	Caribbean Blue	
2122	28 in. (71.1 cm)	Raspberry	
Littmann Classic II Infant Stethoscope			
2114	28 in. (71.1 cm)	Black	
2114R	28 in. (71.1 cm)	Red	
2124	28 in. (71.1 cm)	Caribbean Blue	
2125	28 in. (71.1 cm)	Raspberry	
Littmann Select Stethoscope			
2290	28 in. (71.1 cm)	Black	
2293	28 in. (71.1 cm)	Burgundy	
2294	28 in. (71.1 cm)	Purple	
2296	28 in. (71.1 cm)	Raspberry	
2298	28 in. (71.1 cm)	Royal Blue	
2291	28 in. (71.1 cm)	Caribbean Blue	
Littmann Lightweight II S.E. Stethoscope			
2450	28 in. (71.1 cm)	Black	
2451	28 in. (71.1 cm)	Burgundy	
2452	28 in. (71.1 cm)	Caribbean Blue	
2453	28 in. (71.1 cm)	Lilac	
2454	28 in. (71.1 cm)	Ceil Blue	
2455	28 in. (71.1 cm)	Seafoam Green	
Catalog #	Length/Qty	Tube/Tag Color	
Littmann Teaching Stethoscopes			
Littmann Master Classic II Stethoscope (two headsets)			
2139	40 in. (101.6 cm)	Black	
Littmann Classic II S.E. Stethoscope (two headsets)			
2138	40 in. (101.6 cm)	Black	
Littmann Stethoscope Identification Tags			
2170	5/case	Black	
2171	5/case	Gray	



3M™ Littmann® Stethoscopes are renowned worldwide for their unsurpassed quality. To learn more about them, contact your distributor or visit www.3M.com/Littmann.



Health Care

3M Center, Building 275-4W-02
St. Paul, MN 55144-1000
USA
1 800 228-3957
www.3M.com/Littmann

3M Canada

Post Office Box 5757
London, Ontario N6A 4T1
Canada
1 800 563-2921



*Minimum 10%
Post-Consumer Fiber*

H.I. 4509
Printed in U.S.A.
© 3M 2003 All Rights Reserved
70-2009-0485-5