



Attest™

Sterilization Solutions



3M™ Attest™ Mini Auto-reader 490M System* Quick Reference Guide

BI Preparation and Use for Tabletop Sterilizers



1 Assemble PCD
Create a Process Challenge Device (PCD) by placing a BI** and a chemical integrator in representative packaging.



2 Place
Place PCD in a full load in the location recommended by sterilizer manufacturer.



3 Process
Run sterilization cycle.



4 Cool
After completion of cycle, retrieve PCD, remove BI and allow BI to cool 10 minutes prior to activation.

BI Activation and Incubation



1 Crush
Squeeze BI in activator to close cap and crush media ampoule.



2 Flick
Remove BI and flick wrist to ensure media flows to growth chamber.




3 Inspect
Verify media has transferred to BI growth chamber at the vial bottom.




4 Incubate
Place activated BI in any well. Display shows remaining minutes of incubation or test results.
DO NOT REMOVE BI until the (+) or (-) symbol indicates the test is complete.

Interpretation of 3M™ Sterilization Assurance Core Products


3M™ Attest™ Biological Indicator 1491 and 1492V




3M™ Comply™ Indicator Tape 1322




3M™ Comply™ Indicator Tape 1355



3M™ Attest™ Steam Chemical Integrator 1243



3M™ Attest™ Mini Auto-reader 490M



3M Company
3M Center, Building 275-4E-01
St. Paul, MN 55144 USA

© 2020 3M. All rights reserved. 3M, Attest and Comply are marks and/or registered marks of 3M. Unauthorized use prohibited. 70-2011-7923-4 Rev 9/20

*Refer to 3M™ Attest™ Mini Auto-reader 490M Operator Manuals and 3M™ Attest™ Biological Indicator 1491 and 1492V Instructions for Use for complete instructions and other important information. Refer to the Instructions for Use for complete instructions and other important information related to the use of all other 3M products included in this Guide.

**Refer to the indications for use provided in the biological indicator Instructions for Use. The indications for use should align with the sterilization cycle to be monitored.

