The fluid within the ampule is non-sterile. The applicator is sterile if package is intact.

Product Description

3M™ Cavilon™ Advanced Skin Protectant is a polymeric-cyanoacrylate solution intended for the cover and protection of intact or damaged skin. Upon application to skin, the liquid dries rapidly to form a primary long-lasting waterproof, highly durable film barrier. It is elastomeric, adhering to the contours of the skin and providing a uniform film. The film is transparent and possesses good oxygen and moisture vapor permeability. The polymer-cyanoacrylate is dispersed in a non-stinging solvent. The film is colorless, non-cytotoxic and has a low dermatitis potential. The film adheres to dry, moist or wet skin surfaces (i.e., superficial, partial thickness skin loss) and remains intact during conditions of continuous or repeated exposure to moisture or caustic irritants. It will wear off the skin and does not require removal.

Ingredients: Hexamethyldisiloxane, Acrylic Tetrapolymer, 2-Octyl Cyanoacrylate

Indications for Use

Cavilon Advanced Skin Protectant forms a film barrier intended to cover and protect intact or damaged skin. It is effective in conditions where wet and/or dry skin is frequently or continuously exposed to moisture and caustic irritants such as feces, digestive fluids, wound drainage and urine. The protective film barrier reduces pain associated with Incontinence Associated Dermatitis (IAD) and prevents, stops, and/or reverses the effects of IAD. Cavilon Advanced Skin Protectant also can be used in areas exposed to friction and shear from bedding, clothing, shoes or any other material that would rub against the skin allowing/enabling the skin to heal.

*The protective barrier creates a favorable environment for healing.

Intended Users

Intended for use by medical professionals, or under the guidance of a medical professional.

Patient target groups are those patients with damaged skin, or those at risk of developing skin damage (to include but not limited to: Incontinence-associated Dermatitis, peristomal/perifistula skin damage, periwound skin damage).

Contraindications

Cavilon Advanced Skin Protectant is NOT to be used:

- as a wound dressing for full thickness wounds
- in or around the eyes

Warnings

1. **DANGER! HIGHLY FLAMMABLE!**
2. Cavilon Advanced Skin Protectant is highly flammable until it has completely dried on the skin.
3. The product should only be applied when no ignition sources or heat-producing devices are in use. Allow a **minimum of 90 seconds** for product to dry and vapors to dissipate before using any device that could provide a source of heat or ignition. Pooling of product will significantly extend dry time.
4. Do not use Cavilon Advanced Skin Protectant near fire or flame, heat, sparks and sources of static discharge.
5. Use the product only in well ventilated areas.
6. Avoid use on individuals who are allergic to any of the ingredients.
7. The product is individually packaged for single use only. Reuse could result in increased risk of infection, or inadequate product performance.
8. The product is not intended for applications requiring sterile product (e.g. infusion catheter site protection and care, or surgical site protection).
9. Keep out of the reach of children.

Precautions

- Skin absorption and the effectiveness of topical medications (including: antimicrobials, antifungals, and analgesics) may be reduced or prevented by the presence of the Cavilon Advanced Skin Protectant.
- Use of other barrier products, ointments, creams or lotions may significantly reduce the effectiveness of the product.

- The product can increase the adhesion of some adhesive products (e.g. tape), particularly in the first few days of use. When removing an adhesive product, it is important to exercise care and follow the Directions for Use (Step 10, below)

Directions for Use

Figure 1

Figure 2

Figure 3

- 1) Cleanse skin before applying Cavilon Advanced Skin Protectant. Gently dry areas of intact skin as needed. If areas of erosion (denudement) are present, excess serous or sero-sanguinous drainage may be blotted with a gauze pad if necessary.
- 2) Open package and grasp applicator. See Figure 1. Place thumb at the end of the lever. Aim the sponge end of the applicator **downward** and firmly depress the lever to break the internal ampule. A snapping or popping noise will be noted as the ampule breaks. It is not necessary to continue to press the lever.
- 3) Keep the applicator pointed in a downward position for approximately **10 seconds** to allow fluid to flow into the foam sponge; continued pressure on the lever is not required. Note-fluid will not completely saturate the sponge to the edges.
- 4) Gently press the foam against the skin and then using an even, sweeping motion wipe the foam sponge across the skin. See Figure 2. Apply the product with minimal overlap. Firm pressure and/or downward pressure on the lever or the applicator is **not** needed and may result in excessive wetness or pooling of fluid. Continue application until the entire area of concern has been covered. The applicator may be used in any direction.

Apply the product to *all* skin that is exposed to stool and/or urine, or irritating body fluid. Allow the area to dry for at least 30 seconds. (Figure 3) **or for a minimum of 90 seconds** before using any device that could provide a source of heat or ignition. Refer to Warning Statement #3.

If an area is missed, wait until the fluid has dried completely before applying additional product.

- 5) If Cavilon Advanced Skin Protectant is applied within a skin fold or other area of skin-to-skin contact, make sure that the skin surfaces are separated to allow the fluid to dry completely (*at least 30 seconds*) before allowing skin to return to the normal position.
- 6) **When used under adhesive tapes, dressings, or devices, allow the product to dry for approximately 1 minute before covering with adhesive products.** Refer to Precaution Statement # 3.
- 7) Cleanse affected area as needed. Cavilon Advanced Skin Protectant is waterproof and is not removed by cleansing.
- 8) Reapply one to three times per week. Frequent application may result in buildup of the product.
- 9) The film will wear off the skin and does not require removal. If removal is desired, the film should be removed with an adhesive remover containing hexamethyldisiloxane (HMDS).
- 10) Use an adhesive remover containing HMDS to remove an adhesive product that has been placed over the film. This is especially important for patients with fragile skin.

Storage/ Shelf Life/Disposal

For best results, the product should be stored in a cool and dry environment. Avoid excessive heat.

For shelf life, refer to the expiration date on each package. Refer to facility policy for disposal.

Please report a serious incident occurring in relation to the device to 3M and the local competent authority (EU) or local regulatory authority.

How supplied

Applicators are individually packaged for single patient, one-time use only. The fluid within the ampule is non-sterile. The applicator is sterile if package is intact. Do not use the product if the package is damaged or opened.

Cat. No.	Size
5050G	2.7 ml applicator
5051G	0.7 ml applicator

For further information, please contact your local 3M representative or contact us at 3M.com and select your country.

Symbol Glossary

Symbol Title	Symbol	Description and Reference
Manufacturer		Indicates the medical device manufacturer as defined in Medical Device Regulation (EU) 2017/745 formerly EU Directive 93/42/EEC. Source: ISO 15223, 5.1.1
Authorized Representative in European Community		Indicates the authorized representative in the European Community. ISO 15223, 5.1.2
Date of Manufacture		Indicates the date when the medical device was manufactured. Source: ISO 15223, 5.1.3
Use-by date		Indicates the date after which the medical device is not to be used. Source: ISO 15223, 5.1.4
Batch code		Indicates the manufacturer's batch code so that the batch or lot can be identified. Source: ISO 15223, 5.1.5
Catalogue number		Indicates the manufacturer's catalogue number so that the medical device can be identified. ISO 15223, 5.1.6
Sterilized using ethylene oxide Applicator is sterile if package is intact		Indicates a medical device that has been sterilized using ethylene oxide. Source: ISO 15223, 5.2.3
Non Sterile solution		Indicates a medical device that has not been subjected to a sterilization process. Source: ISO 15223, 5.2.7
Do not use if package is damaged or open		Indicates a medical device that should not be used if the package has been damaged or opened. Source: ISO 15223, 5.2.8
Do not re-use		Indicates a medical device that is intended for one use or for use on a single patient during a single procedure. Source: ISO 15223, 5.4.2
Consult instructions for use or consult electronic instructions for use		Indicates the need for the user to consult the instructions for use. Source: ISO 15223, 5.4.3
Does not contain natural rubber latex		Indicates there is no presence of dry natural rubber nor natural rubber latex as a material of construction within the medical device or the packaging of a medical device. Source: ISO 15223, 5.4.5 and Annex B
Medical Device		Indicates the item is a medical device
Single sterile barrier system		To identify a single layer that makes up the sterile barrier system (SBS) and differentiate it from layers of protective packaging designed to prevent damage to the sterile barrier system and its contents
Unique device identifier		Indicates bar code to scan product information into patient electronic health record
Importer		Indicates the entity importing the medical device into the locale. Source: ISO 15223, 5.1.8

Symbol Title	Symbol	Description and Reference
UKCA mark		Indicates conformity to all applicable regulations and/or directives in the United Kingdom (UK), for products placed on the market in Great Britain (GB) with UK Approved Body involvement. (0086 to be replaced with UK Approved Body ID number)
Swiss Authorised Representative		Indicates the authorised representative in Switzerland. Source: Swissmedic.ch
CE Mark 2797		Indicates conformity to all applicable European Union Regulations and Directives with notified body involvement.
Flammable		Product is classified as flammable. Source: Globally Harmonized Symbols for Classifying Hazardous Chemicals

For more information, see HCBRegulatory.3M.com

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3M Australia Pty. Limited
ABN 90 000 100096
Building A, 1 Rivett Road
North Ryde, NSW 2113
Phone 1300 363 878
Product info: 136 136
Emergency: 1800 097 146

	3M Deutschland GmbH Health Care Business Carl-Schurz-Str. 1 41453 Neuss, Germany
	3M Schweiz GmbH Eggstrasse 91 8803 Rüschlikon Switzerland

UK Responsible Person:
3M United Kingdom PLC
3M Centre, Cain Road
Bracknell, RG12 8HT, United Kingdom

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3M and Cavilon are trademarks of 3M. Used under license in Canada. © 2022, 3M. Unauthorized use prohibited. All rights reserved. 3M et Cavilon sont des marques de commerce de 3M. Utilisées sous licence au Canada. © 2022, 3M. Toute utilisation non autorisée est interdite. Tous droits réservés.