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

3M Deutschland GmbH
Health Care Business
Carl-Schurz-Str. 1
41453 Neuss, Germany

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INSTRUCTIONS FOR USE

Indications for Use

The 3M™ Bair Hugger™ Temperature Management System is intended to prevent and treat hypothermia. In addition, the temperature management system can be used to provide patient thermal comfort when conditions exist that may cause patients to feel too warm or too cold. The temperature management system can be used with adult and pediatric patients.

- NOT STERILE.

Contraindications, Warnings and Cautions

Explanation of Signal Word Consequences

WARNING:

Indicates a hazardous situation which, if not avoided, could result in death or serious injury.

CAUTION:

Indicates a hazardous situation which, if not avoided, could result in minor or moderate injury.

CONTRAINDICATIONS: To reduce the risk of thermal injury:

- Do not apply heat to lower extremities during aortic cross-clamping. Thermal injury may occur if heat is applied to ischemic limbs.

WARNING: To reduce the risk of thermal injury:

- Do not treat patients with the Bair Hugger warming unit hose alone. Always attach the hose to a Bair Hugger warming blanket before providing warming therapy.
- Do not allow the patient to lie on the warming unit hose.
- Do not allow the warming unit hose to directly contact the patient's skin during warming therapy.
- Do not leave neonates, infants, children and other vulnerable patient populations unattended during warming therapy.
- Do not leave patients with poor perfusion unmonitored during prolonged warming therapy.
- Do not place the non-perforated side of the warming blanket on the patient. Always place the perforated side (with the small holes) directly on top of the patient in contact with the patient's skin.
- In the operating room, do not use this warming blanket with any device other than a Bair Hugger 500 series, 700 series or 675 warming unit.

- Do not use a Bair Hugger 200 series warming unit in the operating room.

- Do not use a Bair Hugger 800 series patient adjustable warming unit with any Bair Hugger warming blanket.

- Do not continue warming therapy if the red Over-temp indicator light illuminates and the alarm sounds. Unplug the warming unit and contact a qualified service technician.


- Do not place patient securement device (i.e. safety strap or tape) over the warming blanket.

Internal Use

Directions for Use

1. Place the perforated side of the 3M™ Bair Hugger™ Multi-position upper body blanket (the side with small holes) directly on top of the patient in contact with the patient's skin (Figure A).


2. Place the Multi-position upper body blanket on and/or around patient to cover the greatest amount of skin surface area. Blanket can be safely curved, bent, or conformed to accommodate a range of surgical positions and procedures. Position of blanket may vary based on procedure and needed surgical access. Warming blanket should be placed with adhesive strip toward surgical site. Examples of possible configurations are shown below (Figure B).


WARNING: Do not place patient securement device (i.e. safety strap or tape) over the warming blanket.

WARNING: Do not place the warming blanket directly over a dispersive electrode pad.


WARNING: Do not use a warming blanket over transdermal medication patches.


3. Remove the backing from the adhesive tape strip and adhere the warming blanket to the patient (Figure C). This, together with good surgical draping practice, helps to prevent warmed air from flowing towards the surgical site.

4. Pull the tie-strip tabs, centered along the upper and lower edges of the warming blanket. Tie these strips as needed to secure the blanket in the desired position (Figure D).


5. Two hose ports are provided for clinician preference. If indicated, place the removable hose port card in the hose port that is not being used during warming therapy (Figure E).


6. If the patient is mechanically ventilated (via endotracheal tube or laryngeal mask airway), unfold the head drape and place over the patient's head and neck (Figure F); otherwise, leave the head drape folded.

WARNING: Do not allow the warming blanket or head drape to cover the patient's head or airway when the patient is not mechanically ventilated.

7. Insert the end of the Bair Hugger warming unit hose in the hose port (Figure G). Use a twisting motion to ensure a snug fit. A visual marker is located around the mid-section of the hose end to guide the depth of hose insertion. Support the hose to ensure secure attachment.

WARNING: Do not treat patients with the warming unit hose alone. Always attach the hose to a Bair Hugger warming blanket before providing warming therapy.

8. Select the desired temperature setting on the warming unit to initiate warming therapy. (See the Operator's Manual for your specific Warming Unit Model)

CAUTION: Patient Monitoring Recommendations:

- 3M recommends continuously monitoring core temperature. In the absence of continuous monitoring, monitor the temperature of patients who are incapable of reacting, communicating and/or who cannot sense temperature a minimum of every 15 minutes or according to institutional protocol.
- Monitor cutaneous responses of patients who are incapable of reacting, communicating and/or who cannot sense temperature a minimum of every 15 minutes or according to institutional protocol.
- Adjust air temperature or discontinue therapy when the therapeutic goal is reached, if elevated temperatures are recorded or if there is an adverse cutaneous response in the warmed area.

9. Based on the warming unit model utilized, turn the unit off or to standby mode to discontinue warming therapy. Disconnect the warming unit hose from the warming blanket and discard the blanket per hospital policy.

Please report a serious incident occurring in relation to the device to 3M and the local competent authority (EU) or local regulatory authority.

Symbol Glossary

Symbol Title	Symbol	Description and Reference
Authorized Representative in European Community	EC REP	Indicates the authorized representative in the European Community. ISO 15223, 5.1.2
Batch code	LOT	Indicates the manufacturer's batch code so that the batch or lot can be identified. ISO 15223, 5.1.5
Catalogue number	REF	Indicates the manufacturer's catalogue number so that the medical device can be identified. ISO 15223, 5.1.6
Caution	!	Indicates the need for the user to consult the instructions for use for important cautionary information such as warnings and precautions that cannot, for a variety of reasons, be presented on the medical device itself. Source: ISO 15223, 5.4.4
CE Mark	CE	Indicates conformity to European Union Medical Device Regulation or Directive.
Date of Manufacture	Factory icon	Indicates the date when the medical device was manufactured. Source: ISO 15223, 5.1.3
Do not re-use	Recycling symbol with slash	Indicates a medical device that is intended for one use or for use on a single patient during a single procedure. Source: ISO 15223, 5.4.2

Do not use if package is damaged or open	Box with slash	Indicates a medical device that should not be used if the package has been damaged or opened. Source: ISO 15223, 5.2.8
Green Dot	Green Dot logo	Indicates a financial contribution to national packaging recovery company per European Directive No. 94/62 and corresponding national law. Packaging Recovery Organization Europe.
Manufacturer	Factory icon	Indicates the medical device manufacturer as defined in Medical Device Regulation (EU) 2017/745 formerly EU Directive 93/42/EEC. Source: ISO 15223, 5.1.1
Medical Device	MD	Indicates the item is a medical device.

Natural rubber latex is not present	LATEX	Indicates natural rubber or dry natural rubber latex is not present as a material of construction within the medical device or the packaging of a medical device. Source: ISO 15223, 5.4.5 and Annex B
Non Sterile	NON STERILE	Indicates a medical device that has not been subjected to a sterilization process. Source: ISO 15223, 5.2.7
Unique device identifier	UDI	Indicates bar code to scan product information into patient electronic health record
Use-by date	Hourglass	Indicates the date after which the medical device is not to be used. Source: ISO 15223, 5.1.4

For more information see, HCBGRegulatory.3M.com

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3M, Bair Hugger and the Bair Hugger logo are trademarks of 3M. Used under license in Canada. 3M, Bair Hugger et Bair Hugger logo sont des marques de commerce de 3M, utilisées sous licence au Canada. © 2019, 3M. All rights reserved./Tous droits réservés.