

Replacement of an insufficient composite filling: 1-step impression using 3M™ Impregum™ Super Quick Medium/Ultra Light Body Polyether Impression Materials

Clinical dentistry and photography by: Dr. Gunnar Reich, Munich, Germany
Lab work by: MDT Nicole Freund, Tutzing, Germany



Initial situation: Insufficient composite filling on upper first molar.



Shade selection for 3M™ Lava™ Esthetic Fluorescent Full-Contour Zirconia inlay.



Impression of opposing arch using 3M™ Imprint™ 4 Preliminary VPS Preliminary Impression Material.




Minimally invasive preparation design of molar prior to final impression taking.



Filling of 3M™ Impression Tray with 3M™ Impregum™ Super Quick Medium Body Polyether Material out of the 3M™ Garant™ Cartridge.



Syringing of 3M™ Impregum™ Super Quick Ultra Light Body Polyether Impression Material around the preparation.



Final full arch impression made with 3M™ Impregum™ Super Quick Polyether Material.



Close-up of the final impression with crisp details of the preparation.



For temporary treatment of the preparation 3M™ Protemp™ Plus Temporization Material was used.



Highly esthetic final 3M™ Lava™ Esthetic Fluorescent Full-Contour Zirconia inlay.

Please refer to the instructions for use for complete product information.



3M Oral Care
2510 Conway Avenue
St. Paul, MN 55144-1000
USA
1-800-634-2249

3M Canada
Post Office Box 5757
London, Ontario N6A 4T1
Canada
1-888-363-3685

3M, Garant, Impregum, Imprint, Lava and Protemp are trademarks of 3M Company or 3M Deutschland GmbH. Used under license in Canada. All other trademarks are owned by other companies. © 3M 2020. All rights reserved.

Dr. Gunnar Reich has received an honorarium from 3M Oral Care.

3M.com/Dental