
PICC/CVC Securement Device

- en Securement Device
- fr Dispositif de fixation
- de Fixiersystem
- it Dispositivo di fissaggio
- es Dispositivo de fijación
- nl Fixatiehulpmiddel
- sv Fästanordning
- da Fikseringsanordning
- no Fikseringsenhet
- fi Kiinnitysväline
- pt Dispositivo de fixação
- el Συσκευή Συγκράτησης
- pl Zabezpieczenie
- hu Rögzítőeszköz
- cs Fixační prostředek
- sk Fixačná pomôcka
- sl Varnostni pripomoček
- et Kinnitusseadis
- lv Fiksācijas ierīce
- lt Tvirtinimo įtaisas
- ro Dispozitiv de fixaree
- uk Засіб фіксації
- hr Uređaj za osiguravanje
- bg Изделие за фиксиране
- sr Sredstvo za fiksiranje
- tr Sabitleme Cihazı
- kk Бекіту құрылғысы
- th อุปกรณ์ยึดตำแหน่ง
- sq Pajisja e sigurisë
- mk Уред за прицврстување

Made in USA with globally sourced materials

3M Company
2510 Conway Ave.
St. Paul, MN 55144 USA
1-800-228-3957 (USA only)
3M.com

3M and Tegaderm are trademarks of 3M.
Used under license in Canada.
© 2023, 3M. Unauthorized use prohibited. All rights reserved.
Issue Date: 2023-01, rev. 3
34-8728-3432-9

EC REP 3M Deutschland GmbH
Health Care Business
Carl-Schurz-Str. 1
41453 Neuss, Germany

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en

Product Description

The 3M™ PICC/CVC Securement Device is designed to secure the majority of vascular access devices to the skin. The molded plastic securement device is integrated onto a breathable base with a gentle silicone adhesive.

Indications for Use

The 3M™ PICC/CVC Securement Device is designed to secure the majority of vascular access devices to the skin.

Warnings

- Do not use 3M™ PICC/CVC Securement Device where loss of adhesion could occur, such as with a confused patient, diaphoretic or non-adherent skin, or when the access device is not monitored daily.
- Failure to follow the manufacturer’s instructions for use may result in complications including inadequate securement and skin irritation.
- This is a single use securement device. Reuse and/or repackaging may create a risk of patient or user infection, compromise the structural integrity and/or essential material and design characteristics of the securement device, which may lead to securement device failure.
- This securement device is intended to be used with a securement dressing.

Precautions

- Observe universal blood and body fluid precautions and infection control procedures, during application and removal of the securement device.
- Any active bleeding at the insertion site should be stabilized before applying the securement device.
- The skin should be dry and free of detergent residue to prevent skin irritation and to ensure good adhesion. Allow all preps and protectants to dry completely before applying the securement device.
- The securement device may be used with sutures, if deemed necessary.
- Do not use if securement device restricts flow.
- Do not use on patients with known tape or adhesive allergies.
- Avoid securement device contact with alcohol or acetone: both can weaken bonding of components and securement device adherence.
- Minimize catheter manipulation during application and removal of the securement device.
- Orient the securement device so arrows point toward the catheter insertion site.
- Securement device and catheter position should be routinely inspected.
- Securement device should be replaced at least every 7 days.
- Do not re-sterilize.

Directions for Use

Failure to follow the manufacturer’s instructions for use may result in complications including inadequate securement and skin irritation.

This product is to be used by health care providers.

3M™ PICC/CVC Securement Device: The securement device is designed to secure the majority of vascular access devices to skin.

Site Preparation: Prepare the site according to institution protocol. Clipping of hair where the securement device and dressing will be placed may improve adhesion. Shaving is not recommended. The skin should be clean, dry and free of detergent residue. Allow all preps and skin protectants to dry completely before applying the securement device to prevent skin irritation and to ensure good adhesion.

Any active bleeding at the insertion site should be stabilized before applying the securement device.


 9300 Winnetka Ave. N Brooklyn Park, MN 55445 (763) 468-5700				CUST ID#	PRINT METHOD	SGS#	ROUND	
						7228418-1	3	
Production Art	Description: 3M™ PICC/CVC Securement Device	Item Spec#: 34-8728-3432-9	DATE				CHANGE LOG	INT
	Base File Name: PA: 34872834329.indd GA: Software: InDesign CC 2020 SGS Contact: Email: Printer: Printer Location: Supplier:	Supersedes#:	01.25.23	Amend as per Annotation				GA
		Structure#: SS-00000	02.08.23	Amend as per Annotation				GA
		Cat/Product#: 1879R-2100	02.08.23	Collect Files for Omni				GA
		Reference:	03.01.23	Review Omni files				SM
		Requester: Frank Boeshart						
	Die # / Doc. Size: 5.5" x 8.5"							
INKS	 PROCESS BLACK	MATCH				SPECIAL		
LIMITATION OF LIABILITY – The supplied materials should be thoroughly inspected immediately upon receipt by the recipient for all content, placement and copy accuracy. It is also the responsibility of the recipient printer to examine and proofread all supplied materials before producing image carriers and printing final product. Should the supplied materials be found to be incorrect please contact an SGS representative immediately. The liability of SGS shall be limited to correcting its own product. SGS will not be responsible for any loss or damage including but not limited to the cost of product, print, image carriers, press delays, or other ancillary cost. UPC NOTE – Because of the differences in scanning equipment, SGS cannot guarantee the scanability of UPC barcodes. The Uniform Code Council recommends a UPC symbol be printed black on white at 100% with no truncation. All UPC codes are test scanned and verified before leaving our facility.								

Securement Device Application:

Minimize catheter manipulation during application of the 3M™ PICC/CVC Securement Device.

Refer to Figure 1 for securement device application instructions.

1. Before placing securement device on the skin, orient the securement device with arrows pointing toward the insertion site as shown in **Figure 1a**.
2. Place catheter into securement device as shown in **Figure 1b** and weave lumen(s) under the single plastic arm.
3. Remove liner of attached tape strip and secure lumen(s) to the securement device base as shown in **Figure 1c**.
4. Position securement device on skin at the desired location. Pull and remove liner from one side of base to expose adhesive while holding the securement device in place. Pull and remove liner from the other side of base.
5. Apply pressure to the securement device base to establish good adhesion to the skin.
6. Apply a securement dressing, such as the Tegaderm™ I.V. Advanced Securement Dressing (catalog number 1837 and 1839) or the Tegaderm™ CHG Chlorhexidine Gluconate I.V. Dressing (catalog number 1877(R) and 1879(R)), following the manufacturer's instructions for use.


Site Care:

1. The site should be observed daily for signs of infection or other complications. If infection is suspected, remove the transparent dressing, inspect the site directly, and determine appropriate medical intervention. Infection may be signaled by fever, pain, redness, swelling, or unusual odor or discharge.
2. Inspect the securement device daily and change as necessary, in accordance with facility protocol. Securement device changes should occur at least every 7 days and may be needed more frequently with highly exudative sites or if integrity of the dressing is compromised.

Removal of Securement Device

Minimize catheter manipulation during removal of the 3M™ PICC/CVC Securement Device.

1. Using a low and slow technique, remove dressing in the direction of the catheter insertion site.
2. As the securement device is exposed secure the catheter hub with your gloved finger until securement device is uncovered. Leave remainder of dressing in place over catheter insertion site. **Figure 2**
3. Remove securement device tape strip from catheter lumen(s). **Figure 3**
4. Use gloved finger to stabilize the catheter hub and gently remove the catheter lumen(s) out from under the plastic arm of securement device.
5. Secure end of catheter lumens with tape to skin.
6. Stabilize catheter with one hand. Use your other hand to remove the securement device from skin and remove the remainder of dressing using a low and slow removal technique. **Figure 4**

Note: Dispose of product per facility protocol.

Storage and Shelf Life

For best results, store in a cool, dry place. For shelf life, refer to the expiration date on the package.

If the sterile package is damaged or unintentionally opened before use, discard the product, do not use.

How supplied / Ordering Information

For additional information visit www.3M.com

	3M™ PICC/CVC Securement Device
Catalog #	Securement Device Size
2100	2" x 2 1/8" (5.1 cm x 5.4 cm)

If you have any questions or comments, in the USA please contact the 3M Health Care Customer Help line at 1-800-228-3957.

In Canada, contact 3M Canada Company, P.O. Box 5757, London, Ontario, N6A 4T1, 1-800-364-3577.

For further information outside the United States, contact your local 3M representative or contact us at www.3M.com and select your country.

Please report a serious incident occurring in relation to the securement device to 3M and the local competent authority (EU) or local regulatory authority.

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Symbol Glossary

Symbol Title	Symbol	Description and Reference
Manufacturer		Indicates the medical device manufacturer. Source: ISO 15223, 5.1.1
Authorized Representative in European Community / European Union		Indicates the authorized representative in the European Community / European Union. Source: ISO 15223, 5.1.2, 2014/35/EU, and/or 2014/30/EU
Date of manufacture		Indicates the date when the medical device was manufactured. Source: ISO 15223, 5.1.3
Use-by date		Indicates the date after which the medical device is not to be used. Source: ISO 15223, 5.1.4
Batch code		Indicates the manufacturer's batch code so that the batch or lot can be identified. Source: ISO 15223, 5.1.5
Catalogue number		Indicates the manufacturer's catalogue number so that the medical device can be identified. Source: ISO 15223, 5.1.6
Sterilized using ethylene oxide		Indicates a medical device that has been sterilized using ethylene oxide. Source: ISO 15223, 5.2.3
Do not resterilize		Indicates a medical device that is not to be resterilized. Source: ISO 15223, 5.2.6
Do not use if package is damaged and consult instructions for use		Indicates a medical device that should not be used if the package has been damaged or opened and that the user should consult the instructions for use for additional information. Source: ISO 15223, 5.2.8
Upper limit of temperature		Indicates the upper limit of temperature to which the medical device can be safely exposed. Source: ISO 15223, 5.3.6
Do not re-use		Indicates a medical device that is intended for one single use only. Source: ISO 15223, 5.4.2
Consult instructions for use or consult electronic instructions for use		Indicates the need for the user to consult the instructions for use. Source: ISO 15223, 5.4.3
Does not contain natural rubber latex		Indicates there is no presence of dry natural rubber nor natural rubber latex as a material of construction within the medical device or the packaging of a medical device. Source: ISO 15223, 5.4.5 and Annex B
Medical Device		Indicates the item is a medical device. Source: ISO 15223, 5.7.7
Contains a medicinal substance		Indicates a medical device that contains or incorporates a medicinal substance. Source: ISO 15223, 5.4.7
Unique device identifier		Indicates a carrier that contains unique device identifier information. Source: ISO 15223, 5.7.10
Importer		Indicates the entity importing the medical device into the locale. Source: ISO 15223, 5.1.8
Single sterile barrier system		Indicates a single sterile barrier system. Source: ISO 15223, 5.2.11
UKCA Mark		Indicates conformity to all applicable regulations and/or directives in the United Kingdom (UK), for products placed on the market in Great Britain (GB).

Swiss Authorised Representative		Indicates the authorised representative in Switzerland. Source: Swissmedic.ch
CE Mark		Indicates conformity to all applicable European Union Regulations and Directives with notified body involvement.
Green Dot		Indicates a financial contribution to national packaging recovery company per European Directive No. 94/62 and corresponding national law. Packaging Recovery Organization Europe.

For more information, see HCBRegulatory.3M.com

- Manufacturer
- Fabricant
- Hersteller
- Produttore
- Fabricante
- Fabrikant
- Tillverkare
- Producent
- Produsent
- Valmistaja
- Fabricante
- Κατασκευαστής
- Wytwórca
- Gyártó
- Výrobce
- Yörobca
- Proizvajalec
- Tootja
- Ražotājs
- Gamintojas
- Producător
- Виробник
- Proizvođač
- Производител
- Proizvođač
- Üretici
- Өндүрүш
- ຜູ້ຜິຕ
- Prodhuesi
- Производитель

- Authorized Representative in European Community / European Union
- Représentant autorisé pour la Communauté européenne / l'Union européenne
- Bevollmächtigter Vertreter in der Europäischen Gemeinschaft/ Europäischen Union
- Rappresentante autorizzato nella Comunità Europea/Unione Europea
- Representante autorizado en la Comunidad Europea/Unión Europea
- Gemachtigde vertegenwoordiger in de Europese Gemeenschap/Europese Unie
- Godkänd representant inom EU
- Bemyndiget i EU
- Autorisert representant i Det europeiske fellesskap / EU
- Valtuutettu edustaja Euroopan yhteisön / Euroopan unionin alueella
- Representante Autorizado na Comunidade Europeia/União Europeia
- Εξουσιοδοτημένος αντιπρόσωπος στην Ευρωπαϊκή Κοινότητα / Ευρωπαϊκή Ένωση
- Autoryzowany przedstawiciel we Wspólnocie Europejskiej / Unii Europejskiej
- Hivatalos képviselő az Európai Közösségben / Európai Unióban
- Autorizovaný zástupce v Evropském společenství / Evropské unii
- Splnomocnený zástupca v Európskom spoločenstve/Európskej únii
- Pooblaščen zastopnik v Evropski skupnosti/Evropski uniji
- Volitatud esindaja Euroopa Ühenduses / Euroopa Liidus
- Pilnvarotais pārstāvis Eiropas Kopienā/Eiropas Savienībā

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