

Competitive Product Comparison

Paradigm™ Nano Hybrid Universal Restorative


Gradia® Direct X™

(Manufactured by GC Corporation)


Summary—Paradigm™ Nano Hybrid Universal Restorative offers these advantages:

- Paradigm Nano Hybrid restorative has a greater compressive strength, diametral tensile strength and flexural strength than Gradia Direct X.
- Paradigm Nano Hybrid restorative fracture toughness was measured to be greater than Gradia Direct X.
- The in-vitro wear rate of Paradigm Nano Hybrid restorative is lower than Gradia Direct X.
- The shrinkage of Paradigm Nano Hybrid restorative is less than Gradia Direct X.


Test Results and Discussion:




Strength: The compressive, diametral tensile and flexural strengths of Paradigm Nano Hybrid restorative were measured to be statistically higher than that of Gradia Direct X.



In-vitro wear: Wear rates were determined using a 3-body, in-vitro test. Dimensional loss during the 200,000 cycle test was measured every 40,000 cycles and the wear rates were calculated using linear regression. Paradigm Nano Hybrid restorative exhibited less wear than Gradia Direct X.



Single Notch Fracture Toughness: The single notch fracture toughness of Paradigm Nano Hybrid restorative was measured to be statistically higher than Gradia Direct X.



Volumetric Shrinkage: The volumetric shrinkage of Paradigm Nano Hybrid restorative was compared to Gradia Direct X using the test methodology described by Watts and Cash.

Competitive Product Comparison

	Paradigm™ Nano Hybrid Universal Restorative from 3M ESPE	Gradia® Direct X™ from GC Corporation
Composition	<p>Resin: BIS-GMA, UDMA, BIS-EMA, PEGDMA and TEGDMA</p> <p>Filler System:</p> <ul style="list-style-type: none"> • Surface modified zirconia/silica with a median particle size of approximately 3 microns or less • Non-agglomerated/non-aggregated 20 nanometer surface modified silica particles • 82% by weight (68% by volume) 	<p>Resin:</p> <ul style="list-style-type: none"> • Urethane dimethacrylate 15–30% • Dimethacrylate 0–5% • Camphorquinone <0.1% <p>Filler System:</p> <ul style="list-style-type: none"> • Silica powder 10–20% • Prepolymerised filler 20–30% • Fluoro Alumino-silicate glass 30–40% • Lanthanide nanofiller (posterior only)
Indications	<ul style="list-style-type: none"> • Direct anterior and posterior restorations including occlusal surfaces • Core build-ups • Splinting • Indirect restorations including inlays, onlays and veneers 	Direct restorative for Class I and II cavities
Manufacturer's Claims	<ul style="list-style-type: none"> • Easy to polish • Excellent strength and wear resistance • Does not stick to instruments when adapting • Easy to shape • Does not slump — holds shape prior to cure • Radiopaque • 12 shades that cover most dentists' daily needs 	<ul style="list-style-type: none"> • Excellent wear resistance — and won't wear away opposing enamel, either • Extremely low polymerisation shrinkage stress for reduced post-operative sensitivity • Activated filler particles strengthen matrix for increased fracture toughness • Smooth handling for easy placement, sculptable and brushable • Improved radiopacity
Shades	A1, A2, A3, A3.5, A4, 0A2, 0A3, B1, B2, B3, C2, D3	<p>Standard Shades: X-XBW, X-BW, X-A1, X-A2, X-A3, X-A3.5, X-B1, X-B2, X-C2, X-D2</p> <p>Inside Special Shade: X-A02</p> <p>Outside Special Shade: X-WT</p>
Delivery Options	3g syringes or 0.20g capsules	5.0g syringes or 0.3g unitip

Information obtained from the GC website and Gradia Direct X Instructions for Use.

Customer Service Center: 1300 363 454 (Australia) www.3MESPE.com.au
0800 80 81 82 (New Zealand) www.3MESPE.co.nz

3M ESPE

3M Health Care

3M Australia Pty Limited
 ABN 90 000 100 096
 Building A, 1 Rivett Rd
 North Ryde NSW 2113
 ph: 1300 363 454
www.3MESPE.com.au

3M Health Care

3M New Zealand Limited
 94 Apollo Drive
 Rosedale
 Auckland 0632
 Ph: 0800 80 81 82
www.3MESPE.co.nz

Please Recycle. Printed in Australia.
 © 3M 2012. All rights reserved.

3M, ESPE and Paradigm are trademarks of 3M or 3M ESPE AG. Direct X, GC and Gradia are not trademarks of 3M.