

Art meets science: preserving structure and function with Glass Ceramic Onlays

A case study by Dr. Giuseppe Marchetti, DDS, Italy featuring 3M™ RelyX™ Ultimate Adhesive Resin Cement and 3M™ Scotchbond™ Universal Adhesive



A patient presented with failing amalgam filling in lower molars. Because there was a good amount of existing tooth structure, the decision was made to use a minimally invasive approach. Glass ceramic onlays were used to achieve function with good esthetics.



Fig. 1: Initial situation with failing amalgam fillings.



Fig. 2: Removal of amalgam and caries.



Fig. 3: Final view of cleaned cavities.




Fig. 4: Application of the 3M™ Scotchbond™ Universal Adhesive in self-etch mode on dentin. Rub for 20 seconds and air thin for 5 seconds.




Fig. 5: Followed by light cure of adhesive for 10 sec with 3M™ Elipar™ DeepCure-S LED Curing Light.



Fig. 6: Application of 3M™ Filtek™ Bulk Fill Flowable Restorative in order to render the cavity suitable for impression.



Fig. 7: Final view after composite light curing and polishing.




Fig. 8: Impression with 3M™ Imprint™ 4 VPS Impression Material (light and heavy body).

Art meets science: preserving structure and function with Glass Ceramic Onlays cont.



Fig. 9: Ceramic onlays on the laboratory cast.




Fig. 10: Ceramic onlays ready for cementation.




Fig. 11: After etching with hydrofluoric acid, apply 3M™ Scotchbond™ Universal Adhesive to the inner surface of the onlay, rub for 20 seconds, then air thin for 5 seconds until no more movement of solvent is seen. Apply Scotchbond Universal Adhesive according to tooth surface per Instructions for Use (IFU).



Fig. 12: Apply 3M™ RelyX™ Ultimate Adhesive Resin Cement to adhesively cement the onlays.



Fig. 13: Isolation of neighboring teeth and cementation of the first onlay. Apply 3M™ Scotchbond™ Universal Adhesive to the surface of flowable material, rub for 20 seconds, then air thin for 5 seconds before applying 3M™ RelyX™ Ultimate Adhesive Resin Cement.



Fig. 14: View after cementation of both onlays and excess cement removal.



Fig. 15: Checking occlusion after light curing or allowing the material to self cure per Instructions For Use (IFU).



Fig. 16: Final restorations in place.

3M

3M Oral Care
2510 Conway Avenue
St. Paul, MN 55144-1000 USA

Phone 1-800-634-2249
Web 3M.com/OralCare

3M.com/DentalCements

Before using the products described, please refer to the instructions for use provided with the product packages.

The featured 3M product may be known with an alternative name in different regions.

3M, Elipar, Filtek, Imprint, RelyX and Scotchbond are trademarks of 3M or 3M Deutschland GmbH. Used under license in Canada. All other trademarks are owned by other companies. Please recycle. Printed in U.S.A.
© 3M 2018. All rights reserved. Dr. Marchetti has received an honorarium from 3M Oral Care. 70-2013-7039-5