

3M™ Tube Securement Device


This strong, easy-to-use solution secures non-critical medical tubes made from silicone, PVC and silicone coated latex spanning 6-24 French. The tube securement device provides additional securement to fluid drainage tubing to help stabilize and prevent movement. It is gentle to skin and designed to support pressure injury and medical adhesive-related skin injury (MARS!) prevention practices.

Uniquely shaped

Designed to lift tubes off skin and help promote consistent application.

Sterile device

In single-use packaging.



Reinforced soft cloth backing

Provides extra strength and flexibility.

Two sizes

Small device for medical tubes and catheters 6-12 Fr and large device for medical tubes and catheters 12-24 Fr.



Skin performance: The feather designates products that deliver the securement power you need while minimizing damage to skin.

Reliable securement

Unique adhesive secures medical tubes and catheters for up to 4 days.¹

Ease of use

Easy to apply and removes cleanly in one piece; allows for standardization of technique.

Skin safe

Non-sensitizing silicone adhesive has low potential to cause dermatitis with minimal skin damage and pain upon removal.¹

Designed to help reduce pressure injury risk

With no hard components to cause pressure, product's unique shape minimizes contact between skin and tube.

Considerations and Indications

When it comes to securing tubing, there are a variety of factors to keep in mind, including securement performance, skin protection and clinician ease of use. The 3M™ Tube Securement Device is indicated for non-critical tube securement or in conjunction with sutures, primary securement dressings and/or urinary catheter balloons per facility protocol. Two sizes accommodate tubes spanning 6-24 French.

Applications



Vascular Access

The tube securement device meets the Infusion Therapy Standards of Practice guidelines as an adhesive securement device and should be used in conjunction with sutures and/or primary securement dressings.

Consider using the securement device as secondary securement for primary tubing and extension tubing on central lines and peripheral lines. It can help organize lines with catheters that have multiple lumens.

Surgical Drains

Provides additional securement in conjunction with sutures to fluid drainage tubing to help prevent movement.

Urinary Catheters

When compared to StatLock® Foley Stabilization Device in a comparative prospective study of 24 healthy human subjects showed that the tube securement device:



- 9x faster removal. The average removal time for the tube securement device was 13 seconds compared to 126 seconds for StatLock® Foley Stabilization Device (P<.0001).*
- No removal agent. On average no wipes were needed to remove the tube securement device, while the StatLock® Foley Stabilization Device requires 12 alcohol wipes for removal.*

Feel secure in your choice. Trust the science behind tube securement device for reliable performance.

To order, contact your 3M Health Care Sales Representative or call the 3M Health Care Customer Helpline at 1-800-228-3957. Visit [3M.com/TubeSecurement](https://www.3m.com/TubeSecurement) for more information.

Refer to product Instructions for Use (IFU) for other important information.

Ordering Information

Product	Product Name	Product Number	Tube & Catheter Size	Product Size	Devices/Box	Boxes/Case
	3M™ Tube Securement Device, Small	1500U	6-12 French	1.7 in x 2.4 in (4,3 cm x 6,1 cm)	25	4
	3M™ Tube Securement Device, Large	1501U	12-24 French	2.2 in x 3.5 in (5,5 cm x 9,0 cm)	25	4

*As tested in a comparative, prospective clinical study on 24 healthy human subjects after 3 days of wear followed by device removal from skin per IFU.



3M Company
2510 Conway Avenue
St. Paul, MN 55144-1000 USA
Phone 1-800-228-3957
Web [3M.com/Medical](https://www.3m.com/Medical)

© 2023 3M. All rights reserved.
3M and the other marks shown are marks and/or registered marks. Unauthorized use prohibited.
3M marks used under license in Canada.
All other marks are property of their respective owners.
GL_ENG_70-2013-1619-0