


# Addendum for Operating Manuals positive pressure Full Face Masks

## **Replacing the Plug Connector**

#### **Maintenance**

The plug connector can be loosened or tightened with the tool for push to connect adapters (material number 10035756). Prior to re-assembly, the O-ring at the thread run out is replaced.



- 1a AutoMaXX connector 3S
- 1b AutoMaXX connector Ultra Elite
- 2 ESA connector [assembly with O-Ring]
- 3 PS connector [assembly with O-Ring]

Adapter	For masks	Torque
AutoMaXX	3S-PF and Ultra Elite PF MaXX	10 -1 Nm
ESA	3S-PF and Ultra Elite PF ESA	7,0 +/- 0,5 Nm
PS	3S-PS and Ultra Elite PS	2,5 +/- 0,25 Nm

With every maintenance, check the O-ring in the plug connector.

A damaged O-ring must be replaced immediately.

If it is not damaged, lubricate it (PS connector only) extremely lightly with Parker-O-Lube (material number D0010239).

Regardless of damage O-rings have to be replaced at certain intervals (see table in the operating or maintenance manual).

MSA recommends to check the fit of the plug connector regularly.



#### Attention!

Using torques higher as those listed above could damage the full face mask and destroy its protective properties.

### **Ordering information**

Description	<b>Part No.</b> 10032072
(1a) Quick-connect coupling AutoMaXX / 3S	
(1b) O-Ring for quick connect adapter, 10 pcs	10145638
(1b) Quick-connect coupling AutoMaXX / Ultra Elite	10145637
(2) Quick-connect coupling ESA [including gaskets]	10037519
(2) Gaskets for coupling [kit of 2]	10037709
(3) Plug connector [complete with O-rings]	D2056727
(3) O-rings for plug connector [set of 2]	D2056736
(3) Parker O-lube	D0010239
Tool for push to connect adapter [all]	10035756

Order No.: 10128090/01