## IMPORTANT SAFETY INSTRUCTIONS DO NOT DESTROY

**DANGER:** TO REDUCE RISKS OF BURNS,  
ELECTRIC SHOCK, FIRE, AND PERSONAL INJURY,  
THIS PRODUCT MUST BE USED IN ACCORDANCE  
WITH THE FOLLOWING INSTRUCTIONS:

1. READ ALL INSTRUCTIONS.
  2. DO NOT USE WHILE SLEEPING.
  3. DO NOT USE ON AN INFANT OR ON AN ANIMAL.
  4. THIS PAD IS NOT TO BE USED ON OR BY AN INVALID, A SLEEPING OR UNCONSCIOUS PERSON, A PERSON WITH DIABETES, OR A PERSON WITH POOR BLOOD CIRCULATION.
  5. DO NOT USE ON AREAS OF INSENSITIVE SKIN.
  6. BURNS CAN OCCUR REGARDLESS OF CONTROL SETTING. CHECK SKIN UNDER PAD FREQUENTLY.
  7. **DO NOT USE IN AN OXYGEN ATMOSPHERE.**
  8. **DO NOT USE PINS OR OTHER METALLIC MEANS TO FASTEN THIS PAD IN PLACE.**
  9. **DO NOT SIT ON OR AGAINST OR CRUSH PAD-AVOID SHARP FOLDS.**
  10. **NEVER PULL THIS PAD BY THE SUPPLY CORD.**
  11. **DO NOT USE THE CORD AS A HANDLE.**
  12. CAREFULLY EXAMINE INNER COVER BEFORE EACH USE. DISCARD THE PAD IF INNER COVERING SHOWS ANY SIGNS OF DETERIORATION.
  13. NEVER USE PAD WITHOUT THE REMOVABLE COVER IN PLACE.
  14. THIS UNIT SHOULD NOT BE USED BY CHILDREN WITHOUT ADULT SUPERVISION.
  15. PLACE PAD ON TOP OF AND NOT UNDER THE PART OF THE BODY NEEDING HEAT.
  16. USE THIS PAD ONLY ON A 110-120 VOLT AC CIRCUIT. UNPLUG WHEN NOT IN USE.
  17. DO NOT TAMPER WITH THIS PAD IN ANY WAY.
  18. THERE ARE NO USER-SERVICEABLE PARTS. IF FOR ANY REASON THIS PAD DOES NOT FUNCTION SATISFACTORILY, SEE WARRANTY FOR CONSUMER SERVICE INFORMATION.
  19. DO NOT USE THIS PAD WITH LINIMENT, SALVE OR OINTMENT PREPARATIONS, ESPECIALLY ONES THAT CONTAIN HEAT- PRODUCING INGREDIENTS. SKIN BURNS COULD RESULT.
  20. DO NOT BEND OR PINCH CORD AND DO NOT USE IF CORD IS DAMAGED.
  21. NEVER LEAVE THE APPLIANCE UNATTENDED, ESPECIALLY IF CHILDREN ARE PRESENT.
  22. LOOP CORD LOOSELY WHEN STORING. TIGHT WRAPPING MAY DAMAGE CORD AND INTERNAL PARTS.
  23. **SAVE THESE INSTRUCTIONS.**
- This appliance has a polarized plug (one blade is wider than the other). This plug fits a polarized outlet only one way; if the plug does not fit fully into the outlet, reverse the plug. If it still does not fit, contact a qualified electrician. **DO NOT ATTEMPT TO DEFEAT THIS SAFETY FEATURE.**

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FIGURE / FIGURA 2

### HOW TO USE YOUR HEATING PAD

1. Lay the cloth cover on a flat surface with the printing face up. Open the snap at the bottom.
2. Insert the pad into the cover completely, fitting one side at a time. (Figure 1) Make sure the cord and control remain on the outside of the cover.
3. The pad will lie flat, filling the entire cover. Close snap.
4. Plug into any 110-120 volt AC electrical outlet and set control to the desired heat.

**WARNING:** Never use the heating pad without the cloth cover.

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### FOR MODELS WITH HOT (RED) GEL PACK:

1. Remove gel pack from heating pad.
- a. Place gel pack in microwave.
- b. Make sure gel is spread evenly across pack for proper heating

**CAUTION: UNEVEN GEL COULD CAUSE HOT SPOTS**

- c. Heat at full power (700 Watts) for one minute. NOTE: Heat time may vary based on microwave wattage.
- d. Remove and check heat guide.

2. If not hot enough, return to microwave and heat in ten second intervals, checking Heat Guide until ready.

3. If Heat Guide reads TOO HOT, let cool before use.

### GEL PACK SAFETY REMINDERS:

- a. Gel Pack can be used without heating pad, but must be enclosed in cover.
- b. NEVER HEAT GEL PACK FOR MORE THAN 90 SECONDS AT ONE TIME OR GEL PACK MAY RUPTURE AND LEAK.
- c. NEVER BOIL GEL PACK.

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