
WEIGHT TRACKING SCALE


User Manual
Model HDM561

Congratulations on your purchase of this Health o meter® weight tracking scale. It is so much more than just a scale! You can monitor your weight, weight changes and compare it to your personal goal. It stores personal information for up to 4 different users. Also, you can save up to seven records for each user, and they can be completely private.

The Health o meter® weight tracking scale can help you in achieving your physical best by providing body weight readings and monitoring your weight over time.

KEYS AND ICONS


1. Weight Readout
2. Goal and difference in weight compared to previously saved weight
3. Unit of Measure: Pounds (LB) or Kilograms (KG)
4. User number and memory indicator
5. Information ready to be saved
6. Time/Date start weight and difference in weight compared to previously saved weight
7. Progress tracking stars
8. SET/SAVE button
9. SCROLL/HISTORY button
10. User buttons
11. Private Icon
12. Privacy code numbers


FIRST TIME SETUP

Back of the Scale

1. Select pounds (LB) or kilograms (KG).
2. Install four (4) AA batteries.
3. Place the scale on a flat surface. A hard surface is best for maximum accuracy.

ADJUSTING TIME AND DATE

1. Hold down the "SET/SAVE" button until the HOUR starts blinking.
2. Press or hold down the "SCROLL/HISTORY" button to choose the desired HOUR.
3. Press the "SET/SAVE" button. MINUTES will start blinking.
4. Press or hold down the "SCROLL/HISTORY" button to choose the desired MINUTES.
5. Press the "SET/SAVE" button. "AM" or "PM" will start blinking.


1

RESETTING AN INDIVIDUAL USER'S HISTORY DATA

This process will erase the records for the selected user.

NOTE: PLEASE MAKE SURE YOUR USER ID APPEARS ON THE UPPER RIGHT CORNER OF THE DISPLAY before you reset the user's data. By following this step, you will avoid erasing other user's data when resetting your user's memory.

1. Select the desired "USER" to be reset.
2. Once your user ID is shown, press and hold the "SET" button and then press the "User 1-2" button (if you are User 1 or User 2) or the "User 3-4" button (if you are User 3 or User 4) at the same time for 3 seconds.

Now the individual memory has been reset and the scale will prompt you to program the goal and privacy code.

3. See the PROGRAMMING/RESETTING YOUR WEIGHT GOAL and PRIVACY CODE sections on pages 2 and 3.


REPLACING THE BATTERIES

This unit includes four (4) AA batteries. When "Lo" appears on the display, replace all batteries.

Removing the batteries will reset the Date and Time but the recorded data will not be affected.


Battery Replacement

1. Locate the battery cover on the back of the scale. Push the locking tab and pull the battery cover up.
2. Use the tip of a pen to loosen each battery. Hold the scale at an angle to allow the batteries to slide out.
3. Insert four new batteries, pressing them firmly into their slot.
4. Replace the battery cover.
5. Set the Time/Date (see ADJUSTING TIME AND DATE).

6

6. Press the "SCROLL/HISTORY" button to choose "AM" or "PM".
7. Press the "SET/SAVE" button. Now your clock is set.

8. Press or hold down the "SCROLL/HISTORY" button to choose the desired MONTH.
9. Press the "SET/SAVE" button.
10. Press or hold down the "SCROLL/HISTORY" button to choose the desired DAY.
11. Press the "SET/SAVE" button.
12. Press or hold down the "SCROLL/HISTORY" button to choose the desired YEAR.
13. Press the "SET/SAVE" button. You have setup time and date.

PROGRAMMING/RESETTING YOUR WEIGHT GOAL

1. Press and hold the desired user button until "YES" is displayed.
2. Press the "SCROLL/HISTORY" button to select "YES" or "NO".


NOTE: You can set a goal (by selecting "YES") or just track your weight (by selecting "NO").

3. Press the "SET/SAVE" button to store your selection.
4. If you select "YES", enter your desired goal weight. Each digit is set separately. Scroll to choose your desired weight and press the "SET/SAVE" button to store it.

NOTE: The display will turn off automatically in 20 seconds if no key is pressed.

5. The display will flash the set goal and "SAVED" twice.

NOTE: Now that you have set your goal, the scale will prompt you to enter a privacy code. "YES" will start blinking along with the key icon. (This is optional).

2


COMPLETELY RESETTING THE SCALE

All data will be lost. Erased data cannot be recovered.

Make sure the LCD display is turned off.

1. Turn the scale over.
2. Using a pin, press and hold the button inside the small hole for 3 seconds. The display will show "C".

The scale's memory has been reset.


TROUBLESHOOTING

If "E" appears on the upper row of the display:

- You may have stepped on the scale before it displayed "0.0".
- The scale cannot detect a stable weight. Stand still while weighing.
- The weight has exceeded the scale's capacity (380 lb/172 kg).
- The privacy code has been entered incorrectly.

NOTE: The display will turn off automatically in 20 seconds if no key is pressed.

IT IS NORMAL FOR YOUR WEIGHT TO VARY DURING THE DAY AND FROM ONE DAY TO ANOTHER.

CARING FOR YOUR SCALE

Periodically wipe the scale surface down with a DRY cloth to remove dust.

Weight exceeding the capacity (380 lb/172 kg) may damage your scale. This product is designed for personal use only. Not for commercial or industrial use. This unit is not waterproof; avoid contact with excessive moisture.

7

PROGRAMMING / RESETTING THE PRIVACY CODE

To set up your privacy code you will be required to enter 3 numbers:

1. **YES** and " " will be shown on the display.
2. Press **SCROLL/HISTORY** button to select yes or no.

NOTE: You can set a privacy code (by selecting "YES") or just track your weight without a privacy code (by selecting "NO").

3. If "YES" is chosen, enter your desired 3 number code and press **SET/SAVE**.

NOTE: Each digit is set separately, until you have programmed your entire privacy code.


4. The privacy code and "SAVED" will flash saved twice, followed by " - - - ", then the scale will calibrate to "0.0" with either "LB" or "KG".

NOTE: If the Privacy Code is forgotten, the recorded data for the selected user will have to be erased.

USING YOUR WEIGHT TRACKING SCALE

A. MONITORING YOUR WEIGHT OVER TIME AND COMPARING IT TO YOUR GOAL.

1. Press one of the desired "USER" buttons to select a user. The display will show TIME and DATE; then it will calibrate " - - - " and will show "0.0" with "LB" or "KG".

3

YOUR HEALTHY LIFESTYLE

Health and fitness isn't just about counting calories and cutting fat out of your diet. It is about a balanced combination of eating right and exercising; neither dieting nor exercising alone will give you the kind of results you get when they are used together. To achieve a healthy lifestyle, you must set realistic fitness goals and reduce the number of calories you take in, regardless of whether they come from foods high in protein, carbohydrates, etc.

When you go on a diet, you may lose weight, but you'll also lose beneficial muscle and bone tissue. Instead, reduce the number of calories you are taking in by eating more healthful foods, and include exercise in your regular routine to lose weight faster by burning fat. Exercise makes you feel better; it helps improve circulation, relieves stress, and speeds up your metabolism. So create a plan of exercise and diet appropriate for your needs, body type, and general health. Soon you'll see how this powerful combination makes it easier to look and feel your best.

How does extra weight affect my health?


Extra weight can put you at higher risk for some health problems such as:

- Type 2 diabetes (high blood sugar)
- High blood pressure
- Heart disease and stroke
- Some types of cancer
- Sleep apnea (when breathing stops for short periods during sleep)

Consult with your health care professional for more information.

Helpful Eating Tips

- Make sure no more than 30% of your daily calories come from fat, with no more than 10% saturated fat.
- Never skip meals! Instead, eat sensible portions at every meal and use small, low calories snacks between meals to help control hunger.
- Learn to read food labels and follow recommended serving sizes.
- Balance your calorie intake with your activity level. Eat less on low activity days than you do on days where you enjoy rigorous exercise.
- Follow the USDA's Food Pyramid as your guide to a healthy, balanced eating plan.

Helpful Exercise Tips

- Set goals. Write down what you want to achieve in the short and long term.
- Make exercise a habit. Choose a time every day - even 30 minutes - and stick to it.
- Vary your workout. Mix up your activities to avoid burn out.
- Recruit a friend! Working out as a team keeps you both motivated.

NOTE: These materials are not intended to replace advice from your doctor or fitness professional. Please consult with your physician before beginning any fitness program or fat or weight reduction program. Individual weight loss will vary. **Health o meter® takes no responsibility for individual results, or any claim made by a third party.**

8

NOTE: If you have set a privacy code, the scale will prompt you to enter your code at this time. If no code is entered or the wrong code is entered 3 times in a row, the scale will automatically turn off.

2. Step on the scale and position your feet in a manner that will distribute your weight evenly over the entire scale platform. Stand still while reading the display. Three screens will be displayed, one at a time.
3. Press the "SET/SAVE" button to save your reading.

B. MONITORING YOUR WEIGHT OVER TIME WITHOUT SETTING YOUR GOAL

1. Press the desired "USER" button. The display will show TIME and DATE; then it will calibrate (" - - ") and will show "0.0" with either "LB" or "KG".

NOTE: If you have set a privacy code, the scale will prompt you to enter your code at this time. If no code is entered, the scale will automatically turn off.



2. Step on the scale and position your feet in a manner that will distribute your weight evenly over the entire scale platform. Stand still while reading the display.
3. Press the "SET/SAVE" button to save your reading.

C. WEIGHING WITHOUT SELECTING A USER BUTTON

1. Press the lower corner of the scale platform. The scale will calibrate " - - - " and then will display "0.0" and either "LB" or "KG".
2. Step on the scale and position your feet in a manner that will distribute your weight evenly over the entire scale platform. The display will flash " - - - " several times for approximately 3 seconds before displaying your weight. Stand still while reading the display.

If you would like to save your reading (Only if your profile has been set):

1. STEP OFF THE SCALE.
2. Press the "USER" button where you want to save the information.
3. Enter your privacy code (if set).
4. Press the "SET/SAVE" button to save it.

4

5 YEAR LIMITED WARRANTY

PLEASE DO NOT RETURN THIS PRODUCT TO ANY OF THE ADDRESSES SHOWN OR TO THE PLACE OF PURCHASE.

Sunbeam Products, Inc. ("Sunbeam") warrants that for a period of five years from the date of purchase, this product will be free from defects in material and workmanship. Sunbeam, at its option, will repair or replace this product or any component of the product found to be defective during the warranty period. Replacement will be made with a new or remanufactured product or component. If the product is no longer available, replacement may be made with a similar product of equal or greater value. This is your exclusive warranty. This warranty is valid for the original retail purchaser from the date of initial retail purchase and is not transferable. Keep the original sales receipt. Proof of purchase is required to obtain warranty performance. Sunbeam dealers, service centers, or retail stores selling Sunbeam products do not have the right to alter, modify or any way change the terms and conditions of this warranty. This warranty does not cover normal wear of parts or damage resulting from any of the following: negligent use or misuse of the product, use on improper voltage or current, use contrary to the operating instructions, disassembly, repair or alteration by anyone other than Sunbeam or an authorized Sunbeam service center. Further, the warranty does not cover: Acts of God, such as fire, flood, hurricanes and tornadoes.

What are the limits on Sunbeam's Liability?

Sunbeam shall not be liable for any incidental or consequential damages caused by the breach of any express, implied or statutory warranty or condition.

Except to the extent prohibited by applicable law, any implied warranty or condition of merchantability or fitness for a particular purpose is limited in duration to the duration of the above warranty.

Sunbeam disclaims all other warranties, conditions or representations, express, implied, statutory or otherwise. Sunbeam shall not be liable for any damages of any kind resulting from the purchase, use or misuse of, or inability to use the product including incidental, special, consequential or similar damages or loss of profits, or for any breach of contract, fundamental or otherwise, or for any claim brought against purchaser by any other party.

Some provinces, states or jurisdictions do not allow the exclusion or limitation of incidental or consequential damages or limitations on how long an implied warranty lasts, so the above limitations or exclusion may not apply to you.

This warranty gives you specific legal rights, and you may also have other rights that vary from province to province, state to state or jurisdiction to jurisdiction.

How to Obtain Warranty Service

If you have any question regarding this warranty or would like to obtain warranty service, please call 1-800-672-5625, or visit us at www.healthometer.com and a convenient service center address will be provided to you.

In the U.S.A. this warranty is offered by Sunbeam Products, Inc. located in Boca Raton, Florida 33431.

PLEASE DO NOT RETURN THIS PRODUCT TO ANY OF THE ADDRESSES SHOWN OR TO THE PLACE OF PURCHASE.


© 2005 Sunbeam Products, Inc. All rights reserved. Health o meter® is a trademark of Sunbeam Products, Inc. Boca Raton, FL 33431.
Printed in China.

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