

CVWorks

Standard Applications


2D FIESTA Cine

The workhorse of functional cardiac imaging, 2D FIESTA Cine can be used as a gated functional sequence or non-gated quick localizer to assess the heart's anatomy and function by providing excellent tissue contrast between blood pool, myocardium and valves.


Clinical benefits:

- Qualitative assessment of valvular structure and anatomy
- Aids in evaluation of Arrhythmogenic Right Ventricular Dysplasia (ARVD) and tumors
- Used to calculate LV and RV function and cardiac output in tools such as cmr⁴² and Arterys™

4-chamber



Short axis



Long axis

