



“Finally!”

AMX 115

A 100mA, 200kHz, easy-to-use X-ray

AMX 115 is a commercial configuration of the Brivo XR115



“Fantastic!”



A 100mA X-ray
that promises **great imaging** across anatomies

It is also 40% faster!

This translates to:

- Shorter breath hold
- Lesser patient movement
- Great image quality
- Better patient comfort


AMX 115 helps you easily manage tough cases, such as when your patient is obese or breathless, and in trauma or pediatric cases where uncontrolled patient movement is common.

*6mAs= 60mA x 0.1 sec
*6mAs= 100mA x 0.06 sec


AMX 115 is a commercial configuration of the Brivo XR115

“Superb!”


An **easy-to-use** X-ray




180° rotation



Light, therefore easy to move



Minimum footprint



Easy sliding for convenience

AMX 115 is a well-designed and low-maintenance equipment:

- A less than 2kgf push force – 50% easier to move than most other models in the market
- A 10% smaller footprint than most other models in the market – convenience in narrow spaces
- A weight of less than 120kg – one of the lightest models in the market
- A unique mechanical design for ease-of-use and quicker exams with better SID accuracy

AMX 115 is a commercial configuration of the Brivo XR115

“Awesome!”

A 200kHz X-ray for great imaging and
efficient dosage


This results in:

- Uniform and useful radiation
- Better image quality
- Low scatter radiation
- Focus on patient and technician/physician dose reduction

The AMX 115 200kHz generator is so fast that the X-ray emission is extremely uniform with an efficient X-ray spectrum!



Ultra-high frequency
Lesser unwanted radiation



200kHz benefits
More efficient X-ray
Consistent IQ

“Wow!”

An X-ray that is **ergonomically designed**
for added safety

Cassette cabinet



Accessories holder



Two-sided brakes



Ergonomic handles

AMX 115

is designed ergonomically for safety, comfort and aesthetics:

- Rounded edges and internal cabling
- Minimal pinch points and a 'clean' look
- Ergonomically designed and positioned user touch points (handles, angulated console, or closed cassette cabinet) for greater comfort
- An accessories holder for pens, strips or tapes



“Thank God!”

An X-ray that brings you complete
peace of mind



© 2020 General Electric Company.

AMX 115 is a commercial configuration of the Brivo XR115

GE Healthcare reserves the right to make changes in specifications and features shown herein, or discontinue the product described at any time without notice or obligation. Contact your GE Healthcare representative for the most current information. GE Medical Systems, Inc., doing business as GE Healthcare, a division of General Electric Company. GE and the GE Monogram are trademarks of General Electric Company.

May 2020