



GE HealthCare

AW Family

Advanced Visualization Applications & Platform Compatibility

The AW Family is GE Healthcare's ecosystem of processing and advanced visualization solutions across platforms and applications. This document is intended to overview the compatibility of select AW Family platforms and applications.



AW Family Advanced Visualization Application & Platform Compatibility Overview

Applications listed below are compatible with the following platforms:

- AW4.7 Ext. 18 and higher
- AW Server 3.2 Ext. 4.0 and higher

The applications below are based on the Volume Viewer foundation, therefore access to Volume Viewer (Version 16.0 Ext 4) is a pre-requisite for accessing the following:

- | | |
|-------------------------------|---------------------------|
| • 3D Suite | • FlightPlan for Liver |
| • AdvantageCTC Pro3D EC | • GenIQ |
| • AutoBone Xpress | • GSI Viewer |
| • Bone VCAR | • Hepatic VCAR |
| • CardEP | • Integrated Registration |
| • CardIQ Flow | • Liver Suite |
| • CardIQ Function Xpress | • Lung VCAR |
| • CardIQ Fusion PET | • MR VessellIQ Xpress |
| • CardIQ Fusion SPECT | • OncoQuant |
| • CardIQ Xpress 2.0 | • PET VCAR |
| • Colon VCAR EC | • PROView |
| • CT Perfusion 4D Multi-Organ | • READY View |
| • CT Perfusion 4D Neuro | • Stroke VCAR |
| • DynamicIQ | • TAVI Analysis |
| • Dynamic Shuttle | • Thoracic VCAR |
| • FastStroke | • VessellIQ Xpress |
| • FlightPlan for Embolization | • Volume Viewer |
| • FlightPlan for EVAR | • Volume Viewer Innova |



GE HealthCare

About GE HealthCare

GE HealthCare is a leading global medical technology, pharmaceutical diagnostics, and digital solutions innovator, dedicated to providing integrated solutions, services, and data analytics to make hospitals more efficient, clinicians more effective, therapies more precise, and patients healthier and happier. Serving patients and providers for more than 100 years, GE HealthCare is advancing personalized, connected, and compassionate care, while simplifying the patient's journey across the care pathway. Together our Imaging, Ultrasound, Patient Care Solutions, and Pharmaceutical Diagnostics businesses help improve patient care from prevention and screening, to diagnosis, treatment, therapy, and monitoring. We are an \$18 billion business with 51,000 employees working to create a world where healthcare has no limits.

Follow us on [Facebook](#), [LinkedIn](#), [Twitter](#), [Instagram](#) and [Insights](#) for the latest news, or visit our website [gehealthcare.com](https://www.gehealthcare.com) for more information.

Products mentioned in the material may be subject to government regulations and may not be available in all countries. Shipment and effective sale can only occur after approval from the regulator. Please check with local GE HealthCare representative for details.

GE and the GE monogram are trademarks of General Electric Company. Used under trademark license.

©2023 GE HealthCare.

January 2023

DOC2674331