

Clinical Accessories

Trusted. Reliable. Simple.



Contents

PATIENT MONITORING

ECG

Monitoring Recording Paper

Trunk Cables

Leadwires

Trunk Cable Accessories

Pulse Oximetry

Sensors & Cables

Temperature

Cables & Probes

Invasive Blood Pressure

Cables

Cuffs

Hoses

DINACLICK

SOFT-CUF

CLASSIC-CUF

DURA-CUF

YELLOW CLASSIC-CUF

RADIAL-CUF

Neonatal

Cardiac Output

Cables

Other Accessories

Adequacy of Anesthesia

Cable

Sensors

Neurology

EEG Leadsets

EEG Electrodes

AEP Accessories

Patient Spirometry

Sensors

Tubes

Kits

Gas Monitoring

Airway Adapter Sets

Gas Exhaust Lines

Gas Sampling Lines

Airway Adapters

Water Traps

4 ANESTHESIA & 6 RESPIRATORY CARE

Absorbents

Absorbers

Canisters

Dust Filter

Vaporizer Bottle Adapters

Filters & Breathing Masks

Reservoir Bags & Gas Scavenging

Anesthesia Breathing Bags

Scavenging Bags, Tubes, & Connectors

Patient Breathing Circuits

Circuit Startup Kits

Circuit Kits

Circuits

Circuit Components

Flow Sensors – Respiratory

Respiratory Flow Sensors and

Expiratory Valve Assemblies

Flow Sensors – Anesthesia

Nebulizer Accessories

Anesthesia O2 Cells

Respiratory Safety Guard

Pressure Catheters

DIAGNOSTIC CARDIOLOGY

MAC™, CASE™ and CardioSoft™

ECG

Cables

Leadwires

Accessories

Recording Paper

Electrodes

Holter Monitors

ECG – CM 4000

Leadwire Sets

Accessories

ECG – SEER® 12

Leadwire Sets

Accessories

ECG – SEER® 1000

Leadwire Sets

Accessories

ECG – SEER® Light

Leadwire and Cable Kits

Accessories

MATERNAL INFANT CARE

Maternal Fetal Care

Mini Telemetry

Recording Thermal Paper

Remote Event Marker

Novii™

Novii™ Patch

Cables

Corometrics™

Cables

Transducer Belts

Ultrasound Transducers

Internal Monitoring

Cables

Electrode

Pad

Infant Care

Giraffe™/Panda™

Air filters

Clip

Cover

Phototherapy Accessories

Pad Covers

Resuscitation

T-Piece

Lullaby™ Warmer

Temperature Probes

Temperature Probe Patches

Shuttle™

Patient Restraint

Welcome

GE Healthcare offers you a simple, reliable one-stop solution for high-quality GE Healthcare-compatible clinical accessories. We provide an expanded inventory of products that have been verified by our engineers. The compatible devices listed may not be available for sale and/or cleared/approved in all countries. Check with your representative for product availability.

PATIENT MONITORING



PATIENT MONITORING

ECG

- Leadwires
- Recording Paper
- Trunk Cable
- Trunk Cables Accessories

Pulse Oximetry

- Cables
- Clip
- Other Accessories
- Probes
- Sensors
- Tape & Wrap

Temperature

- Cables
- Probes

Invasive Blood Pressure

- Cables

Cuffs

- Hoses
- DINAClick
- SOFT-CUFF
- DURA-CUF
- CLASSIC-CUF
- YELLOW CLASSIC-CUF
- RADIAL-CUF
- Neonatal

Cardiac Output

- Cables
- Catheter Kits
- Other Accessories

Adequacy of Anesthesia

- Cable
- Electrodes
- Sensors

Neurology

- AEP Accessories
- EEG Electrodes
- EEG Leadsets

Patient Spirometry

- Kits
- Sensors
- Tubes

Gas Monitoring

- Airway Adapter Sets
- Airway Adapters
- Gas Sampling Lines
- Water Traps

ECG



MONITORING RECORDING PAPER



MS Recording Thermal Paper						
Length	CARESCAPE™ B450/B650 Monitors, CARESCAPE™ B850 Monitor with PRN 50-M+ printer	B20(i)/B40(i) Monitors	B1x5 Monitors	CARESCAPE™ V100 and VC150 Monitors	Part Number	Quantity
57.2 mm x 24.4 m				X	2106825-001	10
50 mm x 25.6 m	X	X	X		2106824-001	20
50 mm x 27.4 m	X	X	X		2106823-001	48

For ECG paper for Diagnostic cardiology and Maternal infant care equipment, please go to pages 58 and 68.

Not available in all countries





2106305-003



2106310-002

ECG Trunk Cables – IEC

Description	Length	Part Number	Quantity
Neonatal DIN 3-lead	3.6 m	2106306-003	1
	1.2 m	2106306-004	1
3/5-lead	3.6 m	2106305-003	1
	1.2 m	2106305-004	1
	6 m	2106311-002	1
3/5-lead w/ESU filter	3.6 m	2106310-002	1
6-lead	3.6 m	2106307-003	1
	1.2 m	2106307-004	1
12-lead	3.6 m	2106308-003	1
	1.2 m	2106308-004	1



2106307-003



2106308-004

Product compatibility

- CARESCAPE™ Monitors B450/B650/B850 with CARESCAPE™ Patient Data Module (PDM) or Patient Side Module E-PSM
- B105/B125 patient monitors
- B20/B40 and B20i/B40i patient monitors
- Dash™ 2000/2500/3000/4000/5000 patient monitors
- Solar™ Monitor with CARESCAPE™ Patient Data Module (PDM)
- ECG trunk cables are also compatible with GE Healthcare ECG leadwire sets
- ApexPro™ CH telemetry

Specifications available on request



ECG Trunk Cables – IEC

Description	Length	Part Number	Quantity
Neonatal trunk cable	3.6 m	2080557-003	1
	1.2 m	2080557-004	1

Not available in all countries

ECG LEADWIRES – GRABBER



ECG Leadwires Grabber – IEC Compliant

Description	Connection	Length	Part Number	Quantity
3-lead	Molded	74 cm*	2106390-003	1
		130 cm	2106390-004	
5-lead	Molded	74 cm*	2106389-003	1
		130 cm	2106389-004	
		mix 74 cm, 130 cm	2106389-006	
	Non-molded	74 cm *	2106391-003	
		130 cm	2106391-004	
		mix 74 cm, 130 cm	2106391-006	
5-lead C2-6	Non-molded	74 cm	2106393-004	1
		130 cm	2106393-002	
6-lead	Non-molded	74 cm*	2106398-001	1
		130 cm	2106398-002	
		mix 74 cm, 130 cm	2106398-003	

Product compatibility

- CARESCAPE™ Monitors B450/B650/B850 with CARESCAPE™ Patient Data Module (PDM) or Patient Side Module (E-PSM)
- B105/B125 patient monitors
- B20/B40 and B20i/B40i patient monitors
- Dash™ 2000/2500/3000/4000/5000 patient monitors
- Solar™ monitors with CARESCAPE™ Patient Data Module (PDM)
- ECG trunk cables are also compatible with GE Healthcare ECG leadwire sets
- ApexPro™ CH telemetry
- CARESCAPE™ ONE Monitor

Specifications available on request


*74cm 3-, 5-, and 6-leadwire sets are compatible with ApexPro™

Not available in all countries





2106385-004



2106383-006



2106381-004



2106394-002

ECG Leadwires Snap – IEC Compliant

Description	Connection	Length	Part Number	Quantity
3-lead	Molded	74 cm*	2106385-003	1
		130 cm	2106385-004	
5-lead	Molded	74 cm*	2106383-003	1
		130 cm	2106383-004	
		mix 74 c m, 130 cm	2106383-006	
	Non-molded	74 cm*	2106381-003	
		130 cm	2106381-004	
5-lead C2-6	Non-molded	130 cm	2106394-004	1
		74 cm	2106394-002	
6-lead	Non-molded	74 cm*	2106396-001	1
		130 cm	2106396-002	

Product compatibility

- CARESCAPE™ Monitors B450/B650/B850 with CARESCAPE™ Patient Data Module (PDM) or Patient Side Module (E-PSM)
- B105/B125 patient monitors
- B20/B40 and B20i/B40i patient monitors
- Dash™ 2000/2500/3000/4000/5000 patient monitors
- Solar™ monitors with CARESCAPE™ Patient Data Module (PDM)
- ECG trunk cables are also compatible with GE Healthcare ECG leadwire sets
- ApexPro™ CH telemetry
- CARESCAPE™ ONE Monitor

Specifications available on request

*74cm 3-, 5-, and 6-leadwire sets are compatible with ApexPro™

Not available in all countries





ECG Leadwires – IEC Compliant				
Description	Connection	Length	Part Number	Quantity
4-lead	Snap	74 cm/29 in	2066468-002	1
		130 cm/51 in	2066468-006	
	Grabber	74 cm/29 in	2066468-004	1
		130 cm/51 in	2066468-008	

Product compatibility

- CARESCAPE™ ONE Monitor

Specifications available on request

Not available in all countries



2106403-003



2106403-001



2106399-001



2106399-002

ECG Trunk Cable Accessories

Description		Part Number	Quantity
DIN Converter	3-Lead	2106403-003	1
	5-Lead	2106403-001	1
Leadwire separator		2106399-001	3
Bedsheet clip		2106399-002	1

Product compatibility

- CARESCAPE™ Monitors B450/B650/B850 with CARESCAPE™ Patient Data Module (PDM) or Patient Side Module (E-PSM)
- B105/B125 patient monitors
- B20/B40 and B20i/B40i patient monitors
- Dash™ 2000/2500/3000/4000/5000 patient monitors
- Solar™ monitors with CARESCAPE™ Patient Data Module (PDM)
- ECG trunk cables are also compatible with GE Healthcare ECG leadwire sets
- ApexPro™ CH telemetry

Specifications available on request

Not available in all countries



Pulse Oximetry



TRUSIGNAL™

GE Healthcare provides compatibility with three major brands of Pulse Oximetry (SpO₂) measurement to meet all of your needs:

- GE Healthcare TruSignal™
- MASIMO
- Nellcor



TS-F2-GE



TS-SA4-GE



TS-E2-GE

TruSignal™ Integrated Reusable Sensors

Description	Connector	Patient	Patient Weight	Length	Part Number	Quantity
Finger	GE	Adult, Pediatric	≥ 20 kg	2 m	TS-F2-GE	1
				4 m	TS-F4-GE	1
	Ohmeda			1 m	TS-F1-H	1
	Datex			4 m	TS-F4-N	1
FingerTip	GE	Adult	> 30 kg	4 m	TS-SA4-GE	1
PediTip		Pediatric	15-30 kg	3 m	TS-SP3-GE	1
Ear	GE	Adult, Pediatric	≥ 10 kg	2 m	TS-E2-GE	1
				4 m	TS-E4-GE	1
	Ohmeda				TS-E4-H	1
	Datex				TS-E4-N	1

Not available in all countries





TS-F-D



TS-SA-D/TS-SP-D



TS-E-D



TS-W-D



TS-SE-3



TruSignal™ Reusable Sensors

Description	Patient	Patient Weight	Length	Part Number	Quantity
Finger	Adult, Pediatric	≥ 20 kg	1 m	TS-F-D	1
PediTip	Pediatric	15-30 kg		TS-SP-D	1
Fingertip	Adult	> 30 kg		TS-SA-D	
Ear	Adult, Pediatric	≥ 10 kg		TS-E-D	1
Wrap	Adult, Pediatric, Infant	≥ 3 kg		TS-W-D	1
Sensitive Skin	Adult, Pediatric, Infant, Neonate	No limit		TS-SE-3	3


Requires Compatible Interconnect Cable



TS-G3

TruSignal™ Interconnect Cables

Connector	Length	Part Number	Quantity
GE	3 m	TS-G3	1
Ohmeda		TS-H3	
Datex		TS-N3	



TruSignal™ Disposable Sensors					
Description	Patient	Patient Weight	Length	Part Number	Quantity
Adult/Pediatric	Pediatric, infant	≥ 20 kg	0.3 m	TS-AP-10	10
				TS-AP-25	25
AllFit	Adult, pediatric, infant, neonate	No limit	0.5 m	TS-AF-10	10
				TS-AF-25	25
Adhesive Wrap	Adult	> 20 kg	0.3 m	TS-AAW-10	10
				TS-AAW-25	25
	Pediatric, infant	3-20 kg		TS-PAW-10	10
				TS-PAW-25	25

Requires Compatible Interconnect Cable

TruSignal™ Sensor Accessories				
Description	Use With	Size	Part Number	Quantity
Replacement wide adhesive tape	TruSignal™ Sensitive Skin Sensor, TruSignal™ Wrap Sensor	–	TS-RTW-100	100
Replacement adhesive tape (blue)	TruSignal™ AllFit Sensor	–	TS-RT-100	100
Replacement adhesive tape (teddy bears)			TS-RTB-100	
Replacement foam wrap	TruSignal™ Sensitive Skin Sensor	Small	TS-RWS-24	24
		Medium	TS-RWM-24	
	TruSignal™ Wrap Sensor	Large	TS-RWL-24	
Infant foam sandal	TruSignal™ Sensitive Skin Sensor	–	TS-SND-3	3
Replacement Headband	TruSignal™ Ear Sensor	–	TS-HB-3	3



Product compatibility

- CARESCAPE™ Monitors B450/B650/B850 with CARESCAPE™ Patient Data Module (PDM) or Patient Side Module (E-PSM)
- B105/B125 patient monitors
- B20/B40 and B20i/B40i patient monitors
- Dash™ 2000/2500/3000/4000/5000 patient monitors
- Solar™ monitors with CARESCAPE™ Patient Data Module (PDM)
- ECG trunk cables are also compatible with GE Healthcare ECG leadwire sets
- ApexPro™ CH telemetry
- CARESCAPE™ ONE Monitor

Specifications available on request

Not available in all countries



For CARESCAPE™ B450 / B650 / B850 monitors with a CARESCAPE™ Patient Data Module (PDM) with Masimo™ SpO2 technology.



2027263-002


LNCS and RD Interconnect Cables			
Product	Length	Part Number	Quantity
Masimo LNC-10 interconnect cable	3 m	2027263-002	1
Masimo RD SET GE cable	3,6 m	MASIMO-RD-4085	1



Masimo-RD-4051

Reusable LNCS and RD Sensors					
Sensors	Manufacturer Reference	Patient	Patient Weight	Part Number	Quantity
Masimo Finger sensor	LNCS DC-I	Adult	> 30 kg	2027258-001	1
	LNCS DC-IP	Pediatric	10-50 kg	2027259-001	
Masimo TipClip Ear sensor	LNCS TC-I	Adult	> 30 kg	2027261-001	1
Masimo Finger sensor	RD SET DC-I	Adult	> 30 kg	MASIMO-RD-4050	1
	RD SET DC-IP	Pediatric	10-50 kg	MASIMO-RD-4051	

Requires Compatible Interconnect Cable



2027253-001



MASIMO-RD 4000



MASIMO-RD 4002



MASIMO-RD 4003



MASIMO-RD 4004


Disposable Sensors					
Sensors	Manufacturer Reference	Patient	Patient Weight	Part Number	Quantity
Adhesive LNCS sensors	LNCS NEO	Neonate	< 3 kg	2051069-001	20
	LNCS NEOPT	Neonate	< 1 kg	2051070-001	
	LNCS Adtx	Adult	> 30 kg	2027253-001	
Digit RD SET sensor	RD SET Adt	Adult	> 30 kg	MASIMO-RD 4000	20
	RD Set Inf	Infant	3-20 Kg	MASIMO-RD 4002	
	RD Set Neo	Neonate/Adult	< 3 or > 40 kg	MASIMO-RD 4003	
	RD SET NeoPt	Neonate	< 1 kg	MASIMO-RD 4004	
Digit RD Rainbow SET Sensor	RD Rainbow Set -2 Adt	Adult	> 30 kg	MASIMO-RD 4026	20
	RD Rainbow Set -2 Pdt	Pediatric	10-50 kg	MASIMO-RD 4027	
	RD Rainbow Set -2 Inf	Infant	3-20 Kg	MASIMO-RD 4028	
	RD Rainbow Set -2 Neo	Neonate/Adult	< 3 or > 40 kg	MASIMO-RD 4029	

Requires Compatible Interconnect Cable

These accessories may only be available from GE Healthcare when purchased with compatible GE Healthcare medical devices. Please consult with your local GE Healthcare representative if purchasing replacement accessories.

Not available in all countries

Adhesive Tape		
Product	Part Number	Quantity
Masimo adhesive tape for SpO2 adhesive sensor (RD SET Neo)	MASIMO-RD-4047	102



Reusable LNCS Sensors					
Sensors	Manufacturer Reference	Patient	Patient Weight	Part Number	Quantity
Masimo Finger sensor	LNCS DC-I	Adult	> 30 kg	2027258-001	1
	LNCS DC-IP	Pediatric	10-50 kg	2027259-001	
Masimo TipClip Ear sensor	LNCS TC-I	Adult	> 30 kg	2027261-001	1



2106257-001

REUSABLE LNCS CABLES				
Cables	Manufacturer Reference	Length	Part Number	Quantity
MASIMO LNCS CABLE	RED 25 LNC-10	3 m	2106257-001	1
	RED 25 LNC-4	1.2 m	2106257-002	
	RED 25 LNC-1	30 cm	2106257-003	



MASIMO-RD-4074



MASIMO-RD-4076



MASIMO-RD-4106



MASIMO-RD-4078


Masimo RD Cable				
Cables	Manufacturer Reference	Length	Part Number	Quantity
Masimo RD Cable	RD Rainbow SET 25R-1.5	0.45 m	MASIMO-RD-4074	1
	RD Rainbow SET 25R-05	1.5 m	MASIMO-RD-4076	
	RD Rainbow SET 25R-08	2.4 m	MASIMO-RD-4106	
	RD Rainbow SET 25R-12	3.6 m	MASIMO-RD-4078	

These accessories may only be available from GE Healthcare when purchased with compatible GE Healthcare medical devices. Please consult with your local GE Healthcare representative if purchasing replacement accessories.

Not available in all countries



For CARESCAPE™ B450 / B650 / B850 monitors with a CARESCAPE™ Patient Data Module (PDM) with Nellcor™ SpO2 technology.



Nellcor™ Interconnect Cables

Description	Connector	Length	Part Number	Quantity
OxiMax™ SpO2 Interconnect Cable	GE	1.2 m	2021406-002	1
		3 m	2021406-001	1
	GE for Corometrics™ only	2.9 m	2025350-001	1
OxiSmart SpO2 Interconnect Cable	GE	2.9 m	2006644-001	1
Interconnect Cables*	DOC-4	1.2 m	2059000-001	1
	DOC-8	2.4 m	2059001-001	1
	DOC-10	3 m	2008773-001	1

* These accessories are sold ONLY with the equipment

414248-001



70124033



70124021



70124035



Nellcor™ Reusable Sensors

Sensors	Patient	Patient Weight	Nellcor™ Reference	Part Number	Quantity
DuraSensor Finger Probe	Adult	< 30 kg	DS100A-1	70124021	1
Multisite SpO2 Sensor	Adult	> 1 kg	D-YS	70124033	1
Two-piece SpO2 sensor	Neonatal / Adult	< 3 kg / > 40 kg	OXI-A/N	70124035	1 (delivered with 24 adhesive strips)
	Pediatric, infant	3-40 kg	OXI-P/I	414248-001	1 (delivered with 24 adhesive strips)

Requires Compatible Interconnect Cable

These accessories are sold ONLY with the equipment

Not available in all countries

For CARESCAPE™ B450 / B650 / B850 monitors with a CARESCAPE™ Patient Data Module (PDM) with Nellcor™ SpO2 technology.



70124027



70124032

Nellcor™ Disposable Sensors

Sensors	Nellcor™ Reference	Patient	Patient Weight	Part Number	Quantity
Adhesive SpO2 sensor	MAX-A	Adult	> 30 kg	70124027	24
Long Finger Adhesive sensor	MAXAL			2028117-001	24
Adhesive SpO2 sensor	MAXP	Pediatric	10-50 kg	70124022	24
	MAXI	Infant	3-20 kg	70124026	24
	MAX-N	Adult	> 40 kg	70124032	24
		Neonate	< 3 kg		

Requires Compatible Interconnect Cable
These accessories are sold ONLY with the equipment



Nellcor™ Fixation Clip

Description	Nellcor™ Reference	For Use With	Part Number	Quantity
SpO2 Ear-Clip	D-YSE	D-YS	70124034	1

These accessories are sold ONLY with the equipment



Nellcor™ Sensor Adhesive Tape

Product	Patient	Nellcor™ Reference	Use With	Part Number	Quantity
Sensor Wrap Replacement Tape	Adult/Neonatal	ADH-A/N	OXI-A/N, D-YS	2016130-001	100
	Pediatric/Infant	ADH-P/I	OXI-P/I, D-YS	2016131-001	100

These accessories are sold ONLY with the equipment

Not available in all countries



Temperature



2104178-001



2104178-002



2104179-001

Temperature Cables

Description	Connection	Length	Part Number	Quantity
400 Series for Reusable Probes	Single	0.5 m	2104178-001	1
	Dual	0.5 m	2104178-002	1
400 Series for disposable Probes	Single	3.6 m	2104179-001	1

Product compatibility

- CARESCAPE™ Patient Data Module (PDM), Patient Side Module E-PSMP and E-PT modules
- B105/125 patient monitors
- B20/B40 and B20i/B40i patient monitors

Not available in all countries



TEMPERATURE



2107176-032



2107176-033

Temperature Probes Reusable					
Description	Connection	Matching Code on the Label	Length	Part Number	Quantity
Esophageal/rectal	Adult	1010132011GE	280 cm	2107176-031	1
	Pediatric	1010132010GE	280 cm	2107176-032	1
Skin	Adult	1010132023GE	280 cm	2107176-033	1



Disposable Temperature Probes			
Product	Patient	Part Number	Quantity
Disposable Temperature Probe – General Purpose, 400 series, 9 Fr	Pediatric	M1024229VY	50
Disposable Temperature Probe – General Purpose, 400 series, 12 Fr	Adult	M1024231VY	
Disposable Temperature Probe – Skin, 400 series	Adult - Infant	M1024222VY	

Product compatibility

- CARESCAPE™ Monitors B450/B650/B850 with CARESCAPE™ Patient Data Module (PDM), Patient Side Module E-PSMP and E-PT modules
- B105/125 patient monitors
- B20/B40 and B20i/B40i patient monitors
- Dash™ 3000/4000/5000 patient monitors
- Solar™ 8000M with TRAM module

Specifications available on request

Not available in all countries



Invasive Blood Pressure



Invasive Blood Pressure Cables

Description	Connection	Length	Part Number	Quantity
Argon™	Single	3.6 m	2104166-001	1
		1.2 m	2104166-002	
	Dual	3.6 m	2104166-003	1
		1.2 m	2104166-004	
Edwards™ LS TruWave™	Single	3.6 m	2104162-001	1
		1.2 m	2104162-002	
	Dual	3.6 m	2104162-003	1
		1.2 m	2104162-004	

Product compatibility


- CARESCAPE™ Monitors B450/B650/B850 with CARESCAPE™ Patient Data Module (PDM), Patient Side Module E-PSMP, E-PP, E-PT, E-COP(SV) and E-PiCCO™ modules
- CARESCAPE ONE Monitor
- B105/B125 patient monitors
- B20/B40 and B20i/B40i patient monitors

Specifications available on request

Not available in all countries



INVASIVE BLOOD PRESSURE CABLES



ICU Medical Transpac-IV® - Dual



Utah/Biotrans™ - Single



Utah/Biotrans™ - Dual



Dual Adapter Cable

Invasive Blood Pressure Cables				
Description	Connection	Length	Part Number	Quantity
ICU Medical Transpac-IV®	Single	3.6 m	2104158-001	1
		1.2 m	2104158-002	
	Dual	3.6 m	2104158-003	1
		1.2 m	2104158-004	
Utah/Biotrans™	Single	3.6 m	2104107-001	1
		1.2 m	2104107-002	
	Dual	3.6 m	2104107-003	1
		1.2 m	2104107-004	
Dual Adapter Cable	–	30 cm	2103966-001	1

Product compatibility

- CARESCAPE™ Monitors B450/B650/B850 with CARESCAPE™ Patient Data Module (PDM), Patient Side Module E-PSMP, E-PP, E-PT, E-COP(SV) and E-PiCCO™ modules
- B105/B125 patient monitors
- B20/B40 and B20i/B40i patient monitors

Specifications available on request

Not available in all countries

Cuffs



Hoses					
Description	Patient	Compatibility	Length	Part Number	Quantity
DINAClick Connector	Adult/ Pediatric	V100/VC150, DINAMAP Pro Series, DINAMAP Procure	3.6 m	2058205-001	1
			7.2 m	2058205-002	
		Block style monitor connector for Coro, PDM and some older legacy monitors	1.2 m	2058203-001	
			3.6 m	2058203-002	
			7.2 m	2058203-003	
		CARESCAPE ONE	1.2 m	2087389-001	
			3.6 m	2087389-002	
			7.2 m	2087389-003	
Neo-Snap Connector	Neonate	CARESCAPE ONE	2.4 m	2089791-002	1
			3.6 m	2089791-003	
		Block style monitor connector for Coro, PDM and some older legacy monitors	2.4 m	2058204-002	
			3.6 m	2058204-003	
		V100/VC150, DINAMAP Pro Series, DINAMAP Procure	3.6 m	2058206-001	



Not available in all countries



CONNECTORS

	PART NUMBER	ACCESSORIES	QUANTITY/PACKAGE
	2059294-001	2T DINAClick FEMALE RECEPT	10/PACK
	2055630-001	1T DINAClick 1/8" TUBE FEMALE RECEPT	10/PACK
	2055630-002	1T CONN 2-1 DINAClick 5/32" TUBE FEMALE RECEPT	10/PACK

WALL MANOMETERS AND ACCESSORIES

PART NUMBER	WALL MANOMETER AND ADAPTERS FOR MANUAL INFLATION SYSTEMS	QUANTITY/PACKAGE
2059289-001	INCLUDES COILED TUBING, ADAPTER WITH DINAClick CONNECTOR, BULB AND VALVE	5/BOX
2059289-002	INCLUDES ADAPTER WITH DINAClick CONNECTOR AND 13" TUBING (NO COILED TUBING, NO BULB AND VALVE)	5/PACK

ADAPTERS

	PART NUMBER	DESCRIPTION	QUANTITY/PACKAGE
	2059291-001	Y-ADAPTER: 1T MALE BAYONET TO 2T DINAClick FEMALE RECEPT	10/PACK
	2059295-001	WASH PLUGS – FOR NEO-SNAP FEMALE & DINAClick MALE	50/PACK

Not available in all countries





Our most popular cuff, the CRITIKON™ SOFT-CUF™ is designed with soft, absorbent material with rounded corners for maximizing patient comfort during extended stays at the hospital. These cuffs are ideal for patients with longer term stays or for protocols where the cuff travels with the patient through their total length of time in the hospital. It is economically priced to help support your infection control efforts.



20 cuffs/box

SIZE	COLOR	LIMB CIRCUMFERENCE	2-TUBE DINACLICK	1-TUBE BAYONET CONNECTOR
Infant	Orange	8-13 cm	SFT-P1-2A-5S	2360
Child	Green	12-19 cm	SFT-P2-2A-5S*	2361*
Child Long	Green	12-19 cm	SFT-P2-2A-L-5S	
Small adult	Royal Blue	17-25 cm	SFT-A1-2A-5S*	2362*
Small adult Long	Royal Blue	17-25 cm	SFT-A1-2A-L-5S	
Adult	Navy	23-33 cm	SFT-A2-2A-5S*	2363*
Adult Long	Navy	23-33 cm	SFT-A2-2A-L-5S*	2140*
Large Adult	Pink	31-40 cm	SFT-A3-2A-5S*	2364*
Large Adult Long	Pink	31-40 cm	SFT-A3-2A-L-5S	2141
Thigh	Brown	38-50 cm	SFT-T1-2A-5S	2365

*Part number with an “-INT” suffix, denoting that the Instructions For Use (IFU) is available in multiple languages.



ASSORTMENTS DINACLICK CONNECTOR

PRODUCT	2-TUBE DINACLICK
20 Cuffs : 5 Small Adult, 10 Adult, 5 Large Adult	2059303-501
20 Cuffs : 6 Infant, 10 Child, 4 Small Adult	2059303-502
20 Cuffs : 3 Small Adult, 7 Adult, 3 Adult Long, 3 Large Adult, 4 RADIAL-CUF (forearm)	2059303-503
20 Cuffs : 2 Infant, 2 Child, 2 Small Adult, 6 Adult, 2 Adult Long, 2 Large Adult, 2 RADIAL-CUF (forearm), 2 Thigh	2059303-504
4 cuffs: 1 Small Adult, 1 Adult, 1 Large Adult, 1 RADIAL-CUF	2059303-506

Not available in all countries

CLASSIC-CUF

CRITIKON™ CLASSIC-CUF™ BP cuffs are designed and clinically validated for limited reuse on multiple patients. Made with durable vinyl construction, this time-tested CRITIKON blood pressure cuff is both easy-to-clean and strong enough to withstand repeated inflations – yet economical enough for single patient- use.



5 cuffs/box

SIZE	COLOR	LIMB CIRCUMFERENCE	2-TUBE DINACLICK	1-TUBE BAYONET CONNECTOR
Infant	Orange	8-13 cm	CLA-P1-2A-5S	2300
Child	Green	12-19 cm	CLA-P2-2A-5S	2301
Small Adult	Light Blue	17-25 cm	CLA-A1-2A-5S	2302
Adult	Navy	23-33 cm	CLA-A2-2A-5S	2303
Adult long	Navy	23-33 cm	CLA-A2-2A-L-5S	2306
Large Adult	Pink	31-40 cm	CLA-A3-2A-5S	2304
Large Adult Long	Pink	31-40 cm	CLA-A3-2A-L-5S	2307
Thigh	Brown	38-50 cm	CLA-T1-2A-5S	2305



ASSORTMENTS DINACLICK CONNECTOR

PRODUCT	2-TUBE DINACLICK
3 Cuffs : 1 Small Adult, 1 Adult, 1 Large Adult	2059304-502
20 Cuffs : 5 Small Adult, 10 Adult, 5 Large Adult	2059304-501



Not available in all countries





CRITIKON™ DURA-CUF™ is gentle and designed for durability. Constructed with materials not made with PVC or rubber latex.

Easy-to-use, fast **DINAClick connector** enables nurses to link the cuff to the NIBP hose quicker than a pair of traditional tube connectors. The audible click lets users know there is a proper connection.



5 cuffs/box

SIZE	COLOR	LIMB CIRCUMFERENCE	2-TUBE DINAClick	1-TUBE BAYONET CONNECTOR
Infant	Orange	8-13 cm	DUR-P1-2A-5S	002274
Child	Green	12-19 cm	DUR-P2-2A-5S	002275
Small Adult	Royal Blue	17-25 cm	DUR-A1-2A-5S	002276
Small Adult Long	Royal Blue	17-25 cm	DUR A1-2A-L-5S	
Adult	Navy	23-33 cm	DUR-A2-2A-5S	002277
Adult long	Navy	23-33 cm	DUR-A2-2A-L-5S	002130
Large Adult	Wine	31-40 cm	DUR-A3-2A-5S	002278
Large Adult Long	Wine	31-40 cm	DUR-A3-2A-L-5S	002131
Thigh	Brown	38-50 cm	DUR-T1-2A-5S	002279

ASSORTMENTS DINAClick CONNECTOR

PRODUCT	2-TUBE DINAClick
3 cuffs : 1 Small Adult, 1 Adult, 1 Large Adult	2059301-501
6 cuffs : 2 Infant, 3 Child, 1 Small Adult	2059301-502
6 cuffs : 1 Small Adult, 3 Adult, 1 Adult Long, 1 Large Adult	2059301-503
6 cuffs : 1 Infant, 1 Child, 1 Small Adult, 1 Adult, 1 Large Adult, 1 Thigh	2059301-504

Not available in all countries

YELLOW CLASSIC-CUF

Yellow CRITIKON™ CLASSIC-CUF™ is designed for single patient use within an isolation environment. These vinyl cuffs are priced as disposable-use cuffs and support your infection prevention protocols. They provide the same quality and assurance that you expect from the entire line of CRITIKON cuffs, while tailored for isolation areas.



20 cuffs/box

SIZE	COLOR	LIMB CIRCUMFERENCE	2-TUBE DINACCLICK
Infant	Yellow/Blue	8-13 cm	CLA-P1-2A-X-5S
Child	Yellow/Blue	12-19 cm	CLA-P2-2A-X-5S
Small Adult	Yellow/Blue	17-25 cm	CLA-A1-2A-X-5S
Adult	Yellow/Blue	23-33 cm	CLA-A2-2A-X-5S
Large Adult	Yellow/Blue	31-40 cm	CLA-A3-2A-X-5S
Thigh	Yellow/Blue	38-50 cm	CLA-T1-2A-X-5S



Not available in all countries



GE Healthcare offers the first BP forearm cuff clinically validated for use on obese or overly muscular adult patients. It provides proper fit and comfort to patients, while giving reliable blood pressure readings that support quality clinical care.



SIZE	COLOR	LIMB CIRCUMFERENCE	2-TUBE DINAClick
ADULT	PURPLE	26-36 CM (FOREARM)	SFT-F1-2A-5S

Not available in all countries

NEONATAL



Providing big benefits for your tiniest patients, GE Healthcare's neonatal cuffs are designed with the accuracy you trust and the gentle comfort you expect. We offer a full range of neonatal sizes with the choice of extra soft, absorbent material or easy-to-clean vinyl. The playful, color-coded animal designs help you find the right size at a glance and they are economically priced for single-patient-use.

20 cuffs/box


			2-TUBE NEO-SNAP		1-TUBE NEO-SNAP CONNECTOR	
SIZE	COLOR	LIMB CIRCUMFERENCE	CLASSIC-CUF	SOFT-CUF	CLASSIC-CUF	SOFT-CUF
1	Orange	3-6 cm	CLA-N1-2B	SFT-N1-2B	CLA-N1-1B	SFT-N1-1B
2	LT. Blue	4-8 cm	CLA-N2-2B	SFT-N2-2B	CLA-N2-1B	SFT-N2-1B
3	Green	6-11 cm	CLA-N3-2B	SFT-N3-2B	CLA-N3-1B	SFT-N3-1B
4	Navy	7-13 cm	CLA-N4-2B	SFT-N4-2B	CLA-N4-1B	SFT-N4-1B
5	Pink	8-15 cm	CLA-N5-2B	SFT-N5-2B	CLA-N5-1B	SFT-N5-1B

NEONATAL ASSORTMENT PACK	CUFF TYPE	PART NUMBER	DESCRIPTION
2-Tube Neo-Snap Connectors	SOFT-CUF	2059305-002	2 Neo #1, 3 Neo #2, 5 Neo #3, 5 Neo #4, 5 Neo #5 (20 cuffs/box)
	CLASSIC-CUF	2059306-002	2 Neo #1, 3 Neo #2, 5 Neo #3, 5 Neo #4, 5 Neo #5 (20 cuffs/box)
	CLASSIC-CUF	2059306-003	1 Neo #3, 1 Neo #4, 1 Neo #5 (3 cuffs/pack)
1-Tube Neo-Snap Connectors	SOFT-CUF	2059305-001	2 Neo #1, 3 Neo #2, 5 Neo #3, 5 Neo #4, 5 Neo #5 (20 cuffs/box)
	CLASSIC CUF	2059306-001	2 Neo #1, 3 Neo #2, 5 Neo #3, 5 Neo #4, 5 Neo #5 (20 cuffs/box)

Not available in all countries



Cardiac Output



CARDIAC OUTPUT



Cardiac Output Cables

Description	Length	For use with	Part Number	Quantity
Cardiac Output Cable	3.6 m	CARESCAPE™ Monitors B450, B650 and B850 with the CARESCAPE™ Patient Data Module (PDM)	2104180-001	1
	2.4 m		2104180-002	1

Specifications available on request

Not available in all countries



OTHER ACCESSORIES

Other Accessories				
Description	For use with	Manufacturer Reference	Part Number	Quantity
Catheter Connecting Cable	E-COP, E-COPsv	COM2-CC	16590-HEL	1
Thermodilution Catheter, disposable	E-COP, E-COPsv	131HF7	M1238258	5

These accessories are sold ONLY with the equipment

PiCCO™ CABLES

PiCCO™ Cardiac Output Cables				
Description	Length	For use with	Part Number	Quantity
Continuous Cardiac Output Cable	4.1m	E-PiCCO module	M1154585	1
Injectate Sensor Cable	0.9m		M1182448	1

These accessories are sold ONLY with the equipment



GE Healthcare compatible Cardiac Output probes are available from multiple manufacturers. You can purchase the following product portfolios from our partners directly:

Companies

- Vyair Medical
- Edwards Life sciences
- Getinge
- Argon Medical
- Spectramed

Portfolios

- Injectate Bath Probes
- Flow-through Injectate-Temperature Probes
- PiCCO Accessories
- Temperature Probes

Not available in all countries



Adequacy of Anesthesia



ENTROPY



M1038681

M1174413

Sensors			
Product	Use	Part Number	Quantity
Entropy™ EasyFit Sensor	Disposable	M1174413	25
Entropy™ Sensor		M1038681	

Cable			
Product	Length	Part Number	Quantity
Entropy™ Cable	3.5 m	M1050784	1

Product compatibility

E-ENTROPY module

The GE Healthcare Entropy™ module, E-ENTROPY, and accessories are indicated for adult and pediatric patients older than 2 years within a hospital for monitoring the state of the brain by data acquisition of electroencephalograph (EEG) and frontal electromyograph (FEMG) signals during general anesthesia.

Not available in all countries



Electromyography (EMG)
Requires 5 electrodes



Kinemyography (KMG)
Requires 2 electrodes



888418/897439



888416

Sensors				
Product	Use	Patient	Part Number	Quantity
NMT MechanoSensor	Reusable	Adult	888418	1
		Pediatric	897439	
NMT ElectroSensor		Adult, Pediatric	888416	1
NMT Regional Block Adapter			888417	1

Accessories				
Product	Use	Patient	Part Number	Quantity
NMT Electrodes	Disposable	Adult, Pediatric	57268-HEL	30


Cable			
Product	Length	Part Number	Quantity
NMT Sensor Cable	3,3 m	888414	1

Product compatibility

E-NMT module

Not available in all countries

Neurology



898050

Reusable Preconfigured EEG Leadsets			
Product	Length	Part Number	Quantity
Basic montage	0.6 m	898050	1
AEP montage		898052	1

Not available in all countries





Disposable EEG Electrodes		
Product	Part Number	Quantity
EEG Clip electrodes Ag/AgCl	M1240794	40 pouches of 25 (1,000 total)
	572686	500
EEG Prewired cup electrodes	2083319-001	46 pouches of 10 (460 total)
EEG Prewired wet gel electrodes	2083320-001	42 pouches of 10 (420 total)

All electrodes within a set-up should be similar

Product compatibility

E-EEG(X) module and N-EEG(X) headbox

AEP Accessories		
Product	Part Number	Quantity
Disposable ear tips for use with earphone set	8001233	20

To use with 898052

Product compatibility

E-EEG(X) module and N-EEG(X) headbox

Not available in all countries

Patient Spirometry



896952



8001948



733910-HEL



73393

Sensors				
Description	Type	Approved with Modules	Part Number	Quantity
D-lite™ +	Single Patient Disposable	E-COV, E-COVX, E-CAiOV, E-CAiOVX, E-sCOV, E-sCOVX, E-sCAiOV, E-sCAiOVE, E-sCAiOVX	896952	50
Pedi-lite+			8001948	
D-lite™	Re-usable		733910-HEL	1
Pedi-lite			73393	

Specifications available on request

Not available in all countries



SPIROMETRY



890031

Spirometry Tubes					
Description	Length	Type	Approved with Modules	Part Number	Quantity
Spirometry Tubes	2m	Single Patient Disposable	E-sCOV, E-sCOVX, E-sCAiOV, E-sCAiOVE, E-sCAiOVX, E-COV, E-COVX, E-CAiOV, E-CAiOVX	890031	5
	3m			884101	
	6m			891191-HEL	



D-lite™ ++ set

Kits (include gas sample line, sensor and spirometry tube)					
Description	Length	Type	Approved with Modules	Part Number	Quantity
D-lite™ ++ set	2m	Single Patient Disposable	E-sCOV, E-sCOVX, E-sCAiOV, E-sCAiOVE, E-sCAiOVX	2104297-001	20
	3m			2104297-002	
D-lite™ + kit	2m	Single Patient Disposable	E-sCOV, E-sCOVX, E-sCAiOV, E-sCAiOVE, E-sCAiOVX, E-COV, E-COVX, E-CAiOV, E-CAiOVX	8004381	50
	3m			8004382	
Pedi-lite+ kit	2m	Single Patient Disposable	E-sCOV, E-sCOVX, E-sCAiOV, E-sCAiOVE, E-sCAiOVX, E-COV, E-COVX, E-CAiOV, E-CAiOVX	8002718	50
	3m			M1032634	

Specifications available on request



Pedi-lite+ kit

Not available in all countries



Gas Monitoring



ICU AIRWAY ADAPTER SETS

Airway Gas Sampling Sets					
Description	Type	Approved with Modules	Length	Part Number	Quantity
Adult/Pediatric	Single Patient Disposable	E-miniC, E-sCO, E-sCOV, E-sCOVX	3 m	8090121531	20 pcs/box
			2 m	8090121521	
Infant		E-sCO, E-sCOV, E-sCOVX	2 m	8090121522	20 pcs/box



GAS EXHAUST LINES

Gas Exhaust Lines					
Description	Type	Approved with Modules	Length	Part Number	Quantity
Gas Exhaust Line, white conical and Colder	Single Patient Disposable	E-sCO, E-sCOV, E-sCOVX, E-sCAiO, E-sCAiOE, E-sCAiOV, E-sCAiOVX, E-sCAiOVE, N-CAiO and N-CAiO with new gas exhaust connector; serial number prefix "SUA"	1 m	5514182	1
			0.18 m	5514183	
Gas Exhaust Line, white conical		E-miniC with new gas exhaust connector; serial number prefix "SUB"	2 m	5514184	
Gas Exhaust Line, Luer and Colder		E-sCO, E-sCOV, EsCOVX, E-sCAiO, E-sCAiOV, E-sCAiOVX, E-sCAiOE, E-sCAiOVE and N-CAiO with luer type connector; serial number prefix "SGV"	0.18 m	8004462	
Gas Exhaust Line, Luer and Colder			1 m	8004463	
Gas Exhaust Line, Luer		E-sCO, E-sCOV, EsCOVX, E-sCAiO, E-sCAiOV, E-sCAiOVX, E-sCAiOE, E-sCAiOVE and N-CAiO with luer type connector; serial number prefix "SGV"	2 m	M1031274	5
		E-miniC with luer type connector			



5514182



5514183



5514184



8004462 & 8004463



M1031274



GAS SAMPLING LINES AND ADAPTERS



Gas Sampling Lines					
Description	Type	Approved with Modules	Length	Part Number	Quantity
PVC/PE gas sampling line	Single Patient Disposable	E-miniC, E-CO, E-COV, E-COVX, E-CAiO, E-CAiOV, E-CAiOVX, E-sCO, E-sCOV, E-sCOVX, E-sCAiO, E-sCAiOE, E-sCAiOV, E-sCAiOVX, E-sCAiOVE, N-CAiO	2 m	2097307-001	10
			3 m	2097307-002	
			6 m	2097307-003	5



3090121002

Gas Monitoring Airway Adapter				
Description	Type	Approved with Modules	Part Number	Quantity
Airway Adapter with Luer connector	Single Patient Disposable	E-miniC, E-sCO, E-sCOV, E-sCOVX, E-sCAiO, E-sCAiOE, E-sCAiOV, E-sCAiOVX, E-sCAiOVE, N-CAiO	3090121002	25

WATER TRAPS



881319-HEL



M1200227



876446-HEL



M1182629



8002174

Water Traps – ICU

Description	Type	Approved with Modules	Part Number	Quantity
D-fend™ +	Single Patient Disposable	E-CO, E-COV, E-COVX	881319-HEL	10
D-fend™ Pro+		E-sCO, E-sCOV, E-sCOVX	M1200227	
Mini D-fend™		E-miniC	8002174	

Specifications available on request

Water Traps – Anesthesia

Description	Type	Approved with Modules	Part Number	Quantity
D-fend™	Multi-Patient	E-CAiO, E-CAiOV, E-CAiOVX	876446-HEL	10
D-fend™ Pro		E-sCAiO, E-sCAiOE, E-sCAiOV, E-sCAiOVX, E-sCAiOVE	M1182629	
Mini D-fend™		E-miniC	8002174	

Specifications available on request

Not available in all countries



ANESTHESIA & RESPIRATORY CARE



ANESTHESIA & RESPIRATORY CARE

Absorbents

- Absorbers
- Canisters
- Dust Filter

Vaporizer Bottle Adapters

Filters & Breathing Masks

Reservoir Bags & Gas Scavenging

- Anesthesia Breathing Bags
- Scavenging Bags, Tubes, & Connectors

Patient Breathing Circuits

- Circuits
- Circuit Components
- Circuit Kits
- Circuit Startup Kits

Flow Sensors - Respiratory

- Cables
- Heater
- Respiratory Flow Sensors and Expiratory Valve Assemblies

Flow Sensors – Anesthesia

Nebulizer Accessories

Anesthesia O2 Cells

Respiratory Safety Guard

Pressure Catheters



Anesthesia & Respiratory care



ABSORBENTS

Disposable CO₂ Absorbers

Description	Compatibility	Part Number	Quantity
AMSORB® Plus G-Can Absorber, 1.0L	All Aisys™, Avance™, Aespire™, 9100c and ADU variants	2105489-003	8
AMSORB® Plus Prefilled Cartridge, 1.2L	All Aestiva™ variants	2105489-004	12
AMSORB® Plus CareCan Absorber, 1.4L	Carestation™ 600 series & 750 series*	2105489-006	8
AMSORB® Plus Bubble-Can Universal 1.2L	Apollo™, Cicero™, Fabius CE/GS/Trio™, Julian™, Pallas™, Perseus™, Primus™, Zeus™	2105489-011	6

Reusable CO₂ Absorbent Components

Description	Compatibility	Part Number	Quantity
AMSORB® Plus Jerican Loosefill Absorbent, 5.0L	GE Healthcare Reusable Canisters	2105489-007	2
ABS Reusable CO ₂ Canister	All Aisys™, Avance™, Aespire™ & 9100c variants	1407-7004-000	1
CBS Reusable CO ₂ Canister	Carestation™ 600 series & 750 series*	2071165-001-S	1
Foam Dust Filter	ABS Reusable Canister	2105489-001	40

VAPORIZER BOTTLE ADAPTERS

Easy-Fil Vaporizer Bottle Adapters

Vaporizer Compatibility	Color Indicator	Part Number	Quantity
Isoflurane Easy-Fil	Purple	2106572-001	1
Halothane Easy-Fil	Red	2106572-002	
Sevoflurane Easy-Fil	Yellow	2106572-003	



FILTERS AND BREATHING MASKS



HMEFs and Filters – Disposable

Description	Size	Sampling Port	Part Number	Quantity
HMEF	Neonatal	Included	2106570-006	50
	Pediatric		2106570-009	
	Adult		2106570-010	
Bacterial and Viral Filter	Adult	Included	2106570-007	50
HEPA Filter	Adult		2106570-008	



Anesthesia Breathing Masks – Disposable

Mask Connection	Mask Size	Part Number	Quantity
15mm	0 – Neonatal	2106570-012	50
	1 – Small Child	2106570-013	
22mm	2 – Child	2106570-014	50
	3 – Large Child	2106570-015	
	4 – Adult	2106570-016	50
	5 – Large Adult	2106570-017	
	6 – Extra Large Adult	2106570-018	



Anesthesia Breathing Masks – Reusable

Mask Connection	Mask Size	Part Number	Quantity
15mm	0 – Neonatal	2106570-011	1
		2106570-029	10
	1 – Small Child	2106570-019	1
		2106570-030	10
22mm	2 – Child	2106570-020	1
		2106570-031	10
	3 – Large Child	2106570-021	1
		2106570-032	10
	4 – Adult	2106570-022	1
		2106570-033	10
	5 – Large Adult	2106570-023	1
		2106570-034	10

Not available in all countries



RESERVOIR BAGS AND GAS SCAVENGING

Anesthesia Manual Breathing Bags – Reusable

Bag Connection Size	Material	Bag Size	Part Number	Quantity
22mm	Clear Silicone	0.5L	2106570-001	1
			2106570-025	5
		1.0 L	2106570-002	1
			2106570-026	5
		2.0 L	2106570-003	1
			2106570-027	5
		3.0 L	2106570-004	1
			2106570-028	5



Passive Scavenging Tubes – Disposable

Connection	Length	Type	Part Number	Quantity
30M-30M	5 m	Tube	2105488-006	1
			2105488-022	10



Gas Scavenging Bag – Reusable

Size	Compatability	Connection	Part Number	Quantity
3.0L	Aisys, Avance, Aespire™ (ABS Breathing Systems)	30M	2106570-024	1



Not available in all countries



PATIENT BREATHING CIRCUIT KITS – START UP

All startup kits include patient breathing circuit, reservoir bag tube, reservoir bag, gas sampling lines and gas sampling exhaust lines. These are primarily intended for initial anesthesia machine setup.



Patient Breathing Circuit Kits – Start Up – Adult

Type	Gas Sampling Line	Gas Sampling Exhaust Line	Limb Circuit	Breathing Circuit	Breathing Bag		Part Number	Quantity
				Length	Manual Tube 1.5 M	Size		
Disposable	3m PVC/PE	0.18 m & 1 m	Dual	Expandable up to 3 m	Straight	3L	2096533-001	1
						2L	2096533-003	1
			Single Split Lumen	1.8 m	Angled*	2L	2096536-002	10
					Straight	2L	2096536-003	10
					Angled*	2L	2105488-019	1
Reusable			Dual	1.5 m	Angled*	2L	2105488-024	1
					Straight	2L	2096534-009	1
					Angled*	2L	2096536-001	1

*For use with bag support arm

Not available in all countries




PATIENT BREATHING CIRCUITS

Our wide range of patient breathing circuits and components offer robust and high quality options for both disposable and reusable applications.

Patient Breathing Circuit Kits									
Type	Size	3m PVC/PE Gas Sampling Line	Circuit Type	Breathing Circuit	Breathing Bag		Includes	Part Number	Qty
				Length	Manual Tube	Size			
Disposable	Adult	Yes	Dual Limb	Expandable up to 3 m	Not included	3L	-	2096535-006	20
				1.8 m	Angled*, 1 M		-	2105488-013	20
					Straight, 1.5 M		-	2105488-014	20
		No	Single Split Lumen	1.8 m	Straight	3L	-	2105488-010	20
			Dual Limb	Expandable up to 1.8 m	Not included		-	2105488-011	20
				1.5 m	Not included	-	2x Water Traps	2096535-001	20
					Straight, 0.8 M	-	1x Water Trap	2105488-017	20
			Not Included			Angled*, 2 M	2L		2105488-012
	Neonatal	Yes	Dual Limb	1.5 m	Not included	1L	HME & Filter	2105488-001	20
		No						2105488-021	1
Reusable	Adult	No	Dual Limb	1.5 m	Not included		2x Water Traps	2096534-005	1
	Pediatric							2096534-006	1

*For use with bag support arm



PATIENT BREATHING CIRCUITS



2096534-007
Reusable, Pediatric



2105488-020
Disposable, Pediatric




2105488-015
Disposable, Adult



2105488-009
Disposable, Adult, Inteliflo® Single Split Lumen

Patient Breathing Circuits							
Type	Size	Circuit Type	Breathing Circuit Length	Patient Elbow	Gas Sampling Port Location	Part Number	Quantity
Reusable	Adult	Dual Limb	1.5 m	Double Swivel, Removable	Elbow	2096534-007	1
	Pediatric					2096534-008	1
Disposable	Adult	Single Split Lumen	1.8 m	Fixed, Removable	Elbow	2105488-009	20
		Dual Limb	1.5 m	Double Swivel, Removable	Patient Wye	2096535-002	20
				Fixed, Removable	Elbow	2105488-015	20
	Pediatric			Double Swivel, Removable	Patient Wye	2105488-020	20



Non-Recirculating Breathing Circuits – Jackson Rees

Type	Breathing Circuit Length	Bag Size	Gas Sampling Port Location	Part Number	Quantity
Disposable	0.8 m	1 L	Elbow	2105488-008	20

Not available in all countries



Reusable Patient Breathing Circuit Tubes					
Size	Material	Length	Connections	Part Number	Quantity
Adult	Hytrel®	0.35 m	22F-22F (90 Deg)	2096534-001	1
		0.75 m	22F-22F	2096534-002	1
		0.9 m		2096534-003	1
		1.5 m		2096534-004	1
	Silicone	1 m		2096534-101	1
		1.5 m		2096534-102	1
Pediatric	Silicone	1.5 m	22F-10F	2096534-103	1
	Hytrel®	1.1 m	22F-15F	2105488-007	1
				2105488-023	10



Not available in all countries

PATIENT BREATHING CIRCUITS



Patient Breathing Circuit Components				
Description	Size	Reprocessable	Part Number	Quantity
Y-Piece Connector Kit, 15M-10F	Neonatal	Reusable	2096534-104	10
Patient Angled Y-Piece			2105488-002	10
Patient Straight Y-Piece - Swivel			2105488-004	10
Straight Connector - CBS Bag Arm, 22M-22M	Adult	Reusable	2096533-004	1
ABS Bag Arm Connector, 22M-22M	All		2105488-005	10
Patient Elbow - Swivel	Adult	Reusable	2105488-003	10
Patient Elbow		Disposable	2105488-018	50
Gas Return Adapter, 22M-22F	Adult	Disposable	2105488-016	10

Not available in all countries





Respiratory Flow Sensors and Expiratory Valve Assemblies (EVA)				
Description	Reprocessable	Patient Type	Part Number	Quantity
Flow Sensor Only (Neonatal)	Reusable	Neonatal	1505-3272-000	1
Flow Sensor Only	Single Patient Use	All	1505-3231-001	1
Whole Assembly (housing and flow sensor)			1505-3848-001	1
Housing Only (includes water trap)			1505-8568-001	1
Flow Sensor Only			1505-3231-010	10
Whole Assembly (housing and flow sensor)			1505-3848-010	10
Housing Only (includes water trap)			1505-8568-010	10

Product compatibility

CARESCAPE™ R860 and Engström Ventilators

Not available in all countries



FLOW SENSORS – RESPIRATORY



M1200693



M1188723

Respiratory Flow Sensor and Expiratory Valve Components

Description	Patient Type	Part Number	Quantity
Flow Sensor Cable	Neonatal	1505-5604-000	1
Exhalation Valve Heater	All	M1200693	1
Exhalation Valve Heater Cable		M1188723	1


Product compatibility

CARESCAPE™ R860 and Engström Ventilators

FLOW SENSORS – ANESTHESIA



2087640-001-S



2089610-001-S



2069358-001-S

Anesthesia Flow Sensors

Machine Compatibility	Autoclavable	Part Number	Quantity
Aisys™, Avance™, Aespire™, Aestiva™	Yes	2087640-001-S	1
	No	2089610-001-S	
Carestation™ 600 series & 750 series*	Yes	2069358-001-S	1

*Not cleared approved by the US FDA. Not for sale in the United States.

Not available in all countries





Aerogen® Nebulizer Accessories			
Description	Patient Type	Part Number	Quantity
Aerogen® Pro Nebulizer kit with adult and ped tee and cable	Adult and Pediatric	1505-3846-000	1
Nebulizer head with filler cap	All	AG-AP1000	
Nebulizer Cable		1505-5602-000	
Nebulizer tee with plug	Adult	AG-AP1010	5
	Pediatric	AG-AP1020	
	Neonatal	AG-AP1035	
Nebulizer replacement kit	Pediatric	AG-AP1200	1
		AG-AP1100	1
Aerogen® Solo Nebulizer	Adult	AG-AS3100-S	5
		AG-AS3200-S	10
Nebulizer tee	Adult	AG-AS3010	10
	Pediatric	AG-AS3020	
	Neonatal	AG-AS3035	

Product compatibility

CARESCAPE™ R860 and Engström Ventilators

These accessories are sold ONLY with the equipment.

Not available in all countries



ANESTHESIA O2 CELLS

Anesthesia O2 Cells			
Description	Machine Compatability	Part Number	Quantity
O2 Sensor ONLY	Aisys, Avance, Aespire™	M-10	1
O2 Sensor and cable		2082006-001	
O2 Sensor with Adapter and O-Ring	Carestation™ 600 series & 750 series*	2071018-001-S	

*Not cleared approved by the US FDA. Not for sale in the United States.



2071018-001-S

RESPIRATORY INSPIRATORY SAFETY GUARD

Respiratory Inspiratory Safety Guard			
Machine Compatability	Patient Type	Part Number	Quantity
CARESCAPE™ R860	All	5830446	1
		5514015	10



PRESSURE CATHETERS

Respiratory Pressure Catheters			
Description	Machine Compatability	Part Number	Quantity
Intratracheal Pressure Catheter	CARESCAPE™ R860 and Engström Ventilators	M1045564	10



Not available in all countries



DIAGNOSTIC CARDIOLOGY



DIAGNOSTIC CARDIOLOGY

MAC™, CASE™ and CardioSoft™

ECG

- Accessories
- Cables
- Electrodes
- Leadwires
- Recording Paper

Holter Monitors

ECG – CM 4000

- Leadwire Sets
- Accessories

ECG – SEER® 12

- Leadwire Sets
- Accessories

ECG – SEER® 1000

- Leadwire Sets
- Accessories

ECG – SEER® Light

- Leadwire and Cable Kits
- Accessories

MAC™, CASE™ and CardioSoft™



CABLES AND LEADWIRES – IEC



ECG Trunk Cables – IEC

Description	Length	Part Number	MAC™								CASE™	CardioSoft™
			5	7	VU360	5500	3500	2000	800	600		
Trunk Cable 10-lead IEC	220 cm	2104722-001	X	X			X	X	X	X		

ECG Host Cables – IEC

Description	Length	Part Number	MAC™								CASE™	CardioSoft™
			5	7	VU360	5500	3500	2000	800	600		
Host Patient Cable MAC5500	170 cm	2104403-001				X						
Host Patient Cable CASE	550 cm	2104403-002									X	
Host Patient Cable MAC5500ST	540 cm	2104403-003				X						

Integrated Cable and Leadwire Set – IEC

Description	Length	Part Number	MAC™								CASE™	CardioSoft™
			5	7	VU360	5500	3500	2000	800	600		
ECG Value Cable 10-lead Cable / Leadwire IEC	Chest leads: 80 cm Limb leads: 125 cm Trunk cable: 220cm	2104726-001	X	X			X	X	X	X		



Not available in all countries



LEADWIRES – IEC

ECG Leadwire Set – Banana – IEC												
Description	Length	Part Number	MAC™								CASE™	CardioSoft™
			VU360	5	7	5500	3500	2000	800	600		
Leadwire Set, Base 10, banana, IEC	Limb leads: 125 cm	2104724-002	X	X	X	X	X	X	X	X	X	X
Leadwire Set, 4 Add-on, banana, IEC-HEIM	Chest leads: 80 cm	2104750-001	X			X					X	X
Leadwire Set, Gold, Base 10, banana, IEC	Limb leads: 102 cm Chest leads: 80 cm	2104724-006	X	X	X	X	X	X	X	X	X	X
Leadwire Set, Gold, Base 10, banana, long, IEC	Limb leads: 140 cm Chest leads: 95 cm	2104724-007	X	X	X	X	X	X	X	X	X	X
Leadwire Set, Gold, Add-on, banana, IEC-HEIM	80 cm	2104750-002	X			X					X	X



2104724-002



2104750-001



2104724-006



2104750-002

ECG Leadwire Set – Grabber – IEC												
Description	Length	Part Number	MAC™								CASE™	CardioSoft™
			VU360	5	7	5500	3500	2000	800	600		
Base 10, IEC	Limb lead: 125 cm	2104748-001	X	X	X	X	X	X	X	X	X	X
4 Add-on, IEC-HEIM	Chest lead: 80 cm	2104725-002	X			X					X	X



2104748-001



2104725-002

Not available in all countries





2104729-001



2104767-001

Accessories

Description	Part Number	MAC™								CASE™	CardioSoft™	Quantity
		VU360	5	7	5500	3500	2000	800	600			
Adapter Clip, Wide Mouth Clip	2104729-001	X	X	X	X	X	X	X	X	X	X	14
ECG Leadwire Port Plug	2104767-001	X			X					X	X	4

CARDIOSOFT™



CardioSoft™ Accessories

Description	Part Number	Quantity
Supply BP Bag Belt TONOPORT	2104824-001	1

Not available in all countries




ECG RECORDING PAPER



DCAR ECG Recording Thermal Paper

Folding	White patient data area	Paper size	Red grid wide	Queue	MAC™								CASE™	Part Number	Quantity
					VU360	5	7	5500	3500	2000	800	600			
z-fold	No	A4 *	195 mm	Hole				X	X				X	2104769-001	300 sheets/ pack, 8 packs
			205 mm	Block	X	X	X			X				2104772-001	150 sheets/ pack, 10 packs
		A5	205 mm	Block		X								5684683	150 sheets/ pack, 20 packs
		21.5 cm x 28 cm	200 mm	Hole				X	X				X	2104768-001	300 sheets/ pack, 8 packs
		80 mm x 90 mm	75 mm	Block								X		2104773-001	280 sheets/ pack, 10 packs
		110 mm x 140 mm	100 mm	Block							X			2104775-001	200 sheets/ pack, 10 packs
	Yes	21.5 cm x 28 cm	155 mm	Hole	X			X	X					2104770-001	300 sheets/ pack, 8 packs
				Hole	X	X	X		X	X				2104771-001	150 sheets/ pack, 16 packs

For ECG paper for Patient monitoring, please go to page 6.



Not available in all countries



REUSABLE ECG ELECTRODES

GE Healthcare reusable bulb electrodes with silver-silver chloride (Ag/AgCl) plating are designed to ensure positive contact with the skin and provide maximum signal quality.



Reusable Bulb Electrode		
Description	Part Number	Case Quantity
Reusable Bulb Electrode	2104783-001	6



Bulb Electrode Accessories			
Description	Compatibility	Part Number	Case Quantity
Clamp Limb AHA Value	MAC™ 5, MAC™ 7, MAC™ 600, MAC™ 800, MAC™ 2000 ECG System	2104784-001	4/set
Clamp Limb IEC Value		2104785-001	

Specifications available on request

Not available in all countries

Holter Monitors



CM 4000 LEADWIRE SETS



2105803-005

CM 4000 Leadwire Set – IEC			
Description	Length	Part Number	Quantity
3-leads	80 cm	2105803-001	1
5-leads		2105803-003	1
7-leads		2105803-005	1

Specifications available on request

Not available in all countries

CM 4000 Accessories		
Description	Part Number	Case Quantity
Carry Pouch Reusable	2054933-014	1
Carry Bag Disposable	2054933-015	1
USB Cable	2054933-013	1

Specifications available on request

SEER® 12



2104808-001

SEER® 12 – IEC				
Description	Length	Type	Part Number	Quantity
12-Channel Holter ECG Cable	130 cm	Round	2104808-001	1
	105 cm	Flat	2104809-001	

SEER® 12 ECG cables are compatible with SEER® 12 monitors and CardioMem™ CM 3000-12 BT monitors



2104803-001


SEER® 12 Holter Monitor Accessories

Description	Part Number	Quantity
SEER® 12 lead portable pouch, waist belt, shoulder strap, 1 battery compartment	2104803-001	1
SEER® 12 Carrying pouch, 2 battery compartment	2014606-297	1

2104803-001 is compatible with SEER® 12, CM 3000, CM 3000-12
 2014606-297 must be purchased with service
 Specification available on request

Not available in all countries





SEER® 1000 Leadwire Set – IEC

Description	Length	Part Number	Quantity
5 leadwire, 2 CH	75 cm	2104810-001	1
7 leadwire, 3 CH		2104811-001	
5 leadwire, 2 CH	105 cm	2104812-001	
7 leadwire, 3 CH		2104813-001	
3 leadwire, 3 CH	75 cm	2104818-001	
	105 cm	2104819-001	

Product compatibility

SEER® 1000 ECG cables are compatible with SEER® 1000 24-hour Holter monitor, SEER® 1000 48-hour Holter monitor and SEER® 1000 7-day Holter monitors

Specifications available on request

Not available in all countries



2067634-073

SEER® 1000 Holter Monitor Accessories		
Description	Part Number	Quantity
Pouch with waist belt and shoulder strap	2067634-073	1
Pouch Belt Clip	2067634-166	1

Specification available on request



2067634-166

Not available in all countries

SEER® Light Leadwire Set Kits – IEC			
Description	Length	Part Number	Quantity
Carry Pouch and 2 Channel Patient cable	140 cm	2042686-003	1
Carry Pouch and 3 Channel Patient cable		2042686-004	
Carry Pouch and 2 Channel Patient cable	60 cm	2042686-007	1
Carry Pouch and 3 Channel Patient cable		2042686-008	

Product compatibility

SEER® Light and SEER® Light Extend
Specifications available on request

Not available in all countries

SEER® Light Accessories		
Description	Part Number	Quantity
Compact Flash	2008750-003	1
Connect Kit	2020044-292	1
Controller Kit	2042685-007	1

Product compatibility

SEER® Light and SEER® Light Extend

Specifications available on request



SEER® Light Holter Monitor Accessories		
Description	Part Number	Quantity
SEER® Light, shoulder straps, waist belt, portable pouch	2104827-001	1

Product compatibility

2104827-001 is compatible with SEER® Light and SEER® Light Extend

Specification available on request

Not available in all countries





MATERNAL INFANT CARE

Maternal Fetal Care

Mini Telemetry

- Recording Thermal Paper
- Remote Event Marker

Novii™

- Cables
- Novii™ Patch
- Tape

Corometrics™ and Mini Telemetry

- Cables
- Transducer Belts
- Transducer Button and Loop

Internal Monitoring

- Cables
- Electrode
- Pad

Infant Care

Giraffe™/Panda™

- Air filters
- Clip
- Cover

Phototherapy accessories

- Pad Covers

Resuscitation

- T Piece

Lullaby™ Warmer


- Temperature Probes
- Temperature Probe Patches

Shuttle™

- Patient Restraint

Maternal Fetal Care

- Mini Telemetry
- Novii™
- Corometrics™



MIC Fetal Monitoring Recording Thermal Paper

Scaling	Corometrics™ 170 series	Corometrics™ 250cx series	Part Number	Quantity
30-240 bpm scaling	X	X	2104907-001	40 packs
50-210 bpm scaling	X	X	2104908-001	

For ECG paper for Patient monitoring and Diagnostic cardiology, go to pages 6 and 58.

Specifications available on request

Accessories

Description	Part Number	Quantity
Remote Event Marker	2107349-001	1



2107349-001

Not available in all countries





Novii™ Accessories			
Description		Part Number	Quantity
Novii™ Patch		107-PT-004-10	10
		107-PT-004-50	50
Novii™ Interface Cable	GE FECG	105-PT-102	1
	GE MECG	105-PT-104	
	GE UA	105-PT-106	
3M™ Red Dot™ 2236 Skin Prep Tape		100-PT-007	1

Not available in all countries





MECG Accessories - Vyaire			
Description		Part Number	Quantity
Gray MEGC Cable	3 Lead Set, Model 118	1554AAO	1
	3 Lead Set IEC, Model 118	1554BAO	1
120 Green MEGC Cable	3 Lead Set	1553AAO	1
	3 Lead Set IEC	1553BAO	1
FECG/MECG Y Adapter Cable, 120 Series		1442AAO	1

Sold with equipment only.



2015827-001

Transducer Belts			
Description		Part Number	Quantity
Single patient use	Button-style, blue	2015827-001	10



2108347-001

Ultrasound Transducers		
Description	Part Number	Quantity
Ultrasound Transducer Button and Loop	2108346-001	1
TOCO Transducer Button and Loop	2108347-001	1

Not available in all countries



2109002-001



2109006-001

Fetal Spiral Electrode Accessories			
Description	Type	Part Number	Quantity
Cable	Rectangluar Connector	2109005-001	1
	Round Connector	2109004-001	1
ATTACHMENT PAD	Large, Single Patient Use	2109006-001	50
SPIRAL ELECTRODE	Single Helix, Single Patient Use	2109002-001	50

For ECG Electrodes for Patient Monitoring, go to pages 14 and 15.

Not available in all countries



Infant Care

- Giraffe™/Panda™
- Phototherapy accessories
- Resuscitation



GIRAFFE™/PANDA™



2107356-001



2107357-001

Accessories		
Description	Part Number	Quantity
Air Filters	2107356-001	10
Bag and mask holder clip	2107357-001	1

Not available in all countries





Giraffe™ Hood Covers

Description	Part Number	Quantity
OmniBed™ Cover	6600-0825-800	1
Incubator Cover	6600-0846-800	1

Sold with equipment only

Not available in all countries





Flat cushioned cover with straps
The flat, cushioned BiliSoft™ Cover lets you swaddle baby and phototherapy pad together.



Cushioned cover with developmental positioning foot roll and straps
The BiliSoft™ Nest offers sick babies the boundaries and support they need.

BiliSoft™ Disposable Pad Covers					
Description	BiliSoft™ 1.0	BiliSoft™ 2.0	Size	Part Number	Quantity
Disposable Pad Covers, Flat		X	Small	2105422-001	20
				2105418-001	50
	X		Small	2107201-001	20
				2107200-001	50
	X	X	Large	2105423-001	20
				2105419-001	50
Disposable Nest Pad Covers		X	Small	2105420-001	15
	X			2107202-001	
	X	X	Large	2105421-001	

Product compatibility

The BiliSoft™ Pad Covers and Nest are designed for use with the Bilisoft™ 1.0 and 2.0 Phototherapy System, which is used to treat both premature and full-term infants for hyperbilirubinemia.

Available in small and large sizes to match the size of the Light Pad*

* Note: the BiliSoft™ 2.0 Disposable Pad Cover (Flat and Nest) in Size Large also fits the BiliSoft™ Phototherapy System (Legacy) large pad.

Specifications available on request



Not available in all countries

T-Piece Circuit: Disposable resuscitator designed for neonates, provides accurate reliable ventilation and pulmonary resuscitation.



2109061-001
Disposable T Piece
with Mask

Infant Resuscitator T-Piece Circuit Kit			
Description	Mask Size	Part Number	Quantity
Disposable T-Piece	–	2109059-001	20
Disposable T-Piece with Mask	0	2109060-001	10
	1	2109061-001	10

Product compatibility

With a universal connector, it fits with all commonly used manual resuscitator units and does not require an adapter.

Specifications available on request



Panda™ Warmer



StAr

Not available in all countries



LULLABY™ WARMER

Temperature Probes - Vyaire		
Description	Part Number	Quantity
Disposable Temperature Probe	6600-0196-700	50
	6600-0208-700	10
Reusable Temperature probe – Lullaby™ Warmer	2044384-001	1

Sold with equipment only in most countries.

Temperature Probe Patches - Vyaire		
Description	Part Number	Quantity
Patch Heat Reflecting ROHS	0203-1980-300	50
Heat Reflecting Patches	2044385-001	50

Sold with equipment only in most countries.

SHUTTLE™

Transport		
Description	Part Number	Quantity
Patient restraint	M1169540	5

Not available in all countries





Product may not be available in all countries and regions.
Full product technical specifications are available upon request.
Contact a GE Healthcare representative for more information.

Please visit: www.gehealthcare.com/promotional-locations

Data subject to change.

GE Healthcare provides transformational medical technologies and services to meet the demand for increased access, enhanced quality and more affordable healthcare around the world. GE (NYSE: GE) works on things that matter – great people and technologies taking on tough challenges.

From medical imaging, software & IT, patient monitoring and diagnostics to drug discovery, biopharmaceutical manufacturing technologies and performance improvement solutions, GE Healthcare helps medical professionals deliver great healthcare to their patients.

© GE, 2022. All rights reserved.

Products may not be available in all markets. Consult GE Healthcare for current market availability.

General Electric Company reserves the right to make changes in specifications and features shown herein, or discontinue the product described at any time without notice or obligation. Contact your GE representative for the most current information.

GE, the GE Monogram, D-Lite, Corometrics, Omnibed, Giraffe, Shuttle, Panda, MAC, CASE, CardioSoft, Dash, Lullaby, Novii, CARESCAPE, DINAClick, ApexPro, CRITIKON, SOFTCUF, RADIAL-CUF, CLASSIC-CUF, DURA-CUF, Carestation, Aespire, Aestiva, Aisys Avance, Bilisoft, TruSignal, SEER, Entropy, D-lite and D-fend are trademarks of General Electric Company.

Aerogen is a registered trademark of the Aerogen company. ICU Medical Transpac is a registered trademark of ICU Medical, Inc.; Utah/Biotrans is a registered trademark of Utah Medical Products, Inc. Edwards Lifesciences TruWave is a trademark of Edwards Lifesciences Corporation. Argon is a trademark of Argon Medical Devices. AMSORB and G-Can are trademarks of Armstrong Medical. Hytrel is a registered trademark of E. I. du Pont de Nemours and Company or its affiliates. Inteliflo is a registered trademark of Flexicare Medical Ltd. Nellcor and Oximax are trademarks of a Medtronic Company. Masimo and SET are trademarks of Masimo Corporation. All other trademarks are the property of their respective owners.

Reproduction in any form is forbidden without prior written permission from GE.

JB01042XE Rev3