




GE HealthCare

Mac-Lab Altix XI

Hemodynamic Recording Platform





Streamlined workflow and collaborative documenting at point-of-care

Enhanced insight. Improved clinical functionality. Smooth workflow.

The Altix XI edition of GE Healthcare's Mac-Lab™ recording platform helps improve the Cardiac Cath clinical user experience so clinicians can focus on what matters most: patient care.

What's new in Altix XI?

- **INW in Cloud:** Our platform supports both traditional onsite and cloud-based servers
- **Windows® 11 + Windows Server 2025:** An actively supported OS platform that maximizes hardware efficiency
- **Medis Interface:** Collaboration with Medis QFR®—an innovative AI angiography-derived physiology assessment to support procedural decision-making
- **Cybersecurity:** Enhanced security and Software Of Unknown Pedigree (SOUP) updates



Image courtesy of Medis Medical Imaging

Medis QFR Interface

Mac-Lab interfaces with Medis QFR, enhancing workflow and helping clinicians access additional clinical information. Medis QFR enables physiology-guided lesion assessment across the Percutaneous Coronary Intervention (PCI) workflow without pressure wires or hyperemic agents. The offering helps reduce procedural complexity, minimize invasiveness and support a more streamlined patient and clinician experience.

An enhanced user experience

In Cardiac Cath and the growing structural heart space, Mac-Lab provides customized documentation options and a known workflow, which helps to simplify these complex procedures.

An optimal visual experience

With Altix XI, you have the flexibility to add the EP Cath Command Center—designed to take the clinical user experience to the next level of efficiency. The large display and curved monitor, featuring sophisticated video and layouts, put images and tools at users' fingertips, further enhancing productivity.



Provide more data downstream with DirectConnect

As EMR use continues to increase, Altix XI responds with a focus on interoperability with solutions that fit simply into your existing infrastructure. This allows organizations to easily exchange more data with EMR and other required systems. Altix XI includes DirectConnect, which provides a suite of solutions that increase data flexibility, facilitating the easy exchange of data with your EMR and other required systems.

GE HealthCare is ASC ready

As cardiac care continues to evolve beyond traditional settings, interventional procedures are increasingly emerging in Ambulatory Surgery Centers (ASC). This shift represents one of the most significant transformations in cardiovascular care today—driving the need for scalable, efficient and future-ready technologies. ComboLab™ with Altix XI is built to support this evolution by enabling streamlined workflows, interoperability, operational efficiency and advanced clinical capabilities tailored for ASC demand. Backed by a legacy of innovation, GE HealthCare is helping ASCs confidently prepare for the next generation of cardiac procedures today and beyond.



Built on a strong foundation that seamlessly interfaces security measures within an organization's ecosystem, keeping data private.

For more information, please see our security brochure.



New diagnostic forms allow users to capture more diagnostic event information at the point of care



Provides an option to send a PDF export via DICOM to PACS or other downstream solutions



INW in Cloud: Our platform supports both traditional onsite and cloud-based servers



Windows 11 + Windows Server 2025: An actively supported OS platform that maximizes hardware efficiency



GE HealthCare



Contact your GE HealthCare sales representative if you have any questions.

About GE HealthCare Technologies Inc.

GE HealthCare is a trusted partner and leading global healthcare solutions provider, innovating medical technology, pharmaceutical diagnostics and integrated, cloud-first AI-enabled solutions, services and data analytics. We aim to make hospitals and health systems more efficient, clinicians more effective, therapies more precise and patients healthier and happier. Serving patients and providers for more than 125 years, GE HealthCare is advancing personalized, connected and compassionate care, while simplifying the patient's journey across care pathways. Together, our Imaging, Advanced Visualization Solutions, Patient Care Solutions and Pharmaceutical Diagnostics businesses help improve patient care from screening and diagnosis to therapy and monitoring. We are a \$19.7 billion business with approximately 53,000 colleagues working to create a world where healthcare has no limits.

Follow us on [LinkedIn](#), [X](#), [Facebook](#), [Instagram](#) and [Insights](#) for the latest news, or visit our website <https://www.gehealthcare.com> for more information.

© 2026 GE HealthCare. GE is a trademark of General Electric Company used under trademark license. Subscription offerings, select product features and third-party interfaces may not be available in all countries. Please contact your sales representative for availability. Mac-Lab and ComboLab are trademarks of GE HealthCare.

QFR is a registered trademark of QFR Solutions bv. Microsoft Windows and Microsoft Office are registered trademarks of Microsoft Corporation in the United States and/or other countries. All third-party trademarks remain the property of their respective holders.

September 2026
JB37754XX