




OrthoWorks


High-resolution imaging

2D and 3D MERGE

Alternative orthopedic tissue contrast



3D MERGE
Scaphoid Fracture Staging



3D MERGE



2D MERGE



Routine knee

AIR™ Multi-Purpose Coil, Medium*



FSE PD FatSat
0.4 x 0.5 x 3 mm

FSE PD
0.3 x 0.4 x 3 mm




FSE PD FatSat
0.3 x 0.4 x 3 mm




T1 FSE
0.4 x 0.6 x 3 mm




Routine knee




[FSE](#) PD
0.3 x 0.4 x 2 mm



3D [Cube](#) PD FatSat
0.6 x 0.6 x 1.2 mm




[FSE](#) PD FatSat
0.4 x 0.6 x 3 mm
1:25 min




[FSE](#) PD FatSat
Meniscus
0.6 x 0.6 x 1.0 mm


Osteochondritis dessicans




3D MENSA
0.5 x 0.5 x 1 mm
TR/TE = 18.9 ms / 9.6 ms
FOV = 14 cm
3:31 min



Osteochondritis dessicans



T1 FSE
0.4 x 0.5 x 3 mm
TR/TE = 650 ms / 14 ms
2:52 min



FSE PD FatSat
0.4 x 0.5 x 3 mm
TR/TE = 3168 ms / 33.8 ms
2:39 min

