



OrthoWorks

Isotropic 3D imaging


Isotropic Cube in the knee

10-minute knee exam

Cube PD

[Cube](#) PD and PD FatSat
0.5 x 0.5 x 0.5 mm
total scan time 10:00 min
[HyperSense](#)


Cube PD FatSat




Hip labral tear and osseous cyst

Clinical case:


- Intraosseous acetabular cyst (yellow arrow)
- Anterior labral tear (red arrow)
- Normal labrum (blue arrow)





T1 FSE
0.4 x 0.4 x 3 mm
1:48 min



PROPELLER PD FatSat
0.7 x 0.7 x 3 mm
2:09 min



3D MERGE
0.8 x 0.8 x 0.8 mm
5:11 min



Isotropic Cube in the spine

AIR™ Anterior Array Coil (30ch) + Posterior Array coil

T2 [Cube](#)
0.8 x 0.8 x 0.8 mm
3:42 min
[HyperSense](#)
[HyperCube](#)

