




OrthoWorks

Peripheral nerve imaging

Brachial plexus

AIR™ Anterior Array Coil Coil + Posterior Array + HNU

[Cube](#) STIR
0.7 x 1.3 x 1.2 mm
2:52 min
[HyperSense](#)
[HyperCube](#)




Brachial plexus and radial nerves

AIR™ Anterior Array Coil Coil (30ch) + 48ch Head Coil


[Cube STIR](#)
[HyperSense](#)
[HyperCube](#)



Lower back pain



PROPELLER STIR MB
0.9 x 0.9 x 4.0 mm
3:05 min



DWPEPI MUSE
b750
2.8 x 1.2 x 4.0 mm
4:00 min



DWPEPI MUSE - inverted

