




OrthoWorks


Large FOV and long-bone
imaging with AIR™ Coils

Long-bone imaging – Femur


AIR™ Anterior Array Coil (30ch) + Posterior Array (60ch)




T1 FSE



T1 FSE




FSE STIR



T1 FSE Flex - water


FOV 46cm



FSE STIR
56 cm FOV
2 stations

Long-bone imaging – upper and lower arm

AIR™ Anterior Array Coil (30ch) x 2 + Posterior Array (60ch)



[FSE PD Flex](#) – in-phase, out of phase, water, fat





Long-bone imaging – upper and lower arm

AIR™ Anterior Array Coil (30ch) x 2 + Posterior Array (60ch)




Lower leg

AIR™ Anterior Array Coil (30ch) + Posterior Array (60ch)




Water Image




In Phase


T2 FSE Flex
0.9 x 1.2 x 5 mm
2:54 min



T2 FSE Flex - water
1.2 x 0.8 x 5 mm
2:31 min



DWI with ADC Map



LAVA Flex
1.2 x 0.9 x 2 mm
2:45 min

