



OrthoWorks

Off-isocenter, small FOV
extremity imaging


Fingers



T2 FSE Flex Water
0.6 x 0.5 x 2 mm
1:27 min



T1 FSE
0.3 x 0.4 x 2 mm
50 sec



T2 FSE FatSat
0.4 x 0.5 x 2 mm
1:35 min



T1 FSE FatSat
0.3 x 0.5 x 2 mm
1:27 min

QED T/R Hand/Wrist coil



Elbow


AIR™ Multi-Purpose Coil, Medium*



[Cube](#) PD FatSat
[HyperSense](#) and [HyperCube](#)
0.4 x 0.5 x 0.6 mm
5:01 min




[FSE](#) PD
0.3 x 0.4 x 2.5 mm
3:13 min



[Cube](#) PD FatSat
[HyperSense](#)
0.6 x 0.6 x 0.4 mm
5:13 min



Hand/wrist



FSE PD
0.2 x 0.2 x 1.2 mm



FSE PD
0.2 x 0.2 x 2 mm




FSE PD
0.2 x 0.2 x 2 mm

QED T/R Hand/Wrist coil



Elbow


AIR™ Multi-Purpose Coil, Medium*



FSE PD [Flex](#)
0.4 x 0.6 x 3 mm
3:26 min



T2 FSE
0.3 x 0.3 x 3 mm
4:25 min



FSE PD [Flex](#)
0.4 x 0.6 x 3 mm
3:20 min



