



GE Healthcare

Anesthesia O₂ Cells

Oxygen sensors are critical to Anesthesia machine performance and can be used as back-up if there is no gas module present.

GE Healthcare oxygen sensors are the only approved devices for use with GE Anesthesia machines.

Aisys™, Avance™, Aespire™,
Aestiva™ Series



2082006-001



M-10

Carestation™
600/750* Series



2071018-001-S

Our equipment.
Our accessories.
Simply compatible.

You can trust GE Healthcare accessories and consumables to keep your clinical devices working efficiently.

- Simplified ordering process
- Fast delivery
- One contact point for expert advice
- Quality support at every step

* Carestation 750 not cleared or approved by US FDA. Not for sale in the United States

[gehealthcare.com/products/accessories-and-supplies](https://www.gehealthcare.com/products/accessories-and-supplies)

PRODUCT SPECIFICATIONS

Materials

2071018-001-S

O ₂ Cell Cable:	Polyvinyl Chloride (PVC), Copper
O ₂ Cell Seat:	Acrylonitrile Butadiene Styrene (ABS), Polycarbonate (PC)
Adapter:	Polyphenylsulfone (PPSU)
O-ring:	Silicone

2082006-001 / M-10

O-ring:	Vinyl-methyl-silicone (VMQ)
Housing:	ABS

Materials composition

Flow sensors are not made with BPA, natural rubber latex, PVC, DEHP or other phthalates.

Operating environment

2071018-001-S

Temperature:	10° to 40°C (50° to 104°F)
Humidity:	15% to 95%, non-condensing
Pressure:	618 to 1250 hPa

2082006-001 / M-10

Temperature:	10° to 40°C (50° to 104°F)
Humidity:	15% to 95%, non-condensing
Pressure:	488 to 937 mmHg

Storage

Temperature:	-15° to 50°C (5° to 122°F)
Humidity:	10% to 95% relative humidity
Pressure:	425 to 937 mmHg

Disposal considerations

Handle and dispose of O₂ cells according to site biohazard policies, do not incinerate.

* Carestation 750 not cleared or approved by US FDA. Not for sale in the United States

© 2021 General Electric Company – All rights reserved.

GE Healthcare reserves the right to make changes in specifications and features shown herein, or discontinue the product described at any time without notice or obligation. Contact your GE Healthcare representative for the most current information. GE, the GE Monogram, Aestiva, Aespire, Aisys, Avance and Carestation are trademarks of General Electric Company. GE Healthcare, a division of General Electric Company. GE Medical Systems, Inc., doing business as GE Healthcare.

DOC2465644 1/21

Packaging materials

Primary:	Plastic Clam Shell with Paper
Secondary:	Cardboard box

Usage and Shelf life

2071018-001-S: The O₂ cell and its packaging have an established shelf life of 1 year. They have a useful life of 1,000,000 volume % O₂ hours. Typical use case is approximately 3 years.

2082006-001 / M-10: The O₂ cell and its packaging have an established shelf life of 6 months. They have useful life of more than 750,000 volume % O₂ hours. Typical use case is approximately 3 years.

Regulatory standards

RoHS (Restriction of Hazardous Substances) compliant

EN ISO 80601-2-55:2011(E) - Part 2-55

CE Mark

Product compatibility

Anesthesia O₂ Cells are compatible with the following GE Healthcare equipment:

- 2071018-001-S for use with Carestation 600/750* series
- 2082006-001 and M-10 for use with Aestiva, Aespire, Aisys (CS²) and Avance (CS²) series

PART LIST

Part Number	Description	QTY
2071018-001-S	O ₂ Sensor with adapter and O-ring, CBS	1
2082006-001	O ₂ Sensor and Cable, ABS	1
M-10	O ₂ Sensor, ABS	1

Manufactured by:



IT Gambert GmbH
Hinter dem Chor 21
23966 Wismar, Germany