



Clinical accessories portfolio

CARESCAPE Monitors V100 and VC150



NON-INVASIVE BLOOD PRESSURE

A patient's blood pressure (BP) is often the first step in determining the best clinical course of action. That's why quality data acquisition through your BP cuff is so essential. GE CRITIKON® BP cuffs from GE Healthcare are designed and tested for performance accuracy – from the first touch to the final measurement. By starting with accurate information, you can make well-informed decisions that may ultimately enhance patient outcomes.

For a complete list of cuffs compatible with your equipment, contact your GE Healthcare sales representative.



HOSES

CRITIKON AIR HOSE	PATIENT	LENGTH	PART NUMBER	QUANTITY
2-Tube DINACLICK™ Connector	Adult - Pediatric	3,6 m	2058205-001	1
		7,2 m	2058205-002	
2-Tube Neo-Snap™ Connector	Neonate	3,6 m	2058206-001	1





CUFFS

NEONATAL

Our neonatal BP cuffs are designed with the performance accuracy you trust and the gentle comfort you expect. Each of our 5 sizes is offered in a extra soft absorbent material (SOFT-CUF™ range) or easy-to-clean vinyl (CLASSIC-CUF™ range). Economical enough for pairing with a single patient throughout hospital stay.

20 cuffs/box			2-TUBE NEO-SNAP	
SIZE	COLOR	LIMB CIRCUMFERENCE	CLASSIC-CUF	SOFT-CUF
1	Orange	3-6 cm	CLA-N1-2B	SFT-N1-2B
2	LT. Blue	4-8 cm	CLA-N2-2B	SFT-N2-2B
3	Green	6-11 cm	CLA-N3-2B	SFT-N3-2B
4	Navy	7-13 cm	CLA-N4-2B	SFT-N4-2B
5	Pink	8-15 cm	CLA-N5-2B	SFT-N5-2B

ASSORTMENTS NEO-SNAP CONNECTOR

PRODUCT	2-TUBE NEO-SNAP
20 cuffs SOFT-CUF : 2 Neo #1, 3 Neo #2, 5 Neo #3, 5 Neo #4, 5 Neo #5	2059305-002
20 cuffs CLASSIC-CUF : 2 Neo #1, 3 Neo #2, 5 Neo #3, 5 Neo #4, 5 Neo #5	2059306-002
3 cuffs CLASSIC-CUF : 1 Neo #3, 1 Neo #4, 1 Neo #5	2059306-003



Playful numeric and colorcoded animal designs to help you find the right size at a glance.



Actual illustrations may slightly vary.



CUFFS

DURA-CUF

CRITIKON™ DURA-CUF™ is gentle and designed for durability. Constructed with materials not made with PVC, making them safer for the environment, and natural rubber latex, making them safer for your patient.



Easy-to-use, fast **DINAClick connector** enables nurses to link the cuff to the NIBP hose quicker than a pair of traditional tube connectors. The audible click lets users know there is a proper connection.



5 cuffs/box

SIZE	COLOR	LIMB CIRCUMFERENCE	2-TUBE DINAClick
Infant	Orange	8-13 cm	DUR-P1-2A-5S
Child	Green	12-19 cm	DUR-P2-2A-5S
Small Adult	Royal Blue	17-25 cm	DUR-A1-2A-5S
Small Adult Long	Royal Blue	17-25 cm	DUR A1-2A-L-5S
Adult	Navy	23-33 cm	DUR-A2-2A-5S
Adult long	Navy	23-33 cm	DUR-A2-2A-L-5S
Large Adult	Wine	31-40 cm	DUR-A3-2A-5S
Large Adult Long	Wine	31-40 cm	DUR-A3-2A-L-5S
Thigh	Brown	38-50 cm	DUR-T1-2A-5S



The first GE CRITIKON™ blood pressure cuff, DURA-CUF™, combines comfort and durability.



ASSORTMENTS DINAClick CONNECTOR

PRODUCT

3 cuffs : 1 Small Adult, 1 Adult, 1 Large Adult
6 cuffs : 2 Infant, 3 Child, 1 Small Adult
6 cuffs : 1 Small Adult, 3 Adult, 1 Adult Long, 1 Large Adult
6 cuffs : 1 Infant, 1 Child, 1 Small Adult, 1 Adult, 1 Large Adult, 1 Thigh

2-TUBE
DINAClick

2059301-501
2059301-502
2059301-503
2059301-504



CUFFS

SOFT-CUF



Our most popular cuff, the CRITIKON™ SOFT-CUF™ is designed with soft, absorbent material with rounded corners for maximizing patient comfort during extended stays at the hospital. These cuffs are ideal for patients with longer term stays or for protocols where the cuff travels with the patient through their total length of time in the hospital. It can be cleaned and reused on multiple patients, but is priced economically so it can be disposed of after a single patient use, helping you affordably prevent infections.

20 cuffs/box



SIZE	COLOR	LIMB CIRCUMFERENCE	2-TUBE DINACLICK
Infant	Orange	8-13 cm	SFT-P1-2A-5S
Child	Green	12-19 cm	SFT-P2-2A-5S*
Child Long	Green	12-19 cm	SFT-P2-2A-L-5S
Small adult	Royal Blue	17-25 cm	SFT-A1-2A-5S*
Small adult Long	Royal Blue	17-25 cm	SFT-A1-2A-L-5S
Adult	Navy	23-33 cm	SFT-A2-2A-5S*
Adult Long	Navy	23-33 cm	SFT-A2-2A-L-5S*
Large Adult	Pink	31-40 cm	SFT-A3-2A-5S*
Large Adult Long	Pink	31-40 cm	SFT-A3-2A-L-5S
Thigh	Brown	38-50 cm	SFT-T1-2A-5S

Part number with an "-INT" suffix, denoting that the Instructions For Use (IFU) is available in multiple languages.

ASSORTMENTS DINACLICK CONNECTOR



PRODUCT	2-TUBE DINACLICK
20 Cuffs : 5 Small Adult, 10 Adult, 5 Large Adult	2059303-501
20 Cuffs : 6 Infant, 10 Child, 4 Small Adult	2059303-502
20 Cuffs : 3 Small Adult, 7 Adult, 3 Adult Long, 3 Large Adult, 4 RADIAL-CUF (forearm)	2059303-503
20 Cuffs : 2 Infant, 2 Child, 2 Small Adult, 6 Adult, 2 Adult Long, 2 Large Adult, 2 RADIAL-CUF (forearm), 2 Thigh	2059303-504
4 cuffs: 1 Small Adult, 1 Adult, 1 Large Adult, 1 RADIAL-CUF	2059303-506



Validated quality

CRITIKON BP cuffs are validated with GE Healthcare monitors, and are part of an integrated clinical system.



Better infection control

Source of contamination ignored or underestimated, no other equipment is used more often than a blood pressure cuff without proper cleaning. GE's CRITIKON SOFT-CUF™ cuff provide the optimal solution to protect your patients.



DINACLICK Connector

Easy-to-use, fast DINACLICK connector enables nurses to link the cuff to the monitor faster than other dual-tube connectors.



Sized for comfort and accuracy

Using the wrong size BP cuff can affect determination accuracy. A full range of CRITIKON™ cuff sizes are available, from neonate, pediatric, adult and obese.

NON-INVASIVE BLOOD PRESSURE



PRECISE MEASURES FOR INFORMED CHOICES

GE Healthcare offers **the first BP forearm cuff** clinically validated¹ for use on **obese or overly muscular adult patients**.



By adapting to the anatomy of obese patients, the RADIAL-CUF cuff optimizes the accuracy of the **blood pressure** measurements of this particularly risky patient population.

77% of patients also consider it **more comfortable** than traditional blood pressure cuffs *.



SIZE	COLOR	LIMB CIRCUMFERENCE	2-TUBE DINACLICK
OBESSE ADULT	PURPLE	26-36 CM (FOREARM)	SFT-F1-2A-5S



1. The first clinically validated NIBP Cuff for providing accurate measurements from the forearm, GE Healthcare, 2013.

* GE Healthcare sponsored and compensated Clinimark LLC for managing and executing the clinical portions of this study (March 2008/ December 2011).



SpO₂ SENSORS

Through the legacy of Datex-Ohmeda *, GE SpO₂ technology is a story of innovation that has been going on for more than 20 years. TruSignal™ sensors allow continuous non-invasive monitoring of SpO₂ and heart rate.

We offer a wide range of SpO₂ sensors, designed and tested for their performance and always serving the comfort of your patients.



TECHNOLOGY

GE TRUSIGNAL

For your CARESCAPE V100 and VC150 monitors equipped with TruSignal SpO₂ technology.



REUSABLE INTEGRATED SENSORS					
SENSORS	PATIENT	PATIENT WEIGHT	LENGTH	PART NUMBER	QUANTITY
TruSignal Integrated Finger sensor	Adult, Pediatric	≥ 20 kg	2 m	TS-F2-GE	1
			4 m	TS-F4-GE	
TruSignal Integrated FingerTip	Adult	> 30 kg	4 m	TS-SA4-GE	1
	Pediatric	15 - 30 kg	3 m	TS-SP3-GE	
TruSignal Integrated Ear sensor	Adult, Pediatric	≥ 10 kg	2 m	TS-E2-GE	1
			4 m	TS-E4-GE	




INTERCONNECT CABLES

PRODUCT	LENGHT	PART NUMBER	QUANTITY
TruSignal interconnect cable with GE connector	3 m	TS-G3	1



REUSABLE SENSORS (REQUIRES COMPATIBLE INTERCONNECT CABLE)

SENSORS	PATIENT	PATIENT WEIGHT	LENGHT	PART NUMBER	QUANTITY
TruSignal Finger sensor	Adult, Pediatric	≥ 20 kg	1 m	TS-F-D	1
TruSignal FingerTip sensor	Pediatric	15-30 kg	1 m	TS-SP-D	1
	Adult	> 30 kg	1 m	TS-SA-D	
TruSignal Ear sensor	Adult, Pediatric	≥ 10 kg	1 m	TS-E-D	1
TruSignal Wrap sensor	Adult, Pediatric	≥ 3 kg	1 m	TS-W-D	1
TruSignal Sensitive Skin sensor	Adult, Pediatric, Infant, neonate	No limit	1 m	TS-SE-3	3



DISPOSABLE SENSORS (REQUIRES COMPATIBLE INTERCONNECT CABLE)

SENSORS	PATIENT	PATIENT WEIGHT	LENGHT	PART NUMBER	QUANTITY
TruSignal sensor	Pediatric, infant	≥ 20 kg	0,3 m	TS-AP-10	10
				TS-AP-25	25
TruSignal sensor	Adult, pediatric, infant, neonate	No limit	0,5 m	TS-AF-10	10
				TS-AF-25	25
TruSignal Adhesive Wrap sensor	Adult	> 20 kg	0,3 m	TS-AAW-10	10
				TS-AAW-25	25
TruSignal Adhesive Wrap sensor	Pediatric, infant	3 - 20 kg	0,3 m	TS-PAW-10	10
				TS-PAW-25	25

GE SPO₂ ACCESSORIES

PRODUCT	USE WITH	SIZE	PART NUMBER	QUANTITY
Wide adhesive replacement tape	TS-W-D, TS-SE-3, TS-AF-10, TS-AF-25	-	TS-RTW-100	100
Replacement adhesive tape (blue)	TS-AF-10, TS-AF-25	-	TS-RT-100	100
Replacement adhesive tape (bears)			TS-RTB-100	
Foam wrap replacement	TS-SE-3	Small	TS-RWS-24	24
		Medium	TS-RWM-24	
	TS-W-D	Large	TS-RWL-24	24
Infant foam sandal (2 to 7.7 kg)	TS-SE-3	-	TS-SND-3	3



TECHNOLOGY

NELLCOR™

For your CARESCAPE V100 and VC150 monitors equipped with Nellcor SpO₂ technology.

INTERCONNECT CABLES


PRODUCT	LENGTH	PART NUMBER	QUANTITY
Interconnect cable with GE connector and Nellcor OxiMax™ connector	3 m	2021406-001	1

REUSABLE NELLCOR SENSORS



SENSORS	MANUFACTURER REFERENCE	PATIENT	PATIENT WEIGHT	PART NUMBER	QUANTITY
Nellcor DuraSensor Finger Probe	DS 100A	Adult	< 30 kg	70124021	1
Nellcor SpO ₂ Multisite Sensor	D-YS	Adult	> 1 kg	70124033	1
Nellcor SpO ₂ Sensor	OXI-A/N	Neonate	< 3 kg	70124035	1
		Adult	> 40 kg		
Nellcor SpO ₂ Sensor	OXI-P/I	Pediatric, infant	3-40 kg	414248-001	1 delivered with 50 adhesive strips

DISPOSABLE NELLCOR SENSORS



SENSORS	MANUFACTURER REFERENCE	PATIENT	PATIENT WEIGHT	PART NUMBER	QUANTITY
Nellcor SpO ₂ Finger Sensor	MAX-A	Adult	> 30 kg	70124027	24
Nellcor SpO ₂ Teo Sensor	MAX-N	Adult	> 40 kg	70124032	24
		Neonate	< 3 kg		

FIXATION CLIP



PRODUCT	PART NUMBER	QUANTITY
Nellcor SpO ₂ Ear-Clip D-YSE, for use with 70124033	70124034	1

ADHESIVE TAPE



PRODUCT	PART NUMBER	QUANTITY
Nellcor adhesive tape (ADH-A/N) - Ref: 70124035 (OXI-A/N)	2016130-001	100



TECHNOLOGY

MASIMO™

For your CARESCAPE V100 and VC150 monitors equipped with MASIMO SpO₂ technology.

LNCS AND RD INTERCONNECT CABLES

PRODUCT	DEVICE	LENGTH	PART NUMBER	QUANTITY
Masimo LNC-10 interconnect cable	V100	3 m	2027263-002	1
Masimo LNC-10 interconnect cable	VC150	3 m	2059006-001	1
Masimo RD SET GE cable	V100	3,6 m	MASIMO-RD-4085	1
Masimo RD SET GE cable	VC150	3,6 m	MASIMO-RD-4104	1



REUSABLE LNCS AND RD SENSORS

SENSORS	MANUFACTURER REFERENCE	PATIENT	PATIENT WEIGHT	PART NUMBER	QUANTITY
Masimo Finger sensor	LNCS DC-I	Adult	> 30 kg	2027258-001	1
	LNCS DC-IP	Pediatric	10-50 kg	2027259-001	
Masimo TipClip Ear sensor	LNCS TC-I	Adult	> 30 kg	2027261-001	1
Masimo Finger sensor	RD SET DC-I	Adult	> 30 kg	MASIMO-RD-4050	1
	RD SET DC-IP	Pediatric	10-50 kg	MASIMO-RD-4051	



DISPOSABLE LNCS AND RD SENSORS

SENSORS	MANUFACTURER REFERENCE	PATIENT	PATIENT WEIGHT	PART NUMBER	QUANTITY
Masimo Adhesive Finger Sensor	LNCS Neo	Adult	> 40 kg	2051069-001	20
Masimo Adhesive Teo Sensor		Neonate	< 3 kg		
Masimo Adhesive Finger Sensor	RD SET Adt	Adult	> 30 kg	MASIMO-RD-4000	20
Masimo Adhesive Finger Sensor	RD SET Neo	Adult	> 40 kg	MASIMO-RD-4003	20
Masimo Adhesive Finger or Teo Sensor		Neonate	< 3 kg		
Masimo Adhesive Finger or Teo Sensor	RD SET NeoPt	Premature Baby	< 1 kg	MASIMO-RD-4004	20

ADHESIVE TAPE

PRODUCT	PART NUMBER	QUANTITY
Masimo adhesive tape for SpO ₂ adhesive sensor (RD SET Neo)	MASIMO-RD-4047	102

[illegible]



About GE Healthcare

GE Healthcare is a leading global medical technology and digital solutions innovator, GE Healthcare enables clinicians to make faster, more informed decisions through intelligent devices, data analytics, applications and services, supported by its Edison intelligence platform. With over 100 years of healthcare industry experience and around 50,000 employees globally, the company operates at the center of an ecosystem working toward precision health, digitizing healthcare, helping drive productivity and improve outcomes for patients, providers, health systems and researchers around the world. Follow us on Facebook, LinkedIn, Twitter and Insights , or visit our website www.gehealthcare.com for more information.

The list of accessories displayed in this catalogue is not exhaustive, Contact your GE Healthcare representative to obtain the list of all accessories compatible with your equipment.

www.gehealthcare.co.uk

Product may not be available in all countries and regions. Contact a GE Healthcare Representative for more information.

Data subject to change.

GE, the GE Monogram, CARESCAPE, CRITIKON, DURA-CUF, SOFT-CUF, DINACLICK, TruSignal, D-Fend and D-Lite are trademarks of General Electric

Company.

*All other third-party trademarks are the property of their respective owners.

Reproduction in any form is forbidden without prior written permission from GE.

Nothing in this material should be used to diagnose or treat any disease or condition. Readers must consult a healthcare professional.

© 2021 General Electric Company. All rights reserved.

JB00987XE Rev3