



Neo-Snap

GE Healthcare's neonatal connection system for single tube and dual tube CRITIKON™ blood pressure cuffs

The Neo-Snap™ connector for neonatal blood pressure cuffs is designed to standards that help prevent unintended disconnect. The Neo-Snap connector provides a secure connection with audible click and tactile feedback.

Neo-Snap connector benefits

- **Designed for safety.** Connector meets international standard for improved patient safety
- **Secure connection.** Prevents accidental disconnection
- **Easy-to-use.** Smooth, 360 degree connection that does not have to be aligned to make secure connection – so you can focus more on baby
- **Patient safety.** Not made with natural rubber latex or DEHP



To learn more about GE Healthcare's comprehensive cuff solutions, visit us at [gehealthcare.co.uk/products/patient-monitors/non-invasive-blood-pressure-cuffs](https://www.gehealthcare.co.uk/products/patient-monitors/non-invasive-blood-pressure-cuffs)

NIBP cuff part numbers

Neo-Snap			2-Tube Neo-Snap Connector		1-Tube Neo-Snap Connector	
						
Size	Color	Limb Circumference	SOFT-CUF	CLASSIC-CUF	SOFT-CUF	CLASSIC-CUF
1	Orange	3-6 cm	SFT-N1-2B	CLA-N1-2B	SFT-N1-1B	CLA-N1-1B
2	Light Blue	4-8 cm	SFT-N2-2B	CLA-N2-2B	SFT-N2-1B	CLA-N2-1B
3	Green	6-11 cm	SFT-N3-2B	CLA-N3-2B	SFT-N3-1B	CLA-N3-1B
4	Navy	7-13 cm	SFT-N4-2B	CLA-N4-2B	SFT-N4-1B	CLA-N4-1B
5	Rose	8-15 cm	SFT-N5-2B	CLA-N5-2B	SFT-N5-1B	CLA-N5-1B

Neonatal assortment pack

Neonatal assortment pack	Cuff type	Part number	Description
2-Tube Neo-Snap Connectors	SOFT-CUF	2059305-002	2 Neo #1, 3 Neo #2, 5 Neo #3, 5 Neo #4, 5 Neo #5 (20 cuffs/box)
	CLASSIC-CUF	2059306-002	2 Neo #1, 3 Neo #2, 5 Neo #3, 5 Neo #4, 5 Neo #5 (20 cuffs/box)
	CLASSIC-CUF	2059306-003	1 Neo #3, 1 Neo #4, 1 Neo #5 (3 cuffs/pack)
1-Tube Neo-Snap Connectors	SOFT-CUF	2059305-001	2 Neo #1, 3 Neo #2, 5 Neo #3, 5 Neo #4, 5 Neo #5 (20 cuffs/box)
	CLASSIC CUF	2059306-001	2 Neo #1, 3 Neo #2, 5 Neo #3, 5 Neo #4, 5 Neo #5 (20 cuffs/box)

2-tube Neo-Snap hoses

	Description	Length	Part number	Compatibility
	Care Hose	8 FT/2.4 M	2058204-002	B20, B40, B105, B125, Dash 3000, 4000, 5000, Tram, CARESCAPE Patient Data Module
		12 FT/3.6 M	2058204-003	
	CRITIKON Hose	12 FT/3.6 M	2058206-001	V100/VC150 DINAMAP™ Procure
	CARESCAPE™ ONE	8 FT/2.4 M	2089791-002	NIBP HOSE CARESCAPE ONE NEO-SNAP
	CARESCAPE ONE	12 FT/3.6 M	2089791-003	

Neo-Snap connectors

	Description	Part number	Quantity/package
	Y, 1 Male Barb 1/8" to Neo-Snap Male Plug	2059290-003	10/PACK
	Y, 1 Male Barb 3/32" to Neo-Snap Male Plug	2059290-004	10/PACK
	2T Neo-Snap Male Plug 0.22 Centers Compatibility: B20, B40, B105, B125, Dash 3000, 4000, 5000, Tram, CARESCAPE Patient Data Module and Carescape One	2059293-001	10/PACK
	2T Neo-Snap Male Plug 0.31 Centers Compatibility: V100/VC150 DINAMAP™ Procure	2059293-002	10/PACK
	1T Neo-Snap Male Plug	2059293-003	10/PACK



© 2021 General Electric Company – All rights reserved.

GE Healthcare reserves the right to make changes in specifications and features shown herein, or discontinue the product described at any time without notice or obligation. Contact your GE Healthcare representative for the most current information. GE, the GE Monogram, Neo-Snap, CRITIKON, DINAMAP, and CARESCAPE are trademarks of General Electric Company. GE Healthcare, a division of General Electric Company. GE Medical Systems, Inc., doing business as GE Healthcare.

JB00924XE 7/21