

	Order code	Dispensing system	Viscosity low  high	Impression technique	Recommended max. working time (min:sec)	Intra-oral setting time at 37°C (min:sec)
<b>Super Quick Setting Materials</b>						
3M™ Impregum™ Penta™ Super Quick (Heavy Body)	69406			1-step	0:45	2:00
3M™ Impregum™ Penta™ Super Quick (Medium Body)	69385			Monophase or 1-step	0:45	2:00
3M™ Impregum™ Super Quick (Medium Body)	69417			Monophase or 1-step	0:45	2:00
3M™ Impregum™ Super Quick (Light Body)	69379			1-step	0:45	2:00
3M™ Impregum™ Super Quick (Ultra Light Body)	69418			1-step	0:45	2:00
<b>Quick Setting Materials</b>						
3M™ Impregum™ Penta™ Soft Quick (Heavy Body)*	31769			1-step	1:00	3:00
3M™ Impregum™ Penta™ Soft Quick (Medium Body)	31765			Monophase or 1-step	1:00	3:00
3M™ Impregum™ Garant™ Soft Quick (Medium Body)	31788			Monophase or 1-step	1:00	3:00
3M™ Impregum™ Soft Quick (Light Body)	31767			1-step	1:00	3:00
<b>Regular Setting Materials</b>						
3M™ Impregum™ Penta™ H DuoSoft™ (Heavy Body)	31740			1-step	2:30	3:30
3M™ Impregum™ Penta™ Soft (Medium Body)	31730			Monophase or 1-step	2:45	3:15
3M™ Impregum™ Garant™ Soft (Medium Body)	31785			Monophase or 1-step	1:45	4:15
3M™ Impregum™ Penta™ (Medium Body)	31644			Monophase	2:45	3:15
3M™ Impregum™ Penta™ L DuoSoft™ (Light Body)	31745			1-step	3:15	3:15
3M™ Impregum™ Garant™ L DuoSoft™ (Light Body)	31751			1-step	2:00	3:30

\* 3M™ Impregum™ Penta™ Soft Quick Heavy Body (31769) is not available in New Zealand.

	Order code	Dispensing system	Viscosity low  high	Max. working time (23°C)	Max. intra-oral syringe time (37°C) 	Intra-oral setting time (37°C) 
<b>Super Quick Setting Materials</b>						
3M™ Imprint™ 4 Penta™ Super Quick Heavy	71485			1:15	–	1:15
3M™ Imprint™ 4 Super Quick Heavy	71493			1:15	–	1:15
3M™ Imprint™ 4 Super Quick Ultra-Light	71491			–	0:35	1:15
3M™ Imprint™ 4 Super Quick Light	71490			–	0:35	1:15
3M™ Imprint™ 4 Super Quick Regular	71489			–	0:35	1:15
<b>Regular Setting Materials</b>						
3M™ Imprint™ 4 Penta™ Putty	71486			1:30	–	2:30
3M™ Imprint™ 4 Penta™ Heavy	71484			2:00	–	2:00
3M™ Imprint™ 4 Heavy	71492			2:00	–	2:00
3M™ Imprint™ 4 Light	71488			–	1:00	2:00
3M™ Imprint™ 4 Regular	71487			–	1:00	2:00
<b>Alginate Replacement</b>						
3M™ Imprint™ 4 Penta™ Preliminary Super Quick	71522			1:00	–	1:30
3M™ Imprint™ Preliminary Super Quick	71540			1:00	–	1:30
3M™ Imprint™ 4 Penta™ Preliminary	71521			1:30	–	2:00
3M™ Imprint™ Preliminary*	71539			1:30	–	2:00

\*3M™ Imprint™ Preliminary (71539) is not available in New Zealand.