

**3M** Science.  
Applied to Life.™

## 3M™ Impregum™ Super Quick Polyether Impression Material

# Capture every detail in 2 minutes? Yes, it can.




Now also  
available in  
cartridge  
delivery!



### Quicker set.

With just a two minute setting time, Impregum Super Quick Material lives up to its name and allows for faster impressions in cases like single-unit crowns, implants, inlays/onlays and smaller bridges up to 3 units.



### 100 % Polyether. Trusted precision.

Impregum Super Quick Materials demonstrate excellent performance in moist situations and exhibit great flow behavior, leading to higher levels of accuracy and detailed, high-precision restorations.

The nature of polyether allows you to get the full 45 seconds of working time, avoiding inaccuracies associated with premature setting times common with VPS materials.




### Make your patients feel comfortable.

An improved fresh, minty taste together with a short time in the mouth significantly increases your patients' comfort.

## 3M™ Impregum™ Super Quick Polyether Impression Material

- For monophasic and 1-step impression techniques
- Perfect for impressions for single-unit crowns, implants, inlays/onlays and smaller bridges up to 3 units
- Available in heavy, medium, light and ultra light body consistencies
- Impregum Super Quick Material is scannable with the current generation of scanners for a more efficient workflow.



# The fastest way to polyether precision for smaller cases.

	Dispensing system	Impression technique	Single-unit crown	Small implant (≤ 2 units)	Inlay/onlay (≤ 2 units)	Small bridge (≤ 3 units)	Recommended max. working time (min:sec)	Intra-oral setting time at 37°C (min:sec)
3M™ Impregum™ Penta™ Super Quick (Medium Body)		Monophase or 1-step	✓	✓	✓	✓	0:45	2:00
3M™ Impregum™ Super Quick (Medium Body)		1-step	✓	✓	✓	✓	0:45	2:00
3M™ Impregum™ Penta™ Super Quick (Heavy Body)		1-step	✓	✓	✓	✓	0:45	2:00
3M™ Impregum™ Super Quick (Light Body)		1-step	✓	✓	✓	✓	0:45	2:00
3M™ Impregum™ Super Quick (Ultra Light Body)		1-step	✓	✓	✓	✓	0:45	2:00

## Ordering information



3M™ Impregum™ Super Quick Medium Body/Light Body/Ultra Light Body Introductory Kit (69421)



3M™ Impregum™ Penta™ Super Quick Medium Body/Light Body Introductory Kit (69415)



3M™ Impregum™ Penta™ Super Quick Heavy Body/Light Body Introductory Kit (69407)

### Item no. Product information

#### Introductory Kit for 3M™ Garant™ Dispenser

**69421 3M™ Impregum™ Super Quick Medium Body/Light Body/Ultra Light Body**  
2 3M™ Impregum™ Super Quick MB cartridges – 50 ml each, 1 3M™ Impregum™ Super Quick LB cartridge – 50 ml, 1 3M™ Impregum™ Super Quick ULB cartridge – 50 ml, 5 3M™ Intra-oral Syringes – Purple, 5 3M™ Garant™ Mixing Tips – Purple, 5 3M™ Garant™ Intra-oral Tips – White, 1 3M™ Polyether Tray Adhesive

#### Introductory Kits for 3M™ Pentamix™ 3 Automatic Mixing Unit

**69414 3M™ Impregum™ Penta™ Super Quick Medium Body**  
1 base paste MB – 300 ml, 1 catalyst MB – 60 ml, 1 3M™ Impregum™ Penta™ Super Quick MB cartridge for 3M™ Pentamix™ 3 Automatic Mixing Unit, 10 3M™ Penta™ Mixing Tips – Red, 1 3M™ Penta™ Elastomer syringe, 1 3M™ Polyether Tray Adhesive, 3 3M™ Astringent Retraction Paste capsules, 6 3M™ Dual Arch Trays (3 posterior, 3 anterior trays)

**69415 3M™ Impregum™ Penta™ Super Quick Medium Body/Light Body**  
1 base paste MB – 300 ml, 1 catalyst MB – 60 ml, 2 3M™ Impregum™ Super Quick LB cartridges – 50 ml each, 1 3M™ Impregum™ Penta™ Super Quick MB cartridge for 3M™ Pentamix™ 3 Automatic Mixing Unit, 10 3M™ Penta™ Mixing Tips – Red, 5 3M™ Intra-oral Syringes – Purple, 5 3M™ Garant™ Mixing Tips – Purple, 5 3M™ Garant™ Intra-oral Tips – White, 1 3M™ Polyether Tray Adhesive, 3 3M™ Astringent Retraction Paste capsules, 6 3M™ Dual Arch Trays (3 posterior, 3 anterior trays)

**69407 3M™ Impregum™ Penta™ Super Quick Heavy Body/Light Body**  
1 base paste HB – 300 ml, 1 catalyst HB – 60 ml, 2 3M™ Impregum™ Super Quick LB cartridges – 50 ml each, 1 3M™ Impregum™ Penta™ Super Quick HB cartridge for 3M™ Pentamix™ 3 Automatic Mixing Unit, 10 3M™ Penta™ Mixing Tips – Red, 5 3M™ Intra-oral Syringes – Purple, 5 3M™ Garant™ Mixing Tips – Purple, 5 3M™ Garant™ Intra-oral Tips – White, 1 3M™ Polyether Tray Adhesive, 3 3M™ Astringent Retraction Paste capsules, 6 3M™ Dual Arch Trays (3 posterior, 3 anterior trays)

#### Refills

**69385 3M™ Impregum™ Penta™ Super Quick Medium Body**  
2 base pastes – 300 ml each, 2 catalysts – 60 ml each

**69417 3M™ Impregum™ Super Quick Medium Body**  
4 cartridges – 50 ml each, 10 3M™ Garant™ Mixing Tips purple

**69406 3M™ Impregum™ Penta™ Super Quick Heavy Body**  
2 base pastes – 300 ml each, 2 catalysts – 60 ml each

**69379 3M™ Impregum™ Super Quick Light Body**  
4 cartridges – 50 ml each, 5 3M™ Intra-oral Syringes Purple, 5 3M™ Garant™ Mixing Tips purple

**69418 3M™ Impregum™ Super Quick Ultra Light Body**  
4 cartridges – 50 ml each, 5 3M™ Intra-oral Syringes Purple, 5 3M™ Garant™ Mixing Tips purple

#### Accessories

**69387 3M™ Impregum™ Penta™ Super Quick Medium Body Cartridge**

**69401 3M™ Impregum™ Penta™ Super Quick Heavy Body Cartridge**

**71622 48 3M™ Dual Arch Impression Trays Posterior**

**71623 35 3M™ Dual Arch Impression Trays Anterior**

**30600 3M™ Polyether Tray Adhesive**  
1 bottle – 17 ml



3M Oral Care  
2510 Conway Avenue  
St. Paul, MN 55144-1000  
USA  
1-800-634-2249

3M Canada  
Post Office Box 5757  
London, Ontario N6A 4T1  
Canada  
1-888-363-3685

[3M.com/Dental](http://3M.com/Dental)

3M, Garant, Impregum, Penta and Pentamix are trademarks of 3M Company or 3M Deutschland GmbH. Used under license in Canada. All other trademarks are owned by other companies.

© 3M 2020. All rights reserved. 70-2013-7071-8