


# 3M™ Tegaderm™ +Pad Film Dressing with Non-Adherent Pad 3M™ Medipore™ +Pad Soft Cloth Adhesive Wound Dressing



- All-in-one sterile dressing convenience
- Highly conformable and comfortable
- Unique non-adherent pad



# Trusted Performance Delivered Two Styles

## 3M™ Tegaderm™ +Pad Film Dressing with Non-Adherent Pad

### *A waterproof, absorbent dressing.*

- All-in-one sterile dressing—reduces application time
- A bacterial and viral\* barrier—helps prevent external bacteria and contaminants from entering the site
- A waterproof barrier—patient may bathe or shower without removing dressing
- Conforms to body contours, providing comfort
- Unique, non-adherent absorbent pad

\**In vitro* testing shows that the transparent film dressing provides a viral barrier from viruses 27 nm in diameter or larger while the dressing remains intact without leakage.



Tegaderm™ +Pad on epidural



Tegaderm™ +Pad on hip incision



Tegaderm™ +Pad on arthroscopic incision



Tegaderm™ +Pad

*Full range of sizes to fit every application need.*

### Suggested Applications

#### Primary Dressing Applications

- Post-surgical – especially on areas where swelling or movement is expected
- Acute wounds such as large abrasions, lacerations, and burns
- Light-to-moderately-draining wounds
- Superficial and partial thickness wounds
- I.V. catheter sites
- Fragile skin

in

# of Absorbent Island Dressings

## 3M™ Medipore™ +Pad Soft Cloth Adhesive Wound Dressing

*A gentle, conformable, absorbent dressing.*

- All-in-one sterile dressing—reduces application time
- Conforms to body contours, providing comfort
- Gentle, porous and breathable
- Unique, non-adherent absorbent pad



Medipore™ +Pad on IV catheter site



Medipore™ +Pad on knee



Medipore™ +Pad on appendectomy

### Secondary Dressing Applications

- Over hydrogels, alginates, gauze, or other primary dressings



Medipore™ +Pad  
Full range of sizes to fit every application need.

Catalog No. NHRIC/NDC No.	Dressing Size	Pad Size	Dressings/ Box	Boxes/ Case	HCPCS* Code
<b>3M™ Tegaderm™ +Pad Film Dressing with Non-Adherent Pad</b>					
<b>3582</b> <b>8333-3582-01</b>	2 in. x 2-3/4 in. 5 cm x 7 cm	1 in. x 1-1/2 in. 2,5 cm x 4 cm	50	4	A6203
<b>3584</b> <b>8333-3584-01</b>	2-3/8 in. x 4 in. 6 cm x 10 cm	1 in. x 2-3/8 in. 2,5 cm x 6 cm	50	4	A6203
<b>3586</b> <b>8333-3586-01</b>	3-1/2 in. x 4 in. 9 cm x 10 cm	1-3/4 in. x 2-3/8 in. 4,5 cm x 6 cm	25	4	A6203
<b>3587 Oval</b> <b>8333-3587-01</b>	3-1/2 in. x 4-1/8 in. 9 cm x 10,5 cm	1-3/4 in. x 2-3/8 in. 4,5 cm x 6 cm	25	4	A6203
<b>3588</b> <b>8333-3588-01</b>	6 in. x 6 in. 15 cm x 15 cm	4 in. x 4 in. 10 cm x 10 cm	25	4	A6203
<b>3589</b> <b>8333-3589-01</b>	3-1/2 in. x 6 in. 9 cm x 15 cm	1-3/4 in. x 4 in. 4,5 cm x 10 cm	25	4	A6203
<b>3590</b> <b>8333-3590-01</b>	3-1/2 in. x 8 in. 9 cm x 20 cm	1-3/4 in. x 6 in. 4,5 cm x 15 cm	25	4	A6203
<b>3591</b> <b>8333-3591-01</b>	3-1/2 in. x 10 in. 9 cm x 25 cm	1-3/4 in. x 8 in. 4,5 cm x 20 cm	25	4	A6203
<b>3593</b> <b>8333-3593-01</b>	3-1/2 in. x 13-3/4 in. 9 cm x 35 cm	1-3/4 in. x 11-3/4 in. 4,5 cm x 30 cm	25	4	A6204

Catalog No. NHRIC/NDC No.	Dressing Size	Pad Size	Dressings/ Box	Boxes/ Case	HCPCS* Code
<b>3M™ Medipore™ +Pad Soft Cloth Adhesive Wound Dressing</b>					
<b>3562</b> <b>8333-3562-01</b>	2 in. x 2-3/4 in. 5 cm x 7 cm	1 in. x 1-1/2 in. 2,5 cm x 4 cm	50	4	A6219
<b>3564</b> <b>8333-3564-01</b>	2-3/8 in. x 4 in. 6 cm x 10 cm	1 in. x 2-3/8 in. 2,5 cm x 6 cm	50	4	A6219
<b>3566</b> <b>8333-3566-01</b>	3-1/2 in. x 4 in. 6 cm x 10 cm	1-3/4 in. x 2-3/8 in. 4,5 cm x 6 cm	25	4	A6219
<b>3568</b> <b>8333-3568-01</b>	6 in. x 6 in. 15 cm x 15 cm	4 in. x 4 in. 10 cm x 10 cm	25	4	A6219
<b>3569</b> <b>8333-3569-01</b>	3-1/2 in. x 6 in. 9 cm x 15 cm	1-3/4 in. x 4 in. 4,5 cm x 10 cm	25	4	A6219
<b>3570</b> <b>8333-3570-01</b>	3-1/2 in. x 8 in. 9 cm x 20 cm	1-3/4 in. x 6 in. 4,5 cm x 15 cm	25	4	A6219
<b>3571</b> <b>8333-3571-01</b>	3-1/2 in. x 10 in. 9 cm x 25 cm	1-3/4 in. x 8 in. 4,5 cm x 20 cm	25	4	A6219
<b>3573</b> <b>8333-3573-01</b>	3-1/2 in. x 13-3/4 in. 9 cm x 35 cm	1-3/4 in. x 11-3/4 in. 4,5 cm x 30 cm	25	4	A6220

\* HCPCS codes as referenced in the Local Coverage Determination Surgical Dressing Policy.

**Disclaimer:** HCPCS codes have been provided to assist you in the preparation of insurance claims. Please note, however, that the reimbursement information provided by 3M Health Care and its representatives is intended to provide general information relevant to coverage and coding for 3M products. Insurers' reimbursement policies can vary and the use of the codes discussed here does not guarantee that an insurer will cover or pay at any particular level. Health care providers should exercise independent clinical judgment in choosing the codes which most accurately describe the products provided.



### 3M Critical & Chronic Care Solutions Division

**3M Health Care**  
2510 Conway Avenue  
St. Paul, MN 55144  
USA  
1-800-228-3957  
[www.3M.com/healthcare](http://www.3M.com/healthcare)

**3M Canada**  
Post Office Box 5757  
London, Ontario N6A 4T1  
Canada  
1-800-364-3577  
[www.3M.com/healthcare](http://www.3M.com/healthcare)

3M, Medipore and Tegaderm  
are trademarks of 3M.  
Please recycle.  
Printed in U.S.A.  
© 3M 2008, 2012. All rights reserved.  
70-2008-8925-4 (Bolger)