

3M Science.
Applied to Life.™

C L A R I T Y™ | **ADVANCED**
advanced ceramic brackets

**Advanced
ceramics.
Enhanced
beauty.**

Parts List



Parts List

3M™ Clarity™ ADVANCED Ceramic Brackets

MBT™ Versatile+ Appliance System

MBT™
Versatile+ Appliance System

Maxillary			IN/OUT		0.018 in.		0.022 in.	
Tooth	Torque	Angulation	in.	mm	L	R	L	R
Central	+17	+4	0.026	0.66	006-201	006-202	006-301	006-302
Lateral	+10	+8	0.035	0.89	006-205	006-206	006-305	006-306
Cuspid	0	+8	0.021	0.53	006-207	006-208	006-307	006-308
Cuspid HK	0	+8	0.021	0.53	006-209	006-210	006-309	006-310
Bicuspid Universal	-7	0	0.029	0.74	006-211		006-311	
Bicuspid HK	-7	0	0.029	0.74	006-213	006-212	006-313	006-312
Mandibular			IN/OUT		0.018 in.		0.022 in.	
Tooth	Torque	Angulation	in.	mm	L	R	L	R
Anterior	-6	0	0.045	1.14	006-250		006-350	
Cuspid	0	+3	0.020	0.51	006-251	006-252	006-351	006-352
Cuspid HK	0	+3	0.020	0.51	006-253	006-254	006-353	006-354
1st Bicuspid	-12	+2	0.030	0.76	006-255	006-256	006-355	006-356
1st Bicuspid HK	-12	+2	0.030	0.76	006-257	006-258	006-357	006-358
2nd Bicuspid	-17	+2	0.034	0.86	006-259	006-260	006-359	006-360
2nd Bicuspid HK	-17	+2	0.034	0.86	006-261	006-262	006-361	006-362

Adding Prefix "4" is APC™ Flash-Free Adhesive System – Example: 4006-201

Adding Prefix "3" is APC™ II Adhesive System – Example: 3006-201

Adding Prefix "5" is APC™ PLUS Adhesive System – Example: 5006-201

Single Patient Kits

	0.018 in.	0.022 in.
U/L 5×5 Clarity™ ADVANCED Brackets Cuspid HK	006-100	006-110
U/L 5×5 Clarity ADVANCED Brackets HKs	006-101	006-111
U/L 5×5 Clarity ADVANCED Brackets	006-102	006-112
U/L 3×3 Clarity ADVANCED Brackets Cuspid HK	006-103	006-113
U/L 3×3 Clarity ADVANCED Brackets	006-104	006-114
U 5×5 Clarity ADVANCED Brackets Cuspid HK	006-105	006-115
U 5×5 Clarity ADVANCED Brackets HKs	006-106	006-116
U/L 5×5 Clarity ADVANCED Brackets	006-107	006-117
U 3×3 Clarity ADVANCED Brackets Cuspid HK	006-108	006-118
U/L 3×3 Clarity ADVANCED Brackets	006-109	006-119
Adding Prefix "4" is APC™ Flash-Free Adhesive System – Example: 4006-100		
Adding Prefix "3" is APC™ II Adhesive System – Example: 3006-100		
Adding Prefix "5" is APC™ PLUS Adhesive System – Example: 5006-100		

3M™ Clarity™ ADVANCED Ceramic Brackets


CLARITY™ | ADVANCED
advanced ceramic brackets

Roth* Prescription


Maxillary			IN/OUT		0.018 in.		0.022 in.	
Tooth	Torque	Angulation	in.	mm	L	R	L	R
Central	+12	+5	0.026	0.66	006-215	006-216	006-315	006-316
Lateral	+8	+9	0.035	0.889	006-217	006-218	006-317	006-318
Cuspid HK	0	+11	0.021	0.53	006-221	006-222	006-321	006-322
Bicuspid HK	-7	0	0.029	0.74	006-213	006-212	006-313	006-312
Mandibular			IN/OUT		0.018 in.		0.022 in.	
Tooth	Torque	Angulation	in.	mm	L	R	L	R
Anterior	-1	0	0.045	1.14	006-264		006-364	
Cuspid HK	-11	7	0.020	0.51	006-267	006-268	006-367	006-368
1st Bicuspid HK	-17	0	0.030	0.76	006-271	006-272	006-371	006-372
2nd Bicuspid HK	-22	0	0.034	0.87	006-275	006-276	006-375	006-376
					Adding Prefix "4" is APC™ Flash-Free Adhesive System – Example: 4006-215			
					Adding Prefix "3" is APC™ II Adhesive System – Example: 3006-215			
					Adding Prefix "5" is APC™ PLUS Adhesive System – Example: 5006-215			

Single Patient Kits

	0.018 in.	0.022 in.
U/L 5x5 Clarity™ ADVANCED Brackets HK	006-121	006-131
U/L 3x3 Clarity ADVANCED Brackets Cuspid HK	006-123	006-133
U 5x5 Clarity ADVANCED Brackets Cuspid HK	006-125	006-135
U 5x5 Clarity ADVANCED Brackets HK	006-126	006-136
U 3x3 Clarity ADVANCED Brackets Cuspid HK	006-128	006-138
Adding Prefix "4" is APC™ Flash-Free Adhesive System – Example: 4006-121		
Adding Prefix "3" is APC™ II Adhesive System – Example: 3006-121		
Adding Prefix "5" is APC™ PLUS Adhesive System – Example: 5006-121		



Complement the aesthetic Clarity™ ADVANCED Ceramic Brackets with tooth colored ligature ties.



Unitek™ Ligature Ties –
Preformed Coated
Part No. 390-170
.012 diameter
100 per pack

*3M version of this prescription. No endorsement by the Doctor is implied.



Parts List *continued*

3M™ Clarity™ ADVANCED Ceramic Brackets

CLARITY™ | ADVANCED
advanced ceramic brackets

High Torque Brackets

Maxillary			IN/OUT		0.018 in.		0.022 in.	
Tooth	Torque	Angulation	in.	mm	L	R	L	R
Central	+22	+5	0.026	0.66	006-223	006-224	006-323	006-324
Lateral	+14	+8	0.035	0.89	006-225	006-226	006-325	006-326
Cuspid HK	+7	+5	0.021	0.53	006-229	006-230	006-329	006-330
Bicuspid Universal	0	0	0.029	0.74	006-281*		006-381*	
Mandibular			IN/OUT		0.018 in.		0.022 in.	
Tooth	Torque	Angulation	in.	mm	L	R	L	R
Cuspid HK	+7	+5	0.020	0.51	006-279	006-280	006-379	006-380
Adding Prefix "4" is APC™ Flash-Free Adhesive System – Example: 4006-201								
Adding Prefix "3" is APC™ II Adhesive System – Example: 3006-201								
Adding Prefix "5" is APC™ PLUS Adhesive System – Example: 5006-201								

High Torque Kits

	0.018 in.	0.022 in.
U/L 3x3 Clarity ADVANCED Hi TQ Brackets Cuspid HK	006-143	006-153
U 5x5 Clarity ADVANCED Hi TQ Brackets Cuspid HK	006-170*	006-180*
U 3x3 Clarity ADVANCED Hi TQ Brackets Cuspid HK	006-148	006-158
Adding Prefix "4" is APC™ Flash-Free Adhesive System – Example: 4006-100		
Adding Prefix "3" is APC™ II Adhesive System – Example: 3006-100		
Adding Prefix "5" is APC™ PLUS Adhesive System – Example: 5006-100		

*Not available with APC™ PLUS Adhesive.



3M Oral Care
2510 Conway Avenue
St. Paul, MN 55144-1000 USA
Phone 1-800-423-4588
Web 3M.com/ortho

3M Canada
Health Care Division
300 Tartan Dr.
London, ON N5V 4M9 Canada
Phone 1-800-443-1661

3M, APC, Clarity, MBT, and Unitek are trademarks of 3M.
Used under license in Canada.
© 3M 2016. All rights reserved. 70-2021-5261-0 1612