

**3M** Science.  
Applied to Life.™

**Esthetic Restoration Solution**

**A beautiful smile  
she can't wait to share.**



# Because you want every patient to leave with a beautiful smile.

Patients love the way highly polished restorations look and feel. But creating a high-end, natural-looking gloss can take more of your valuable time. So what do you do ... sacrifice higher esthetics for higher productivity? Or spend the extra time to create a beautiful smile ... because, "that's just the way it is?"

The days of time-consuming esthetic restorations are gone. With our Esthetic Restoration Solution, you'll love how simple it is to create natural-looking restorations with a gorgeous paste-like gloss.

Oh, don't be surprised if word of your great work gets around ... because a beautiful smile is worth sharing.



The truth is, the composite you use matters. You can spend a lot of time polishing most composites, only to find the initial gloss doesn't last long.

Or, you can create smiles your patients will love by using our Esthetic Restoration Solution. Thanks to 3M's nanotechnology, Filtek™ Supreme Ultra Universal Restorative is easy to polish and also offers unsurpassed polish retention.

Plus, nanotechnology makes Filtek Supreme Ultra restorative beautifully strong, giving you the option of using it in any place and for any Class.

Having one restorative that can do it all is a beautiful thing.



Class I



Class II



Class III



Class IV



Class V

## Nanotechnology. For beautifully strong restorations ... anywhere.

# Diamond particles. For a gorgeous, paste-like gloss.

The esthetic appearance of a restoration depends on the shade, shape and the gloss achieved by final polishing.

With the 3M™ ESPE™ Sof-Lex™ Diamond Polishing System, you can get a gorgeous, diamond paste-like gloss—with the convenience of a rubberized system!

The Sof-Lex™ Pre-Polishing Spiral (beige) smooths and removes scratches in restorations that develop during contouring and prepares the surface for high-gloss polishing. The Sof-Lex™ Diamond Polishing Spiral (pink) delivers the final high-gloss polish.

You'll find it reassuring to know that you no longer need to stop and switch between traditional points, cups, discs and brushes. The flexible spirals adapt to all tooth surfaces—anterior and posterior, convex and concave. They can be sterilized and reused, making them an economical choice for your practice.

## A difference you can see!



Filtek™ Supreme Ultra  
Universal Restorative polished  
with Sof-Lex™ Diamond  
Polishing System

TPH Spectra® Universal  
Composite polished with  
Enhance® Finishing System  
and PoGo® Polishing System

Notice a clearer reflection when  
Filtek™ Supreme Ultra Universal Restorative is  
polished with the Sof-Lex™ Diamond Polishing System.

# Now simple esthetics can be as easy as 1, 2, 3 ...

A factor in the success of this solution is 3M's unprecedented history in abrasive technologies. We leveraged that expertise to create a polishing solution that can help you maximize the esthetic benefits of Filtek™ Supreme Ultra Universal Restorative with ease and simplicity.

Filtek™ Supreme Ultra  
Universal Restorative

1



Sof-Lex™  
Pre-Polishing Spiral

2



Sof-Lex™ Diamond  
Polishing Spiral

3




# Ordering Information



## Simple Esthetics

While Filtek™ Supreme Ultra Universal Restorative is known for its wide range of shades and opacities, most single shade restorations are made with one of 8 Body shades.



### Filtek™ Supreme Ultra Universal Restorative

#### Item # Product Information

<b>6029SU</b>	<b>Filtek™ Supreme Ultra Universal Restorative Single Opacity Trial Kit—Capsule</b> Kit includes: 20 – 0.2g capsules: 10 each of Body shades A1 and A2; 1 – 1.5ml Scotchbond™ Universal Adhesive Vial; 1 – 3ml Scotchbond™ Universal Etchant Syringe; 5707SD Restorative Dispenser; Technique Guides; Instructions for Use
<b>6028SU</b>	<b>Filtek™ Supreme Ultra Universal Restorative Single Opacity Trial Kit—Syringe</b> Kit includes: 1 – 4g syringe of Body shades A1 and A2; 1 – 1.5ml Scotchbond™ Universal Adhesive Vial; 1 – 3ml Scotchbond™ Universal Etchant Syringe; Technique Guides; Instructions for Use
<b>6029DU</b>	<b>Filtek™ Supreme Ultra Universal Restorative Dual Opacity Trial Kit—Capsule</b> Kit includes: 40 – 0.2g capsules: 10 each of Dentin shades A3 and A4 and Enamel shades A2 and A3; 1 – 1.5ml Scotchbond™ Universal Adhesive Vial; 1 – 3ml Scotchbond™ Universal Etchant Syringe; 5707SD Restorative Dispenser; Technique Guides; Instructions for Use
<b>6028DU</b>	<b>Filtek™ Supreme Ultra Universal Restorative Dual Opacity Trial Kit—Syringe</b> Kit includes: 1 – 4g syringe of Dentin shades A3 and A4 and Enamel shades A2 and A3; 1 – 1.5ml Scotchbond™ Universal Adhesive Vial; 1 – 3ml Scotchbond™ Universal Etchant Syringe; Technique Guides; Instructions for Use
<b>6029M</b>	<b>Filtek™ Supreme Ultra Universal Restorative Master's Kit—Capsule</b> Kit includes: 132 – 0.2g capsules: 3 each of A1D, A2D, A3D, A4D, B3D, C4D, WD; 6 each of A1B, A2B, A3B, A3.5B, A4B, B1B, B2B, C2B; 3 each of B3B, B5B, C1B, C3B, D2B, D3B, WB, XWB; 6 of A2E; 3 each of A1E, A3E, B1E, B2E, D2E, WE, XWE; 3 each of Translucent shades AT, BT, CT, GT. All in plastic reusable case with 1 Shade Selection Wheel; Technique Guides; Instructions for Use

#### DENTIN SHADE REFILLS

Refill SHADE	Item # SYRINGE	Item # CAPSULE
A1D	6028A1D	6029A1D
A2D	6028A2D	6029A2D
A3D	6028A3D	6029A3D
A4D	6028A4D	6029A4D
B3D	6028B3D	6029B3D
C4D	6028C4D	6029C4D
White	6028WD	6029WD

#### BODY SHADE REFILLS

Refill SHADE	Item # SYRINGE	Item # CAPSULE
A1B	6028A1B	6029A1B
A2B	6028A2B	6029A2B
A3B	6028A3B	6029A3B
A3.5B	6028A3.5B	6029A3.5B
A4B	6028A4B	6029A4B
A6B	6028A6B	6029A6B
B1B	6028B1B	6029B1B
B2B	6028B2B	6029B2B
B3B	6028B3B	6029B3B
B5B	6028B5B	6029B5B
C1B	6028C1B	6029C1B
C2B	6028C2B	6029C2B
C3B	6028C3B	6029C3B
D2B	6028D2B	6029D2B
D3B	6028D3B	6029D3B
White	6028WB	6029WB
Extra White	6028XWB	6029XWB

#### TRANSLUCENT SHADE REFILLS

Refill SHADE	Item # SYRINGE	Item # CAPSULE
Clear	6028CT	6029CT
Blue	6028BT	6029BT
Gray	6028GT	6029GT
Amber	6028AT	6029AT

SYRINGE SHADE REFILLS: Include 1 – 4g syringe; Instructions for Use

CAPSULE SHADE REFILLS: Dentin, Enamel and Translucent Shade Refills include 10 – 0.2g capsules; Body Shade Refills include 20 – 0.2g capsules; Instructions for Use

#### ENAMEL SHADE REFILLS

Refill SHADE	Item # SYRINGE	Item # CAPSULE
A1E	6028A1E	6029A1E
A2E	6028A2E	6029A2E
A3E	6028A3E	6029A3E
B1E	6028B1E	6029B1E
B2E	6028B2E	6029B2E
D2E	6028D2E	6029D2E
White	6028WE	6029WE
Extra White	6028XWE	6029XWE

### Sof-Lex™ Diamond Polishing System

#### Item # Product Information

<b>5092-IM</b>	<b>Sof-Lex™ Diamond Polishing System—Introductory Kit</b> Kit includes: 5 Sof-Lex™ Pre-Polishing Spirals (Beige); 5 Sof-Lex™ Diamond Polishing Spirals (Pink); 1 Sof-Lex™ RA mandrel; Instructions for Use; Technique Card
<b>5090</b>	<b>Sof-Lex™ Pre-Polishing Spirals—Refill</b> Includes: 15 Sof-Lex Pre-Polishing Spirals (Beige)
<b>5091</b>	<b>Sof-Lex™ Diamond Polishing Spirals—Refill</b> Includes: 15 Sof-Lex Diamond Polishing Spirals (Pink)



[www.3M.com/EstheticRestorations](http://www.3M.com/EstheticRestorations)



3M Oral Care  
2510 Conway Avenue  
St. Paul, MN 55144-1000 USA

Phone 1-800-634-2249  
Web 3M.com/dental

3M Canada  
Post Office Box 5757  
London, Ontario N6A 4T1  
Canada

Phone 1-888-363-3685

3M, ESPE, Filtek and Sof-Lex are trademarks of 3M or 3M Deutschland GmbH. Used under license in Canada. Printed in USA. © 3M 2016. All rights reserved. 70-2013-0720-7