



# This is for the passionate, who embrace a busy day.



#### 1 Shade. 1 Universal opacity. It doesn't get any simpler.

Efficiency is key to navigating a tight schedule. That's one reason dentists use a single-shade restoration for about 80% of their cases.\*

Knowing the positive impact it would have on your practice if those restorations were easier, we created 3M™ Filtek™ Universal Restorative.

It starts with simplifying shade selection. Just 8 designer shades plus an Extra White (XW) cover the 19 VITA classical and bleach shades.\* Plus, the shades have a universal opacity that meets most clinical needs.

- NaturalMatch technology (made up of nanofillers, proprietary lowstress monomers and pigments) enables one-shade esthetics to blend more naturally with dentition.
- Fluorescent pigments enhance the lifelike appearance.
- 3M filler technology allows for increased radiopacity.

Now you can keep dark areas from diminishing your esthetic results.

3M's first Pink Opaquer effectively masks metal and stained dentition.

\*3M Internal Data

## It's for patients who take pride in a beautiful smile ...



#### **Extra White** restorations for patients who bleach

The tooth-whitening business continues to expand,1 but dentists tell us it's difficult to find composites that can match the whitest teeth. We answered with an XW, to help you meet patient expectations. XW

#### Shade blending that passes the "friend test"

With the unique shade-matching effect of 3M NaturalMatch technology, even anterior work can be virtually unnoticeableeven by close friends. The material's universal opacity helps it blend more invisibly with surrounding dentition.

1. Expected size of the global teeth whitening products market by 2024. Hexa Research. (2017, August 17). Teeth Whitening Products Market Size Worth USD 7.40 Billion By 2024 [Press release]. Retrieved from www.hexaresearch.com/press-release/global-teeth-whitening-products-market.

#### A glamorous gloss that lasts

The initial polish of any restorative can be impressive, but many dull over time. Thanks to 3M's patented TRUE nanofiller technology, Filtek Universal restorations polish quickly to a high gloss, and polish can even be improved with toothbrush abrasion.\*

100K magnification of nanocluster. SEM photo courtesy of Dr. J. Perdigão, University of Minnesota.

\*3M Internal Data

# and for the bold who say, "Yes!" to life. **Dual "endurance"** technologies

Chewing, brushing, grinding, gritting ... restorations get quite a workout. We raised the bar on durability by combining two of 3M's top technologies. Low-stress monomer technologies pioneered in our bulkfill composites, help protect against shrinkage and shrinkage stress. 3M's TRUE nanotechnology provides high strength and excellent wear resistance.



#### Even the capsule is designed to make restorations easier.

Our capsule has a shear thinning design which allows the material to extrude like a flowable—then become more viscous for shaping. The tip is longer (8mm), and the opening flexes during placement enabling good adaptation in hard to reach areas.



### This is 3M™ Filtek™ **Universal Restorative**

It's for all of you who want your restorations, your composite, your day,





#### **Key features**

- Simplified shade selection—just 8 designer shades to cover most of your patients
- Universal opacity to blend more naturally
- Extra White (XW) for matching bleached teeth
- Fluorescent pigments to give restorations a lifelike appearance
- Pink Opaquer to easily mask dark areas
- 3M's nanotechnology for easy polish that lasts
- Proprietary low-stress monomers for fracture resistance
- Increased radiopacity for easier identification at recall
- Capsule warming—Filtek Universal Restorative capsules can be safely warmed.\*

\*See Instructions for Use

#### **Indications**

- Direct anterior and posterior restorations (including occlusal surfaces)
- Core buildups
- Splinting
- Indirect restorations including inlays, onlays and veneers



#### **Ordering Information**

#### 3M™ Filtek™ Universal Restorative

Item # P	Product Information  3M™ Filtek™ Universal Restorative Essential Kit, Capsules  Kit includes: 80 – 0.2g capsules: 20 each of A1, A2, A3, A3.5; 10 – 0.2g capsules of PO; Technique Guide; Instructions for Use			
	<b>M™ Filtek™ Universal Re</b> it includes: 5 – 4g syring		Kit, Syringes , A3, A3.5 and PO; Technique Guide; Instructions for Use.	
Refill	Item #	Item #		
SHADES	SYRINGE	CAPSULE	\ \ \ \ \ \ \ \ \ \ \ \ \ \ \ \ \ \ \ \	
— A1	6555A1	6550A1 ——		
A2	6555A2	6550A2		
A3	6555A3	6550A3		
A3.5	6555A3.5	6550A3.5	designer shades -	
A4	6555A4	6550A4		
B1	6555B1	6550B1		
B2	6555B2	6550B2		
— D3	6555D3	6550D3 ——		
XW—Extra White	Bleach 6555XW	6550XW	/ "; \	
PO—Pink Opaque	er 6555PO	6550PO		
SHADE REFILL	CAPSULES	—includes 20 – 0.2	2g capsules SYRINGES—includes 1 – 4g syringe	
PINK OPAQUER (PO)	REFILL CAPSULES	—includes 10 – 0.2	g capsules SYRINGES—includes 1 – 4g syringe	

All refills include a Technique Guide and Instructions for Use.

#### www.3M.com/FiltekUniversal



3M Oral Care 2510 Conway Avenue St. Paul, MN 55144-1000 USA

Phone 1-800-634-2249 Web 3M.com/dental 3M Canada Post Office Box 5757 London, Ontario N6A 4T1 Canada

Phone 1-888-363-3685

3M and Filtek are trademarks of 3M or 3M Deutschland GmbH. Used under license in Canada. Printed in USA. © 3M 2019. All rights reserved. 70-2013-7114-6