

8 Designer shades. 1 Universal opacity. It doesn't get any simpler.

Efficiency is key to navigating a tight schedule. That's one reason dentists use a single-shade restoration for about 80% of their cases.* Making those restorations easier is why we created 3M™ Filtek™ Universal Restorative.

It starts with simplifying shade selection.

Just 8 designer shades plus an Extra White (XW) cover the 19 VITA classical and bleach shades.* Plus, the shades have a universal opacity that meets most clinical needs.


- NaturalMatch technology (made up of nanofillers, proprietary low-stress monomers and pigments) enables one-shade esthetics to blend more naturally with dentition.
- Fluorescent pigments enhance the lifelike appearance.
- 3M filler technology allows for increased radiopacity.

Enhancements like an Extra White (XW) and a Pink Opaquer (PO) help you deliver excellent esthetics in unique situations.

Create Extra White restorations to meet the expectations of patients who bleach.



Keep dark areas from diminishing esthetic results. A Pink Opaquer effectively masks metal and stained dentition.





A new class of universal composite

From the moment of cure, restorations are challenged by chewing, brushing and grinding. We raised the bar for durability by combining two proven proprietary dental technologies for the first time in a universal composite.

Low-stress monomer technologies

Pioneered in our bulk-fill composites, they help protect restorations from problems caused by shrinkage and shrinkage stress.



3M's patented TRUE nanotechnology

Gives restorations excellent polish that lasts, as well as high strength and excellent wear resistance.

100K magnification of nanocluster. SEM photo courtesy of Dr. J. Perdigão, University of Minnesota.

Key features

- **Simplified shade selection**—just 8 designer shades to cover most of your patients
- Universal opacity to blend more naturally
- Extra White (XW) for matching bleached teeth
- Fluorescent pigments to give restorations a lifelike appearance
- Pink Opaquer to easily mask dark areas
- 3M's nanotechnology for easy polish that lasts
- Proprietary low-stress monomers for fracture resistance
- Increased radiopacity for easier identification at recall
- Capsule warming—3M[™] Filtek[™] Universal Restorative capsules can be safely warmed.*

*See Instructions for Use

Indications

- Direct anterior and posterior restorations (including occlusal surfaces)
- Core buildups
- Splinting
- Indirect restorations including inlays, onlays and veneers



Ordering Information

3M™ Filtek™ Universal Restorative

Item #	Product Information				
	3M™ Filtek™ Universal Restorative Essential Kit, Capsules Kit includes: 80 – 0.2g capsules: 20 each of A1, A2, A3, A3.5; 10 – 0.2g capsules of PO; Technique Guide; Instructions for Use.				
		Muniversal Restorative Essential Kit, Syringe 5 – 4g syringes: 1 each of A1, A2, A3, A3.5 and PO; Technique Guide; Instructions for Use.			
Refill		Item #	Item #		
SHADES		SYRINGE	CAPSULE		
—A1		6555A1	6550A1	designer shades -	
A2		6555A2	6550A2		
A3		6555A3	6550A3		
A3.5		6555A3.5	6550A3.5		
A4		6555A4	6550A4		
B1		6555B1	6550B1		
B2		6555B2	6550B2		
— D3		6555D3	6550D3 ——	─	
XW—Extra Whi	te	6555XW	6550XW	, , , , ,	
PO—Pink Opaq	uer	6555PO	6550PO		
SHADE REFILLS		CAPSULES-	-includes 20 – 0.2	g capsules SYRINGES—includes 1 – 4g syringe	
PINK OPAQUER (PO) REFILLS		CAPSULES-	-includes 10 - 0.2g	g capsules SYRINGES—includes 1 – 4g syringe	

All refills include a Technique Guide and Instructions for Use.

www.3M.com/FiltekUniversal



3M Oral Care 2510 Conway Avenue St. Paul, MN 55144-1000 USA

Phone 1-800-634-2249 Web 3M.com/dental 3M Canada Post Office Box 5757 London, Ontario N6A 4T1 Canada

Phone 1-888-363-3685

3M and Filtek are trademarks of 3M or 3M Deutschland GmbH. Used under license in Canada. Printed in USA. © 3M 2019. All rights reserved. 70-2013-7115-3