

Impregum™ Soft

Polyether Impression Material



Captures more detail than any VPS
no matter how *ultra* or *extra* it claims to be.

Summary of advantages

- Better initial hydrophilicity than any VPS – provides void-free impressions in moist conditions
- “Snap-set behavior” of *soft* polyether is more forgiving
- Excellent flow properties – captures fine details for precise impressions
- Less-rigid, *soft* polyether is easier to remove – great for dual-arch trays
- Quick-Step material cuts working and setting times by at least 33% over Impregum Soft Polyether Impression Material

“3M ESPE has successfully addressed my major concerns – rigidity, taste and set time – with the new Impregum Soft Quick Step Impression Material”

– Dr. Brian Kinsey, Atlanta, GA

Better than any VPS for capturing detail in a wet environment.*

Unlike VPS material, Impregum™ Soft Polyether Impression Material displaces moisture right from the start of the mix – when it counts – for void-free impressions.



Water drop on Aquasil Impression Material prior to set (left). Initially, Aquasil shows hydrophobic behavior. Water drop on Impregum™ Soft Polyether Impression Material prior to set (right).

Source: Mondon M., Ziegler Chr., Reusch B., “Surface Wettability During the Curing Process of Impression Materials,” AADR, #618, 2001.

Better than any VPS for achieving precise-fitting restorations, due to polyether’s unique snap-set behavior.*

Impregum Soft polyether material will not start setting before the working time ends, and when it does set, it does so immediately. The result is precise-fitting restorations without distortion.

Better than any VPS at capturing fine detail for better-fitting crowns, reducing crown adjustments.*

Polyether material flows better into the toughest spots, yet does not slump – capturing more detail than any VPS. In fact, it offers the best flow behavior during the entire working time.

The Shark Fin Test proves it ...
the taller the fin, the better the flow!



Results shown are at the end of the working time recommended by the manufacturer.*

*Data on file. Go to www.3MESPE.com/GreatImpressions.

3M ESPE

Now you can get the accuracy of polyether with the convenience of handgun dispensing.

Impregum™ Soft Polyether Tray Impression Material is now available in two delivery options:

- Handgun dispenser: Fast, simple and convenient
- Pentamix™ 3 Automatic Mixing Unit: Get homogeneous, void-free mix at the push of a button

It's easier than ever to try Impregum Soft polyether material. See for yourself how using polyether impression material can reduce remakes, retakes and adjustment times, making your practice more efficient and more profitable.



More proven solutions brought to you by 3M ESPE, the worldwide leader in impression materials.

Ordering Information

Item #	Product Information
P31739	Impregum™ Penta™ Soft Heavy Body/Light Body Introductory Kit Kit includes: 1 base paste HB – 300ml; 1 catalyst HB – 60ml; 1 Penta™ cartridge HB; 10 Penta mixing tips, red; 2 cartridges LB – 50ml; 10 cartridge mixing tips, white; 10 intraoral tips, white; 1 cartridge dispenser; 1 polyether tray adhesive – 17ml
P31779	Impregum™ Penta™ Soft Quick Step Heavy Body/Light Body Introductory Kit Kit includes: 1 base paste HB – 300ml; 1 catalyst HB – 60ml; 1 Penta cartridge HB; 10 Penta mixing tips, red; 10 cartridge mixing tips, white; 10 intraoral tips, white; 1 polyether tray adhesive – 17ml; 2 cartridges LB – 50ml
P31735	Impregum™ Penta™ Soft Medium Body Introductory Kit Kit includes: 1 base paste MB – 300ml; 1 catalyst MB – 60ml; 1 Penta cartridge MB; 10 Penta mixing tips, red; 1 Penta elastomer syringe; 1 polyether tray adhesive – 17ml
P31771	Impregum™ Penta™ Soft Quick Step Medium Body Introductory Kit Kit includes: 1 base paste MB – 300ml; 1 catalyst MB – 60ml; 1 Penta cartridge MB; 10 Penta mixing tips, red; 1 Penta elastomer syringe; 1 polyether tray adhesive – 17ml
31786	Impregum™ Soft Tray Introductory Kit Kit includes: 3 Tray cartridges – 50ml; 1 LB cartridge – 50ml; 10 cartridge mixing tips, purple; 10 cartridge mixing tips, white; 10 intraoral tips, white; 1 polyether tray adhesive – 17ml
31789	Impregum™ Soft Quick Step Tray Introductory Kit Kit includes: 3 Tray cartridges – 50ml; 1 LB cartridge – 50ml; 10 cartridge mixing tips, purple; 10 cartridge mixing tips, white; 10 intraoral tips, white; 1 polyether tray adhesive – 17ml

Item #	Product Information
31741	Impregum™ Penta™ Soft Heavy Body Refill Pack 2 base paste – 300ml; 2 catalyst – 60ml
31769	Impregum™ Penta™ Soft Quick Step Heavy Body Refill Pack 2 base paste – 300ml; 2 catalyst – 60ml
31731	Impregum™ Penta™ Soft Medium Body Refill Pack 2 base paste – 300ml; 2 catalyst – 60ml
31765	Impregum™ Penta™ Soft Quick Step Medium Body Refill Pack 2 base paste – 300ml; 2 catalyst – 60ml
31785	Impregum™ Soft Tray Refill Pack 4 cartridges – 50ml; 10 cartridge mixing tips, purple
31785J	Impregum™ Soft Tray Jumbo Refill Pack 30 cartridges – 50ml
31788	Impregum™ Soft Quick Step Tray Refill Pack 4 cartridges – 50ml; 10 cartridge mixing tips, purple
31788J	Impregum™ Soft Quick Step Tray Jumbo Refill Pack 30 cartridges – 50ml
31746	Impregum™ Penta™ Soft Light Body Refill Pack 2 base paste – 300ml; 2 catalyst – 60ml
31752	Impregum™ Soft Light Body Refill Pack 4 cartridges – 50ml; 20 cartridge mixing tips, white
31767	Impregum™ Soft Quick Step Light Body Refill Pack 4 cartridges – 50ml; 20 cartridge mixing tips, white

Customer Care Center: 1-800-634-2249

www.3MESPE.com/GreatImpressions

3M ESPE

Dental Products

3M Center
Building 275-2SE-03
St. Paul, MN 55144-1000
U.S.A.
1-800-634-2249

3M Canada

Post Office Box 5757
London, Ontario N6A 4T1
Canada
1-888-363-3685

Please recycle. Printed in U.S.A.
© 3M 2008. All rights reserved.
70-2013-0155-6

3M, ESPE, Impregum, Penta and Pentamix are trademarks of 3M or 3M ESPE AG. Used under license in Canada. Aquasil and Exafast are not trademarks of 3M.